r-biopharm®

## RIDASCREEN® Tetracyclin

**REF** R3505

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# RIDASCREEN® Tetracyclin

## Brief information

RIDASCREEN® Tetracyclin (Art. No. R3505) is a competitive enzyme immunoassay for the quantitative analysis of tetracyclins in milk, milk powder, cheese, yoghurt, honey, meat, fish, shrimp and whole egg.  
(see chapter 1. Intended use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:           milk, honey: dilution  
  milk powder: reconstitution, dilution  
  cheese: homogenization, extraction  
  yoghurt: incubation, centrifugation, dilution  
  honey: dissolve  
  meat, fish, shrimp, whole egg: homogenization,  
  extraction

Time requirement:           sample preparation (for 10 samples)  
  milk ..... approx. 30 min  
  milk powder ..... approx. 45 min  
  cheese ..... approx. 90 min  
  yoghurt ..... approx. 40 min  
  honey ..... approx. 30 min  
  meat ..... approx. 45 min  
  fish, shrimp ..... approx. 30 min  
  whole egg ..... approx. 50 min  
  test implementation (incubation time) ... 1 h 30 min

Limit of detection:           milk ..... approx. 0.7 µg/L  
(depending on matrix)   milk powder ..... approx. 0.8 µg/kg  
  cheese ..... approx. 1.0 µg/kg  
  yoghurt ..... approx. 0.6 µg/kg  
  honey ..... approx. 2.0 µg/kg  
  meat ..... approx. 0.7 µg/kg  
  fish ..... approx. 1.0 µg/kg  
  shrimp ..... approx. 0.5 µg/kg  
  whole egg ..... approx. 1.2 µg/kg

Recovery rate: (corresponding to the Standard substance)	milk .....	approx. 111 %
	milk powder .....	approx. 120 %
	cheese .....	approx. 70 %
	yoghurt .....	approx. 95 %
	honey.....	approx. 98 %
	meat .....	approx. 78 %
	fish .....	approx. 74 %
	shrimp.....	approx. 91 %
	whole egg .....	approx. 56 %
Specificity:	Tetracycline (standard substance).....	100 %
	Chlortetracycline .....	approx. 70 %
	Rolitetracycline .....	approx. 34 %
	Demeclocycline .....	approx. 26 %
	Oxytetracycline .....	approx. 13 %
	Minocycline.....	approx. 3 %
	Doxycycline .....	approx. 2 %

The specificity of the RIDASCREEN® Tetracyclin test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website <https://food.r-biopharm.com/media/technical-guides/>.

### **Related product and accessories for tetracyclin determination**

RIDA® Tetracyclin Spiking solution (Art. No. R3599)

## 1. Intended use

RIDASCREEN® Tetracyclin is a competitive enzyme immunoassay for the quantitative analysis of tetracyclins in milk, milk powder, cheese, yoghurt, honey, meat, fish, shrimp and whole egg.

## 2. General information

Aureomycin (chlortetracycline) was isolated 1948 by Duggan as a metabolite of the actinomyces species *streptomyces aureofaciens*. This was the first antibiotic substance of the group of tetracyclines. In Germany tetracycline, chlortetracycline and oxytetracycline are authorized for veterinary use. Therefore health risks for the consumer can occur.

In all animal species used for food production, tetracycline residues (sum of parent substance and 4-epimer) are limited in the EU law, commission regulation (EU) No 37/2010 with the following limits: 100 µg/kg in muscle and in milk.

## 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with tetracycline-protein-conjugate. Tetracycline standards or sample solutions and anti-tetracycline antibodies are added. Free tetracycline and immobilized tetracycline compete for the tetracycline antibody binding sites (competitive enzyme immunoassay). Any unbound antibody is then removed in a washing step and enzyme labeled secondary antibody (conjugate), which is directed against the anti-tetracycline antibody, is added. After removing unbound enzyme labeled antibodies by a washing step, substrate/chromogen is added to the wells and incubated. Bound conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the tetracycline concentration in the sample.

## 4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Sample buffer 1	White	Ready to use		60 mL
Sample buffer 2	Brown	Ready to use		60 mL
Standard 1	White	Concentrate 10x	0 µg/L	1.3 mL
Standard 2	White	Concentrate 10x	0.5 µg/L	1.3 mL
Standard 3	White	Concentrate 10x	1.5 µg/L	1.3 mL
Standard 4	White	Concentrate 10x	3 µg/L	1.3 mL
Standard 5	White	Concentrate 10x	6 µg/L	1.3 mL
Standard 6	White	Concentrate 10x	18 µg/L	1.3 mL
Wash buffer salt Tween		Salt for dissolving		
Conjugate	Red	Ready to use		13 mL
Antibody	Black	Ready to use		6 mL
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		13 mL
Stop solution	Yellow	Ready to use		14 mL

## 5. Reagents required but not provided

### 5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Laboratory mixer (Stomacher, ultra turrax; for meat, fish, and shrimp samples)
- Centrifuge (at least 3,000 x g) + centrifugal vials with cap (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (60 °C / 140 °F; for fluctuation range please refer to the instructions of the water bath manufacturer) (for yoghurt samples)
- ultrasonic bath (optional for milk powder and honey samples)
- Graduated pipettes
- Variable 20 - 200 µL and 200 - 1000 µL micropipettes
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET (Art. No. Z9996FF)

## 5.2 Reagents

- Distilled water (dist. water) or deionized water

For honey, meat, fish and shrimp samples:

- 20 mM PBS buffer: 0.55 g  $\text{NaH}_2\text{PO}_4 \times \text{H}_2\text{O}$  + 2.85 g  $\text{Na}_2\text{HPO}_4 \times 2 \text{H}_2\text{O}$  + 9 g NaCl, adjust pH to 7.4 with 1 M NaOH, fill up to 1000 mL with distilled water

For cheese samples:

- methanol 10 % (v/v)

For meat samples:

- n-hexane

For whole egg samples:

- 50 mM succinic acid buffer: dissolve 5.9 g of succinic acid in 500 mL of distilled water, adjust pH to 4.0 with 1 N NaOH, fill up to 1000 mL with distilled water

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

The standards contain tetracycline. Particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

Do not reuse wells of the microtiter strips (coated microtiter plate, see chapter 10.2.). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 0.6 absorbance units ( $A_{450\text{ nm}} < 0.6$ ) for standard 1

## 9. Sample preparation

The samples should be stored in a cool place, protected against light.

### 9.1 Milk

#### Milk with fat content > 1.5 %:

- Centrifuge fatty milk samples: 10 min / 3,000 g / 10 °C (50 °F)  
(if a refrigerated centrifuge is not available, chill sample to 10 °C (50 °F) prior to centrifugation)
- Remove the upper cream layer completely (e.g. with a pasteur pipette) and dilute the resulting skimmed milk 1:10 (1+9) with sample buffer 2 in a new vial (e. g. 50 µL milk + 450 µL sample buffer 2)
- Use 50 µL per well in the assay

#### Milk with fat content ≤ 1.5 %:

- Dilute skimmed milk samples 1:10 (1+9) with sample buffer 2  
(e. g. 50 µL milk + 450 µL sample buffer 2)
- Use 50 µL per well in the assay

## 9.2 Milk powder

- Weigh 10 g of milk powder and add pre-warmed distilled water (60 °C / 140 °F) to a final volume of 100 mL
- Shake / rotate for 10 min or longer until milk powder is dissolved completely optional for milk powders which do not dissolve easily: incubate for 3 min in an ultrasonic bath
- Centrifuge: 10 min / 3,500 g / 10 °C (50 °F)
- Remove the upper cream layer completely (e.g. with a pasteur pipette) and dilute the resulting skimmed milk 1:10 (1+9) with sample buffer 2 in a new vial (e. g. 50 µL milk + 450 µL sample buffer 2)
- Use 50 µL per well in the assay

## 9.3 Cheese

- Add 20 mL of methanol 10 % (v/v) to 5 g of cheese
- Homogenize the sample (use stomacher, mixer or ultra turrax)
- Transfer the homogenized sample into a 50 mL centrifugal vial
- Shake 10 min upside-down
- Centrifuge: 15 min / 3,000 g / 4 °C (39.2 °F)
- Transfer 1 mL of the aqueous (middle) phase into a 1.5 mL vial
- Centrifuge 5 min / 20,000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- Dilute supernatant 1:5 (1+4) in sample buffer 2 (e.g. 100 µL supernatant + 400 µL sample buffer 2)
- Use 50 µL per well in the assay

## 9.4 Yoghurt

- Transfer 5 g of the sample into a centrifugal vial
- Incubate the sample for 15 min at 50 °C (122 °F), e.g. in a waterbath
- Vortex until the sample is homogenized completely
- Centrifuge: 10 min / 4000 g / 10 °C (50 °F)
- Dilute the supernatant 1:10 (1+9) with sample buffer 2 (e. g. 50 µL supernatant + 450 µL sample buffer 2)
- Use 50 µL per well in the assay

## 9.5 Honey

- Weigh 1 g of honey in a centrifugal vial
- Dilute 1:50 (1+49) with 20 mM PBS buffer, pH 7.4 (see chapter 5.2.) (e.g. 1 g honey + 50 ml 20 mM PBS buffer, pH 7.4)
- Optional for honey samples which do not dissolve easily: incubate 5 min in an ultrasonic bath
- Mix intensively for 2 min on a vortex

- Before use in the test shake the vials upside down briefly
- Use 50 µL per well in the assay

### 9.6 Meat (beef, pork, poultry)

- Homogenize the sample (use stomacher, mixer or ultra turrax)
- Transfer 1 g of the homogenized sample into a centrifugal vial and add 9 ml 20 mM PBS buffer pH 7.4 (see chapter 5.2.)
- Vortex sample and buffer
- Shake 10 min for extraction
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- Transfer 1 mL of supernatant into a new vial
- Add 2 mL of hexane
- Vortex for 10 s
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- Use 50 µL of the lower aqueous phase per well in the assay

### 9.7 Fish and shrimp

- Homogenize the sample (use stomacher, mixer or ultra turrax)
- Transfer 1 g of the homogenized sample into a centrifugal screw cap vial and add 9 mL 20 mM PBS buffer pH 7.4 (see 5.2.)
- Vortex sample and buffer
- Shake/rotate the sample for 10 min for extraction
- Centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- Use 50 µL of the upper aqueous phase per well in the assay

### 9.8 Whole egg

- Homogenize egg yolk and egg white of one egg thoroughly
- Transfer 4 g of the mixed egg into a 50 mL polypropylene vial and add 20 ml of 50 mM succinic acid buffer
- Mix 15 min at room temperature (20 - 25 °C / 68 - 77 °F) by shaking / rotating
- Centrifuge: 15 min / 4,000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- Dilute supernatant 1:10 (1+9) with 20 mM PBS buffer, pH 7.4 (e.g. 100 µL of supernatant + 900 µL 20 mM PBS buffer, H 7.4)
- Use 50 µL per well in the assay

## 10. Test procedure

### 10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **standards** are provided as concentrates. In order to produce tetracycline standards ready to use, dilute 50 µL standard concentrate with 450 µL **sample buffer** in each case (1:10 dilution) and mix thoroughly. Use glass vials.

For honey, meat, fish, shrimp and whole egg samples use **sample buffer 1** for dilution of concentrates.

For milk, milk powder, cheese and yoghurt, use **sample buffer 2** for dilution of concentrates.

Standard 1:	50 µL standard concentrate (0 µg/L) + 450 µL buffer	0 µg/L
Standard 2:	50 µL standard concentrate (0.5 µg/L) + 450 µL buffer	0.05 µg/L
Standard 3:	50 µL standard concentrate (1.5 µg/L) + 450 µL buffer	0.15 µg/L
Standard 4:	50 µL standard concentrate (3 µg/L) + 450 µL buffer	0.3 µg/L
Standard 5:	50 µL standard concentrate (6 µg/L) + 450 µL buffer	0.6 µg/L
Standard 6:	50 µL standard concentrate (18 µg/L) + 450 µL buffer	1.8 µg/L

**The standards must be prepared freshly for each test series.**

A PBS-Tween buffer is needed as **wash buffer**, please use the washing buffer salt (pouch) contained in the kit (see chapter 4.). Dissolve the total content of the pouch in one liter of distilled water. The ready to use wash buffer expires after approx. 4 weeks at 2 - 8 °C (36 - 46 °F).

Alternatively: Dissolve the contents of the pouch in 100 mL of distilled water to obtain a 10fold concentrated wash buffer. This solution expires after approx. 8 weeks, stored at room temperature (20 - 25 °C / 68 - 77 °F). Use one part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use wash buffer.

Components should be stored immediately at 2 - 8 °C (35 - 46 °F) when no longer required.

## 10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

It is recommended to pipette the antibody, the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore covering the microtiter plates is recommended.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Pipette 50  $\mu\text{L}$  of each diluted standard or sample (prepared according to chapter 9.) in duplicate to the wells,
3. Add 50  $\mu\text{L}$  of anti-tetracycline antibody to each well, mix gently by shaking the plate manually and incubate for 1 h at room temperature (20 - 25  $^{\circ}\text{C}$  / 68 - 77  $^{\circ}\text{F}$ ).
4. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250  $\mu\text{L}$  diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
5. Add 100  $\mu\text{L}$  of the conjugate to each well, mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25  $^{\circ}\text{C}$  / 68 - 77  $^{\circ}\text{F}$ ).
6. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250  $\mu\text{L}$  diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
7. Add 100  $\mu\text{L}$  of substrate/chromogen to each well, mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25  $^{\circ}\text{C}$  / 68 - 77  $^{\circ}\text{F}$ ) in the dark.
8. Pipette 100  $\mu\text{L}$  of stop solution into each well. Mix gently by shaking the plate manually and measure the extinction at 450 nm. Read within 30 min after addition of stop solution.

## 11. Evaluation

Special software, **RIDASOFT® Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays.

For the evaluation it should be clarified, that all quality criteria are fulfilled for the current test run. The course of the standard curve is shown in the Quality Assurance Certificate (certificate of analysis) enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates semilogarithmic against the tetracycline concentration [ $\mu\text{g/L}$ ].

## 12. Result interpretation

In order to obtain the tetracycline concentration in  $\mu\text{g/L}$  ( $\mu\text{g/kg}$ ) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factors are as follows:

milk, yoghurt, meat, shrimp, fish.....	10
milk powder .....	10
cheese .....	25
honey .....	50
whole egg.....	60

## 13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

## 14. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control. Antibiotic-free and antibiotic containing (spiked) samples should be used.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- To do spike experiments to ensure an accurate and correct test procedure.
- To contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de) if automates (e.g. ThunderBolt® / Bolt™) are used.

**For further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Version overview

Version number	Chapter and title
2011-05-23	Release version
2015-10-30	General revision
2021-12-16	General revision Changes made: <ul style="list-style-type: none"><li>– Revalidation of the test: the matrices sausage, butter, curd, yoghurt with fruits, kefir, cream and sour cream were not considered</li><li>– 4. Reagents provided: amount of conjugate and substrate/chromogen was increased from 10 mL to 13 mL</li><li>– Additions in chapters 6 and 7</li><li>– New: Chapters 13. Limitations of the method and 14. Recommendation</li></ul>
2023-01-09	Current version Changes made: <ul style="list-style-type: none"><li>– Mistake correction in chapter 10.2: conjugate is ready to use</li></ul>

## Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.