r-biopharm®

RIDASCREEN® Tetracycline Plus

REF R3506

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RIDASCREEN® Tetracycline Plus

Brief information

RIDASCREEN® Tetracycline Plus (Art. No. R3506) is a competitive enzyme immunoassay for the quantitative analysis of tetracycline and its analogues in milk, meat, fish, shrimps, honey and eggs (see chapter 1. Intended use).

Apart from distilled water or deionized water and methanol, all reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) ... approx. 30 min
test implementation (incubation time) 45 min

Limit of detection: (depending on matrix)

Milk	approx. 1.4 µg/L (ppb)
Meat	approx. 1.7 µg/kg (ppb)
Shrimps, Fish	approx. 1.7 µg/kg (ppb)
Honey	approx. 0.8 µg/kg (ppb)
Egg	approx. 0.9 µg/kg (ppb)

Detection capability (CC β):

Milk	approx. 5 µg/kg (ppb)
Meat	approx. 5 µg/kg (ppb)
Shrimps, Fish	approx. 5 µg/kg (ppb)
Honey	approx. 5 µg/kg (ppb)
Egg	approx. 5 µg/kg (ppb)

Recovery rate: (corresponding to the Standard substance)

Milk	approx. 110 %
Meat	approx. 95 %
Shrimps, Fish	approx. 95 %
Honey	approx. 95 %
Egg	approx. 90 %

Specificity: Tetracycline (Standard substance) 100 %

4-epi-Tetracyclin.....	approx. 70 %
Oxytetracyclin	approx. 59 %
4-epi-Oxytetracyclin.....	approx. 42 %
Chlortetracyclin.....	approx. 57 %
4-epi-Chlortetracyclin.....	approx. 33 %
Doxycyclin	approx. 24 %
Demeclocyclin.....	approx. 50 %
Minocyclin	<1 %

Further information is contained in the validation report.

The specificity of the RIDASCREEN® Tetracycline Plus test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection (LoD) and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website <https://food.r-biopharm.com/media/technical-guides/>.

Related product and accessories for tetracycline determination

R3598 RIDA® Tetracycline Spiking solution

1. Intended use

RIDASCREEN® Tetracycline Plus (Art. No. R3506) is a competitive enzyme immunoassay for the quantitative analysis of tetracyclines in milk, meat, fish, shrimps, honey and egg.

For detailed results and further information on validation data with other food matrices, please refer to the validation report.

2. General information

In all animal species used for food production, tetracycline residues (sum of parent substance and 4-epimer) are limited in the EU law, Commission Regulation (EU) No. 37/2010 with the following limits:

Muscle	100 µg/kg
Milk	100 µg/kg
Egg	200 µg/kg

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against tetracycline.

Tetracycline standards or sample solutions and tetracycline enzyme conjugate are added. Free tetracycline and tetracycline enzyme conjugate compete for the tetracycline antibody binding sites (competitive enzyme immunoassay). Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the tetracycline concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Cap Color	Format	Volume
Microtiter plate	-	ready to use	1 plate 96 wells
Sample dilution buffer	transparent	4-fold concentrate	1 bottle 22 mL
Standard	black	lyophilized	2 vials, 3 µg/L 2 mL (after reconstitution)
Wash buffer salt Tween	-	salt for dissolving	1 pouch 1 L (after dissolving)
Conjugate	blue	lyophilized	1 vial 6 mL (after reconstitution)
Stop solution	yellow	ready to use	1 bottle 14 mL
Substrate/Chromogen Red Chromogen Pro	brown	ready to use	1 bottle 13 mL
Extraction buffer	transparent	ready to use	1 bottle 90 mL

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Fume hood
- Vortex mixer
- Centrifuge (at least 4,000 x g) + centrifugal vials with cap (e.g. 15 mL centrifuge tubes from Greiner)
- Shaker
- Graduated pipettes
- Variable 20 - 200 µL and 200 - 1000 µL micropipettes
- If available: 8-channel micropipette 100 - 300 µL
- If available: microtiter plate shaker
- Multipipette with 2.5 ml combi tips
- Aluminum foil or parafilm
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)

5.2 Reagents

- Distilled water (dist. water) or deionized water
- Methanol

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate, see chapter 10.2 Test procedure). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (36 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 1

9. Sample preparation

9.1 Meat, fish, shrimps and egg

- Weigh 1 gram of homogenized meat, fish, shrimps or egg sample.
- Add 2 mL of methanol and vortex for 1 minute.
- Add 1 mL of extraction buffer, shortly vortex and mix head over head for 5 minutes.
- Centrifuge for 5 minutes at 4000 g (20 - 25 °C / 68 - 77 °F).
- Dilute the upper methanol layer 1:10 (1 + 9) in the sample dilution buffer (e. g. 30 µL upper layer + 270 µL sample dilution buffer)
- Mix and use 50 µL per well in the assay.

9.2 Honey

- Weigh 1 gram of the sample and add 1 mL of water, mix well until dissolved.
- Add 2 mL of methanol and vortex for 1 minute.
- Add 1 mL of extraction buffer, shortly vortex and mix head over head for 5 minutes.
- Centrifuge for 5 minutes at 4000 g (20 - 25 °C / 68 - 77 °F).
- Dilute the upper methanol layer 1:10 (1 + 9) in the sample dilution buffer (e. g. 30 µL upper layer + 270 µL sample dilution buffer)
- Mix and use 50 µL per well in the assay.

9.3 Milk*

- Pipette 1 mL of milk.
- Add 2 mL of methanol and vortex for 1 minute.
- Add 1 mL of extraction buffer, shortly vortex and mix head over head for 5 minutes.
- Centrifuge for 5 minutes at 4000 g (20 - 25 °C / 68 - 77 °F).
- Dilute the upper methanol layer 1:10 (1 + 9) in the sample dilution buffer (e. g. 30 µL upper layer + 270 µL sample dilution buffer)

- Mix and use 50 µL per well in the assay.

* In the case of milk powder, prepare liquid milk according to the manufacturer's instructions. If the instructions are not known, dilute the milk powder 10-fold in water (e.g. 10 g up to a volume of 100 ml). When calculating the concentration per gram of the sample, account for the additional dilution factor used in preparing the milk from powder.

10. Test procedure

10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The sample dilution buffer is provided as a 4-fold concentrate. The concentrate is diluted 4 times with distilled water (e.g. 5 mL sample dilution buffer concentrate + 15 mL distilled water). The diluted sample dilution buffer should be prepared fresh before the analysis.

The standard dilution buffer is needed to dilute the standards. This buffer contains 10 % of methanol. Prepare standard dilution buffer fresh before use (e.g. 1 mL methanol + 9 mL sample dilution buffer). Mix well.

The conjugate is provided lyophilized. Add 6 mL of the sample dilution buffer to the conjugate vial and mix. The reconstituted conjugate is stable at - 20 °C (- 4 °F) for at least 3 months.

The standards are provided lyophilized. Prepare a dilution range of tetracycline standards. The reconstituted standards are stable at - 20 °C (- 4 °F) for at least 3 months. Add 2 mL of the standard dilution buffer to the tetracycline standard vial and mix. This solution contains tetracycline concentration at 3 µg/L. Pipette 0.5 mL of this solution into a clean tube and add 0.5 mL of the standard dilution buffer. Continue to make a dilution range of 1.5, 0.75, 0.38 and 0.19 µg/L. Use standard dilution buffer directly as 0 µg/L standard.

A PBS-Tween buffer is required as wash buffer, please use the wash buffer salt (pouch) contained in the kit (see chapter 4. Reagents provided). Dissolve the content of the pouch in one liter of distilled water. The ready-to-use wash buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternatively: Dissolve the content of the pouch in 100 mL of distilled water to obtain a 10-fold concentrated wash buffer. This solution expires after approx. 8 - 12 weeks, stored at room temperature (20 - 25 °C / 68 - 77 °F). Use one part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

Components should be stored immediately at 2 - 8 °C (36 - 46 °F) when no longer required.

10.2 Test procedure

Standards and samples must be tested on one microtiter plate at one time. Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore, cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µL of each standard (0 – 0.19 – 0.38 – 0.75 – 1.5 – 3 µg/L) or prepared sample in duplicate to the wells; use a new pipette tip for each standard or sample.
3. Add 50 µL of conjugate to each well (multi-dispenser or 8-channel pipette) and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Using an 8-channel pipette, fill the wells with ready to use wash buffer (250 µL per well) (see chapter 10.1 Test preparation). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
5. Add 100 µL of the substrate/chromogen to each well (multi-dispenser or 8-channel pipette) and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
6. Add 100 µL of stop solution to each well (multi-dispenser or 8-channel pipette) and measure the absorbance at 450 nm within 5 min after addition of stop solution.

11. Evaluation

A special software, RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF), is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter function.

For the evaluation it should be clarified that all quality criteria are fulfilled for the current test run. The course of the standard curve is shown in the certificate of analysis (CoA) enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates semilogarithmic against the tetracycline concentration [$\mu\text{g}/\text{kg}$].

12. Result interpretation

In order to obtain the tetracycline concentration in $\mu\text{g}/\text{L}$ ($\mu\text{g}/\text{kg}$) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factors are as follows:

Milk	20
Meat	20
Fish	20
Shrimps	20
Honey	20
Egg	20

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

14. Recommendation

In order to ensure a high analytical performance, we recommend analyzing each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

To ensure a high analytical performance, we recommend:

- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- To do spike experiments to ensure an accurate and correct test procedure.
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt® / Bolt™) are used.

15. Further application notes

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2024-12-10	Release version

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

1. In conformance with the German Civil Code (“BGB”) R-Biopharm AG provides a limited warranty (“Gewährleistung”) against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm’s own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG’s liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG’s sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser’s expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements (“Gewährleistung”) implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

R-Biopharm AG
Postanschrift / Postal Address:
An der neuen Bergstraße 17
64297 Darmstadt, Germany
Sitz / Corporate Seat: Pfungstadt
Tel.: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-40
E-mail: info@r-biopharm.de
www.r-biopharm.com

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are registered trademarks of R-Biopharm AG.