r-biopharm®

RIDASCREEN® Nitrofurantoin (SEM)

REF R3724

Enzymimmunoassay zur quantitativen Bestimmung
von Semicarbazid (SEM)

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RIDASCREEN® Nitrofuran (SEM)

Brief information

RIDASCREEN® Nitrofuran (SEM) (Art. No. R3724) is a competitive enzyme immunoassay for the quantitative analysis of semicarbazide (SEM) in shrimp, meat (bovine, porcine and chicken) and fish (see chapter 1. Intended use).

All reagents required for performing the enzyme immunoassay, including the standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 duplicate determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: Homogenization, derivatization, extraction, centrifugation, evaporation, defatting
Depending on the workflow in the laboratory, a 3-hour or an over-night derivatization may be performed. For highest precision, the over-night method is recommended.

Time requirement: Method I (over-night)
Sample preparation approx. 3 h
Derivatization approx. 16 h at 37 °C

Method II (3 hour)
Sample preparation approx. 3 h
Derivatization approx. 3 h at 50 °C

Test implementation (incubation time) 45 min

Limit of detection: Shrimp approx. 230 ng/kg (ppt)
(corresponding to the Meat approx. 200 ng/kg
standard substance) Fish approx. 200 ng/kg

Detection capability (cc β): Shrimp 250 ng/kg
(corresponding to the Meat 350 ng/kg
standard substance) Fish 350 ng/kg

Recovery rate:	Shrimp	approx. 100 %
(corresponding to the standard substance)	Meat	approx. 116 %
	Fish.....	approx. 85 %

Specificity:	Nitrophenyl-(NP) SEM (Standard substance)	
	100 %
	Nitrofurazon	1.37 %
	Semicarbazide (SEM).....	0.02 %
	Furaltadon, Furazolidon, Nitrofurantoin..	< 0.01 %
	Chloramphenicol, AMOZ, AOZ, AHD.....	< 0.01 %
	2-NP-AMOZ, 2-NP-AOZ, 2-NP-AHD.....	< 0.01 %

The specificity of the RIDASCREEN® Nitrofurantoin (SEM) test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from the values determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice. Such lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products and accessories for nitrofurantoin determination

- RIDA® Nitrofurantoin (SEM) Spiking Solution (Art. No. R3797)
- RIDASCREEN® Nitrofurantoin (AOZ) (Art. No. R3703)
- RIDASCREEN® Nitrofurantoin (AHD) (Art. No. R3713)
- RIDASCREEN® Nitrofurantoin (AMOZ) (Art. No. R3722)
- RIDA® Nitrofurantoin (AOZ) Spiking Solution (Art. No. R3798)
- RIDA® Nitrofurantoin (AHD) Spiking Solution (Art. No. R3796)
- RIDA® Nitrofurantoin (AMOZ) Spiking Solution (Art. No. R3799)

1. Intended use

RIDASCREEN® Nitrofurantoin (SEM) (R3724) is a competitive enzyme immunoassay for the quantitative analysis of semicarbazide (SEM) in shrimp, fish, meat (bovine, porcine and chicken).

2. General

Nitrofurans are synthetic broad-spectrum antibiotics, which are frequently used in animal production due to their excellent antibacterial and pharmacokinetic properties. They have also been used as growth promoters during the production of shrimps, poultry and pigs. Long term animal experiments have shown that the parent compounds and their metabolites have carcinogenic and mutagenic characteristics. This led to the prohibition of nitrofurans for the treatment of animals used for food production. In 1993, the EU banned the nitrofurans furaltadone, nitrofurantoin and nitrofurazone for use in animals used as sources of food, and in 1995 the use of furazolidone was also prohibited. The analysis of nitrofurans are based on the detection of the tissue bound metabolites of nitrofurans. The parent compounds are difficult to detect accurately since they are metabolized very rapidly after treatment. The tissue bound nitrofurantoin metabolites however are present for a long time after administration and they are used to detect nitrofurantoin abuse.

Parent compound	Metabolite	Abbreviation
Nitrofurantoin	1-Aminohydantoin	AHD
Furaltadone	3-Amino-5-morpholinomethyl-2-oxazolidinone	AMTZ
Furazolidone	3-Amino-2-oxazolidinone	AOZ
Nitrofurazone	Semicarbazide	SEM

Prior to analysis, the metabolites have to be derivatized by incubation with 2-nitrobenzaldehyde (2-NBA) into NP-AHD, NP-AMTZ, NP-AOZ and NP-SEM.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies and anti-SEM antibodies. Standards or sample and SEM conjugate are added into the wells. Free and enzyme conjugated SEM compete for the antibody binding sites (competitive enzyme immunoassay). Any unbound conjugate is then removed in a washing step. The

conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the SEM concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 duplicate determinations (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume	
Microtiter plate	-	Ready-to-use	96 wells	
2-Nitrobenzaldehyde		Solid	100 mg	
Standard 1	White	Ready-to-use	0 ng/L	1.3 mL
Standard 2	White	Ready-to-use	100 ng/L	1.3 mL
Standard 3	White	Ready-to-use	300 ng/L	1.3 mL
Standard 4	White	Ready-to-use	900 ng/L	1.3 mL
Standard 5	White	Ready-to-use	2700 ng/L	1.3 mL
Standard 6	White	Ready-to-use	8100 ng/L	1.3 mL
Wash buffer salt Tween	-	Salt for dissolving		
Conjugate	Red	Ready-to-use	7.5 mL	
Substrate/ Chromogen Red Chromogen Pro	Brown	Ready-to-use	13 mL	
Stop solution	Yellow	Ready-to-use	14 mL	

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Scale
- Centrifuge + centrifugal vials with cap (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker (Vortex)
- Mixer (Stomacher, Ultraturrax)
- Rotary evaporator
- Magnetic stirrer
- Pasteur pipettes
- Graduated pipettes
- Variable 20 - 200 μ L and 200 - 1000 μ L micropipettes
- If necessary: Multistepper or 8-channel pipette for 50 μ L and 100 μ L
- Microtiter plate spectrophotometer (450 nm)

- Optional: RIDASOFT® Win.NET (Art. No. Z9996FF)

5.2 Reagents

- Distilled or deionized water
- 1 M HCl
- 1 M NaOH
- 0.1 M K₂HPO₄
- Ethyl acetate
- n-Hexane
- Dimethyl sulfoxide (DMSO)

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate, see chapter 10.2.). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the red-stained substrate/chromogen prior to test implementation
- A value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 1

9. Preparation of Samples

The samples should be stored in a cool place, protected against light.

For the derivatization step, a solution of 10 mM 2-nitrobenzaldehyde (see chapter 4.) in dimethyl sulfoxide (DMSO) has to be prepared directly before use (i.e. dissolve 7.6 mg 2-nitrobenzaldehyde in 5 mL DMSO).

9.1 Shrimps

- Homogenize the sample, use stomacher or mixer
- Prepare 10 mM 2-nitrobenzaldehyde in dimethyl sulfoxide (DMSO) (see chapter 5.2.)
- Mix 1 g of the homogenized sample with 4 mL distilled water, 0.5 mL 1 M HCl and 100 μ L 10 mM 2-nitrobenzaldehyde in DMSO (see chapter 5.2) by vortexing

Derivatization

- Incubate for 3 h at 50 °C (122 °F) or over-night (approx. 16 h) at 37 °C (98.6 °F)
- Add 5 mL 0.1 M K_2HPO_4 , 0.4 mL 1 M NaOH and 10 mL ethyl acetate, vortex for 30 s
- Centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- Transfer 5 mL of the ethyl acetate layer (upper layer) into a new vial and evaporate to dryness at 60 °C (140 °F)

Note: If the sample has a high fat content, oily residues that cannot be evaporated may occur. In this case, the evaporation step is completed and sample preparation can be continued.

- Dissolve the residue in 1 mL n-hexane and vortex with 1 mL wash buffer for 30 s
- Centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- Remove 750 μ L from the lower, aqueous phase and use 50 μ L per well in the assay

9.2 Meat and poultry

- Homogenize the sample, use stomacher or mixer
- Prepare 10 mM 2-nitrobenzaldehyde in dimethyl sulfoxide (DMSO) (see chapter 5.2.)
- Mix 1 g of the homogenized sample with 4 mL distilled water, 0.5 mL 1 M HCl and 300 μ L 10 mM 2-nitrobenzaldehyde in DMSO (see chapter 5.2) by vortexing

Derivatization

- Incubate for 3 h at 50 °C (122 °F) or over-night (approx. 16 h) at 37 °C (98.6 °F)
- Add 5 mL 0.1 M K_2HPO_4 , 0.4 mL 1 M NaOH and 10 ml ethyl acetate, vortex for 30 s
- Centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- Transfer 5 mL of the ethyl acetate layer (upper layer) into a new vial and evaporate to dryness at 60 °C (140 °F)

Note:

If the sample has a high fat content, oily residues that cannot be evaporated may occur. In this case, the evaporation step is completed and sample preparation can be continued.

- Dissolve the residue in 1 mL n-hexane and vortex with 1 mL wash buffer for 30 s
- Centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- Remove 750 μ L from the lower, aqueous phase and use 50 μ L per well in the assay

9.3 Fish

- Homogenize the sample, use stomacher or mixer
- Prepare 10 mM 2-nitrobenzaldehyde in dimethyl sulfoxide (DMSO) (see chapter 5.2.)
- Mix 1 g of the homogenized sample with 4 mL distilled water, 0.5 mL 1 M HCl and 200 μ L 10 mM 2-nitrobenzaldehyde in DMSO (see chapter 5.2) by vortexing

Derivatization

- Incubate for 3 h at 50 °C (122 °F) or over-night (approx. 16 h) at 37 °C (98.6 °F)
- Add 5 mL 0.1 M K₂HPO₄, 0.4 mL 1 M NaOH and 10 mL ethyl acetate, vortex for 30 s
- Centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- Transfer 5 mL of the ethyl acetate layer (upper layer) into a new vial and evaporate to dryness at 60 °C (140 °F)

Note: If the sample has a high fat content, oily residues that cannot be evaporated may occur. In this case, the evaporation step is completed and sample preparation can be continued.

- Dissolve the residue in 1 mL n-hexane and vortex with 1 mL wash buffer for 30 s
- Centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- Remove 750 µL from the lower, aqueous phase and use 50 µL per well in the assay

10. Test implementation

10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

As wash buffer, a PBS Tween buffer is needed; please use the enclosed buffer salt (see chapter 4. Reagents provided). To prepare the buffer, dissolve the entire contents of the pouch in 1 L distilled water. The dissolved wash buffer can be stored for approximately 4 weeks at 2 - 8 °C.

Alternative: Dissolve the contents of the pouch in 100 ml distilled water (10-fold concentrate). The solution can be stored for approximately 8 weeks at room temperature (20 - 25 °C). To prepare the ready-to-use solution, mix 1 part of the 10-fold concentrate with 9 parts of distilled water.

Return all reagents to 2 - 8 °C immediately after use.

10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µL of each standard or prepared sample to separate duplicate wells.
3. Add 50 µL of the conjugate to the bottom of each well, mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C).
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL washing buffer (see chapter 10.1.) and pour out the liquid again. Repeat two more times.
5. Add 100 µL of substrate/chromogen solution to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C) in the dark.
6. Add 100 µL of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 15 minutes after addition of stop solution.

11. Results

Special software, **RIDASOFT® Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays.

For the evaluation it should be clarified, that all quality criteria are fulfilled for the current test run. The course of the standard curve is shown in the Quality Assurance Certificate (certificate of analysis) enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a semi-logarithmic plot against the SEM concentration [ng/L].

12. Result interpretation

In order to obtain the SEM concentration in ng/kg actually contained in a sample, the concentration read from the standard curve must be further multiplied by the corresponding dilution factor. When working according to the instructions, the dilution factors are as follows:

Shrimp:	2
Meat:	2
Fish:	2

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

14. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt® / Bolt™) are used.

15. Further application notes

Application notes are available on request.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2021-03-23	Release version
2022-10-27	Current version General revision Changes made: <ul style="list-style-type: none">– “Brief information”, chapter 6, 7 and 10 updated– Chapters 13, 14, 15 added

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.