## **RIDASCREEN® Nitrofuran (DNSH)**

**REF R3725**

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# RIDASCREEN® Nitrofurazone (DNSH)

## Brief information

RIDASCREEN® Nitrofurazone (DNSH) (Art. No. R3725) is a competitive enzyme immunoassay for the quantitative analysis of nitrofurazone metabolite 3,5-dinitrosalicylic acid hydrazide (DNSH) in shrimps (see chapter 1. Intended use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: Homogenization and extraction

Time requirement: Sample preparation (for 10 samples) ..... approx. 3 h  
Derivatisation ..... approx. 3 h at 50 °C  
Test implementation (incubation time)..... 45 min

Detection capability (CC $\beta$ ): Shrimps ..... 250 ng/kg  
(corresponding to the standard substance)

Recovery rate: Shrimps ..... approx. 82 %  
(corresponding to the standard substance)

Specificity:

NP-DNSH (standard substance).....	100 %
Nitrofurazone.....	122 %
DNSH .....	10 %
DNSA .....	0.5 %
DNSHA.....	55 %
SEM .....	< 0.1 %
AMOZ .....	< 0.1 %
AOZ .....	< 0.1 %
AHD .....	< 0.1 %
Nitrofurazone .....	< 0.1 %

Furaltadon .....	< 0.1 %
Furazolidon .....	< 0.1 %
Chloramphenicol .....	< 0.1 %
Nitrofurantoin .....	< 0.1 %

The specificity of the RIDASCREEN® Nitrofurantoin (DNSH) test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection (LoD) and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website <https://food.r-biopharm.com/media/technical-guides/>.

### **Related product and accessories for nitrofurantoin determination**

RIDA® Nitrofurantoin (DNSH) Spiking Solution	(Art. No. R3795)
RIDASCREEN® Nitrofurantoin (AHD)	(Art. No. R3713)
RIDASCREEN® Nitrofurantoin (AMOZ)	(Art. No. R3722)
RIDASCREEN® Nitrofurantoin (AOZ)	(Art. No. R3703)
RIDASCREEN® Nitrofurantoin (SEM)	(Art. No. R3724)
RIDA® Nitrofurantoin (AHD) Spiking Solution	(Art. No. R3796)
RIDA® Nitrofurantoin (AMOZ) Spiking Solution	(Art. No. R3799)
RIDA® Nitrofurantoin (AOZ) Spiking Solution	(Art. No. R3798)
RIDA® Nitrofurantoin (SEM) Spiking Solution	(Art. No. R3797)

## 1. Intended use

RIDASCREEN® Nitrofurantoin (DNSH) (Art. No. R3725) is a competitive enzyme immunoassay for the quantitative analysis of nitrofurantoin metabolite 3,5-dinitrosalicylic acid hydrazide (DNSH) in shrimps.

## 2. General information

Nitrofurantoin is a nitrofurantoin antibiotic banned as a feed additive in the European Union and other countries. Nitrofurantoin is metabolized to 3,5-dinitrosalicylic acid hydrazide (DNSH) in living organisms. DNSH is a marker for the detection of illegal use of nitrofurantoin in animal husbandry. In accordance with Commission Regulation (EU) 2019/1871 a reference point for action (RPA) of 0.5 µg/kg for DNSH and other nitrofurantoin metabolites applies from 28 November 2022.

Parent compound	Metabolite	Abbreviation
Nitrofurantoin	3,5-dinitrosalicylic acid hydrazide	DNSH
Nitrofurantoin	1-Aminohydantoin	AHD
Furaltadone	3-Amino-5-morpholinomethyl-2-oxazolidinone	AMOZ
Furazolidone	3-Amino-2-oxazolidinone	AOZ
Nitrofurazone	Semicarbazide	SEM

Prior to analysis, the metabolites have to be derivatized by incubation with 2-Nitrobenzaldehyde (2-NBA) into NP-DNSH, NP-AHD, NP-AMOZ, NP-AOZ, and NP-SEM.

## 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies and specific anti-DNSH antibodies. Standards or samples and enzyme-conjugated DNSH are added into the wells. Free and enzyme-conjugated DNSH compete for the anti-DNSH antibody binding sites. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells and bound enzyme-conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the DNSH concentration in the sample.

## 4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
2-Nitrobenzaldehyd	-	solid		100 mg
Standard 1	White	Ready-to-use	0 ng/L	1.3 mL
Standard 2	White	Ready-to-use	100 ng/L	1.3 mL
Standard 3	White	Ready-to-use	250 ng/L	1.3 mL
Standard 4	White	Ready-to-use	625 ng/L	1.3 mL
Standard 5	White	Ready-to-use	1563 ng/L	1.3 mL
Standard 6	White	Ready-to-use	3906 ng/L	1.3 mL
Wash buffer salt Tween		Salt for dissolving		
Conjugate	Red	Ready-to-use		7.5 mL
Substrate/Chromogen Red Chromogen Pro	Brown	Ready-to-use		13 mL
Stop solution	Yellow	Ready to use		14 mL

## 5. Reagents required but not provided

### 5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of  $\pm 0.01$  g) and weighing vessels
- Vortex mixer
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 4,000 x g) + centrifugal vials with cap (e.g. 50 mL centrifuge tubes from Greiner)
- Rotary evaporator
- Magnetic stirrer
- Graduated pipettes
- Variable 20 - 200  $\mu$ L and 200 - 1000  $\mu$ L micropipettes
- If available: 8-channel micropipettes for 50  $\mu$ L, 100  $\mu$ L, and 250  $\mu$ L
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET (Art. No. Z9996FF)

## 5.2 Reagents

- Distilled water (dist. water) or deionized water
- 1 M HCl
- 1 M NaOH
- 0.1 M K<sub>2</sub>HPO<sub>4</sub>
- Ethyl acetate
- n-Hexane
- Dimethyl sulfoxide (DMSO)

For the derivatization step, a solution of 10 mM 2-nitrobenzaldehyde (see chapter 4.) in dimethyl sulfoxide (DMSO) has to be prepared directly before use (i.e. dissolve 7.6 mg 2-nitrobenzaldehyde in 5 mL DMSO).

## 6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

The standards contain nitrofuran metabolite (NP-DNSH). Particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

Do not reuse wells of the microtiter strips (coated microtiter plate, see chapter 10.2.). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided, and further store at 2 - 8 °C (36 - 46 °F).

The substrate/chromogen and conjugate are light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- Bluish coloration of the substrate/chromogen prior to test implementation
- Value of less than 0.8 absorbance units ( $A_{450\text{ nm}} < 0.8$ ) for the zero standard

## 9. Sample preparation

The samples should be stored in a cool place, protected against light.

For the derivatization step, a solution of 10 mM 2-nitrobenzaldehyde (see chapter 4.) in dimethyl sulfoxide (DMSO) has to be prepared directly before use (i.e. dissolve 7.6 mg 2-nitrobenzaldehyde in 5 mL DMSO).

### 9.1 Shrimps

- Homogenize the sample, use stomacher or mixer
- Prepare 10 mM 2-nitrobenzaldehyde in dimethyl sulfoxide (DMSO) (see chapter 5.2.)
- Mix 1 g of the homogenized sample with 4 mL distilled water, 0.5 mL 1 M HCl and 100 µL 10 mM 2-nitrobenzaldehyde (in DMSO) by vortexing.

### Derivatization

- Incubate for 3 h at 50 °C (122 °F)
- Add 5 mL 0.1 M  $K_2HPO_4$ , 0.4 mL 1 M NaOH and 10 mL ethyl acetate
- Vortex for 30 s
- Centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- Transfer 5 mL of the ethyl acetate layer (upper layer) into a new vial and evaporate to dryness at 60 °C (140 °F)

Note: If the sample has a high fat content, oily residues that cannot be evaporated may occur. In this case, the evaporation step is completed and sample preparation can be continued.

- Dissolve the residue in 1 mL n-hexane and vortex with 1 mL wash buffer (see chapter 10.1) for 30 s

- Centrifuge: 10 min / 4000 g / at room temperature (20 - 25 °C / 68 - 77 °F)
- Remove 750 µL from the lower, aqueous phase and use 50 µL per well in the assay

## 10. Test procedure

### 10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

As **wash buffer** a PBS tween buffer is needed. Please use the wash buffer salt (see chapter 4.) contained in the kit. Dissolve the entire buffer salt in one liter of distilled water. The ready to use washing buffer expires after approx. 4 weeks at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the pouch in 100 mL of distilled water to obtain a 10-fold concentrated washing buffer. This solution expires after approx. 8 weeks, store at room temperature (20 - 25 °C). Use 1 part 10-fold concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

Components should be stored immediately at 2 - 8 °C (36 - 46 °F) when no longer required.

### 10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µL of each standard or prepared sample extract to separate duplicate wells.
3. Add 50 µL of the conjugate (red cap) to the bottom of each well, mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C) in the dark.

4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL washing buffer (see chapter 10.1.) and pour out the liquid again. Repeat two more times.
5. Add 100 µL of substrate/chromogen (brown cap) solution to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C) in the dark.
6. Add 100 µL of the stop solution (yellow cap) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 15 minutes after addition of stop solution.

## 11. Evaluation

Special software, **RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays.

For the evaluation it should be clarified, that all quality criteria are fulfilled for the current test run. The course of the standard curve is shown in the Quality Assurance Certificate (certificate of analysis) enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates semi-logarithmic against the DNSH concentration [ng/L].

## 12. Result interpretation

In order to obtain the DNSH concentration in ng/L (ng/kg) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the **dilution factor is 2**.

### 13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

### 14. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- To do spike experiments in matrix to ensure an accurate and correct test procedure.
- To contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de) if automates (e.g. ThunderBolt® / Bolt™) are used.

**For further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

### Version overview

Version number	Chapter and title
2023-11-28	Release version

## Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.