r-biopharm®

RIDASCREEN® DNSH

REF R3740

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RIDASCREEN® DNSH

Brief information

RIDASCREEN® DNSH (Art. No. R3740) is a competitive enzyme immunoassay for the quantitative analysis of nifursol metabolite 3,5-dinitrosalicylic acid hydrazide (DNSH) in meat and seafood (see chapter 1. Intended use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)..... approx. 4.5 h
test implementation (incubation time)..... 60 min

Limit of detection: meat..... approx. 0.13 µg/kg (ppb)
(depending on matrix) fish and shellfish..... approx. 0.15 µg/kg

Detection capability (CC β):
(depending on matrix) meat..... approx. 0.25 µg/kg
fish and shellfish..... approx. 0.25 µg/kg

Recovery rate: meat (Ø)..... approx. 92 %
(corresponding to the fish and shellfish (Ø) approx. 95 %
standard substance)

Specificity: DNSH (standard substance)..... 100 %
Nifursol..... 119 %
AHD..... < 0.1 %
AMOZ..... < 0.1 %
AOZ < 0.1 %
SEM < 0.1 %

The specificity of the RIDASCREEN® DNSH test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

Related product and accessories for nitrofurantoin determination

RIDASCREEN® Nitrofurantoin (AHD)	(Art. No. R3713)
RIDASCREEN® Nitrofurantoin (AMOZ)	(Art. No. R3722)
RIDASCREEN® Nitrofurantoin (AOZ)	(Art. No. R3703)
RIDASCREEN® Nitrofurantoin (SEM)	(Art. No. R3724)
RIDA® Nitrofurantoin (AHD) Spiking Solution	(Art. No. R3796)
RIDA® Nitrofurantoin (AMOZ) Spiking Solution	(Art. No. R3799)
RIDA® Nitrofurantoin (AOZ) Spiking Solution	(Art. No. R3798)
RIDA® Nitrofurantoin (SEM) Spiking Solution	(Art. No. R3797)

1. Intended use

RIDASCREEN® DNSH (Art. No. R3740) is a competitive enzyme immunoassay for the quantitative analysis of nifursol metabolite 3,5-dinitrosalicylic acid hydrazide (DNSH) in meat (poultry, pork, beef) and seafood (fish, shrimp, and prawn).

2. General information

Nifursol is a nitrofurantoin antibiotic banned as a feed additive in the European Union and other countries. Nifursol is metabolized to 3,5-dinitrosalicylic acid hydrazide (DNSH) in living organisms. DNSH is a marker for the detection of illegal use of nifursol in animal husbandry. The detection of DNSH by LC-MS/MS requires derivatization of this metabolite with 2-nitrobenzaldehyde to NP-DNSH, like other nitrofurantoin metabolites such as SEM, AHD, AMOZ and AOZ. The antibody used in this test detects DNSH without derivatization with 2-nitrobenzaldehyde. In accordance with Commission Regulation (EU) 2019/1871 a reference point for action (RPA) of 0.5 µg/kg for DNSH and other nitrofurantoin metabolites applies from 28 November 2022.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against DNSH. In the first step of the assay standards and samples are added. DNSH present in the standard solutions or in the samples bind to the immobilized antibody. After incubation, the content of the wells is emptied. Then DNSH enzyme conjugate is added to each well and it binds to any available antibody. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the DNSH concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Format		Volume/Amount
Microtiter plate	Ready to use		96 wells
Extraction buffer	Ready to use		2 x 50 mL
Dilution buffer	Concentrate	10x	30 mL
Standard 1 (meat)	Lyophilized	4.4 µg/kg	2 x
Standard 2 (seafood)	Lyophilized	3.6 µg/kg	2 x
Wash buffer	Concentrate	20x	30 mL
Conjugate	Concentrate	100x	200 µL
Protease enzyme	Concentrate	40x	2 x 1.3 mL
Substrate/Chromogen	Ready to use		12 mL
Stop solution	Ready to use		12 mL

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g) and weighing vessels
- Fume hood
- Vortex mixer
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 4,000 x g) + centrifugal vials with cap (e.g. 15 mL centrifuge tubes from Greiner)
- If available: 8-channel micropipette 100 - 300 µL
- If available: microtiter plate shaker
- Variable 20 - 200 µL and 200 - 1000 µL micropipettes
- Multipipette with 2.5 ml combitips
- Aluminium foil or parafilm
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET (Art. No. Z9996FF)

5.2 Reagents

- Distilled water (dist. water) or deionized water
- Acetonitrile
- 1 M hydrochloric acid (HCl)

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

The standards contain DNSH. Particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate, see chapter 10.2.). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (36 - 46 °F).

The substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the substrate/chromogen prior to test implementation
- Value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for the zero standard

9. Sample preparation

The samples should be stored in a cool place, protected against light.

- Add 1 mL of diluted protease enzyme (see chapter 10.1) to 1 g of homogenized sample and vortex for 10 s
- Incubate the sample for 1 h at 50 °C (122 °F)
- Add 100 µL of 1 M hydrochloric acid (HCl) to the sample, vortex for 10 s
- Incubate the sample for 3 h at 50 °C (122 °F)
- Add 1 mL of acetonitrile to the samples and vortex for 10 s
- Add 1 mL of extraction buffer to the sample and vortex/shake for 30 s
- Centrifuge the sample for 5 min at 4000 g (20 - 25 °C / 68 - 77 °F)
- Dilute the upper acetonitrile layer 1:8 (1 + 7) in the dilution buffer (see chapter 10.1) (e. g. 50 µL upper layer + 350 µL dilution buffer)
- Mix well and use 100 µL directly in the test

10. Test procedure

10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **dilution buffer** is provided as a 10-fold concentrate. Before use, mix the concentrated buffer and dilute 1:10 (1 + 9) with distilled water (e.g. 10 mL buffer concentrate + 90 mL dist. water). Prepare fresh before use.

The **wash buffer** is provided as a 20-fold concentrate. Before use, the buffer has to be diluted 1:20 (1 + 19) with distilled water (e.g. 1 mL buffer concentrate + 19 mL dist. water). For each strip 20 mL of diluted rinsing buffer is required. Prepare dilutions freshly before use.

The **conjugate** is provided as a concentrate (100-fold). Since the diluted conjugate has a limited stability, only the amount which actually is needed should be diluted. Spin down the conjugate in the vial by a short centrifugation step (1 min, 1000 x g). The conjugate concentrate has to be diluted 1:100 (1 + 99) with dilution buffer (e.g. 20 µL concentrate + 1980 µL dilution buffer). For two strips 1.6 mL of diluted conjugate is required. Store unused concentrated conjugate at 2 - 8 °C.

The **protease enzyme** is provided as a concentrate (40-fold). Before pipetting, the conjugate concentrate should be shaken carefully and mixed well. The protease enzyme has to be diluted 1:40 (1 + 39) with dist. water (e.g. 0.5 mL

enzyme + 19.5 mL dist. water). For each sample 1 mL of diluted protease enzyme is required.

The **standard dilution buffer** is needed to dilute the standards. It is the dilution buffer containing 12.5 % of acetonitrile. Prepare standard dilution buffer fresh before use (e. g. 1.25 mL acetonitrile + 8.75 mL dilution buffer). Mix well.

The **standards** are provided lyophilized. Prepare a dilution range of DNSH standards. For optimum calibration select the standard corresponding to the samples used. The reconstituted standards are stable at - 20 °C (- 4 °F) for 4 weeks.

Meat:

Add 2 mL of the **standard dilution buffer** to DNSH “meat” standard and mix. This solution contains DNSH concentration corresponding to 4.4 µg/kg of DNSH in a sample. Pipette 0.4 mL of this solution into a clean tube and add 0.6 mL of the standard dilution buffer. Continue to make a dilution range of 1.76, 0.70, 0.28 and 0.11 µg/kg.

Seafood:

Add 2 mL of the **standard dilution buffer** to the DNSH “seafood” standard and mix. This solution contains DNSH concentration corresponding to 3.6 µg/kg of DNSH in a sample. Pipette 0.4 mL of this solution into a clean tube and add 0.6 mL of the standard dilution buffer. Continue to make a dilution range of 1.44, 0.58, 0.23 and 0.09 µg/kg.

Components should be stored immediately at 2 - 8 °C (36 - 46 °F) when no longer required.

10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Pipette 100 μL of the standard dilution buffer in duplicate as **blank** (e. g. wells G1, G2) and 100 μL of the standard dilution buffer in duplicate to the wells as **maximum signal** (e. g. wells A1, A2).
3. Pipette 100 μL of each standard or sample (prepared according to chapters 9. and 10.1) in duplicate to the remaining wells, mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 $^{\circ}\text{C}$ / 68 - 77 $^{\circ}\text{F}$) (in the dark).
4. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. **Note:** do not wash.
5. Add 100 μL of the diluted conjugate (see chapter 10.1) to each well, except blanks (e. g. wells G1, G2), mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 $^{\circ}\text{C}$ / 68 - 77 $^{\circ}\text{F}$) (in the dark).
6. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 300 μL diluted washing buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
7. Add 100 μL of substrate/chromogen to each well, mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 $^{\circ}\text{C}$ / 68 - 77 $^{\circ}\text{F}$) in the dark.
8. Pipette 100 μL of stop solution into each well. Mix gently by shaking the plate manually and measure the extinction immediately at 450 nm.

11. Evaluation

Special software, **RIDASOFT[®] Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN[®] enzyme immunoassays.

For the evaluation it should be clarified, that all quality criteria are fulfilled for the current test run. The course of the standard curve is shown in the Quality Assurance Certificate (certificate of analysis) enclosed in the test kit.

Remark for the calculation without software:

Subtract the mean absorbance values of the blanks (G1, G2) from the individual ones of the wells containing the standards and the samples.

The absorbance values of the six standards and the samples (mean values of the duplicates) are divided by the mean absorbance value of the zero standard (maximum signal) and multiplied by 100.

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = B/B_0 (\%)$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates semilogarithmic against the DNSH concentration [$\mu\text{g}/\text{kg}$]. The DNSH concentration in meat and seafood expressed in $\mu\text{g}/\text{kg}$ can be read directly from the calibration curve.

12. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

13. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control. Antibiotic-free and antibiotic containing (spiked) samples can be used.
- To do spike experiments to ensure an accurate and correct test procedure.
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt® / Bolt™) are used.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2022-10-11	Release version
2022-12-19	Current version Editorial changes and volume information error correction in chapter 10.1 (only german part)

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.