RIDASCREEN[®]FAST Casein

REF R4612

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RIDASCREEN®FAST Casein

Brief information

RIDASCREEN®FAST Casein (Art. No. R4612) is a sandwich enzyme immunoassay for the quantitative analysis of casein in food validated for the method (see chapter 1).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization and extraction
Time requirement:	sample preparation: for 10 samples..... approx. 20 min for 10 heated samples..... approx. 30 min test implementation (incubation time)..... 30 min
Standard material:	Casein from SIGMA Aldrich (Art. No. C5890)
Limit of detection: (depending on matrix)	Allergen extraction buffer 0.12 mg/kg (ppm) Casein* 0.19 mg/kg (ppm) Casein RIDA® Extractor 2..... 0.71 mg/kg (ppm) Casein* 0.41 - 0.95 mg/kg (ppm) Casein *mean value
Limit of quantification:	Allergen extraction buffer 0.5 mg/kg (ppm) Casein RIDA® Extractor 2..... 2.5 mg/kg (ppm) Casein
Specificity:	The antibodies used specifically detect caseins of cow's milk. There is a cross-reactivity to sheep's, goat's and buffalo's milk. Any further cross-reactivity is not known. Especially, a cross-reactivity to β -lactoglobulin does not exist. Further information regarding cross-reactivity is contained in the validation report.

Cross-reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross-reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments (see chapter 13).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

Related products for milk determination

RIDA® Extractor 2 (Art. No. R4613)

RIDASCREEN®FAST Milk (Art. No. R4652)

RIDASCREEN®FAST β -Lactoglobulin (Art. No. R4912)

RIDASCREEN® β -Lactoglobulin (Art. No. R4901)

Bioavid Lateral Flow Milch / Milk (Art. No. BL623-15)

Bioavid Lateral Flow Casein (Art. No. BLH714-15)

1. Intended use

RIDASCREEN®FAST Casein (Art. No. R4612) is a sandwich enzyme immunoassay test for the quantitative analysis of contaminations by casein in food. Due to the large number of different food products, the following samples were examined as representative for different food product categories within the scope of the test development: ice cream, wine, chocolate, sausage and cookies.

The test kit should generally be able to provide usable results for other foods within the methodological limits (see chapter 13). In individual cases, the user must check the suitability of the test kit by carrying out appropriate experiments before use.

For detailed results and further information on validation data with other food matrices, please refer to the validation report. Further applications are regularly validated in our laboratories, which we make available in our application notes (see chapter 15).

2. General information

Casein refers to a family of proteins commonly found in mammalian milk. The term “casein” is derived from the Latin word “*caseus*”, which means “cheese”. Cow’s milk contains 3.2 % proteins, which consist of approx. 10 % β -lactoglobulin (leading protein of whey), approx. 80 % caseins and approx. 10 % other proteins. The most important allergen for children is β -lactoglobulin while the caseins become the dominant allergen later in adults. Caseins are also heat-stable allergens.

Cow's milk is one of the most important allergenic components of food, especially for children. That is why the labeling of casein or milk is mandatory in many countries worldwide. Although there are no legal limits for casein, food producers are strongly advised to test for very low casein concentrations to protect people with allergies and to avoid allergen-related recalls.

Furthermore, Casein is a protein that is often used as a fining agent in wine production. In the fining process, casein is added to bind the suspended particles that settle out of the wine.

Casein is a rough flaked curdling protein, which forms micelles in the milk and precipitates under acidic conditions. The group of caseins consists of the α -caseins, β -casein, κ -casein and γ -casein (proteolytic protein fragment of β -casein by the milk protease plasmin). κ -casein can be cleaved into a hydrophobic (para κ -casein) and into a water soluble polar component (macropeptid) by proteolysis, e.g. by using lab ferment. Caseins can be present as an ingredient or as a contaminant in raw and processed food products. Skim milk powder or caseins / caseinates are often added to food products (e.g. in sausages) to increase the protein content or as thickening agent. According to the **regulation (EU) No. 1169/2011**, milk and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against casein. By adding the standard or sample solution to the wells, casein present in the sample will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed

in a washing step. Following the washing step, a solution containing antibodies conjugated to peroxidase is added. This conjugate binds to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product. A stop solution is added which results in a color change from blue to yellow. The absorbance of the solution which is proportional to the casein concentration in the sample is measured photometrically at 450 nm and expressed as mg/kg casein.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Extractor 2	Blue	Concentrate	2x	30 mL
Allergen extraction buffer	Green	Concentrate	10x	100 mL
Additive 1	Blue			2 g
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 mL
Standard 2*	Transparent	Ready to use	0.5 mg/kg	1.3 mL
Standard 3*	Transparent	Ready to use	1.5 mg/kg	1.3 mL
Standard 4*	Transparent	Ready to use	4.5 mg/kg	1.3 mL
Standard 5*	Transparent	Ready to use	13.5 mg/kg	1.3 mL
Wash buffer	Brown	Concentrate	10x	100 mL
Conjugate	Red	Concentrate	11x	0.7 mL
Conjugate buffer	Black	Ready to use		7 mL
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		13 mL
Stop solution	Yellow	Ready to use		14 mL

* The concentration values of the standards already take into account the dilution factor 20, which results from the sample preparation according to chapter 9.2. Thus, the casein concentrations of the samples can directly be read from the standard curve.

If sample preparation with Extractor 2 is used (see chapter 9.1), a dilution factor of 100 is obtained. The values read from the standard curve must be multiplied by factor 5 in addition.

5. Reagents required but not provided

5.1. Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (60 °C / 140 °F and 100 °C / 212 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12 μ m)
- Graduated pipettes
- Graduated cylinders
- Variable 20 - 200 μ L and 200 - 1000 μ L micropipettes
- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100 μ L
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)

5.2. Reagents

- Distilled water (dist. water) or deionized water
- 1 M sodium hydroxide (NaOH)
- 1 M hydrochloric acid (HCl)
- If necessary: bovine serum albumin (BSA, e.g. Serva, fraction V, protease free, Art. No. 11926)

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate and pre-plate, if necessary (see chapter 10.2). Use separate pipette tips for each standard and each sample extract to avoid cross-contamination.

The Extractor 2 is harmful to health. It contains mercaptoethanol. It should be worked **under a chemical hood** and skin contact should be avoided (use gloves!).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation.
- Value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 1.2$) for standard 5.

9. Preparation of Samples

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation.

- Carry out the sample preparation in a room isolated from the ELISA procedure.

The samples should be stored in a cool place and protected against light before extraction.

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The Allergen extraction buffer is provided as a 10-fold concentrate and must be diluted prior use. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that, dilute the heated buffer concentrate 1:10 (1+9) with dist. water (e.g. 900 mL dist. water + 100 mL buffer concentrate). The **diluted Allergen extraction buffer (AEB)** is stable at 20 - 25 °C (68 - 77 °F) for approx. 4 weeks or at 2 – 8 °C (35 - 46 °F) for 12 weeks.

The AEB is required for the dilution of the extracts for use in ELISA. For the dilution of matrices containing pine, sunflower or pumpkin seeds, the diluted AEB must be additionally mixed with BSA (final concentration: 2.5 %; e.g. 10 mL diluted AEB + 0.25 g BSA (= BSA-AEB)).

For the preparation of the final **Allergen extraction buffer containing Additive 1 (A-AEB)**, which is required for the extraction according to chapter 9.1, weigh 1.35 g of Additive 1 in a glass beaker and add 15 mL 1 M NaOH. Stir until the Additive 1 is solved. Then, fill 700 mL diluted AEB (see above) in a measuring cylinder. Add the 15 mL Additive 1 solution by stirring constantly. Transfer any residues of the Additive 1 solution into the measuring cylinder by rinsing with diluted AEB. Adjust the Additive 1 containing Allergen extraction buffer (A-AEB) to pH 9 with 1 M HCl and fill up to 750 mL with diluted AEB.

750 mL A-AEB are sufficient for 45 samples. The buffer can be used for approx. 3 weeks at room temperature 20 - 25 °C (68 - 77 °F) (do **not** store in the refrigerator). Discard the buffer if crystals are present. Crystal formation may begin after 10 days. We therefore recommend preparing a fresh buffer solution every week. Alternatively, a reduced volume of the A-AEB of e.g. 375 mL can be prepared. When preparing a smaller volume of A-AEB, the composition described above must be adhered to accordingly. Use clean bottles when preparing the buffer. Particles of dust can initiate crystallization.

Make sure to heat the A-AEB in time in the 60 °C (140 °F) water bath (the 100 °C / 212 °F water bath is used for the extraction of the samples).

The **Extractor 2** is provided as 2fold concentrate and has to be diluted 1:2 (1+1) with dist. water (e.g. 30 mL Extractor 2 + 30 mL dist. water). The complete diluted Extractor 2 is sufficient for 15 samples and can be used for approx.

3 month at room temperature (20 - 25 °C / 68 - 77 °F). Additional Extractor 2 can be ordered at R-Biopharm with the article number R4613.

In the following section, the following abbreviations are used:

- AEB: final diluted Allergen extraction buffer
- A-AEB: AEB with addition of Additive 1

9.1. Extraction with Extractor 2 and A-AEB as universal sample extraction

Heat the A-AEB to 60 °C (140 °F) before sample extraction.

Homogenize a representative, adequately big amount of a solid sample (e.g. 50 g; grind it thoroughly to powder and mix well) or mix well a sample in case of liquid foods.

- Transfer 1 g (or in case of liquid samples 1 mL) of homogenized sample to a new vial, add 4 mL prepared Extractor 2 (see chapter 9), close the vial and mix well.
- Cook it for 10 min at 100 °C (212 °F) in a water bath (if a lump forms, it must be dissolved in the still warm liquid by shaking the vial).
- Let the sample cool down shortly on the lab bench (1 - 2 min).
- Add 16 mL (or 15 mL in case of liquid samples) pre-heated A-AEB (see chapter 9) to the cooked sample.
- Mix vigorously (e.g. vortexer).
- Extract for 10 min at 60 °C (140 °F) in a water bath.
- Let the sample cool down shortly in ice water (3 - 5 min).
- Filter sample or centrifuge for 10 min at > 2,500 x g; if possible at 4 °C (39 °F).

Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge.

- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The extract (supernatant of centrifugation step or filtrate) can be stored undiluted in a well-sealed container at 2 - 8 °C until used in the test (shelf life approx. 3 days). Unused extracts can also be stored undiluted at -20 °C (- 4 °F) for several months.
- For use in the test, extracts prepared with Extractor 2 (supernatant of centrifugation step or filtrate) must be diluted with AEB (for dilution see chapter 10.2). The diluted sample extracts have a limited shelf life and must be used in the test within 30 minutes.

- If further dilutions are required (samples with absorbance values ($A_{450\text{ nm}}$) > Standard 5), these extract dilutions must be made before dilution with AEB for test use (see chapter 10.2).
- For these dilutions, the following buffer should be used to keep the composition of the extract constant:

A-AEB 16 mL

Extractor 2 2 mL

Dest. water 2 mL

This is followed by the normal dilution with AEB (see chapter 10.2).

Remark

The extracts prepared with Extractor 2 and A-AEB according to chapter 9.1 can be also used with the RIDASCREEN®FAST β -Lactoglobulin (Art. No. R4912) and RIDASCREEN®FAST Milk (Art. No. R4652) assays.

9.2. Extraction with AEB for food like non-hydrolyzed infant formula, ice cream, wine, chocolate, beverages, sausage, rice crispies

The following extraction procedure is only suited for the mentioned food since they can also be prepared without Extractor 2, Additive 1 and cooking step.

Heat the AEB to 60 °C (140 °F) before sample extraction

Homogenize a representative, adequately big amount of a solid sample (e.g. 50 g; grind it thoroughly to powder and mix well) or mix well a sample in case of liquid foods.

- Transfer 1 g (or in case of liquid samples 1 mL) homogenized sample to a new vial and add 20 mL (or 19 mL in case of liquid samples) pre-heated AEB (see chapter 9).

1 ml wine can also be extracted with 9 mL diluted and heated AEB. Doing this increases the method's sensitivity (LOD: 0.12 mg/L; LOQ: 0.25 mg/L).

- Mix vigorously (e.g. vortexer).
- Extract for 10 min at 60 °C (140 °F) in a water bath by shaking occasionally.
- Let the sample cool down shortly in ice water (3 - 5 min).
- Filter sample or centrifuge for 10 min at > 2,500 x g; if possible at 4 °C (39 °F).

Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed ($> 10,000 \times g$) for 10 min in a microcentrifuge.

- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- Use the sample extracts (supernatant of centrifugation step or filtrate) immediately (within 30 minutes) in the assay. A longer time period may influence the recovery. Alternatively, the extracts can be stored in a well-sealed vial at 2 - 8 °C until it is used (shelf life approx. 1 day). Unused extracts can be stored undiluted at -20 °C (-4 °F) for several months.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in conjugate buffer (e.g. 2 mL conjugate buffer + 200 µL conjugate, sufficient for 2 microtiter strips).

The **wash buffer** is provided as a 10-fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with dist. water (i.e. 900 mL dist. water + 100 mL buffer concentrate). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

Kit components should be stored immediately at 2 - 8 °C (35 - 46 °F) when no longer required.

10.2. Test procedure

The extracts prepared according to chapter 9.1 (with Extractor 2 and A-AEB) must be diluted 1:5 (1+4) with AEB (see chapter 9) before use (e.g. 400 µL AEB + 100 µL sample).

Individual samples extracted with Extractor 2 and A-AEB (e.g. corn or corn products, pine, sunflower or pumpkin seeds) may produce non-specific matrix effects in the ELISA. To avoid these, such sample extracts can be diluted 1:5 with BSA-AEB (see chapter 9) before use in the test

(e.g. 400 µL BSA-AEB + 100 µL sample). For this purpose, the diluted AEB is mixed with BSA (final concentration 2.5 %): e.g. 10 mL diluted AEB + 0.25 g BSA.

Use the diluted sample extracts immediately (within 30 minutes) in the assay. A longer time period may influence the recovery.

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1) should be used as a pre-plate to avoid a time shift over the microtiter plate due to pipetting. All standards and samples are pipetted into the uncoated plate (at least 150 µL per well) and then exactly 100 µL are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore, cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µL of each standard or sample extract (prepared according to chapter 9 and - if necessary - diluted according to chapter 10.2) in duplicate to the wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
4. Add 100 µL of the diluted conjugate (see chapter 10.1) to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three cycles).

6. Add 100 µL of substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µL of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

11. Evaluation

A special software, **RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter or cubic spline function. Once an evaluation method has been selected, it should be retained and not be switched between the two functions. Due to its mathematical background, the cubic spline function cannot calculate concentrations outside the measuring range (< standard 2 and > standard 5). The 4-parameter function enables the calculation of values between standard 1 and standard 2 (refer to chapter 13).

It must be clarified that all quality criteria are met for the current test run. The course of the standard curve can be taken from the quality assurance certificate (certificate of analysis, CoA), which is available via the QR code on the test kit label. As the absorbance values in the laboratory may differ from those stated on the certificate, it is recommended to compare the ratios of the standards to each other with those on the certificate. For this purpose, the $B/B_{\max}$ values (the ratio of the absorbance values of the standards to the highest standard) are compared with each other. In the current test run, these should be comparable to the ratios of the standards on the certificate.

A further sample dilution and new determination is recommended for samples with an absorption ($A_{450\text{nm}}$) > standard 5. In case of a further dilution, the additional dilution factor must be taken into account when calculating the allergen concentration.

The assay can be also evaluated when running in single determinations. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not consistent with standards like EN 15633-1 and EN 15842. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, when performed as a single determination, a

higher degree of inaccuracy can be expected, and when performing repeated tests, a higher degree of variation in the results can be expected.

Extraction with Extractor 2 and A-AEB (see chapter 9.1):

When working according with this extraction, the sample dilution factor is 100. A dilution factor of 20 is already taken into account with the standard concentrations (see chapter 4*). Hence, values read from the standard curve need to be multiplied with factor 5 only.

Extraction with AEB (see chapter 9.2):

When working in accordance with this extraction, the sample dilution factor is 20. Hence, a dilution factor of 20 is already taken into account with the standard concentrations (see chapter 4*), the casein concentration can directly be read from the standard curve.

If wine is extracted with 9 mL AEB only, a lower dilution of 1:10 needs to be respected when calculating the casein concentration. Values read from the standard curve need to be multiplied with factor 0.5 in this case.

12. Result interpretation

The test result is calculated in mg casein per kg food and thus indicates the protein concentration.

Higher absorbance values ($A_{450\text{nm}}$) of the calibration curve, especially for the zero standard, as mentioned on the certificate may be a result of insufficient washing or a casein contamination of reagents.

Results between LoD and LoQ indicate a low allergen concentration in the sample. Depending on the matrix tested, values below the LoQ can still be determined with sufficient precision ($\text{CV} < 30\%$). However, values in this range are generally subject to greater uncertainty due to the higher fluctuation range of the test. If the precision of the test has not been validated with a specific sample matrix, results below the measuring range should therefore not be reported with a quantitative value, but qualitatively "< LoQ". Further information on this can be found in the validation report.

A result below the LoD does not exclude an allergen contamination below the detection limit of the assay, or that other allergen components, such as lipids, may be present in a sample. The result should be reported accordingly.

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LOD, LOQ, upper limit of measurement range).

Raw fish strongly binds casein. Therefore, the recovery may only be 10 %. This is also true for cooked fish if casein is added before cooking.

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). Sufficient accuracy is given with a fluctuation of max ± 1 %.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

In general, matrix effects in immunological test methods can be suppressed by adding foreign protein (depending on the test, e.g., BSA, gelatin, skim milk powder) during extraction. The influence of the addition on recovery should be checked by spiking experiments.

The recommendation for the use of BSA mentioned under chapter 10.2 should only be applied if a matrix effect is actually present. It is also not limited to the samples mentioned, but is to be understood as a basic recommendation for samples with matrix effect.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

Cross-reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. In single cases, cross-reactivities can be lost or only become apparent through the processing of foods.

For evaluation of the cross-reactivity only one representative sample was analyzed. Other samples may show a different result. All analyzed cross-reactivities are described in the validation report.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

14. Recommendation

In order to ensure a high analytical performance, we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in standards like EN 15633-1 and EN 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control and to ensure an accurate and correct test procedure. Allergen-free and allergen-containing (spiked) samples should be used. An example of a spike experiment is given in the validation report.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction may be necessary.
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt® / Bolt™) are used.

15. Further application notes

- Sample extraction for heavily liquid absorbing matrices using RIDA® Extractor 2 (Art. No. R4613).
- More sensitive allergen analysis for wine (dilution 1:10).
- RIDASCREEN®FAST Allergen - Swabbing method for the qualitative analysis of allergens in the production line or for laboratory equipment.

- Preparation of food samples with the RIDA® Extraction Solution (colorless) (Art. No. R7098).

This sample preparation shows similar results as the method with A-AEB and Extractor 2. However, the so prepared extracts cannot be used with the RIDASCREEN®FAST β -Lactoglobulin (Art. No. R4912) and in the RIDASCREEN®FAST Milk (Art. No. R4652).

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2019-05-03	Release version
2021-06-16	Previous version
2022-05-06	Previous version
2025-10-20	Current version Changes made: – General linguistic and formal revision – Addition and revision of related products for milk determination – Addition of intended use in chapter 6 – Addition of crystal formation and preparation of A-AEB in chapter 9 – Revision of handling single determinations in chapter 11 – Revision and additions in chapter 12, 13 and 14 – New disclaimer – Revised Board of Management

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

1. In conformance with the German Civil Code ("BGB") R-Biopharm AG provides a limited warranty ("Gewährleistung") against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm's own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG's liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG's sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser's expenses, losses, or damages from loss of goodwill, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements ("Gewährleistung") implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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are registered trademarks of R-Biopharm AG.