RIDASCREEN®FAST Milk

REF R4652

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RIDASCREEN®FAST Milk

Brief information

RIDASCREEN®FAST Milk (Art. No. R4652) is a sandwich enzyme immunoassay for the quantitative determination of milk protein in food validated for the method (see chapter 1).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization and extraction
Time requirement:	sample preparation (for 10 samples) approx. 45 min test implementation (incubation time).....30 min
Standard material:	The RIDASCREEN® standard material is calibrated to the reference material NIST RM 1549a whole milk powder
Limit of detection: (matrix-dependent)	0.7 mg/kg (ppm) milk protein
Limit of quantification:	2.5 mg/kg (ppm) milk protein corresponds to 7.1 mg/kg (ppm) skim milk powder
Specificity:	The antibodies used specifically detect caseins and β -lactoglobulin of cow's milk. There is cross-reactivity to sheep's, goat's and buffalo's milk. Further information are contained in the validation report.

Cross-reactivities of the used antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross-reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments (see chapter 13).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

Related products for milk protein determination

RIDASCREEN®FAST Casein (Art. No. R4612)

RIDASCREEN®FAST β -Lactoglobulin (Art. No. R4912)

RIDASCREEN® β -Lactoglobulin (Art. No. 4901)

Bioavid Lateral Flow Milch / Milk (Art. No. BL623-15)

Bioavid Lateral Flow Casein incl. Hook Line (Art. No. BLH714-15)

RIDA® Extractor 2 (Art. No. R4613)

1. Intended use

RIDASCREEN®FAST Milk (Art. No. R4652) is a sandwich enzyme immunoassay test developed for the quantitative analysis of milk proteins in food, which may contain milk or milk powder. Due to the large number of different food products, the following samples were examined representative for different food product categories within the scope of the test development: cookies, infant formula, chocolate, ice cream and sausage.

The test kit should generally be able to provide usable results for other foods within the methodological limits (see chapter 13). In individual cases, the user must check the suitability of the test kit by carrying out appropriate experiments before use.

For detailed results and further information on validation data with other food matrices, please refer to our brochure validation report. Further applications are regularly validated in our laboratories, which we make available in our application notes (see chapter 15).

2. General information

Food allergy is defined as an adverse immunologic response to a dietary protein like peanut or egg. Allergies are IgE-mediated dysfunctions of the immune system. In contrast, food intolerance like lactose intolerance is caused by lack of certain enzymes (= metabolic disease). Antibodies are not involved

in the disease's development. The allergy to milk is to be distinguished from the lactose intolerance: In the case of intolerance, the symptoms are caused by a certain type of sugar (milk sugar, lactose), in the case of allergy, however, by certain proteins in the milk. Milk allergy is an allergy in the classic sense, e.g. the immune system fights a substance that is actually harmless to the body and allergic symptoms arise. Milk is consumed directly, as a drink, or as a processed product such as cream, butter, yogurt or cheese. Cow's milk powder is widely used in baked goods, sauces and dessert.

Cow's milk contains 3.2 % proteins which consist of approx. 10 % β -lactoglobulin (leading protein of whey) and 80 % caseins. Casein is heat-stable, while β -lactoglobulin is heat-sensitive. Most milk allergy sufferers are allergic to both protein fractions in milk. The most important allergen for children is β -lactoglobulin while the caseins become to be dominant later in adults. The allergen can be present as an ingredient or as a contamination in raw and cooked / processed products. According to the **regulation (EU) No. 1169/2011** milk and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

Milk contains approx. 3.2 % milk protein. Thus, a sample which contains 1 mg/kg milk protein corresponds to a milk content of approx. 31 mg/kg milk.

3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against casein-proteins and β -lactoglobulin-proteins. By adding the standard or sample solution to the wells, milk proteins present in the sample will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Following the washing step, a solution containing casein- and β -lactoglobulin antibodies conjugated to peroxidase is added. This conjugate binds to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product. A stop solution is added which results in a color change from blue to yellow. The absorbance of the solution which is proportional to the milk proteins concentration in the sample is measured photometrically at 450 nm and expressed as mg/kg milk proteins.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
3 x Extractor 2	Blue	Concentrate	2x	30 mL
Allergen extraction buffer	Green	Concentrate	10x	100 mL
Additive 1	Blue			2 g
Standard 1*	Transparent	Ready to use	0 mg/kg	1.3 mL
Standard 2*	Transparent	Ready to use	2.5 mg/kg	1.3 mL
Standard 3*	Transparent	Ready to use	7.5 mg/kg	1.3 mL
Standard 4*	Transparent	Ready to use	22.5 mg/kg	1.3 mL
Standard 5*	Transparent	Ready to use	67.5 mg/kg	1.3 mL
Wash buffer	Brown	Concentrate	10x	100 mL
Conjugate	Red	Concentrate	11x	0.7 mL
Conjugate buffer	Black	Ready to use		7 mL
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 mL
Stop solution	Yellow	Ready to use		14 mL

* The dilution factor 100, which results after sample preparation in chapter 9, has already been considered for the standard concentrations. Therefore, the milk protein concentrations of the samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1. Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (60 °C / 140 °F and 100 °C / 212 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12 μ m)
- Graduated pipettes
- Graduated cylinder
- Variable 20 - 200 μ L and 200 - 1000 μ L micropipettes

- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100 µL
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)

5.2. Reagents

- Distilled water (dist. water) or deionized water
- 1 M sodium hydroxide (NaOH)
- 1 M hydrochloric acid (HCl)
- If necessary: bovine serum albumin (BSA, e.g. Serva, fraction V, protease free, Art. No. 11926)

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate and pre-plate, if necessary, see chapter 10.2). Use separate pipette tips for each standard and each sample extract to avoid cross-contamination.

The Extractor 2 is harmful to health. It contains mercaptoethanol. It should be worked **under a chemical hood** and skin contact should be avoided (use gloves!).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation.
- Value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 1.2$) for standard 5.

9. Preparation of Samples

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment may lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation.
- Carry out the sample preparation in a room isolated from the ELISA procedure.

The samples should be stored in a cool place and protected against light before extraction.

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The Allergen extraction buffer is provided as a 10-fold concentrate and must be diluted prior use. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that, dilute the heated buffer concentrate 1:10 (1+9) with dist. water (e.g. 900 mL dist. water + 100 mL buffer concentrate). The **diluted Allergen extraction buffer (AEB)** is stable at 2 - 8 °C (35 - 46 °F) for 12 weeks. Current data show

that the AEB can alternatively be stored at room temperature (20 - 25 °C / 68 - 77 °F) for up to four weeks (not part of the AOAC *Performance Tested Methods*SM (PTM) approvals).

The AEB is required for the dilution of the extracts. For the dilution of matrices containing pine, sunflower or pumpkin seeds, the diluted AEB must be additionally mixed with BSA (final concentration: 2.5 %; e.g. 10 mL diluted AEB + 0.25 g BSA (= **BSA-AEB**)).

For the preparation of the final **Allergen extraction buffer containing Additive 1 (A-AEB)**, weigh 1.35 g of Additive 1 in a glass beaker and add 15 mL 1 M NaOH. Stir until the Additive 1 is solved. Then, fill 700 mL diluted AEB (see above) in a measuring cylinder. Add the 15 mL Additive 1 solution by stirring constantly. Transfer any residues of the Additive 1 solution into the measuring cylinder by rinsing with diluted AEB. Adjust the Additive 1 containing Allergen extraction buffer (A-AEB) to pH 9 with 1 M HCl and fill up to 750 mL with diluted AEB.

750 mL A-AEB are sufficient for 45 samples. The buffer can be used for approx. 3 weeks at room temperature 20 - 25 °C (68 - 77 °F) (do **not** store in the refrigerator!). Discard the buffer if crystals are present. Crystal formation may begin after 10 days. We therefore recommend preparing a fresh buffer solution every week. Alternatively, a reduced volume of the A-AEB of e.g. 375 mL can be prepared. When preparing a smaller volume of A-AEB, the composition described above must be adhered to accordingly. Use clean bottles when preparing the buffer. Particles of dust can initiate crystallization.

Make sure to heat the A-AEB in time in the 60 °C (140 °F) water bath (the 100 °C / 212 °F water bath is used for the extraction of the samples).

The **Extractor 2** is provided as 2fold concentrate and has to be diluted 1:2 (1+1) with dist. water (e.g. 30 mL Extractor 2 + 30 mL dist. water). The complete diluted Extractor 2 is sufficient for 15 samples and can be used for approx. 3 month at room temperature (20 - 25 °C / 68 - 77 °F). Additional Extractor 2 can be ordered at R-Biopharm with the article number R4613.

In the following section, the following abbreviations are used:

- AEB: final diluted Allergen extraction buffer
- A-AEB: AEB with addition of Additive 1

9.1. Extraction

Heat the A-AEB to 60 °C (140 °F) before sample extraction.

Homogenize a representative, adequately big amount of a solid sample (e.g. 50 g; grind it thoroughly to powder and mix well) or mix well a sample in case of liquid foods.

- Transfer 1 g (or in case of liquid samples 1 mL) of homogenized sample to a new vial, add 4 mL prepared Extractor 2 (see chapter 9), close the vial and mix well.
- Cook it for 10 min at 100 °C (212 °F) in a water bath (if a lump forms, it must be dissolved in the still warm liquid by shaking the vial).
- Let the sample cool down shortly on the lab bench (1 - 2 min).
- Add 16 mL (or 15 mL in case of liquid samples) pre-heated A-AEB (see chapter 9.1) to the cooked sample.
- Mix vigorously (e.g. vortexer).
- Extract for 10 min at 60 °C (140 °F) in a water bath.
- Let the sample cool down shortly in ice water (3 - 5 min).
- Filter sample or centrifuge for 10 min at $> 2,500 \times g$; if possible at 4 °C (39 °F).

Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed ($> 10,000 \times g$) for 10 min in a microcentrifuge.

- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The extract (supernatant of centrifugation step or filtrate) can be stored undiluted in a well-sealed container at 2 - 8 °C until used in the test (on the same day). Unused extracts can also be stored undiluted at -20 °C (-4 °F) for one month.
- For use in the test, the extracts (supernatant from centrifugation step or filtrate) must be diluted with AEB (for dilution see chapter 10.2). The diluted sample extracts have a limited shelf life and must be used in the test within 30 minutes.
- If further dilutions are required (samples with absorbance values ($A_{450 \text{ nm}}$) $>$ Standard 5), these extract dilutions must be made before dilution with AEB or BSA-AEB for test use (see chapter 10.2).
- For these dilutions, the following buffer should be used to keep the composition of the extract constant:

A-AEB	16 mL
Extractor 2	2 mL
Dest. water	2 mL

This is followed by the normal dilution with AEB (see chapter 10.2).

Remark

Sample extracts with Extractor 2 and A-AEB can also be used both in the RIDASCREEN®FAST β -Lactoglobulin (Art. No. R4912) and in the RIDASCREEN®FAST Casein (Art. No. R4612) assays.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in conjugate buffer (e.g. 2 mL conjugate buffer + 200 μ L conjugate, sufficient for 2 microtiter strips).

The **wash buffer** is provided as a 10-fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with dist. water (e.g. 900 mL dist. water + 100 mL buffer concentrate). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

Kit components should be stored immediately at 2 - 8 °C (35 - 46 °F) when no longer required.

10.2. Test procedure

The extracts prepared with Extractor 2 and A-AEB must be diluted 1:5 (1+4) with AEB (see chapter 9) before usage (e.g. 400 μ L AEB + 100 μ L sample).

Individual samples extracted with Extractor 2 and A-AEB (e.g. corn or corn products, pine, sunflower or pumpkin seeds) may produce non-specific matrix effects in the ELISA. To avoid these, such sample extracts can be diluted 1:5 with BSA-AEB (see chapter 9) before use in the test (e.g. 400 μ L BSA-AEB + 100 μ L sample). For this purpose, the diluted AEB is mixed with BSA (final concentration 2.5 %): e.g. 10 mL diluted AEB + 0.25 g BSA.

Use the diluted sample extracts immediately (within 30 minutes) in the assay. A longer time period may influence the recovery.

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µL per well) and then exactly 100 µL are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore, cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µL of each standard or sample extract (prepared according to chapter 9 and diluted according to chapter 10.2) in duplicate to the wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat three more times (a total of four wash cycles).
4. Add 100 µL of the diluted conjugate (see chapter 10.1) to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat three more times (a total of four wash cycles).
6. Add 100 µL of the reddish substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µL of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 min¹ after addition of stop solution.

¹ Current data show that measurement is possible up to 30 min after addition of the stop solution (not part of the AOAC *Performance Tested Methods*SM (PTM) approval). For the AOAC-PTM approval, a measurement within 10 min was specified.

11. Evaluation

A special software, **RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter or cubic spline function. Once an evaluation method has been selected, it should be retained and not be switched between the two functions. Due to its mathematical background, the cubic spline function cannot calculate concentrations outside the measuring range (< standard 2 and > standard 5). The 4-parameter function enables the calculation of values between standard 1 and standard 2 (refer to chapter 13). The 4-parameter evaluation is not part of the AOAC-PTM approval. For AOAC-PTM approval, the cubic spline function was used.

The test result is calculated in mg milk-protein per kg of food and thus indicates the protein concentration.

When working in accordance with this extraction, the sample dilution factor is 100. Hence, a dilution factor of 100 is already taken into account with the standard concentrations (see chapter 4*), the milk-protein concentration can directly be read from the standard curve.

A further sample dilution and new determination is recommended for samples with an absorption ($A_{450\text{nm}}$) > standard 5. In case of further dilution, the additional dilution factor must be taken into account when calculating the result.

It must be clarified that all quality criteria are met for the current test run. The course of the standard curve can be taken from the quality assurance certificate (certificate of analysis, CoA), which is available via the QR code on the test kit label. As the absorbance values in the laboratory may differ from those stated on the certificate, it is recommended to compare the ratios of the standards to each other with those on the certificate. For this purpose, the B/B_{max} values (the ratio of the absorbance values of the standards to the highest standard) are compared with each other. In the current test run, these should be comparable to the ratios of the standards on the certificate.

The assay can be also evaluated when running in single determinations. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not

consistent with standards like EN 15633-1 and EN 15842 and also does not correspond to the procedure for the AOAC Approval. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, when performed as a single determination, a higher degree of inaccuracy can be expected, and when performing repeated tests, a higher degree of variation in the results can be expected.

12. Result interpretation

Higher absorbance values ($A_{450\text{nm}}$) of the calibration curve, especially for the zero standard, as mentioned on the certificate may be a result of insufficient washing or a hazelnut contamination of reagents.

Results between LoD and LoQ indicate a low allergen concentration in the sample. Depending on the matrix tested, values below the LoQ can still be determined with sufficient precision ($\text{CV} < 30\%$). However, values in this range are generally subject to greater uncertainty due to the higher fluctuation range of the test. If the precision of the test has not been validated with a specific sample matrix, results below the measuring range should therefore not be reported with a quantitative value, but qualitatively "< LoQ". Further information on this can be found in the validation report.

A result below the LoD does not exclude an allergen contamination below the detection limit of the assay, or that other allergen components, such as lipids, may be present in a sample. The result should be reported accordingly.

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LOD, LOQ, upper limit of measurement range).

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max $\pm 1\%$.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

In general, matrix effects in immunological test methods can be suppressed by adding foreign protein (depending on the test, e.g., BSA, gelatin, skim milk powder) during extraction. The influence of the addition on recovery should be checked by spiking experiments.

The recommendation for the use of BSA mentioned under chapter 10.2 should only be applied if a matrix effect is actually present. It is also not limited to the samples mentioned, but is to be understood as a basic recommendation for samples with matrix effect.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

Cross-reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. In single cases, cross-reactivities can be lost or only become apparent through the processing of foods.

For evaluation of the cross-reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross-reactivities are described in the validation report.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

Remarks:

Different milk components like milk powder, whey (β -lactoglobulin), caseinates, etc. are added to food products. RIDASCREEN®FAST Milk is calibrated to skim milk powder. Therefore, if the analyzed food contains:

- Skim milk powder, whole milk powder or milk, then the ELISA results are quantitative.
- Whey (β -lactoglobulin), then the ELISA results may be underestimated.
- Caseinates, then the ELISA results may be underestimated.

If a quantitative result for caseins and / or β -lactoglobulin is requested, the sample should be tested in the RIDASCREEN®FAST Casein (Art. No. R4612) and / or the RIDASCREEN®FAST β -Lactoglobulin (Art. No. R4912).

If the result is below the detection limit, further milk constituents (e.g. sugars like lactose) may be present in the food sample.

14. Recommendation

In order to ensure a high analytical performance, we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in standards like EN 15633-1 and EN 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control and control and ensurance of an accurate and correct test procedure. Allergen-free and allergen-containing (spiked) samples should be used. An example of a spike experiment is given in the validation report.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 to 7.5) to neutral prior to extraction may be necessary.
- The recovery for fragmented proteins is reduced in sandwich ELISA, such samples should be analyzed with a competitive ELISA test systems like the RIDASCREEN® β -Lactoglobulin (Art. No. R4901).
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt® / Bolt™) are used.

15. Further application notes

- Sample extraction for heavily liquid absorbing matrices using RIDA® Extractor 2 (Art. No. R4613).

- RIDASCREEN®FAST Allergen - Swabbing method for the qualitative analysis of allergens in the production line or for laboratory equipment.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2015-07-06	Release version
2021-07-01	Previous version
2021-11-15	Previous version
2025-10-20	Previous version
2026-01-14	Current version Changes made: – Chapter 11: Note regarding evaluation with 4-parameter function

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

1. In conformance with the German Civil Code ("BGB") R-Biopharm AG provides a limited warranty ("Gewährleistung") against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm's own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG's liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG's sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser's expenses, losses, or damages from loss of goodwill, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements ("Gewährleistung") implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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are registered trademarks of R-Biopharm AG.