RIDA® QUICK Aflatoxin RQS

REF R5208

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RIDA®QUICK Aflatoxin RQS

Brief information

RIDA®QUICK Aflatoxin RQS (R5208) is an immunochromatographic test for the quantitative determination of aflatoxin in corn. Aflatoxin extraction is performed with 70 % methanol or 50 % ethanol.

The test kit contains 20 test strips, each of which can be used for one analysis.

Results are evaluated with

- the RIDA®SMART APP software (ZRSAM1000) and an approved smartphone or
- the RIDA®SMART APP software (ZRSAM1000) installed on an Android smartphone and the RIDA®SMART BOX (ZRSA-SB) or
- the RIDA®QUICK SCAN reader.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) approx. 10 min
test implementation (incubation time) 3 min

Limit of detection (LOD): < 2 µg/kg (ppb)

Limit of quantification (LOQ): 2.6 µg/kg (70 % methanol extraction)
..... 2.9 µg/kg (50 % ethanol extraction)

Measuring range: 2 - 75 µg/kg and 50 - 300 µg/kg
Quantitative evaluation with the RIDA®SMART APP

Specificity: The RIDA®QUICK Aflatoxin RQS test detects aflatoxin in corn. Further information is contained in the validation report.

The specificity of the RIDA®QUICK Aflatoxin RQS test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample

matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products

RIDASCREEN® Aflatoxin Total	(Art No. R4701)
RIDASCREEN®FAST Aflatoxin	(Art No. R5202)
RIDA®QUICK Aflatoxin RQS ECO	(Art No. R5209)
RIDA®QUICK Mycotoxin ECO Extractor	(Art. No. R5000)
RIDA® Aflatoxin column	(Art No. R5001 / R5002)
RIDA®SMART APP Software	(Art No. ZRSAM1000)
RIDA®SMART BOX	(Art No. ZRSA-SB)
Trilogy® certified reference material and standards for Aflatoxin (ISO 17034)	

1. Intended use

RIDA®QUICK Aflatoxin RQS is an immunochromatographic test in strip format for the quantitative determination of aflatoxin in corn. Aflatoxin extraction is performed using 70 % methanol or 50 % ethanol.

2. General information

Aflatoxins are secondary metabolites of the fungi species *Aspergillus flavus*, *parasiticus* and *nomius*. These fungal species are found in humid, tropical areas and can contaminate plant-based food and feed. This contamination usually takes place in the countries of cultivation. Aflatoxins belong to the strongest natural occurring carcinogenic substances.

3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-aflatoxin antibody detects aflatoxin in the sample. During incubation of the test strip, a band pattern (test line / control line) forms that is used to determine the concentration of aflatoxin. The intensity of the test line depends on the aflatoxin concentration of the sample. It increases as the aflatoxin concentration increases. It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the aflatoxin concentration in the sample increases.

The test strip is evaluated using the RIDA®SMART APP software (ZRSAM1000) or the RIDA®QUICK SCAN reader.


Fig. 1: RIDA®QUICK Aflatoxin RQS test strip.

4. Reagents provided

Each kit contains sufficient materials for max. 20 determinations. Each test kit contains:

Component	Cap color	Format	Volume
20 x Test strip	-	ready to use	One determination per test strip (separately packed)
1 x Mobile solvent	white	ready to use	4.5 mL
1 x Dilution buffer	transparent	ready to use	20 mL
1 x RIDA®SMART APP cover	-	ready to use	Evaluation cover for use with RIDA®SMART APP

5. Equipment and reagents required but not provided

5.1 Equipment

- Gloves
- Laboratory grinder
- Weighing scale
- Graduated cylinder
- Shaker
- Centrifuge (≥ 2000 g)
- 100 μ l pipette
- Latest version of the RIDA®SMART APP software (ZRSAM1000) activated and installed on a compatible smartphone.

5.2 Reagents

- 70 % methanol (e.g. 70 mL methanol (100 %) / 30 mL dist. water)
- 50 % ethanol (e.g. 50 mL ethanol (denatured) / 50 mL dist. water)

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should be carried out only by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 47 °F). Do not freeze the test strips or any test kit components.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. Do not remove the test strips from the test strip packaging until immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature. Do not allow methanol or ethanol to evaporate.

The sample should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

8.1 Measurement range 2 - 75 µg/kg

- Weigh 10 g of the ground and homogenized sample and add 20 mL methanol (70 %) or 20 mL ethanol (50 %)
- Close the tube and shake the sample for 5 min (using a shaker)
- Centrifuge (1 min, 2000 g) until a particle-free solution can be obtained
- Add 100 µL of the particle-free supernatant to 200 µL of mobile solvent and mix carefully (= Dilution I)
- Use 100 µL of Dilution I in the test

Note:

Please avoid methanol evaporation during sample preparation.

If necessary, the sample weight can be increased, e.g. 20 g sample in 40 mL 70 % methanol or 40 mL ethanol (50 %).

8.2 Measurement range 50 - 300 µg/kg

For an extended measuring range of 50 - 300 µg/kg, Dilution I from sample preparation 8.1. must be diluted further.

- Dilute 100 µL of Dilution I with 900 µL dilution buffer and mix carefully (= Dilution II)
- Use 100 µL of Dilution II in the test

9. Test procedure

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

- Remove a test strip from its foil pouch
- Apply 100 µL of prepared sample solution (Dilution I or Dilution II, see 8. Sample preparation) to the application area of the test strip (see Fig. 1)
- Incubate the test strip at room temperature
- After exactly 3 min of incubation, evaluate the test strip with the RIDA®SMART APP or the RIDA®QUICK SCAN

10. Evaluation

The left line in the reaction area is the test band (test line, see Fig. 1). Its appearance and intensity depends on the aflatoxin concentration of the sample. The sample is contaminated with aflatoxin, if the control band (control line) is visible and the test band (test line) is also visible.

The right line in the reaction area is a control band (control line) and must be present after each test procedure. If this line is not visible, the test was not performed correctly or the reagents were not in order. In this case, the test should be repeated with a new test strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample was prepared according to 8.1 and is highly contaminated with aflatoxin, the control line will be only faintly visible. In this case the test result might be reported as „>75 ppb or invalid result“ (RIDA®SMART APP) or “>75 ppb or invalid“ (RIDA®QUICK SCAN reader). In this instance, dilute the sample as described in 8.2 and repeat the test with a new test strip.

To obtain comparable results, the test strips must always be evaluated after 3 min.

10.1 Evaluation with RIDA®SMART APP

First, please read the instructions for use of the RIDA®SMART APP software attentively. The test strip evaluation has to be performed by using the RIDA®SMART APP software according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on our website www.r-biopharm.com.

For the test strip evaluation the RIDA®SMART APP cover is needed which is included in the test kit.

RIDA®SMART APP software applications for evaluation:

Matrix	Measuring range	RIDA®SMART APP application
Corn	70 % Methanol	Methanol 2 - 75 ppb
		Methanol 50 - 300 ppb
Corn	50 % Ethanol	Ethanol 2 - 75 ppb
		Ethanol 50 - 300 ppb

10.2 Evaluation with RIDA®QUICK SCAN

First, carefully read the user guide for the RIDA®QUICK SCAN and proceed as instructed. Point 3 contains the description for evaluating the test strips.

To obtain optimal results every test kit contains lot-specific standard curve parameters. These parameters can be found on the certificate of analysis (CoA) and have to be used only with this specific lot. The CoA contains parameters for several standard curves. Depending on the extraction method used, select the corresponding parameters.

To enter the parameters manually into the RIDA®QUICK SCAN reader use the plain text information on the CoA. To enter the data with the barcode scanner use the 2D barcodes. In both cases, the checksum shown on the RIDA®QUICK SCAN screen needs to be identical to the respective checksum on the certificate (number below the 2D barcode).

RIDA®QUICK SCAN evaluation methods:

Measuring range	RIDA®QUICK SCAN method
2 - 75 µg/kg	Afla ECO Low
50 - 300 µg/kg	Afla ECO High

11. Recommendation

In order to ensure a high analytical performance we recommend:

- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- To do spike experiments to ensure an accurate and correct test procedure.

12. Sensitivity

This test determines aflatoxin contaminations in the measurement range of 2 - 300 µg/kg in combination with the RIDA®SMART APP software and in the measurement range of 2 - 300 µg/kg in combination with the RIDA®QUICK SCAN reader.

13. Further application notes

Further application notes are available on request.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2019-11-25	Release version
2020-12-10	General revision
2024-02-16	Current version General revision Changes made: <ul style="list-style-type: none">- 1. Intended use: products added- 3. Test principle: Figure added- 4. Reagents provided: universal Smart App cover- 6. Warnings and precautions for the use completed- 11./12./13. Additions- Imprint updated

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.