RIDA® QUICK Ochratoxin ECO

REF R5404

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RIDA®QUICK Ochratoxin ECO

Brief information

RIDA®QUICK Ochratoxin ECO (Art. No. R5404) is an immunochromatographic test for the quantitative determination of ochratoxin in corn and wheat.

The test kit contains 20 test strips, each of which can be used for one analysis. All reagents required for the assay are contained in the test kit. The results are evaluated quantitatively with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000).

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)... approx. 10 min
test implementation (incubation time) approx. 3 - 5 min

Limit of detection (LOD): 2 µg/kg (ppb)

Limit of quantification (LOQ): 4 µg/kg

Measuring range: 2 - 30 µg/kg and 30 - 150 µg/kg
Quantitative evaluation with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000)

Specificity: The RIDA®QUICK Ochratoxin ECO test reacts with ochratoxin in corn and wheat.

Further information is contained in the validation report.

The specificity of the RIDA®QUICK Ochratoxin ECO test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products and accessories for ochratoxin determination

RIDASCREEN® Ochratoxin A 30/15	(Art. No. R1312)
RIDASCREEN®FAST Ochratoxin A	(Art. No. R5402)
RIDA® Ochratoxin A column	(Art. No. R1303)
RIDA®SMART Mycotoxin ECO Extractor	(Art. No. R5000)
RIDA®SMART APP software	(Art. No. ZRSAM1000)
RIDA®SMART BOX	(Art. No. ZRSA-SB)
Validated Android smartphone	(see product catalogue)

1. Intended use

RIDA®QUICK Ochratoxin ECO is a quantitative immunochromatographic test in strip format for the quantitative determination of ochratoxin in corn and wheat.

2. General information

The mycotoxin ochratoxin A is formed by fungi of the genus *Aspergillus* and *Penicillium*. In addition to its pronounced nephrotoxicity, ochratoxin A exhibits hepatotoxic, teratogenic, carcinogenic and immunosuppressive properties. The ingestion of contaminated food and feed of plant origin, but also of animal origin, poses a particular health risk to humans. Ochratoxin A has already been detected in pig blood and kidneys and in human blood and breast milk.

3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-ochratoxin antibody detects ochratoxin in the sample. During incubation of the test strip, a band pattern (test line / control line) forms that is used to determine the concentration of ochratoxin. The intensity of the test line depends on the ochratoxin concentration of the sample. It increases as the ochratoxin concentration increases. It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the ochratoxin concentration in the sample increases.

The results are evaluated quantitatively with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000).


Figure 1: RIDA®QUICK Ochratoxin ECO test strip.

4. Reagents provided

Each kit contains sufficient materials for max. 20 determinations. Each test kit contains:

Component	Cap color	Format	Volume
20 x Test strip	-	Ready-to-use	One determination per test strip (separately packed)
1 x Extraction solution	Transparent	10x concentrate	110 mL
1 x RIDA®SMART APP cover	-	Ready-to-use	Evaluation cover for use with RIDA®SMART APP

5. Reagents required but not provided

5.1 Equipment

- Laboratory grinder
- Balance
- Graduated cylinder
- Optional: Shaker (e.g. Heidolph Multi-Reax or equivalent)
- Centrifuge (≥ 2000 g)
- Variable 20 - 200 μ L and 200 - 1000 μ L micropipettes
- RIDA®SMART BOX (Art. No. ZRSA-SB) in combination with a smartphone and the latest version of the RIDA®SMART APP software (Art. No. ZRSAM1000)
- A compatible Android smartphone (see product catalogue) and the latest version of the RIDA®SMART APP software (Art. No. ZRSAM1000)

5.2 Reagents

- Distilled or deionized water

6. Warnings and precautions for the users

The product / test is only suitable for use within the scope of its intended use.

This test should be carried out only by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze the test strips or any test kit components.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. Do not remove the test strips from the test strip packaging until immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

The sample should be stored in a cool place, protected against light.

8.1 Extraction solution

A **ready-to-use extraction solution** is required for the extraction. Extraction solution is available as a **10-fold concentrate** and must be diluted 1:10 (1 + 9) with deionized or distilled water before use (e.g. 100 mL concentrate + 900 mL dist. water).

The diluted extraction solution has a shelf life of one week at 2 - 8 °C (36 - 46 °F). If turbidity occurs in the diluted extraction buffer (e.g. caused by contamination), it must be discarded.

8.2 Extraction for corn and wheat

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- Weigh 10 g of the crushed and mixed sample and add 50 mL (for corn) or 30 mL (for wheat) of ready-to-use extraction solution.
- Close tube and shake the sample vigorously for 5 min (manually or with shaker / vortex).
- Allow mixture to sediment (≥ 3 min).
- Use 1 - 2 mL of supernatant and centrifuge (1 min at 2000 g).
- Take 100 μ L of the clear supernatant for use (see chapter 9. Test procedure).
 - Corn samples are stable for 48 hours at 4 - 23 °C (39 - 73 °F) after extraction.
 - Wheat samples are stable for 30 min at 4 - 23 °C (39 - 73 °F) after extraction.

9. Test procedure

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

9.1 Corn

Measuring range: 2 - 30 µg/kg

- Remove a test strip from its foil pouch.
- Apply 100 µL of the sample solution (see chapter 8. Sample preparation) on the application area of the test strip (see fig. 1).
- Incubate the test strip at room temperature for 3 min.
- Evaluate the result after exactly 3 min incubation time with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000).

Measuring range: 30 - 150 µg/kg

- Mix 100 µL of clear supernatant with 500 µL of extraction solution (= dilution 1).
- Remove a test strip from its foil pouch.
- Apply 100 µL of dilution 1 on the application area of the test strip (see fig. 1).
- Incubate the test strip at room temperature for 3 min.
- Evaluate the result after exactly 3 min incubation time with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000).

9.2 Wheat

Measuring range: 2 - 30 µg/kg

- Remove a test strip from its foil pouch.
- Apply 100 µL of the sample solution (see chapter 8. Sample preparation) on the application area of the test strip (see fig. 1).
- Incubate the test strip at room temperature for 5 min.
- Evaluate the result after exactly 5 min incubation time with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000).

Measuring range: 30 - 150 µg/kg

- Mix 100 µL of clear supernatant with 1000 µL of extraction solution (= dilution 1).
- Remove a test strip from its foil pouch.
- Apply 100 µL of dilution 1 on the application area of the test strip (see fig. 1).
- Incubate the test strip at room temperature for 5 min.
- Evaluate the result after exactly 5 min incubation time with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000).

10. Evaluation

The left line in the reaction area is the test band (test line, see fig. 1). Its appearance and intensity depends on the ochratoxin concentration of the sample. The sample is contaminated with ochratoxin, if the control band (control line) is visible and the test band (test line) is also visible.

The right line in the reaction area is a control band (control line) and must be present after each test procedure. If this line is not visible, the test was not performed correctly or the reagents were not in order. In this case, the test should be repeated with a new test strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample is highly contaminated with ochratoxin, the control line will be only faintly visible. In this case the test result might be reported as “invalid result“. In this instance, dilute the sample as described in chapter 9. Test procedure and repeat the test with a new test strip.

To obtain comparable results, the test strips must always be evaluated after 3 min for corn or 5 min for wheat, respectively.

Evaluation with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000)

First, please read the instructions for use of the RIDA®SMART APP software (Art. No. ZRSAM1000) and the RIDA®SMART BOX (Art. No. ZRSA-SB) attentively. The test strip evaluation has to be performed by using the

RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000) according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on R-Biopharm RIDA®SMART APP product page <https://food.r-biopharm.com/products/rida-smart-app/>.

For the test strip evaluation with a validated Android smartphone the RIDA®SMART APP cover is needed which is included in the test kit. The RIDA®SMART APP cover must be used only with the appropriate lot of the test kit. For further information please check the Certificate of Analysis (CoA) which can be downloaded from <https://eifu.r-biopharm.com/>. Please make sure the correct RIDA®SMART APP evaluation cover is present prior to starting the test procedure and test strip evaluation.

RIDA®SMART APP software applications for evaluation:

Matrix	Measuring range	RIDA®SMART APP application
Corn	2 - 30 µg/kg	Corn 2 - 30 ppb
	30 - 150 µg/kg	Corn 30 - 150 ppb
Wheat	2 - 30 µg/kg	Wheat 2 - 30 ppb
	30 - 150 µg/kg	Wheat 30 - 150 ppb

11. Recommendation

In order to ensure a high analytical performance we recommend:

- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- To do spike experiments to ensure an accurate and correct test procedure.

12. Sensitivity

This test determines ochratoxin contaminations in the measurement range of 2 - 150 µg/kg with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone or a validated Android smartphone in combination with the RIDA®SMART APP software (Art. No. ZRSAM1000).

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2023-07-21	Release version

Explanation of symbols

General symbols:

-  Follow the instructions for use
-  Batch number
-  Expiry date (YYYY-MM)
-  Storage temperature
-  Article number
-  Number of test determinations
-  Manufacturing date (YYYY-MM)
-  Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.