
r-biopharm®

RIDA® QUICK Zearalenon RQS

REF R5504

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RIDA[®]QUICK Zearalenon RQS

Brief information

RIDA[®]QUICK Zearalenon RQS (Art. No. R5504) is an immuno-chromatographic test for the quantitative determination of zearalenone in corn (see chapter 1. Intended use).

The test kit contains 20 test strips, each of which can be used for one analysis. All reagents required for the assay are contained in the test kit.

Results are evaluated with

- the RIDA[®]SMART APP software (ZRSAM) and an approved smartphone or
- the RIDA[®]SMART APP software (ZRSAM) installed on an Android smartphone and the RIDA[®]SMART BOX (ZRSA-SB).

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)... approx. 10 min
test implementation (incubation time)..... 5 min

Limit of detection: approx. 50 µg/kg

Measuring range: 50 - 1000 µg/kg

Specificity: The RIDA[®]QUICK Zearalenon RQS test reacts with zearalenone in corn samples.

Further information is contained in the validation report.

The specificity of the RIDA[®]QUICK Zearalenon RQS test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products and accessories for zearalenone determination

RIDASCREEN® Zearalenon	(Art. No. R1401)
RIDASCREEN®FAST Zearalenon	(Art. No. R5502)
RIDASCREEN®FAST Zearalenon SC	(Art. No. R5505)
RIDA®SMART APP software	(Art. No. ZRSAM)
RIDA®SMART BOX	(Art. No. ZRSA-SB)
Trilogy® certified reference material and standards for Zearalenon (ISO 17034)	

1. Intended use

RIDA®QUICK Zearalenon RQS (Art. No. R5504) is an immuno-chromatographic test in strip format for the quantitative determination of zearalenone in corn.

2. General information

The mycotoxin zearalenone is formed by fungi of the genus *Fusarium*. Zearalenone is a phytohormone which displays, apart from its anabolic properties, mainly estrogenic effects. Because of its estrogenic properties, zearalenone may induce fertility disorders in animals with clinical signs of hypoestrogenism - an aspect of a disease which although reported mainly in hogs, is described in other species such as cow, horse and sheep, too.

The potential health risk for humans induced by this mycotoxin, which is taken up with foods of vegetable or animal origin, is extensively discussed.

3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-zearalenone antibody detects zearalenone in the sample. During incubation of the test strip, a continuous band pattern (test line and control line) needs to be formed that is used to determine the concentration of zearalenone. The intensity of the test line depends on the zearalenone concentration of the sample. It increases as the zearalenone concentration increases. A weak test band can also occur with a negative sample (contaminated below the detection limit). It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the zearalenone concentration in the sample increases.

The test strip is evaluated using the RIDA®SMART APP software (Art. No. ZRSAM).


Fig. 1: RIDA®QUICK Zearalenon RQS test strip.

4. Reagents provided

Each kit contains sufficient materials for max. 20 determinations. Each test kit contains:

Component	Cap color	Format	Volume
20 x Test strip	-	Ready to use	One determination per test strip (separately packed)
1 x Mobile solvent	White	Ready to use	13 mL
1 x RIDA®SMART APP cover	-	Ready to use	Evaluation cover for use with RIDA®SMART APP

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Laboratory grinder
- Balance
- Graduated cylinder
- Optional: Shaker (e.g. Rock-it 360, Trilogy® Analytical Laboratory Inc. or equivalent)
- Optional: Centrifuge ($\geq 2,000$ g)
- 100 μ L and 600 μ L micropipette
- Latest version of the RIDA®SMART APP software (Art. No. ZRSAM) activated and installed on a compatible smartphone
- Optionally: RIDA®SMART BOX (Art. No. ZRSA-SB)

5.2 Reagents

- Distilled water (dist. Water) or deionized water
- 70 % methanol (e.g. 70 mL 100 % methanol + 30 mL distilled water), exact adjustment necessary

6. Warnings and precautions for the users

The product/test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze the test strips or any test kit components.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. Do not remove the test strips from the test strip packaging until immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature. Do not allow methanol to evaporate.

The sample should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- Weigh 10 g of ground and homogenized sample into a suitable container and add 20 mL 70 % methanol.
- Close the tube and shake the sample vigorously approx. 5 min (manually or with shaker / vortex).
- Let the solution settle for 3 min or centrifuge.
- Mix 600 µL of the temperate mobile solvent with 100 µL of clear supernatant.
- Apply 100 µL of the sample solution on the application area of the test strip (see fig. 1).

Note

Please avoid methanol evaporation during sample preparation.

If necessary, the sample weight can be increased, e.g. 20 g sample in 40 mL 70 % methanol. In this case the volume of methanol has to be adapted.

9. Test procedure

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

- Remove a test strip from its foil pouch.
- Apply 100 µL of the sample solution (see chapter 8. Sample preparation) on the application area of the test strip (see Fig. 1).
- Incubate the test strip at room temperature.
- Evaluate the result after 5 min (+ max. 15 seconds) of incubation time with the RIDA®SMART APP.

10. Evaluation

The left line in the reaction area is the test band (test line, see fig. 1). Its appearance and intensity depend on the zearalenone concentration of the sample. The sample is contaminated with zearalenone, if the control band (control line) is visible and the test band (test line) is also visible.

The right line in the reaction area is a control band (control line) and must be present after each test procedure. If this line is not visible, the test was not performed correctly, or the reagents were not in order. In this case, the test should be repeated with a new test strip. If the right band is missing again, please inform your local distributor.

Please note:

- Fragmented bands:

Results will be invalid and test strips are inoperative if (testing or control) bands appear fragmented.

- Band detection:

The result field in the RIDA[®]SMART APP must always be visually checked. The correct test- and control-bands (red bands) have to be framed as shown in Fig. 2. Results are invalid if another band is framed (see, e.g. Fig. 3).


Fig. 2: Band detection valid


Fig. 3: Band detection invalid

- Reported result:

If the sample is highly contaminated with zearalenone, the control line will be only faintly visible. In this case the test result might be reported as “> 1000 ppb or invalid result“ (RIDA[®]SMART APP).

To obtain comparable results, the test strips must always be evaluated after 5 min.

10.1 Evaluation with RIDA[®]SMART APP software in combination with an approved smartphone

First, please read the instructions for use of the RIDA[®]SMART APP software (Art. No. ZRSAM) attentively. The test strip evaluation has to be performed by using the RIDA[®]SMART APP software according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA[®]SMART APP on our website www.r-biopharm.de. To analyze this test with RIDA[®]SMART APP software you must update to the latest software version. On the Certificate of analysis (CoA) you find the latest version by scanning the QR Code.

For the test strip evaluation the RIDA®SMART APP cover is needed which is included in the test kit.

RIDA®SMART APP software applications for evaluation:

Matrix	Measuring range	RIDA®SMART APP application
Corn	50 - 1000 µg/kg	50 - 1000 ppb

10.2 Evaluation with RIDA®SMART APP software in combination with the RIDA®SMART BOX

First, please read the instructions for use of the RIDA®SMART APP software and of the RIDA®SMART BOX attentively. The test strip evaluation with RIDA®SMART APP software in combination with the RIDA®SMART BOX has to be performed according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on our website www.r-biopharm.com. To analyze this test with RIDA®SMART APP software you must update to the latest software version. On the Certificate of analysis (CoA) you find the latest version by scanning the QR Code.

The RIDA®SMART software applications that need to be selected can be found in chapter 10.1.

11. Recommendation

In order to ensure a high analytical performance, we recommend:

- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- In order to ensure accurate and correct test procedure, please analyze artificially contaminated samples

12. Sensitivity

This test determines zearalenone contaminations in the measurement range of 50 - 1000 µg/kg in combination with the RIDA®SMART APP software.

13. Further application notes

Further application notes are available on request.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2011-10-24	Release version
2016-06-30	General revision
2021-12-07	General revision
2024-02-28	General revision
2025-02-10	General revision
2025-06-24	Current version Changes made: Art. No. RIDA®SMART APP changed

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

1. In conformance with the German Civil Code (“BGB”) R-Biopharm AG provides a limited warranty (“Gewährleistung”) against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm’s own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG’s liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG’s sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser’s expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements (“Gewährleistung”) implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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are registered trademarks of R-Biopharm AG.