RIDA® QUICK Zearalenon ECO

REF R5506

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RIDA®QUICK Zearalenon ECO

Brief information

RIDA®QUICK Zearalenon ECO (Art. No. R5506) is an immuno-chromatographic test for the quantitative determination of zearalenone in corn and wheat (see chapter 1. Intended use).

The test kit contains 20 test strips, each of which can be used for one analysis. All reagents required for the assay are contained in the test kit.

Results are evaluated with

- the RIDA®SMART APP software (Art. No. ZRSAM) and an approved smartphone or
- the RIDA®SMART APP software (Art. No. ZRSAM) installed on an Android smartphone and the RIDA®SMART BOX (Art. No. ZRSA-SB).

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)... approx. 10 min
test implementation (incubation time)..... 5 min

Limit of detection (LoD): 15 µg/kg (corn)
22.5 µg/kg (wheat)

Limit of quantification (LoQ): 40 µg/kg

Measuring range: 15 – 500 µg/kg (corn) or 22.5 – 500 µg/kg (wheat)
and 500 - 4000 µg/kg (ppb)
Quantitative evaluation with the RIDA® SMART APP
software (Art. No. ZRSAM)

Specificity: The RIDA®QUICK Zearalenon ECO test reacts with zearalenone in corn and wheat samples.

Further information is contained in the validation report.

The specificity of the RIDA®QUICK Zearalenon ECO test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer

system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products and accessories for zearalenone determination

RIDASCREEN® Zearalenon ECO	(Art. No. R1403)
RIDASCREEN®FAST Zearalenon	(Art. No. R5502)
RIDASCREEN®FAST Zearalenon SC	(Art. No. R5505)
RIDA®SMART APP software	(Art. No. ZRSAM)
RIDA®SMART BOX	(Art. No. ZRSA-SB)
Trilogy® certified reference material and standards for Zearalenon (ISO 17034)	

1. Intended use

RIDA®QUICK Zearalenon ECO (Art. No. R5506) is an immuno-chromatographic test in strip format for the quantitative determination of zearalenone in corn and wheat.

2. General information

The mycotoxin zearalenone is formed by fungi of the genus *Fusarium*. Zearalenone is a phytohormone which displays, apart from its anabolic properties, mainly estrogenic effects. Because of its estrogenic properties, zearalenone may induce fertility disorders in animals with clinical signs of hypoestrogenism - an aspect of a disease which although reported mainly in hogs, is described in other species such as cow, horse and sheep, too.

The potential health risk for humans induced by this mycotoxin, which is taken up with foods of vegetable or animal origin, is extensively discussed.

3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-zearalenone antibody detects zearalenone in the sample. During incubation of the test strip, a continuous band pattern (test line and control line) needs to be formed that is used to determine the concentration of zearalenone. The intensity of the test line depends on the zearalenone concentration of the sample. It increases as the zearalenone concentration increases. A weak test band can also occur with a negative sample (contaminated below the detection limit). It must be possible

to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the zearalenone concentration in the sample increases.

The test strip is evaluated using the RIDA®SMART APP software (Art. No. ZRSAM).


Fig. 1: RIDA®QUICK Zearalenon ECO test strip.

4. Reagents provided

Each kit contains sufficient materials for max. 20 determinations. Each test kit contains:

Component	Cap color	Format	Volume
20 x Test strip	-	Ready to use	One determination per test strip (separately packed)
1x Extraction solution	transparent	10x concentrate	40mL
1 x RIDA®SMART APP cover	-	Ready to use	Evaluation cover for use with RIDA®SMART APP

5. Reagents required but not provided

5.1 Equipment

- Protective equipment (Lab coat, lab gloves and lab goggles)
- Laboratory grinder
- Balance
- Graduated cylinder and beaker glass or screw top bottle
- Optionally: Shaker (e.g. Rock-it 360, Trilogy® Analytical Laboratory Inc. or equivalent)
- Centrifuge ($\geq 2,000$ g) or filter paper (e.g. Whatman No. 1 or equivalent)
- 100 μ L and 800 μ L micropipette
- Latest version of the RIDA®SMART APP software (Art. No. ZRSAM) activated and installed on a compatible smartphone
- Optionally: RIDA®SMART BOX (Art. No. ZRSA-SB)

5.2 Reagents

- Distilled water (dist. Water) or deionized water

6. Warnings and precautions for the users

The product/test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze the test strips or any test kit components.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. Do not remove the test strips from the test strip packaging until immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

The sample should be stored in a cool place, protected against light.

A **ready-to-use extraction solution** is required for the extraction. Extraction solution is available as a **10-fold concentrate** and must be diluted 1:10 (1 + 9)

with deionized or distilled water before use (e.g. 100 mL concentrate + 900 mL dist. water).

The diluted extraction solution has a shelf life of one week at 2 – 8 °C (36- 46 °F). If turbidity occurs in the diluted extraction buffer (e.g. caused by contamination), it must be discarded. A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- Weigh 5 g of ground and homogenized sample into a suitable container and add 15 mL of ready-to-use extraction solution.
- Close the tube and shake the sample vigorously approx. 5 min (manually or with shaker / vortex).
- Allow mixture to sediment (≥ 2 min) until 2 mL of supernatant can be taken
- Use 1 – 2 mL of supernatant and centrifuge (1 min at 2000g) or filtrate.

Extraction stability

Corn samples are stable for 24 hours at 23 °C (73°F) after extraction.

Wheat samples are stable for 30 minutes at 23 °C (73°F) after extraction.

8.1 Measuring range 15 – 500 µg/kg (corn) /22.5 – 500 µg/kg (wheat)

- Apply 100 µL of the particle- free supernatant on the application area of the test strip (see fig. 1).

8.2 Measuring range 500 – 4000 µg/kg

For the measuring range of 500 – 4000 µg/kg the particle- free supernatant from sample preparation 8. must be further diluted.

- Dilute 100 µL of the particle- free supernatant in 800 µL of the ready-to-use extraction solution
- Use 100 µL of the dilution in the test.

Note

If necessary, the sample weight can be increased . In this case the volume of extraction buffer has to be adapted. e.g. 10 g sample in 30 mL ready-to-use extraction solution.

9. Test procedure

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

- Remove a test strip from its foil pouch.
- Apply 100 µL of the sample solution (see chapter 8. Sample preparation) on the application area of the test strip (see Fig. 1).
- Incubate the test strip at room temperature.
- Evaluate the result after 5 min (+ max. 15 seconds) of incubation time with the RIDA®SMART APP.

10. Evaluation

The left line in the reaction area is the test band (test line, see fig. 1). Its appearance and intensity depend on the zearalenone concentration of the sample. The sample is contaminated with zearalenone, if the control band (control line) is visible and the test band (test line) is also visible.

The right line in the reaction area is a control band (control line) and must be present after each test procedure. If this line is not visible, the test was not performed correctly, or the reagents were not in order. In this case, the test should be repeated with a new test strip. If the right band is missing again, please inform your local distributor.

Please note:

- Fragmented bands:


Results will be invalid and test strips are inoperative if (testing or control) bands appear fragmented.

- Band detection:

The result field in the RIDA®SMART APP must always be visually checked. The correct test- and control-bands (red bands) have to be framed as shown in Fig. 2. Results are invalid if another band is framed (see, e.g. Fig. 3).

■ Test line ■ Control line

■ Test line ■ Control line

Fig. 2: Band detection valid

Fig. 3: Band detection invalid

- Reported result:

If the sample is highly contaminated with zearalenone, the control line will be only faintly visible. In this case the test result might be reported as “> 500 ppb or invalid result” (RIDA®SMART APP). In this instance, dilute the sample as described. To obtain comparable results, the test strips must always be evaluated after 5 min.

10.1 Evaluation with RIDA®SMART APP software in combination with an approved smartphone

First, please read the instructions for use of the RIDA®SMART APP software (Art. No. ZRSAM) attentively. The test strip evaluation has to be performed by using the RIDA®SMART APP software according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on our website www.r-biopharm.de. To analyze this test with RIDA®SMART APP software you must update to the latest software version. On the Certificate of analysis (CoA) you find the latest version by scanning the QR Code.

For the test strip evaluation the RIDA®SMART APP cover is needed which is included in the test kit.

RIDA®SMART APP software applications for evaluation:

Matrix	Measuring range	RIDA®SMART APP application
corn	15 - 500 µg/kg	Corn 15 - 500 ppb
	500 - 4000 µg/kg	Corn 500 - 4000 ppb
wheat	22.5 - 500 µg/kg	Wheat 22.5 - 500 ppb
	500 - 4000 µg/kg	Wheat 500 - 4000 ppb

10.2 Evaluation with RIDA®SMART APP software in combination with the RIDA®SMART BOX

First, please read the instructions for use of the RIDA®SMART APP software and of the RIDA®SMART BOX attentively. The test strip evaluation with RIDA®SMART APP software in combination with the RIDA®SMART BOX has to be performed according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on our website www.r-biopharm.com. To analyze this test with RIDA®SMART APP software you must update to the latest software version. On the Certificate of analysis (CoA) you find the latest version by scanning the QR Code.

The RIDA®SMART software applications that need to be selected can be found in chapter 10.1.

11. Recommendation

In order to ensure a high analytical performance, we recommend:

- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- In order to ensure accurate and correct test procedure, please analyze artificially contaminated samples

12. Sensitivity

This test determines zearalenone contaminations in the measurement range of 15 - 4000 µg/kg (corn) or 22.5 - 4000 µg/kg (wheat) in combination with the RIDA®SMART APP software.

13. Further application notes

Further application notes are available on request.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2026-07-27	Release version

Explanation of symbols

General symbols:

-  Follow the instructions for use
-  Batch number
-  Expiry date (YYYY-MM-DD)
-  Storage temperature
-  Article number
-  Number of test determinations
-  Manufacturing date (YYYY-MM-DD)
-  Manufacturer + address

Disclaimer

1. R-Biopharm AG provides a limited statutory warranty ("Gewährleistung") as mandated and pursuant to German law against defects in title and material. For products with a defined shelf life, the warranty applies for the stated shelf life. For use cycle limited products, it applies until the stated cycle limit is reached, but not beyond the general warranty period. For all other products, the general warranty period is twelve (12) months from the date the risk transfers.
Such warranty is subject to timely notice of the alleged defect, giving sufficient detail to assess the claim. A product's defectiveness is assessed based on the specifically agreed quality, its suitability for the contractually stated purpose, and whether all agreed-upon accessories, documentation, and instructions were provided.
2. R-Biopharm AG is not liable for issues arising from:
 - a) The failure to read, understand, or follow product instructions and information for use or safety
 - b) The use of unqualified personnel
 - c) The failure to follow industry standards such as Good Laboratory Practices
 - d) The use of unsuitable controls, sample matrices, or procedures
 - e) Other misuse
 - f) unilateral product modifications
 - g) Improper storage by others
 - h) External chemical, electromagnetic, mechanical, or electrolytic influences that exceed the defined parameters.
 - i) Harm caused by external events beyond R-Biopharm's control, including but not limited to burglary, theft, lightning, fire, water, or acts of nature.
3. R-Biopharm AG's liability is limited to foreseeable, contract-typical losses in case of negligent breaches of obligations that are essential to achieving the agreement's purpose and that the other party could reasonably rely on (material contract obligations). R-Biopharm AG is not liable for losses resulting from non-material contract obligations.
4. ALL OTHER WARRANTIES OR GUARANTEES OF ANY KIND, EXPRESS OR IMPLIED, ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, COURSE OF DEALINGS BETWEEN THE PARTIES OR OTHERWISE. R-Biopharm AG assumes no liability for consequential damages, in particular loss of profit, interruptions, or other indirect damages.
5. The foregoing is not intended to limit liability to the extent such liability cannot be excluded or limited under applicable German law or convention, including cases of intent, willful misconduct, or gross negligence, personal injury or death, express warranties or assumptions of procurement risks, strict liability, or other non-waivable legal mandates.
6. The foregoing liability provisions apply equally to R-Biopharm AG, its directors, officers, employees, and agents.

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are registered trademarks of R-Biopharm AG.