RIDASCREEN®FAST Fumonisin ECO

REF R5603

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RIDASCREEN®FAST Fumonisin ECO

Brief information

RIDASCREEN®FAST Fumonisin ECO (Art. No. R5603) is a competitive enzyme immunoassay for the quantitative analysis of fumonisin in corn and feed (see chapter 1. Intended use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for maximum of 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: Extraction, centrifugation/filtration, dilution

Time requirement: Sample preparation (for 10 samples)
Corn and feed approx. 10 min
Test implementation (incubation time)...approx. 8 min

Limit of detection: Corn < 0.25 mg/kg (ppm)
(corresponding to the Standard substance) Feed..... < 0.25 mg/kg (ppm)

Recovery rate: Corn approx. 103 %
(naturally contaminated samples - Trilogy®) Feed approx. 113 %

Specificity: Fumonisin B1 approx. 100 %
Fumonisin B2 approx. 62 %
Fumonisin B3 approx. 70 %

The specificity of the RIDASCREEN®FAST Fumonisin ECO test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice manual. It lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA analyses with them. The

manual can be retrieved, printed and downloaded from www.r-biopharm.com/products/food-feed-analysis.

Related product and accessories for fumonisin determination

RIDASCREEN® Fumonisin ECO (Art. No. R3411)

RIDA® QUICK Fumonisin RQS ECO (Art. No. R5606)

TRILOGY® Liquid Standard Fumonisin B1, B2 (Art. No. TAS-MM18LZ1-2)

TRILOGY® Liquid Standard Fumonisin B1 (Art. No. TAS-M18LZ1-2)

TRILOGY® Liquid Standard Fumonisin B2 (Art. No. TAS-M19LZ1-2)

TRILOGY® Dried Standard Fumonisin B1, B2 (Art. No. TAS-MM18DZ1-2)

1. Intended use

RIDASCREEN®FAST Fumonisin ECO is a competitive enzyme immunoassay for the quantitative determination of fumonisin in corn and feed.

2. General information

Fumonisin are carcinogenic, neuro-, hepato- and pneumotoxic metabolites of *Fusarium moniliforme*, a mould fungi, which grows host-specific on corn. The toxic concentrations of fumonisin differ significantly depending on the animal species. A concentration of approx. 5 - 10 mg/kg fumonisin in feed induces neuro-toxic effects in horses. In pigs, the ingestion of 4 - 16 mg/kg body weight may result in liver cirrhosis and more than 16 mg/kg bw. may lead to pulmonary edema. Chickens tolerate higher concentrations of fumonisin in feed, up to 75 mg/kg. Cattle seem to be insensitive to high fumonisin concentrations.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-fumonisin antibodies.

Fumonisin standards or sample solutions, fumonisin enzyme conjugate and anti-fumonisin antibodies are added. Free fumonisin and fumonisin enzyme conjugate compete for the fumonisin antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-fumonisin antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue

product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the fumonisin concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for maximum of 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate M	-	Ready to use		48 wells
ECO extractor	Transparent	Concentrate	10x	1x 120 mL
Dilution buffer	Transparent	Ready to use		60 mL
Standard 1*	White	Ready to use	0 mg/L	1.3 mL
Standard 2*	White	Ready to use	0.25 mg/L	1.3 mL
Standard 3*	White	Ready to use	0.50 mg/L	1.3 mL
Standard 4*	White	Ready to use	1 mg/L	1.3 mL
Standard 5*	White	Ready to use	3 mg/L	1.3 mL
Standard 6*	White	Ready to use	9 mg/L	1.3 mL
Wash buffer salt Tween		Dissolve the salt		
Conjugate	Red	Ready to use		3 mL
Antibody	Black	Ready to use		3 mL
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 mL
Stop Solution	Yellow	Ready to use		14 mL

*) The dilution factor 50 for the sample has already been considered. Therefore, the fumonisin concentrations (mg/kg or mg/L) of samples can be read directly from the standard curve.

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Grinder (mill), laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Graduated cylinder: (plastic or glass) 250 mL
- Centrifuge (at least 3,500 x g) + centrifugal vials with cap (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Alternative option: Filter paper: Whatman No. 1 or equivalent
- (Horizontal) shaker

- Vortex mixer
- Variable 20 - 200 µL and 200 - 1000 µL micropipettes
- Multi stepper or 8-channel pipette for 50, 100 and 300 µL
- Microtiter plate spectrophotometer (450 nm)

5.2 Reagents

- Distilled water (dist. water) or deionized water

6. Warnings and precautions for the users

This test should be carried out only by trained laboratory employees. The instruction for use must be strictly followed.

The standards contain fumonisin, particular care should be taken. Avoid contact of the reagent with the skin (use gloves).

Decontamination of the glassware and fumonisin solutions is best carried out using a sodium hypochlorite solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate (see chapter 10.2.)). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (36 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers. Exceptions are checked and explicitly listed in the instructions for use.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to addition in the wells
- Extinction less than 0.8 ($E_{450\text{ nm}} < 0.8$) for zero standard

9. Sample preparation

The samples should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure (recommended particle size: 500 μm).

9.1. Extraction buffer

For extraction, the **ready-to-use ECO extractor** is needed. To obtain the ready-to-use ECO extractor, dilute the ECO Extractor (10x concentrate) 1:10 (1 + 9) with distilled or deionized water. The ready-to-use ECO Extractor is stable for one week at 2 - 8 °C (36 - 46 °F). If turbidity occurs in the ready-to-use ECO extractor (possibly caused by contamination), discard or do not use it. The ECO extractor is interchangeable between RIDASCREEN® Mycotoxin ELISA lots and products.

9.2. Extraction of corn and feed

Bring all reagents and samples to room temperature (20 - 25 °C / 68 - 77 °F) before use and perform the sample preparation at room temperature.

- Weigh 5 g of ground and homogenized sample into a suitable container (e.g. 50 mL tube) and add 25 mL of ready-to-use ECO extractor. (Alternative possibility: Weigh 10 g of the ground and homogenized sample and 50 mL of diluted ECO Extractor.)
- Vortex the sample until completely soaked.

- Shake the sample vigorously for 5 minutes (manually or with shaker at 420 rpm).
- Centrifuge for 5 min at 3,500 g and room temperature (20 - 25 °C / 36 - 46 °F) or filtrate.
- Dilute supernatant 1:10 (1 + 9) with dilution buffer (e.g. 100 µL supernatant + 900 µL dilution buffer).
- Add 50 µL of the diluted supernatant per well in the assay.

Note: Samples measured above of standard 6 (> 9 mg/kg (ppm)) are recommended to be further diluted 1:10 to at least 1:64 with dilution buffer. The resulting dilution factor must be taken into account in the calculation.

Remark: The extracted analyte is stable for at least 6 h at room temperature before and after dilution in dilution buffer.

10. Test implementation

10.1 Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use and perform the test at room temperature.

As **wash buffer** a PBS tween buffer is needed. Please use the wash buffer salt contained in the kit (see chapter 4.). Dissolve the entire buffer salt in one liter of distilled water. The ready to use washing buffer expires after approx. 4 weeks at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the envelope in 100 mL of distilled water to obtain a 10-fold concentrated washing buffer. This 10-fold concentrate expires after approx. 8 weeks when stored at room temperature (20 - 25 °C / 68 - 77 °F). Use 1 part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use wash buffer.

Components should be stored immediately at 2 - 8 °C (36 - 46 °F) when no longer required.

10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow wells to dry between work steps. Accurate washing is very important. Do not allow microwells to dry up totally and avoid prolonged intervals between the working steps. Reproducibility in

any ELISA is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the ELISA test procedure.

It is recommended to pipette the conjugate, the antibody, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipet 50 μ L of standard or prepared sample into separate wells; use a new pipette tip for each standard or sample.
3. Add 50 μ L of conjugate to each well.
4. Add 50 μ L of antibody to each well. Mix gently by shaking the plate manually and incubate for 5 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Wash each well with 250 μ L wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat the washing step two more times (a total of three wash cycles).
6. Add 100 μ L of substrate/chromogen to each well, mix gently by shaking the plate manually and incubate for 3 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Pipette 100 μ L of stop solution into each well. Mix gently by shaking the plate manually and measure the extinction at 450 nm. Read within 10 min after addition of stop solution.

11. Evaluation

Special software, **RIDASOFT® Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using Cubic Spline function.

For the evaluation it should be clarified, that all quality criteria are fulfilled for the current test run. The course of the standard curve is shown in the Quality Assurance Certificate (certificate of analysis) enclosed in the test kit.

Please note that the following method must be used for the evaluation: RidawinFF.NET \ FOOD \ Mycotoxins FAST \ R5603 FAST Fumonisin ECO.met.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates semilogarithmic against the fumonisin concentration [mg/kg]. The fumonisin concentration in mg/kg corresponding to the extinction of each sample can be read from the calibration curve.

12. Result interpretation

Results between LoD and LoQ indicate a low mycotoxin concentration in the sample. Calculated result show a high uncertainty in this area due to the method's high variation below LoQ. Therefore, such results should not be reported with a quantitative value, but qualitative as "< LoQ".

A result below the LoD does not exclude a mycotoxin contamination below the detection limit of the assay. The result should be reported accordingly.

A further dilution and new detection of samples is recommended for absorbance values ($A_{450 \text{ nm}}$) > standard 6. In case of a further dilution, the additional dilution factor must be taken into account when calculating the mycotoxin concentration.

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure, and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max ± 1 %.

14. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis.

This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- Pre-flush pipette tips prior to pipetting.
- Carry along test controls for quality control. Mycotoxin-free and mycotoxin containing samples should be used.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2022-10-14	Release version

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.