r-biopharm®

RIDASCREEN®FAST DON SC

REF R5905

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RIDASCREEN®FAST DON SC

Brief information

RIDASCREEN®FAST DON SC (Art. No. R5905) is a competitive enzyme immunoassay for the quantitative determination of deoxynivalenol (DON) in cereals, malt and feed (see chapter 1. Intended use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: Extraction and filtration

Time requirement: Sample preparation (for 10 samples) .. approx. 10 min
 Test implementation (incubation time)..... 8 min

Limit of detection: 0.074 mg/kg (ppm)
(depending on matrix)

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

Related product and accessories for DON determination

RIDASCREEN® DON (R5906)

RIDASCREEN®FAST DON (R5901/R5902)

RIDA®QUICK DON RQS ECO (R5911)

TRILOGY® dried standards, liquid standards, certified reference material and QC material available (see product catalogue)

1. Intended use

RIDASCREEN®FAST DON SC is a competitive enzyme immunoassay for the quantitative determination of DON in cereals, malt and feed.

2. General information

DON belongs to the trichothecene group of mycotoxins and is formed by fungi of the genus *Fusarium*. DON often occurs in plant products particularly in cereals. Of the trichothecene mycotoxins DON, 3-acetyl- and 15-acetyl-deoxynivalenol are the most frequently occurring toxins in Europe and Northern America. The toxin concentrations found in wheat, corn or rice are often in the ppm range. Due to their high cytotoxic and immunosuppressive properties these toxins pose a risk to human and animal health.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-DON antibodies. DON standards or sample solutions, DON enzyme conjugate and anti-DON antibodies are added. Free DON and DON enzyme conjugate compete for the DON antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-DON antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the DON concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume
Microtiter plate M	-	Ready to use	48 wells
Standard 1*	white	Ready to use	0 mg/L 1.3 mL
Wash buffer salt Tween	-	Dissolve the salt	-
Conjugate	red	Ready to use	3 mL
Antibody	black	Ready to use	3 mL
Substrate/Chromogen Red Chromogen Pro	brown	Ready to use	10 mL
Stop solution	yellow	Ready to use	14 mL

* Only standard 1 (zero standard, 0 mg/L) is included in the test kit. The standard curve (B/B_0) is provided with the certificate of analysis (CoA) of the test kit. For the calculation of results see chapter 11. Evaluation.

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Grinder (mill), laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Graduated cylinder: (plastic or glass) 100 mL
- Filter paper: Whatman No. 1 or equivalent
- (Horizontal) shaker
- Variable 20 - 200 μ L and 200 - 1000 μ L micropipettes
- Multi stepper or 8-channel pipette for 50, 100 and 250 μ L
- Microtiter plate spectrophotometer (450 nm)

5.2 Reagents

- Distilled water (dist. water) or deionized water

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

Decontamination of the glassware and toxin solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate, see chapter 10.2 Test procedure). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for standard 1 (zero standard)

9. Sample preparation

The samples should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

- Weigh 5 g of ground sample into a suitable container and add 100 mL of distilled water*
- Blend the sample by ultra-turrax (or equivalent) for two minutes or shake vigorously for three minutes (manually or with shaker)
- Filter the extract through Whatman No. 1 filter paper (or equivalent)
- Use 50 µL of the filtrate per well in the test

* Sample size may be increased if required, but the volume of water must be adapted accordingly, e.g. 25 g in 500 mL of distilled water or 50 g in 1000 mL of distilled water.

10. Test procedure

10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The specific reaction starts with the addition of the specific antibody. Not more than three strips should be applied to the test, when a single step pipette is used. More strips (up to 6) can be applied, when a multistep pipette is used.

Components should be stored immediately at 2 - 8 °C (35 - 46 °F) when no longer required.

The **DON standard 1** (zero standard, 0 mg/L) is provided ready to use. B/B₀-values of DON standards 2 to 8 (0.074, 0.222, 0.666, 1.333, 2, 4, and 6 mg/L) are reported on the certificate of analysis (CoA) of the test. The standard curve is calculated with the RIDASOFT® Win.NET (see chapter 11. Evaluation) according to those values. The dilution factor 20 for the sample has been considered in this calculation.

As **wash buffer** a PBS tween buffer is needed. Please use the wash buffer salt contained in the kit (see chapter 4. Reagents provided). Dissolve the entire buffer salt in one liter of distilled water. The ready to use washing buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the pouch in 100 mL of distilled water to obtain a 10fold concentrated washing buffer. Use 1 part of this concentrate and dilute with 9 parts of distilled water to obtain the ready to use washing buffer. The 10fold concentrate expires after approx. 8 - 12 weeks, store at room temperature (20 - 25 °C / 68 - 77 °F).

10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

It is recommended to pipette the conjugate, the antibody, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for the standard and samples to be run. Record standard and sample positions.
2. Pipet 50 µL of standard or prepared sample into separate wells. Use a new pipette tip for the standard or each sample.
3. Add 50 µL of conjugate to the bottom of each well.
4. Add 50 µL of antibody to each well. Mix gently by shaking the plate manually and incubate for 5 min (± 1) at room temperature (20 - 25 °C / 68 - 77 °F).
5. Dump the liquid out of the wells into a sink. Tap the microwell holder upside down onto a clean filter towel (three times in a row) to remove all remaining liquid from the wells. Using a wash bottle or multichannel pipette, fill the wells (250 µL per well) with wash buffer (see chapter 10.1. Test preparation). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.
6. Add 100 µL of substrate/chromogen to each well. Mix gently by shaking the plate manually and incubate for 3 min (± 0.5) at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µL of stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Evaluation

Special software, **RIDASOFT® Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using logit-log function.

For the evaluation it should be clarified, that all quality criteria are fulfilled for the current test run. The course of the standard curve is shown in the certificate of analysis (CoA) enclosed in the test kit.

For calculation of results you need to transfer the measured absorbance for standard 1 (0 ppm) as well as the B/B₀ values for standards 2 to 8 (0.074, 0.222, 0.666, 1.333, 2, 4, and 6 ppm) provided with the kit certificate, into the RIDASOFT® Win.NET. From this the software program calculates the standard curve and the content of DON in the samples.

12. Result interpretation

Results between LOD and LOQ indicate a low mycotoxin concentration in the sample. Calculated result show a high uncertainty in this area due to the method's high variation below LOQ. Therefore, such results should not be reported with a quantitative value, but qualitative as "< LOQ".

A result below the LOD does not exclude a mycotoxin contamination below the detection limit of the assay. The result should be reported accordingly.

A further dilution and new detection of samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 8 (see chapter 11. Evaluation). In case of a further dilution, the additional dilution factor must be taken into account when calculating the mycotoxin concentration.

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result.

14. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that

this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- Pre-flush pipette tips prior to pipetting.
- Carry along test controls for quality control. Mycotoxin-free and mycotoxin containing samples should be used.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2016-09-28	Release version
2023-09-27	Current version General revision, changes made: <ul style="list-style-type: none">– Removal of FGIS (2014-052) certificate– Changes in chapter 5 + 6 + 7 + 11– New chapter 12 + 13 + 14 + version overview + explanation of symbols + new disclaimer

Explanation of symbols

General symbols:

-  Follow the instructions for use
-  Batch number
-  Expiry date (YYYY-MM)
-  Storage temperature
-  Article number
-  Number of test determinations
-  Manufacturing date (YYYY-MM)
-  Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.