r-biopharm®

RIDA® QUICK DON RQS ECO

REF R5911

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RIDA[®]QUICK DON RQS ECO

Brief information

RIDA[®]QUICK DON RQS ECO (Art. No. R5911) is an immunochromatographic test for the quantitative determination of deoxynivalenol (DON) in grain (wheat, corn, barley, oats) (see chapter 1. Intended use).

The test kit contains 20 test strips, each of which can be used for one analysis. All reagents required for the assay are contained in the test kit. Results are evaluated with

- the RIDA[®]SMART APP software (ZRSAM1000) and an approved smartphone or
- the RIDA[®]SMART APP software (ZRSAM1000) installed on an Android smartphone and the RIDA[®]SMART BOX (ZRSA-SB) or
- the RIDA[®]QUICK SCAN reader.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) ... approx. 10 min
test implementation (incubation time) 3 min

Limit of detection: 0.15 mg/kg (ppm)

Limit of quantification: 0.25 mg/kg

Measuring range: 0.25 - 7.5 mg/kg und 7.5 - 50 mg/kg
Quantitative evaluation with the RIDA[®]SMART APP

Specificity: The RIDA[®]QUICK DON RQS ECO test reacts with DON in wheat, corn, barley and oats samples.

Further information is contained in the validation report.

The specificity of the RIDA[®]QUICK DON RQS ECO test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances,

the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products and accessories for DON determination

RIDA®QUICK Mycotoxin ECO Extractor	(Art. No. R5000)
RIDASCREEN® DON	(Art. No. R5906)
RIDASCREEN®FAST DON	(Art. No. R5901/R5902)
RIDASCREEN®FAST DON SC	(Art. No. R5905)
RIDA®SMART APP Software	(Art. No. ZRSAM1000)
RIDA®SMART BOX	(Art. No. ZRSA-SB)
Trilogy® certified reference material and standards for DON (ISO 17034)	

1. Intended use

RIDA®QUICK DON RQS ECO (Art. No. R5911) is a quantitative immunochromatographic test in strip format for the quantitative determination of DON in grains (wheat, corn, barley, oats). The DON extraction is carried out with water.

2. General information

DON is a trichothecene mycotoxin. It is produced by fungi in the *Fusarium* genus. DON can be found in plant products and grains in particular. DON is the most prevalent toxin of the more than 100 known trichothecenes in Europe and North America.

3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-DON antibody detects DON in the sample. During incubation of the test strip, a band pattern (test line / control line) forms that is used to determine the concentration of DON. The intensity of the test line depends on the DON concentration of the sample. It increases as the DON concentration increases. It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the DON concentration in the sample increases. The test strip is evaluated using the RIDA®SMART APP software (Art. No. ZRSAM1000) or the RIDA®QUICK SCAN reader.


Fig. 1: RIDA®QUICK DON RQS ECO test strip.

4. Reagents provided

Each kit contains sufficient materials for max. 20 determinations. Each test kit contains:

Component	Cap color	Format	Volume
20 x Test strip	-	Ready to use	One determination per test strip (separately packed)
1 x Mobile solvent	White	Ready to use	25 mL
1 x RIDA®SMART APP cover	-	Ready to use	Cover for use with RIDA®SMART APP

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Laboratory grinder
- Weighing scale
- Graduated cylinder
- Optional: Shaker
- Optional: Centrifuge (≥ 2000 g)
- Optional: Filter paper (e.g. Whatman No. 1 or equivalent)
- 100 - 1000 μ L pipette
- Latest version of the RIDA®SMART APP software (Art. No. ZRSAM1000) activated and installed on a compatible smartphone
- Optional: RIDA®SMART BOX (Art. No. ZRSA-SB)

5.2 Reagents

- Distilled water (dist. Water) or deionized water

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze the test strips or any test kit components.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. Do not remove the test strips from the test strip packaging until immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

The sample should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

8.1. Measurement range 0.25 - 7.5 mg/kg

- Weigh 5 g of the ground and homogenized sample into a suitable container (sealable tube) and add 25 mL of distilled or deionized water
- Close the tube and shake the sample vigorously for 30 s (manually or with shaker / vortex)
- Let the solution settle (≥ 3 min), filtrate or centrifuge (e.g. 1 min, 2000 g) to get a particle-free supernatant

- Dilution I: Dilute 100 µL of the particle-free supernatant with 500 µL of the mobile solvent and mix thoroughly (e.g. invert tube five times)
- Use 100 µL of Dilution I in the test

8.2. Measurement range 7.5 - 50 mg/kg

- Dilution II: Dilute 100 µL of Dilution I (see chapter 8.1) with 700 µL of the mobile solvent and thoroughly mix (e.g. invert tube five times)
- Use 100 µL of Dilution II in the test

Note:

If necessary, the sample weight can be increased, e.g. 10 g sample in 50 mL dist. Water. In this case the volume of water has to be adapted.

9. Test procedure

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

- Remove a test strip from its foil pouch.
- Apply 100 µL of the sample solution (see chapter 8. Sample preparation) on the application area of the test strip (see Fig. 1).
- Incubate the test strip at room temperature.
- Evaluate the result after 3 min incubation time with the RIDA®SMART APP software or RIDA®QUICK SCAN reader.

10. Evaluation

The left line in the reaction area is the test band (test line, see fig. 1). Its appearance and intensity depends on the DON concentration of the sample. The sample is contaminated with DON, if the control band (control line) is visible and the test band (test line) is also visible.

The right line in the reaction area is a control band (control line) and must be present after each test procedure. If this line is not visible, the test was not performed correctly or the reagents were not in order. In this case, the test should be repeated with a new test strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample is highly contaminated with DON, the control line will be only faintly visible. In this case the test result might be reported as „>7,5

ppm or invalid result“ (RIDA®SMART APP) or “>7,5 ppm or invalid“ (RIDA®QUICK SCAN reader). In this instance, dilute the sample as described in chapter 8.2 and repeat the test with a new test strip.

To obtain comparable results, the test strips must always be evaluated after 3 min.

10.1 Evaluation with RIDA®SMART APP

First, please read the instructions for use of the RIDA®SMART APP software (Art. No. ZRSAM1000) attentively. The test strip evaluation has to be performed by using the RIDA®SMART APP software according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on our website www.r-biopharm.de.

For the test strip evaluation the RIDA®SMART APP cover is needed which is included in the test kit-

RIDA®SMART APP software applications for evaluation:

Matrix	Measuring range	RIDA®SMART APP application
Wheat	0.25 - 7.5 mg/kg	Wheat 0,25 - 7,5 ppm
	7.5 - 50 mg/kg	Wheat 7,5 - 50 ppm
Corn	0.25 - 7.5 mg/kg	Corn 0,25 - 7,5 ppm
	7.5 - 50 mg/kg	Corn 7,5 - 50 ppm
Barley/oats	0.25 - 7.5 mg/kg	Barley_ Oats 0,25 - 7,5 ppm
	7.5 - 50 mg/kg	Barley_ Oats 7,5 - 50 ppm

10.2 Evaluation with RIDA®QUICK SCAN

First, carefully read the user guide for the RIDA®QUICK SCAN and proceed as instructed. Chapter 3. contains the description for evaluating the test strips.

To obtain optimal results every test kit contains lot-specific standard curve parameters. These parameters can be found on the Certificate of analysis (CoA) and have to be used only with this specific lot. The CoA contains parameters for several standard curves. Depending on the extraction method used, select the corresponding parameters.

To enter the parameters manually into the RIDA®QUICK SCAN reader use the plain text information on the CoA. To enter the data with the barcode scanner use the 2D barcodes. In both cases, the checksum shown on the RIDA®QUICK SCAN screen needs to be identical to the respective checksum on the certificate (number below the 2D barcode).

RIDA®QUICK SCAN evaluation methods:

Measuring range	RIDA®QUICK SCAN method
0.25 - 7.5 mg/kg	DON Low
7.5 - 50 mg/kg	DON High

11. Sensitivity

This test determines DON contaminations in the measurement range of 0.25 - 50 mg/kg (0.25 - 7.5 mg/kg and 7.5 - 50 mg/kg) in combination with the RIDA®SMART APP software with the RIDA®QUICK SCAN reader.

12. Further application notes

Further application notes are available on request. Please contact your local distributor.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2017-11-27	Release version
2019-02-27	General revision
2021-10-14	General revision
2024-02-20	Current version General revision Changes made: <ul style="list-style-type: none">– 1. Intended use: products added– 4. Reagents provided: Smart App cover simplified– 6. Warnings and precautions for the users completed– Imprint updated

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.