RIDASCREEN® Egg

REF R6411

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RIDASCREEN® Egg

Brief information

RIDASCREEN® Egg (Art. No. R6411) is a sandwich enzyme immunoassay for the quantitative analysis of native and processed hen's egg in food validated for the method (see chapter 1).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization and extraction
Time requirement:	sample preparation with AEP (for 10 samples)..... approx. 20 min sample preparation with A-AEP for processed samples (for 10 samples)..... approx. 30 min test implementation (incubation time) 50 min
Standard material:	NIST reference material 8445 (whole egg powder)
Limit of detection: (matrix-depended)	0.13 mg/kg (ppm) whole egg powder* 0.04 - 0.20 mg/kg (ppm) whole egg powder *mean value
Limit of quantification:	0.25 mg/kg (ppm) whole egg powder corresponds to 0.12 mg/kg whole egg protein
Specificity:	The used antibodies specifically detect the egg white proteins ovalbumin and ovomucoid from hen's egg. Cross-reactivity is given to egg from geese, quails and ostrich. Further information are contained in the validation report.

The cross-reactivities of the antibodies used were determined for pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross-reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments (see chapter 13).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

Related products for the detection of egg

RIDASCREEN®FAST Egg (Art. No. R6402)

RIDASCREEN®FAST Lysozym (Art. No. R6452)

bioavid Lateral Flow Egg incl. Hook Line (Art. No. BLH708-15)

1. Intended use

RIDASCREEN® Egg (Art. No. R6411) is a sandwich enzyme immunoassay for the quantitative analysis of native and processed hen's egg in foods. Due to the large number of different food products, the following samples were examined representatively for different product groups within the scope of the test development: ice cream, noodles, salad dressing, wine, chocolate, and cookies.

The test kit should generally be able to provide usable results for other foods within the methodological limits (see chapter 13). In individual cases, the user must check the suitability of the test kit by carrying out appropriate experiments before use.

For detailed results and further information on validation data with other food matrices, please refer to the validation report. Further applications are regularly validated in our laboratories, which we make available in our Application Notes.

2. General information

Egg allergy is very common in children, although it often disappears with age. Egg protein can be present as an ingredient or as contamination in raw or heated foods.

Egg white contains approx. 10 % - 11 % protein. Allergologically of importance are four main allergens that make up 80 % of egg white protein content. The main allergens include ovomucoid (11 %), ovalbumin (54 %), ovotransferrin (12 %), and lysozyme (3.5 %). In contrast, the proteins in the egg yolk have only moderate allergenicity. According to the **regulation (EU) No. 1169/2011**, egg and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia, China and New Zealand.

3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with antibodies against ovalbumin and ovomucoid. By adding the standard or sample solution to the wells, egg protein (ovalbumin and ovomucoid) present in the sample will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Following the washing step, a solution containing antibody conjugated to peroxidase is added. This conjugate is bound to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue end product. A stop solution is added which results in a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance of the solution which is proportional to the ovalbumin and ovomucoid concentration in the sample is measured photometrically at 450 nm and expressed as mg/kg whole egg powder.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Lid color	Format	Concentration	Contents
Microtiter plate		Ready for use		96 wells
Allergen extraction buffer	Green	Concentrate	10x	100 mL
Extractor 4	Blue	Ready to use		100 mL
Additive 1	Blue			2 g
Additive ECO	Green			2.5 g
Standard 1	Transparent	Ready to use	0 mg/kg*	1.3 mL
Standard 2	Transparent	Ready to use	0.25 mg/kg*	1.3 mL
Standard 3	Transparent	Ready to use	0.5 mg/kg*	1.3 mL
Standard 4	Transparent	Ready to use	1 mg/kg*	1.3 mL
Standard 5	Transparent	Ready to use	2 mg/kg*	1.3 mL
Wash buffer	Brown	Concentrate	10x	100 mL
Conjugate	Red	Ready to use		11 mL
Substrate/Chromogen	Brown	Ready to use		13 mL
Stop solution	Yellow	Ready to use		14 mL

* The concentration values of the standards already take into account the dilution factor 20, which results from the sample preparation according to chapter 9.1 and 9.2. Thus, the whole egg powder concentration of the samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1 Equipment:

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- pH meter
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (60 °C / 140 °F and 100 °C / 212 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12 μ m)
- Graduated pipettes
- Graduated cylinder
- Variable 20 - 200 μ L and 200 - 1000 μ L micropipettes

- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100 µL
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)

5.2. Reagents:

- Distilled water (dist. water) or deionized water
- 5 M hydrochloric acid (HCl)
- 1 M sodium hydroxide (NaOH)
- Skim milk powder (SMP) (food quality; egg-free)

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the plate (coated plate and pre-plate, if necessary, see chapter 10.2). Use separate pipet tips for each standard and each sample extract to avoid cross-contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the reagents at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation.
- Value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 1.2$) for standard 5.

9. Preparation of samples

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation.
- Carry out the sample preparation in a room isolated from the ELISA procedure.

The samples should be stored in a cool place, protected against light.

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The Allergen extraction buffer is provided as a 10fold concentrate and must be diluted prior use. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that, dilute the heated buffer concentrate 1:10 (1+9) with dist. water (e.g. 900 mL dist. water + 100 mL buffer concentrate). The **diluted Allergen extraction buffer (AEB)** is either stable at 20 - 25 °C (68 - 77 °F) for approx. 4 weeks or at 2 - 8 °C (35 - 46 °F) for 12 weeks.

For extraction according chapter 9.1, as an alternative to weighing the skim milk powder (SMP) to the sample, the SMP can be added to the AEB (**SMP-AEB**). Add 2.5 g of skim milk powder per 100 mL to the diluted AEB. Only

prepare as much SMP-AEB as needed. The shelf life of the SMP-AEB is one day.

Extractor 4 is ready to use and solid at 2 - 8 °C (35 - 46 °F). Before use, the Extractor 4 must be completely dissolved at room temperature (20 - 25 °C / 68 - 77 °F) or by heating it in a water bath at 37 °C (99 °F).

The **Extractor Egg** must be made fresh on the day of the extraction. Mix 0.5 g Additive ECO with 25 mL Extractor 4 and shake upside down until all crystals have dissolved. In order to dissolve the Additive ECO better in the Extractor 4, the mixture can be heated in a water bath at 37 °C (99 °F). A volume of 25 mL Extractor Egg will be enough for the extraction of about 25 samples.

For the preparation of the final **Allergen extraction buffer containing Additive 1 (A-AEB)**, which is required for the extraction according to chapter 9.2, weigh 1.35 g of Additive 1 in a beaker and add 15 mL 1 M NaOH. Stir until the Additive 1 is solved. Then, fill 700 mL diluted AEB (see above) in a measuring cylinder. Add the 15 mL Additive 1 solution by stirring constantly. Transfer any residues of the Additive 1 solution into the measuring cylinder by rinsing with diluted AEB. Adjust the Additive 1 containing Allergen extraction buffer (A-AEB) to pH 9 with 1 M HCl and fill up to 750 mL with diluted AEB.

750 mL A-AEB are sufficient for 39 samples. The buffer can be used for approx. 3 weeks at room temperature 20 - 25 °C (68 - 77 °F) (do **not** store in the refrigerator). Discard the buffer if crystals are present. Crystal formation may begin after 10 days. We therefore recommend preparing a fresh buffer solution every week. Alternatively, a reduced volume of the A-AEB of e.g. 375 mL can be prepared. When preparing a smaller volume of A-AEB, the composition described above must be adhered to accordingly. Use clean bottles when preparing the buffer. Particles of dust can initiate crystallization.

For extraction according chapter 9.2, as an alternative to weighing the skim milk powder (SMP) to the sample, the SMP can be added to the A-AEB (**SMP-A-AEB**). Add 0.5 g of skim milk powder per 19 mL A-AEB. Only prepare as much SMP-A-AEB as needed. The shelf life of the SMP-A-AEB is one day.

In the following section, the following abbreviations are used:

- AEB: final diluted Allergen extraction buffer
- SMP-AEB final diluted Allergen extraction buffer with addition of SMP
- A-AEB: AEB with addition of Additive 1
- SMP-A-AEB A-AEB with addition of SMP

9.1 Sample extraction using AEB for **unprocessed / unheated** foods (e.g. ice cream, wine, chocolate, salad dressing)

Heat the AEB (or SMP-AEB) to 60 °C (140 °F) before sample extraction.

Homogenize a representative, adequately big amount of a solid sample (e.g. 50 g; grind it thoroughly to powder and mix well) or mix well a sample in case of liquid foods. Melt or grate chocolate.

- Transfer 1 g (or in case of liquid samples 1 mL) of homogenized sample to a new vial, add 0.5 g SMP and 20 mL (or 19 mL in case of liquid samples) pre-heated AEB (see chapter 9). Alternatively, add 20 mL (or 19 mL for liquid samples) of pre-heated SMP-AEB.
- Mix thoroughly until the sample is completely suspended (e.g. vortexer).
- Incubate for 10 minutes at 60 °C (140 °F), shaking occasionally.
- Let the sample cool down (e.g. in a water bath with ice).
- Filter sample or centrifuge for 10 min at > 2,500 x g; if possible at 4 °C (39 °F).

Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge.

- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The extract (supernatant of centrifugation step or filtrate) can be stored up to one day in a well-closed container at 2 - 8 °C (35 - 46 °F). Unused extracts can be stored at -20 °C (-4 °F) for three months.
- If further dilutions are required (samples with absorbance values ($A_{450\text{nm}}$) > Standard 5), these extract dilutions must be done with SMP-AEB (see chapter 9).

9.2 Sample extraction with Extractor Egg and A-AEB for **processed / heated** samples (e.g. noodles and cookies)

Heat the A-AEB (or SMP-A-AEB) to 60 °C (140 °F) before sample extraction.

Homogenize a representative, adequately big amount of a solid sample (e.g. 50 g; grind it thoroughly to powder and mix well) or mix well a sample in case of liquid foods. Melt or grate chocolate. Melt or grate chocolate.

- Transfer 1 g (or in case of liquid samples 1 mL) of homogenized sample to a new vial, add 0.5 g SMP and mix with 1 mL Extractor Egg and 19 mL (or

18 mL in case of liquid samples) pre-heated A-AEB (see chapter 9). Alternatively, add 19 mL (or 18 mL for liquid samples) of pre-heated SMP-A-AEB.

- Mix thoroughly (e.g. vortexer).
- Incubate for 10 minutes at 60 °C (140 °F), shaking occasionally.
- Let the sample cool down (e.g. in a water bath with ice).
- Filter sample or centrifuge for 10 min at $> 2,500 \times g$; if possible at 4 °C (39 °F).

Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed ($> 10,000 \times g$) for 10 min in a microcentrifuge.

- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The extract (supernatant of centrifugation step or filtrate) can be stored up to 4 weeks in a well-closed container at 2 - 8 °C (35 - 46 °F). Unused extracts can be stored at -20 °C (-4 °F) for three months.
- If further dilutions are required (samples with absorbance values ($A_{450\text{nm}}$) $>$ Standard 5), these extract dilutions must be done in a mixture of 1 mL Extractor Egg and 19 mL SMP-A-AEB (see chapter 9).

10. Test procedure

10.1 Test preparation

Bring all the reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10-fold concentrate. Before use, the buffer has to be diluted 1:10 (1 + 9) with dist. water (e.g. 900 mL dist. water + 100 mL buffer). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 μL per well) and then

exactly 100 µL are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

Avoid direct sunlight during all incubations. Therefore, cover the microtiter plates.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µL of each standard or sample (prepared according to chapter 9.1 or 9.2) in duplicate to the wells and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat four more times (a total of five wash cycles).
4. Add 100 µL conjugate each into the respective wells and incubate for 20 minutes at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat four more times (a total of five wash cycles).
6. Add 100 µL substrate/chromogen into the wells, and incubate for 10 minutes at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µL of stop solution into each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 min after addition of stop solution.

11. Evaluation

A special software, **RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter function.

When working in accordance with this extraction, the sample dilution factor is 20. Hence, a dilution factor of 20 is already taken into account with the standard concentrations (see chapter 4*), the egg concentration can directly be read from the standard curve.

It must be clarified that all quality criteria are met for the current test run. The course of the standard curve can be taken from the quality assurance certificate (certificate of analysis, CoA), which is available via the QR code on the test kit label. As the absorbance values in the laboratory may differ from those stated on the certificate, it is recommended to compare the ratios of the standards to each other with those on the certificate. For this purpose, the $B/B_{\max}$ values (the ratio of the absorbance values of the standards to the highest standard) are compared with each other. In the current test run, these should be comparable to the ratios of the standards on the certificate.

A further dilution and new detection of samples is recommended for absorbance values ($A_{450\text{nm}}$) > standard 5. Further dilutions should be made with the respective extraction buffer (see chapter 9.1 and 9.2). In case of a further dilution, the additional dilution factor must be taken into account when calculating the egg concentration.

The assay can be also evaluated when running in single determinations. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not consistent with standards like EN 15633-1 and EN 15842. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, when performed as a single determination, a higher degree of inaccuracy can be expected, and when performing repeated tests, a higher degree of variation in the results can be expected.

12. Result interpretation

The test is calibrated against the NIST reference material 8445 (whole egg powder). The result is expressed as mg/kg whole egg powder. Reference material 8445 contains 49 % +/- 1 % total protein. Multiply by a factor of 0.49 to convert the result into mg/kg total protein; multiply by a factor of 0.263 for the result in mg/kg egg white protein.

Higher absorbance values ($A_{450\text{nm}}$) of the calibration curve, especially for the zero standard, as mentioned on the certificate may be a result of insufficient washing or a egg contamination of reagents.

Results between LoD and LoQ indicate a low allergen concentration in the sample. Depending on the matrix tested, values below the LoQ can still be determined with sufficient precision ($CV < 30\%$). However, values in this range are generally subject to greater uncertainty due to the higher fluctuation range of the test. If the precision of the test has not been validated with a specific sample matrix, results below the measuring range should therefore not be reported with a quantitative value, but qualitatively "< LoQ". Further information on this can be found in the validation report.

A result below the LoD does not exclude an egg contamination below the detection limit of the assay, or that other allergen components, such as lipids, may be present in a sample. The result should be reported accordingly.

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LoD, LoQ, upper limit of measurement range).

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max. $\pm 1\%$.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive or increased results, but also reduce or suppress a

correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

In general, matrix effects in immunological test methods can be suppressed by adding foreign protein (depending on the test, e.g., BSA, gelatin, skim milk powder) during extraction. The influence of the addition on recovery should be checked by spiking experiments.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on recovery and assay results.

Cross-reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. In single cases, cross-reactivities can be lost or only become apparent through the processing of foods.

For evaluation of the cross-reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross-reactivities are described in the validation report.

For detailed results and further information about other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

14. Recommendation

In order to ensure a high analytical performance, we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control and ensurance of an accurate and correct test procedure. Allergen-free and allergen-containing (spiked) samples should be used. An example of a spike experiment is given in the validation report.

- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- To perform PCR (e.g. SureFood®) for confirmation of the result.
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt® / Bolt™) are used.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2019-08-14	Release version
2024-05-13	Previous version
2025-10-20	Current version Changes made: – Addition of the registered trademark RIDASOFT® – General linguistic revision – Addition of crystal formation and preparation of A-AEB in chapter 9 – Transfer of the note for single determinations to chapter 11 – Revision and additions in chapter 12, 13 and 14 – New disclaimer

Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM-DD)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM-DD)
	Manufacturer + address

Patent Marking:

The extraction means in this product contains sulfite. Food inspection methods using a sulfite-containing extractant as in this product and/or corresponding detection kits are subject to the following patents of MORINAGA & Co., Ltd.: European Patent EP 2 224 239 B1, Australian Patent AU 2008 330 507 B2, United States Patent US 8 859 212 B2, Japanese Patent JP 5 451 854 B2. The patent holder has granted R-Biopharm AG a license to use, and sell products that employ, said protected technology in the above-mentioned territories.

Disclaimer

1. In conformance with the German Civil Code ("BGB") R-Biopharm AG provides a limited warranty ("Gewährleistung") against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm's own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG's liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG's sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser's expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements ("Gewährleistung") implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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are registered trademarks of R-Biopharm AG