r-biopharm®

## RIDASCREEN® Walnut

**REF** R6601

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# RIDASCREEN® Walnut

## Brief information

RIDASCREEN® Walnut (Art. No. R6601) is a sandwich enzyme immunoassay for the quantitative analysis of raw and roasted walnut or walnut protein in food validated for the method (see chapter 1).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)... approx. 20 min  
test implementation (incubation time)..... 50 min

Standard material: walnut paste  
(protein content:  $15.8 \pm 0.32$  g / 100 g)

Limit of detection: 0.12 mg/kg (ppm) walnut (mean value)  
(depending on matrix) (corresponds to 0.019 mg/kg walnut protein)  
0.00 - 0.30 mg/kg walnut  
(corresponds to 0.00 - 0.047 mg/kg walnut protein)

Limit of quantification: 1 mg/kg (ppm) walnut  
(corresponds to 0.158 mg/kg (ppm) walnut protein)

Specificity: The antibodies specifically detect raw and roasted walnut proteins. There is a slight cross-reactivity to pecan. No cross-reactivity was found against a further 98 tested foods. Further information is contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments (see chapter 13. Limits of the method).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

### **Related product and accessories for walnut determination**

Bioavid lateral flow Walnut incl. hook line (Art. No. BLH707-15)

SureFood® ALLERGEN Walnut (Art. No. S3607)

SureFood® ALLERGEN 4plex Peanut / Hazelnut / Walnut + IAC  
(Art. No. S3402)

## 1. Intended use

RIDASCREEN® Walnut (Art. No. R6601) is a sandwich enzyme immunoassay for the quantitative analysis of raw and roasted walnut or walnut protein in food. Due to the large number of different food products, the following samples were examined as representative for different food product categories within the scope of the test development: cookies, chocolate, ice cream and breakfast cereals. It can be assumed that the assay is also suitable for the analysis of other foods. However, this must be verified by the user before applying the test kit to that food category.

For detailed results and further information on validation data with other food matrices, please refer to the validation report.

## 2. General information

The walnut (*Juglans regia*) is one of the most popular tree nuts in the world. Walnuts are very common ingredients in chocolate products and are also consumed as snacks and in baked goods such as cookies. Furthermore, they are used as an additive for ice cream, dressings and salads.

Walnuts are one of the eight foods that are the most common cause of food allergies. Eating less than one milligram walnut protein can trigger allergic reactions in sensitized individuals. The main danger for sensitized individuals is the contamination of walnuts in food (e.g. traces of walnuts in chocolate, which should not contain walnuts). Therefore, it is important to test foods that may have been contaminated during the production process for the presence of unintentional walnut traces.

According to Regulation (EU) No. 1169/2011, walnuts and their products must be declared as an allergen on food labels. There are similar regulations in the USA, Canada, Australia and New Zealand, for example.

## 3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against walnut proteins. By adding the standard or sample solution to the wells, walnut proteins present will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Following the washing step, a solution containing antibody conjugated to peroxidase is added. This conjugate is bound to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed.

Any unbound conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue end product. A stop solution is added which results in a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance of the solution which is proportional to the walnut protein concentration in the sample is measured photometrically at 450 nm and expressed as mg/kg (ppm) walnut or walnut protein.

#### 4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
<b>Microtiter plate</b>	-	Ready to use		96 wells
<b>Allergen extraction buffer</b>	Green	<b>Concentrate</b>	<b>10x</b>	100 mL
<b>Standard 1*</b>	Transparent	Ready to use	0 mg/kg	1.3 mL
<b>Standard 2*</b>	Transparent	Ready to use	1 mg/kg	1.3 mL
<b>Standard 3*</b>	Transparent	Ready to use	3 mg/kg	1.3 mL
<b>Standard 4*</b>	Transparent	Ready to use	9 mg/kg	1.3 mL
<b>Standard 5*</b>	Transparent	Ready to use	18 mg/kg	1.3 mL
<b>Wash buffer</b>	Brown	<b>Concentrate</b>	<b>10x</b>	100 mL
<b>Conjugate</b>	Red	Ready to use		11 mL
<b>Substrate/Chromogen</b> Red Chromogen Pro	Brown	Ready to use		13 mL
<b>Stop solution</b>	Yellow	Ready to use		14 mL

\* The dilution factor 20, which results after sample preparation, has already been considered for the standard concentrations. Therefore, the walnut concentration of the samples can directly be read from the standard curve.

#### 5. Reagents required but not provided

##### 5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of  $\pm 0.01$  g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials with cap (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (60 °C / 140 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)

- Fluted filter (pore size 8 - 12 µm)
- Graduated pipettes
- Graduated cylinders
- Variable 20 - 200 µL and 200 - 1000 µL micropipettes
- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100 µL
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)

## 5.2 Reagents

- Distilled water (dist. water) or deionized water
- Walnut-free skim milk powder (SMP; food quality)

## 6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

Do not reuse wells of the microtiter strips (coated microtiter plate and pre-plate if necessary, see chapter 10.2). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (36 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 1.2 absorbance units ( $A_{450\text{ nm}} < 1.2$ ) for standard 5

## 9. Sample preparation

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation.
- Carry out the sample preparation in a room isolated from the ELISA procedure.

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The Allergen extraction buffer is provided as a 10-fold concentrate and must be diluted prior use. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that, dilute the heated buffer concentrate 1:10 (9 + 1) with dist. water (e.g. 900 mL dist. water + 100 mL buffer concentrate). The **diluted Allergen extraction buffer (AEB)** is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks or at 2 - 8 °C (36 - 46 °F) for 12 weeks.

Make sure that the AEB is heated in time in the 60 °C (140 °F) water bath.

In the following chapter, the following abbreviations are used:

- AEB: final diluted Allergen extraction buffer
- SMP: skim milk powder

## 9.1 Extraction with AEB and SMP

Homogenize a representative, adequately big amount of a solid sample (e.g. 50 g; grind it thoroughly to powder and mix well) or mix well a sample in case of liquid foods.

- Transfer 1 g (or in case of liquid samples 1 mL) of homogenized sample to a new vial and add 1 g SMP.
- Add 20 mL (or in case of liquid samples 19 mL) pre-heated AEB (see chapter 9).
- Mix vigorously (e.g. vortex).
- Extract for 10 min at 60 °C (140 °F) in a water bath by shaking occasionally.
- Let the sample cool down (e.g. 3 - 5 min in ice water).
- Filter or centrifuge sample for 10 min at > 2,500 x g; if possible at 4 °C (39 °F).

(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge.)

- Transfer the supernatant to a new vial.
- If the supernatant is not free of particles after centrifugation, filter the extract.
- The extract (supernatant of centrifugation step or filtrate) should be tested by ELISA on the same day if stored at room temperature (20 - 25 °C / 68 - 77 °F). It can be stored for up to two days in a well-closed container at 2 - 8 °C (36 - 46 °F). Unused extracts can be stored at - 20 °C (- 4 °F) for one month.

### Remark:

As an alternative to weighing the SMP to the sample, the SMP can be added to the AEB (1 g SMP per 20 mL AEB). The shelf life of the SMP-AEB is only 1 day.

In case of samples with concentrations above the measuring range, further dilutions should be made with the SMP-AEB (1 g SMP per 20 mL AEB).

## 10. Test procedure

### 10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10-fold concentrate. Before use, the buffer has to be diluted 1:10 (9 + 1) with dist. water (e.g. 900 mL dist. water + 100 mL buffer concentrate). Prior to dilution, dissolve eventually formed crystals by

incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

Components should be stored immediately at 2 - 8 °C (36 - 46 °F) when no longer required.

## 10.2 Test procedure

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µL per well) and then exactly 100 µL are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µL of each standard or sample (prepared according to chapter 9) in duplicate to the wells and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat three more times (a total of four wash cycles).
4. Add 100 µL of the conjugate to each well and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat three more times (a total of four wash cycles).
6. Add 100 µL of the reddish substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.

7. Add 100 µL of stop solution into each well. Mix gently by shaking the plate manually and measure the extinction at 450 nm. Read within 10 min after addition of stop solution.

## 11. Evaluation

A special software, **RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter function.

It must be clarified that all quality criteria are met for the current test run. The course of the standard curve can be taken from the quality assurance certificate (certificate of analysis, CoA), which is available via the QR code on the test kit. As the absorbance values in the laboratory may differ from those stated on the certificate, it is recommended to compare the ratios of the standards to each other with those on the certificate. For this purpose, the  $B/B_{\max}$  values (the ratio of the absorbance values of the standards to the highest standard) are compared with each other. In the current test run, these should be similar to the ratios of the standards on the certificate.

When working according to these instructions for use, the samples are diluted 1:20 during extraction. The sample dilution factor of 20 is already considered in the concentration data of the standards (see chapter 4\*). The concentration of walnut in the sample can therefore be read directly from the standard curve.

The assay can be also evaluated when running in single well per extract. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not consistent with standards like EN 15633-1 and EN 15842. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single wells per extract.

## 12. Result interpretation

The test result is calculated in mg walnut per kg food. Since the test is calibrated against walnut paste with a protein content of 15.8 %, the result must be multiplied by 0.158 to obtain mg/kg walnut protein.

Results between LoD (Limit of Detection) and LoQ (Limit of Quantification) indicate a low allergen concentration in the sample and such results should not be reported with a quantitative value, but qualitative as “< LoQ”.

A result below the LoD does not exclude an allergen contamination below the detection limit of the assay, or that other allergen components, such as lipids, may be present in a sample. The result should be reported accordingly.

A further dilution and new detection of samples is recommended for absorbance values ( $A_{450\text{ nm}}$ ) > standard 5. In case of a further dilution, the additional dilution factor must be considered when calculating the walnut concentration. Further dilutions should be made with the SMP-AEB (1 g per 20 mL AEB).

Compared to the certificate, higher absorbance values ( $A_{450\text{ nm}}$ ) for the standard curve, especially for the zero standard, may be a result of insufficient washing or walnut contamination of single kit components.

### **13. Limits of the method**

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

The protein content and protein composition may vary in different walnut varieties. Therefore, different walnut varieties may give different results, as the test was calibrated against exemplary walnut varieties used with the standard material.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LoD, LoQ, upper limit of measurement range).

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max.  $\pm 1$  %.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and interlaboratory comparisons.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When

analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

Cross reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

## 14. Recommendation

To ensure a high analytical performance, we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control. Walnut-free and walnut-containing (spiked) samples should be used. An example of a spiking experiment is given in the validation report.
- To do spike experiments to ensure an accurate and correct test procedure.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- To perform PCR (e.g. SureFood®) for confirmation of the result.
- To contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de) if automates (e.g. ThunderBolt® / Bolt™) are used.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

## Version overview

Version number	Chapter and title
2024-08-30	Release version
2024-12-03	Current version Changes made: – Correction: no cross-reactivity to anis

## Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

## Disclaimer

1. In conformance with the German Civil Code (“BGB”) R-Biopharm AG provides a limited warranty (“Gewährleistung”) against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.  
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm’s own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
  - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
  - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
  - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
  - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG’s liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG’s sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser’s expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements (“Gewährleistung”) implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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are registered trademarks of R-Biopharm AG.