RIDASCREEN® Risk Material

REF R6701

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RIDASCREEN® Risk Material

Brief Overview

RIDASCREEN® Risk Material (Art. No.: R6701) is a sandwich enzyme immunoassay for the detection of specified risk materials in processed (heated) meat and sausages, as defined in Annex V of Regulation (EC) No. 999/2001, for the prevention, control, and eradication of certain transmissible spongiform encephalopathies (see Chapter 1: Intended Use). These include central nervous system (CNS) tissues - such as brain, spinal cord, eyes, and tonsils (hereinafter referred to as “CNS”). For detection of CNS in or on raw meat and sausage products, as well as on contaminated surfaces, RIDASCREEN® Risk Material 10/5 (Art. No.: R6703) is recommended.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization, extraction, centrifugation and dilution
Time requirement:	sample preparation (10 samples) approx. 1 h test implementation (incubation time) 1 h
Limit of Detection (LoD):	≤ 0.2 % for CNS tissue* in processed (heated) meat and sausages (Matrix-dependent)
Limit of Quantification (LoQ):	0.2 % for CNS tissue* in processed (heated) meat and sausages (Matrix-dependent)
Specificity:	The RIDASCREEN® Risk material test allows detection of CNS tissue from cattle (incl calves), sheep, goat, horse, poultry and pig. The test does not allow identification of tissue type (e.g. brain, spinal cord or ocular).

*The test is calibrated specifically to bovine CNS (brain) tissue

We additionally refer you to our ELISA Manual to enhance the quality of your ELISA testing process. It outlines minimum standards regarding conditions to be observed when using and performing ELISA analysis with R-Biopharm AG tests. The manual can be accessed, printed and downloaded at:

<https://food.r-biopharm.com/media/technical-guides/>.

Additional products and accessories for CNS tissue detection

RIDASCREEN® Risk Material 10/5 R6703

Swabs and sample tubes Z00101

1.Intended use

RIDASCREEN® Risk Material is a sandwich enzyme immunoassay for the detection of risk material (CNS) in processed (heated) meat and sausages.

2. General information

BSE (bovine spongiform encephalopathy) was first diagnosed in the United Kingdom in 1986. The mechanism of infection is not yet fully understood. However, it is undisputed that the pathologically altered prion protein (prP^{Sc}) plays a major role in the transmission and development of this neurodegenerative disease.

Feeding inadequately heat-treated meat-and-bone meal likely led to the spread of the disease, resulting in the culling of over 180,000 animals by the end of 2000. Experimental infection using infectious brain material successfully induced disease in various animal species. Similar disease patterns have been described in sheep, goats, cats, and other mammals.

It is considered highly likely that the new variant of Creutzfeldt-Jakob Disease (vCJD) is caused by the consumption of infectious meat and meat products. The mortality rate of this disease, which had claimed 111 lives in the UK by May of 2002, is 100%. As a protective measure, the European Union banned the feed of meat-and-bone meal (2000/766/EC) and the processing of specified risk material in food (2000/418/EC and 2001/999/EC).

3. Test principle

The assay is based on the antigen-antibody reaction. Detection of specified risk material is achieved through the quantification of Glial Fibrillary Acidic Protein (GFAP), a cellular marker found in particular high concentrations in CNS tissue. The wells of the microtiter strips are coated with specific antibodies against GFAP. If CNS tissue is present, the glial fibrillary acidic protein contained then binds to the specific capture antibodies. A peroxidase-conjugated antibody against GFAP (enzyme conjugate) is then used to detect the bound antigen. Unbound enzyme conjugate is removed in a subsequent washing step. Detection is carried out by adding substrate/chromogen. The bound enzyme conjugate converts the chromogen into a blue end product. The addition of the stop reagent causes a color change from blue to yellow. It is measured photometrically at 450nm. The absorbance of the solution is proportional to the concentration of specified risk material in the sample.

4. Package (Kit) Content

Each Kit contains sufficient reagents for max. 96 measurements (including standard analyses). Each package contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 Cavities
Sample buffer	White	Concentrate	10 x	10 mL
Extraction buffer	White	Concentrate	10 x	100 mL
Standard 1*	White	Ready to use	0 %	1.3 mL
Standard 2*	White	Ready to use	0.2 %	1.3 mL
Standard 3*	White	Ready to use	1.0 %	1.3 mL
Standard 4*	White	Ready to use	2.0 %	1.3 mL
Wash buffer salt Tween	-	Salt for dissolution	-	-
Conjugate	Red	Ready to use	-	6 mL
Substrate/Chromogen	Brown	Ready to use	-	10 mL
Stop solution	Yellow	Ready to use	-	14 mL

*) The concentration values already account for the dilution factor resulting from sample preparation. This allows the CNS tissue concentrations (%) of the samples to be read directly from the standard curve.

5. Materials required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal glass vials
- Impact mill, mortar, ultra-turrax or homogenizer
- water bath 100 °C (212 °F)
- graduated pipettes
- variable 20 µL - 200 µL and 200 - 1000 µL micropipettes
- Laboratory gloves
- If necessary, a scale (measuring range at least up to 2 g, accuracy of ± 0.01 g)
- Graduated cylinder (plastic or glass) 25mL
- Vortexer If necessary, Multistep or 8-channel pipette for 50 µL, 100 µL and 250 µL
- Optional: RIDA® ABSORBANCE 96 (Art. Nr. ZRA96FF)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. Nr. Z9996FF)

5.2. Reagents:

- Distilled or deionized water

6. Precautions

This test is intended exclusively for use within its specified purpose.

It must be performed only by trained laboratory personnel. The instructions for use must be followed precisely.

Safety information regarding its components can be found in the product Safety Data Sheets (SDS), available at:

<https://food.r-biopharm.com/products/ridascreen-risk-material/>

The microtiter strip wells (coated microtiter plate included in a Kit) must not be reused. Use separate pipette tips for each standard and each sample extract to prevent cross-contamination.

All reagents and other materials must be disposed of responsibly after use to protect human health and the environment. Observe all applicable national laws and regulations when disposing of materials.

7. Reagent Storage and Handling

Store all reagents between 2 °C to 8 °C (35 °F to 46° F). Never freeze any Test Kit components.

To remove microtiter strips, only open the foil pouch once it has reached room temperature (20 °C to 25 °C) (68 °F to 77 °F) to prevent condensation from forming in the wells.

Any unused wells and desiccant should be stored together in the sealed foil pouch and kept refrigerated between 2 °C to 8 °C (35 °F to 46 °F).

Protect the reddish-colored substrate/chromogen from direct light exposure, as it is sensitive to light.

Product quality can no longer be guaranteed once the expiration date has passed (refer to the Test Kit's outer label for the expiration date).

Individual reagents are not interchangeable with Kits from different Lot numbers.

8. Signs of Reagent Deterioration

- A bluish discoloration of the reddish substrate/chromogen solution prior to dispensing into the wells
- An extinction value below 0.6 ($E_{450\text{ nm}} < 0.6$) for Standard 4 (2.0%)

9. Sample Preparation

Keep samples refrigerated until use.

Before starting, allow all reagents and samples to reach room temperature (20 °C to 25 °C) (68 °F to 77 °F) and perform all sample preparation steps at room temperature.

9.1. Component Preparation

Before use, dilute the sample buffer, which is supplied as a tenfold concentrate, at a ratio of 1:10 with distilled water (e.g. 10 mL concentrate plus 90 mL distilled water). We recommend preparing the dilution in a single step (i.e. mix 10 mL of buffer concentrate with 90 mL of distilled water). If crystals are detected within the concentrate, they must be dissolved beforehand by warming the solution in a water bath at 37 °C (98 °F).

The diluted sample buffer should be stored at 2 °C to 8 °C (35 °F to 46 °F) and may be used until the expiration date indicated on the buffer label.

The **extraction buffer** is provided as a 10fold concentrate. Therefore the buffer concentrate has to be diluted 1:10 (1+9) with dist. water before use. Dissolve eventually occurring crystals in a water bath at 37 °C (98 °F) before dilution. We recommend dilution of the buffer concentrate in only one step (i.e. 100 mL concentrate + 900 mL dist. water).

The diluted buffer can be used up to the expiry date printed on the buffer label. To avoid new crystallization, the buffer can be stored at room temperature (20 °C– 25 °C / 68 °F- 77 °F).

9.2. Meat and Meat product samples

- Before analysis, thoroughly mix the sample to ensure reliable test results. CNS tissue contamination may be distributed unevenly within a sample, which can lead to inconsistent results.
- mix 2 g of the homogenized sample with 10 mL diluted extraction buffer (see 9.1.) in a centrifugal glass vial and shake vigorously for suspending (vortex)
- heat the sample for 15 min at 100 °C (212 °F)
- chill down the sample and centrifuge: 10 min / 3500 g / 10 °C (50 °F)
- remove the upper fat layer completely (draw off, e.g. by using a pasteur pipette)
- dilute the fat free supernatant 1:20 (1+19) with diluted sample dilution buffer (see 9.1.), e.g. 50 µL of supernatant + 950 µL of diluted buffer
- use 50 µL per well in the assay

10. Test Execution

10.1. Test Preparation

Bring all reagents to room temperature (20 °C to 25 °C) (68 °F to 77 °F) before use.

A PBS-Tween buffer is needed as **washing buffer**, please use the washing buffer salt (pouch) contained in the kit (see 4.). Dissolve the total content of the pouch in one liter of distilled water. The ready to use washing buffer expires after approx. 4 - 6 weeks at 2 °C- 8 °C (36 °F- 46 °F).

Alternative:

- Dissolve the contents of the sachet in 100 mL of distilled water to create a 10x concentrate. To prepare the ready-to-use solution, mix 1 part of the 10x concentrate with 9 parts of distilled water.
- The 10x concentrate can be stored at room temperature (20 °C to 25 °C) (68°F to 77 °F) for approximately 8 to 12 weeks.

Immediately return any unused reagents to storage from 2 °C to 8 °C (36 °F- 46 °F).

10.2. Test Execution

Avoid drying out the wells between steps. Consistent and thorough washing is critical for reliable test results. Always follow the specified washing sequences to ensure uniform treatment across all wells.

Use a multichannel or multistepper pipette when dispensing the conjugate, substrate/chromogen, and stop solution. This helps prevent timing delays across the microtiter plate.

During all incubation steps, keep the microtiter plate covered and protect it from direct sunlight.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
2. Pipet 50 µL of standard solution or prepared sample to separate wells.
3. Add 50 µL of enzyme conjugate (red cap) to the bottom of each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 °C- 25 °C / 68 °F- 77 °F)

4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times (total of three wash cycles).
5. Add 100 µL of substrate/chromogen solution (brown cap) to each well and incubate for 30 min at room temperature (20 °C- 25 °C, 68 °F- 77 °F) in the dark.
6. Add 100 µL of stop solution (yellow cap) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 minutes after addition of stop solution.

11. Interpretation

Color intensity of selected Standards serves as interpretive reference.

R-Biopharm offers an optional interpretation software tool specifically developed for the RIDASCREEN® enzyme immunoassays:

RIDASOFT® Win.NET Food & Feed (Art. Nr. Z9996FF).

Use 4-Parameter function to evaluate the results.

Verify that the test run meets all quality criteria from the certificate of analysis before evaluating. The standard curve shape is also provided in the certificate.

12. Interpretation of results

If a result is below the LoD, report it as “<Detection limit (%)” since results below the LoD cannot categorically rule out all presence of CNS tissue.

13. Method Limitations

Follow all instructions (especially for sample prep, test execution, and storage) to achieve the specified detection limit.

Detection results and limits are influenced by several factors, including the sample matrix, laboratory conditions, level of processing, and extraction method.

This test is designed to detect central nervous system (CNS) tissue by identifying the Glial Fibrillary Acidic (GFAP) marker protein. It is important to note that the test is

calibrated specifically to bovine CNS (brain) tissue. When interpreting results, consider that marker protein concentrations may differ across species.

Marker protein concentrations may naturally vary within the same tissue type and between different types of CNS tissue (e.g. brain, spinal cord, eye), even among animals of the same species. Such variability can be significant. It may lead to differing levels of marker protein in the sample, which in turn can result in positive test outcomes at varying concentrations of CNS tissue.

To ensure accurate interpretation of results, always account for both species-specific and tissue-specific variability in marker protein content.

14. Recommendation

For reliable results:

- Pre-rinse pipette tips with the standard solution or sample extract before use
- Include control tests in routine quality checks to validate results and detect deviations.

15. Additional Applications

For additional details regarding the product and its uses, please contact:
info@r-biopharm.de.

Version overview

Version number	Chapter and title
2010-07-01	Release version
2025-11-21	Current version Changes made: Comprehensive linguistic and formal revision Added and corrected product listings related to the detection of central nervous system (CNS) tissue

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

1. R-Biopharm AG provides a warranty that its products meet the specified requirements and are free from material and workmanship defects for 12 months following the transfer of risk. For products with limited shelf life, the warranty remains valid until the stated expiration date or, for limited-use items, until the stated number of uses has been reached. Warranty claims must be promptly and correctly submitted by the customer. A product is considered defective if it does not have the agreed features, is not suitable for its intended contractual use, or is not provided with the required accessories and instructions (collectively "Subjective Requirements").

2. R-Biopharm AG assumes no liability for any consequences arising from the following:

- (a) Failure to read, comprehend, or adhere to product usage or safety instructions;
- (b) Employment of individuals who do not possess the necessary training or qualifications to use the product;
- (c) Non-compliance with relevant industry standards and practices, including Good Laboratory Practices;
- (d) Utilization of control materials, samples, sample matrices, or processing that are not appropriate for use with the product;
- (e) Any other improper use of the product;
- (f) Modifications or alterations of the product;
- (g) Improper storage by the customer or others;
- (h) Exposure to chemical, electromagnetic, mechanical, or electrolytic influences beyond standard ranges documented by R-Biopharm AG;
- (i) Damage or malfunctions resulting from external events beyond the control of R-Biopharm AG, such as burglary, theft, lightning, fire, water damage, or Force Majeure.

3. R-Biopharm AG is liable in cases of fraud, gross negligence, or intent, for injury to life, body, or health, for the assumption of a guarantee or procurement risk under Section 276 of the German Civil Code (BGB), or where liability is mandated by law.

For slight negligence, R-Biopharm AG is only liable for breaches of material contractual obligations necessary for the contract's purpose and relied on by the customer. In these cases, liability is limited to typical, foreseeable damage. R-Biopharm AG assumes no liability for other minor negligent breaches.

ALL OTHER EXPRESS OR IMPLIED WARRANTIES OR GUARANTEES OF ANY KIND ARE EXCLUDED, REGARDLESS OF WHETHER THEY ARE ASSUMED THROUGH CUSTOM, PRACTICE, COURSE OF DEALING, OR OTHERWISE.

R-Biopharm AG accepts no liability for consequential damages, in particular lost profits, production delays, or other indirect losses.

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© 2025 R-Biopharm AG
are registered trademarks of R-Biopharm AG.