RIDASCREEN® Risk Material 10/5

REF R6703

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RIDASCREEN® Risk Material 10/5

Brief Overview

RIDASCREEN® Risk Material 10/5 (Art. No.: R6703) is a sandwich enzyme immunoassay for the detection of specified risk materials in or on raw meat and sausage products, as well as on contaminated surfaces, as defined in Annex V of Regulation (EC) No. 999/2001, for the prevention, control, and eradication of certain transmissible spongiform encephalopathies (see Chapter 1: Intended Use). These include central nervous system (CNS) tissues - such as brain, spinal cord, eyes, and tonsils (hereinafter referred to as “CNS”).

This test is not suitable for the analysis of heat-treated meat or sausage products. For this purpose, RIDASCREEN® Risk Material (Art. No.: R6701) is recommended.

All reagents required for performing the enzyme immunoassay, including internal standards, are included in the Test Kit. The Test Kit is sufficient for a maximum of 96 determinations (including standard determinations). A microtiter plate photometer is required for evaluation.

Sample preparation:	Transfer the swab-collected or weighed sample into sample buffer. If necessary, vortex and centrifuge
Time requirement:	sample preparation (10 meat samples) approx. 10 min test implementation (incubation time) 15 min
Limit of Detection (LoD):	≤ 0.1 % for CNS tissue* in raw meat and meat products (Matrix-dependent)
Limit of Quantification (LoQ):	0.1 % for CNS tissue* in raw meat and meat products (Matrix-dependent)

*The test is calibrated specifically to bovine CNS (brain) tissue.

Specificity: The RIDASCREEN® Risk material 10/5 test allows detection of CNS tissue from cattle (incl. calves), sheep, goat, horse, poultry and pig. The test does not allow identification of tissue type (e.g. brain, spinal cord or ocular).

We additionally refer you to our ELISA Manual to enhance the quality of your ELISA testing process. It outlines minimum standards regarding conditions to be observed when using and performing ELISA analysis with R-Biopharm AG tests. The manual can be accessed, printed and downloaded at: <https://food.r-biopharm.com/media/technical-guides/>.

Additional products and accessories for CNS tissue detection

RIDASCREEN® Risk Material R6701
Swabs and sample tubes Z0010

1. Intended Use

RIDASCREEN® Risk Material 10/5 is a sandwich enzyme immunoassay for the detection of risk material (CNS) in or on raw meat, meat products, as well as on contaminated surfaces.

2. General information

BSE (bovine spongiform encephalopathy) was first diagnosed in the United Kingdom in 1986. The mechanism of infection is not yet fully understood. However, it is undisputed that the pathologically altered prion protein (prP^{Sc}) plays a major role in the transmission and development of this neurodegenerative disease.

Feeding inadequately heat-treated meat-and-bone meal likely led to the spread of the disease, resulting in the culling of over 180,000 animals by the end of 2000. Experimental infection using infectious brain material successfully induced disease in various animal species. Similar disease patterns have been described in sheep, goats, cats, and other mammals.

It is considered highly likely that the new variant of Creutzfeldt-Jakob Disease (vCJD) is caused by the consumption of infectious meat and meat products. The mortality rate of this disease, which had claimed 111 lives in the UK by May of 2002, is 100%.

As a protective measure, the European Union banned the feed of meat-and-bone meal (2000/766/EC) and the processing of specified risk material in food (2000/418/EC and 2001/999/EC).

3. Testing Principles

The assay is based on the antigen-antibody reaction. Detection of specified risk material is achieved through the quantification of Glial Fibrillary Acidic Protein (GFAP), a cellular marker found in particular high concentrations in

CNS tissue. The wells of the microtiter strips are coated with specific antibodies against GFAP. If CNS tissue is present, the glial fibrillary acidic protein contained then binds to the specific capture antibodies. A peroxidase-conjugated antibody against GFAP (enzyme conjugate) is then used to detect the bound antigen. Unbound enzyme conjugate is removed in a subsequent washing step. Detection is carried out by adding substrate/chromogen. The bound enzyme conjugate converts the chromogen into a blue end product. The addition of the stop reagent causes a color change from blue to yellow. It is measured photometrically at 450nm. The absorbance of the solution is proportional to the concentration of specified risk material in the sample.

4. Package (Kit) Content

Each Kit contains sufficient reagents for 96 measurements (including standard analyses). Each package contains:

Component	Cap color	Format		Volume
Microtiterplate	-	Ready to use		96 Cavities
Sample buffer	White	Concentrate	10 x	10 mL
Standard 1*	White	Ready to use	0 %	1.3 mL
Standard 2*	White	Ready to use	0.1 %	1.3 mL
Standard 3*	White	Ready to use	0.2 %	1.3 mL
Standard 4*	White	Ready to use	0.4 %	1.3 mL
Wash buffer salt Tween	-	Salt for dissolution	-	-
Conjugate	Red	Ready to use	-	6 mL
Substrate/Chromogen	Brown	Ready to use	-	10 mL
Stop solution	Yellow	Ready to use	-	14 mL
2x Pasteur Pipette	-	Ready to use		1 mL

*) The concentration values already account for the dilution factor resulting from sample preparation. This allows the CNS tissue concentrations (%) of the samples to be read directly from the standard curve.

5. Additional Reagents Required – Required Accessories

- Variable 20 to 200µl and 200 to 1000µl micropipettes
- Mikrotiter Plate Photometer (450 nm)
- Laboratory gloves
- If necessary, a scale (measuring range at least up to 50mg, accuracy of $\pm 0.001\text{g}$)
- Impact mill, mortar, ultra-turrax or homogenizer
- Graduated cylinder (plastic or glass) 25mL
- Centrifuge (min. 3,500 x g) & centrifuge test tubes with cap
- Vortexer If necessary, Multistep or 8-channel pipette for 50µL, 100µL and 250µL
- Optional: RIDA® ABSORBANCE 96 (Art. Nr. ZRA96FF)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. Nr. Z9996FF)

As sampling accessories under Item No.: Z0010 from R-Biopharm

- Swabs
- Sample Tubes

5.1 Reagents

- Distilled or deionized water

6. Precautions

This test is intended exclusively for use within its specified purpose.

It must be performed only by trained laboratory personnel. The instructions for use must be followed precisely.

Safety information regarding its components can be found in the product Safety Data Sheets (SDS), available at:

<https://food.r-biopharm.com/de/produkte/ridascreen-risk-material-105/>

The microtiter strip wells (coated microtiter plate included in a Kit) must not be reused. Use separate pipette tips for each standard and each sample extract to prevent cross-contamination.

All reagents and other materials must be disposed of responsibly after use to protect human health and the environment. Observe all applicable national laws and regulations when disposing of materials.

7. Reagent Storage and Handling

Store all reagents between 2°C to 8 °C (35° to 46° F). Never freeze any Test Kit components.

To remove microtiter strips, only open the foil pouch once it has reached room temperature (20°C to 25 °C) (68° to 77°F) to prevent condensation from forming in the wells.

Any unused wells and desiccant should be stored together in the sealed foil pouch and kept refrigerated between 2°C to 8 °C (35° to 46°F).

Protect the reddish-colored substrate/chromogen from direct light exposure, as it is sensitive to light.

Product quality can no longer be guaranteed once the expiration date has passed (refer to the Test Kit's outer label for the expiration date).

Individual reagents are not interchangeable with Kits from different Lot numbers.

8. Signs of Reagent Deterioration

- A bluish discoloration of the reddish substrate/chromogen solution prior to dispensing into the wells
- An extinction value below 0.6 ($E_{450nm} < 0.6$) for Standard 4 (0.4%)

9. Sample Preparation

Keep samples refrigerated until use.

Before starting, allow all reagents and samples to reach room temperature (20°C to 25°C) (68° to 77°F) and perform all sample preparation steps at room temperature.

9.1 Component Preparation

Before use, dilute the sample buffer, which is supplied as a tenfold concentrate, at a ratio of 1:10 with distilled water (e.g. 10ml concentrate plus 90ml distilled water). We recommend preparing the dilution in a single step (i.e. mix 10 mL of buffer concentrate with 90 mL of distilled water). If crystals are detected within the concentrate, they must be dissolved beforehand by warming the solution in a water bath at 37°C (98°F).

The diluted sample buffer should be stored at 2°C to 8°C (35° to 46° F) and may be used until the expiration date indicated on the buffer label.

For each determination, prepare 1 mL of diluted buffer in a sample tube using one of the graduated pipettes provided in the Test Kit.

9.2 Raw meat and meat product (sausage) samples

- Before analysis, thoroughly mix the sample to ensure reliable test results. CNS tissue contamination may be distributed unevenly within a sample, which can lead to inconsistent results.
- For homogenized samples, dip the dry swab several times into the raw sample mass. For pieces of solid meat, wipe the surface thoroughly with the swab, rotating it to ensure full contact between the entire swab surface and the sample.
- Transfer the sample by pressing the swab several times against the inner wall of the tube while rotating it in 1 mL of diluted sample buffer (see section 9.1).
- (Alternatively, weigh 50 mg of sample and add 1 mL of diluted sample buffer [see section 9.1]).
- Vortex the sample thoroughly.
- If necessary, centrifuge the sample: 5 minutes at $\geq 3500 \times g$ at room temperature (20°C to 25°C) (68° to 77° F).
- Use 50 µL of the supernatant per well in the test.

9.3 Surfaces and Devices

- Moisten a swab with diluted sample buffer (see section 9.1).
- Using a crisscross motion, thoroughly wipe a surface area or relevant device area of approximately 10×10 cm (4 x 4 in).
- Transfer the sample into 1 mL of diluted sample buffer (see section 9.1) by pressing and rotating the swab against the inner wall of the tube.
- Vortex the sample thoroughly.
- Apply a 50 µL sample for each well in the test.

9.4 Sample Storage:

Store samples suspended in sample buffers from 2°C to 8°C (35° to 46° F). for up to 2 days.

Before testing, bring the samples to room temperature (20°C to 25°C) (68 to 77° F) and briefly mix by shaking.

10. Test Execution

10.1 Test Preparation

Bring all reagents to room temperature (20°C to 25°C) (68° to 77 °F) before use.

A PBS-Tween buffer is required as the **wash buffer**. Please use the enclosed buffer envelope (See Section 4). To prepare the buffer, dissolve the entire contents of the envelope in 1 liter of distilled water. Once prepared, the wash buffer can be stored at 2°C to 8°C (35° to 46° F).

for approximately 4 to 6 weeks.

Alternative:

- Dissolve the contents of the sachet in 100 ml of distilled water to create a tenfold concentrate. To prepare the ready-to-use solution, mix 1 part of the tenfold concentrate with 9 parts of distilled water.
- The tenfold concentrate can be stored at room temperature (20°C to 25°C) (68° to 77° F) for approximately 8 to 12 weeks.

Immediately return any unused reagents to storage from 2°C to 8°C (35° to 46°F)

10.2 Test Execution

Avoid drying out the wells between steps. Consistent and thorough washing is critical for reliable test results. Always follow the specified washing sequences to ensure uniform treatment across all wells.

Use a multichannel or multistepper pipette when dispensing the conjugate, substrate/chromogen, and stop solution. This helps prevent timing delays across the microtiter plate.

During all incubation steps, keep the microtiter plate covered and protect it from direct sunlight.

1. Insert as many wells into the holder frame as needed for all standards and samples. Record the positions of the standards and samples. No more than 4 strips of 8 wells each (i.e. 28 samples plus 4 standards) should be processed at the same time.
2. Pipette 50µl of each standard solution or prepared sample into the designated well.

3. Add 50 µL of conjugate (Red Cap) to each well. Mix gently by shaking the plate manually and incubate for 10 minutes (± 1 minute) at room temperature (20° to 25°C (68° to 77° F)).
4. Empty the wells by firmly tapping out the liquid, then remove residual fluid by striking the plate three times onto absorbent laboratory paper. Fill each well with 250 µL of diluted wash buffer (see Section 10.1), then empty as described above. Repeat this step two more times (total of three wash cycles).
5. Add 100 µL of substrate/chromogen solution (Brown Cap) to each well and incubate in the dark for 5 minutes (± 1 minute) at room temperature (20° to 25°C (68° to 77° F)).
6. Add 100 µL of stop reagent (Yellow Cap) to each well and mix gently by shaking the plate manually. Measure absorbance at 450 nm within 15 minutes after adding the stop reagent.

11. Interpretation

Color intensity of selected Standards serves as interpretive reference.

R-Biopharm offers an optional interpretation software tool specifically developed for the RIDASCREEN® enzyme immunoassays:

RIDASOFT® Win.NET Food & Feed (Art. Nr. Z9996FF).

Use linear regression to evaluate the results.

Verify that the test run meets all quality criteria from the certificate of analysis before evaluating. The standard curve shape is also provided in the certificate.

Graphical Interpretation:

Extinction values for the Standards are plotted on a standard graph paper coordinate system, using extinction values on the y-axis and risk-tissue concentrations [%] on the x-axis. A calibration line through the standard data points is then applied. The CNS tissue concentration (%) for each sample can then be determined by locating its extinction value on the calibration curve and finding the corresponding concentration.

12. Interpretation of results

If a result is below the LoD, report it as “<Detection limit (%)” since results below the LoD cannot categorically rule out all presence of CNS tissue.

13. Method Limitations

Follow all instructions (especially for sample prep, test execution, and storage) to achieve the specified detection limit.

Detection results and limits are influenced by several factors, including the sample matrix, laboratory conditions, level of processing, and extraction method.

This test is designed to detect central nervous system (CNS) tissue by identifying the Glial Fibrillary Acidic (GFAP) marker protein.

It is important to note that the test is calibrated specifically to bovine CNS (brain) tissue. When interpreting results, consider that marker protein concentrations may differ across species. Marker protein concentrations may naturally vary within the same tissue type and between different types of CNS tissue (e.g. brain, spinal cord, eye), even among animals of the same species. Such variability can be significant.

It may lead to differing levels of marker protein in the sample, which in turn can result in positive test outcomes at varying concentrations of CNS tissue.

To ensure accurate interpretation of results, always account for both species-specific and tissue-specific variability in marker protein content.

14. Recommendation

For reliable results:

- Pre-rinse pipette tips with the standard solution or sample extract before use
- Include control tests in routine quality checks to validate results and detect deviations.

15. Additional Applications

For additional details regarding the product and its uses, please contact: info@r-biopharm.de.

Version overview

Version number	Chapter and title
2010-07-01	Release version
2025-11-21	Current version Changes made: Comprehensive linguistic and formal revision Added and corrected product listings related to the detection of central nervous system (CNS) tissue

Explanation of symbols

General symbols:

-  Follow the instructions for use
-  Batch number
-  Expiry date (YYYY-MM-DD)
-  Storage temperature
-  Article number
-  Number of test determinations
-  Manufacturing date (YYYY-MM-DD)
-  Manufacturer + address

7. Disclaimer

1. R-Biopharm AG provides a warranty that its products meet the specified requirements and are free from material and workmanship defects for 12 months following the transfer of risk. For products with limited shelf life, the warranty remains valid until the stated expiration date or, for limited-use items, until the stated number of uses has been reached. Warranty claims must be promptly and correctly submitted by the customer. A product is considered defective if it does not have the agreed features, is not suitable for its intended contractual use, or is not provided with the required accessories and instructions (collectively "Subjective Requirements").
2. R-Biopharm AG assumes no liability for any consequences arising from the following:
 - (a) Failure to read, comprehend, or adhere to product usage or safety instructions;
 - (b) Employment of individuals who do not possess the necessary training or qualifications to use the product;
 - (c) Non-compliance with relevant industry standards and practices, including Good Laboratory Practices;
 - (d) Utilization of control materials, samples, sample matrices, or processing that are not appropriate for use with the product;
 - (e) Any other improper use of the product;
 - (f) Modifications or alterations of the product;
 - (g) Improper storage by the customer or others;
 - (h) Exposure to chemical, electromagnetic, mechanical, or electrolytic influences beyond standard ranges documented by R-Biopharm AG;
 - (i) Damage or malfunctions resulting from external events beyond the control of R-Biopharm AG, such as burglary, theft, lightning, fire, water damage, or Force Majeure.
3. R-Biopharm AG is liable in cases of fraud, gross negligence, or intent, for injury to life, body, or health, for the assumption of a guarantee or procurement risk under Section 276 of the German Civil Code (BGB), or where liability is mandated by law.

For slight negligence, R-Biopharm AG is only liable for breaches of material contractual obligations necessary for the contract's purpose and relied on by the customer. In these cases, liability is limited to typical, foreseeable damage. R-Biopharm AG assumes no liability for other minor negligent breaches.

ALL OTHER EXPRESS OR IMPLIED WARRANTIES OR GUARANTEES OF ANY KIND ARE EXCLUDED, REGARDLESS OF WHETHER THEY ARE ASSUMED THROUGH CUSTOM, PRACTICE, COURSE OF DEALING, OR OTHERWISE.

R-Biopharm AG accepts no liability for consequential damages, in particular lost profits, production delays, or other indirect losses.

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© 2025 R-Biopharm AG
are registered trademarks of R-Biopharm AG.