RIDASCREEN®FAST Gliadin

REF R7002

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RIDASCREEN®FAST Gliadin

Brief information

RIDASCREEN®FAST Gliadin (Art. No. R7002) is a R5 sandwich enzyme immunoassay for the quantitative analysis of prolamins from wheat (gliadin), rye (secalin) and barley (hordein) in food validated for the method (see chapter 1).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction and centrifugation

Time requirement: sample preparation
Cocktail (patented) (for 10 samples)2 h
Cocktail ECO (for 10 samples) approx. 35 min
test implementation (incubation time) 30 h

Standard material: The RIDASCREEN® standard material is the standard of the Prolamin Working Group (PWG-Gliadin, <http://www.wgpat.com/handling.html>).

Limit of detection: 0.5 mg/kg (ppm) gliadin or to 1 mg/kg (ppm) gluten
(matrix-dependent)

Limit of quantification: 5 mg/kg (ppm) gliadin or to 10 mg/kg (ppm) gluten

Specificity: The monoclonal antibody R5 reacts with the gliadin-fractions from wheat and corresponding prolamins from rye and barley. Further information are contained in the validation report.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice manual. It lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA-analysis. The brochure can be retrieved, printed and downloaded from the website:

<https://food.r-biopharm.com/media/technical-guides/>

Related products and accessories for gluten/gliadin determination

RIDASCREEN® Gliadin (Art. No. R7001)
RIDASCREEN® Gliadin competitive (Art. No. R7021)
RIDASCREEN® Total Gluten (Art. No. R7041)
RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051)
RIDA®QUICK Gliadin (Art. No. R7003 / R7004 / R7005)
Cocktail (patented) (Art. No. R7006 / R7016)
Cocktail ECO (Art. No. R7080)
RIDA® Extraction Solution (colorless) (Art. No. R7098)
Set of 3 processed Gliadin Assay Controls (Art. No. R7012)
SureFood® ALLERGEN Gluten (Art. No. S3606)
SureFood® QUANTARD Allergen 40 (Art. No. S3301)
SureFood® ALLERGEN 4plex Cereals (Art. No. S7006)

1. Intended use

RIDASCREEN®FAST Gliadin (Art. No. R7002) is a R5 sandwich enzyme immunoassay for the quantitative determination of contaminations in gluten-free declared foods by prolamins from wheat (gliadin), rye (secalin) and barley (hordein). Due to the large number of different food products, the following samples were examined representatively for different product groups within the scope of the test development: buckwheat flour, butter cookie and rice flour. It can be assumed that the test is also suitable for the analysis of other foodstuffs; this is to be verified by the user himself.

For detailed results and further information on validation data with other food matrices please refer to the Validation Report. Further applications are regularly validated in our laboratories, which we make available in our Application Notes (see chapter 15).

The sample preparation using the Cocktail (patented) (Art. No. R7006 / R7016) is the official R5-Mendéz method according to the Codex Alimentarius and the AOAC.

The faster sample preparation using the environmental-friendlier Cocktail ECO (Art. No. R7080) is convenient for the screening of samples. The Cocktail ECO has an extraction efficiency of approx. 70 - 110% compared to Cocktail (patented).

2. General

The use of wheat flour and gluten in foodstuffs is extremely common because of their heat stability and useful effects on e.g. texture, moisture retention and flavour. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley. Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (CODEX STAN 118/1979) two categories for labeling of food according to the gluten content now exist:

- 1) Food products which contain less than 20 mg/kg can be labeled as "**gluten-free**".
- 2) Food products labeled as "**very low gluten**" can have a gluten content above 20 and up to 100 mg/kg.

The threshold of 20 mg/kg has been adopted by many national legislations in many countries. The prolamin content (e.g. gliadin) of gluten is per definition generally assumed to be 50 % (CODEX STAN 118-1979).

3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific R5 antibodies against gliadins. By adding the standard or sample solution to the wells, gliadin present in the solutions will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Following the washing step, a solution containing R5 antibodies conjugated to peroxidase is added. This conjugate binds to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product.

The addition of the stop solution leads to a color change from blue to yellow. The absorbance of the solution, which is proportional to the gliadin concentration in the sample, is measured photometrically at 450 nm and expressed as ng/mL (ppb) gliadin.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		48 wells
Buffer	White	Concentrate	5x	60 mL
Standard 1	Transparent	Ready to use	0 ng/mL	1.3 mL
Standard 2	Transparent	Ready to use	10.0 ng/mL	1.3 mL
Standard 3	Transparent	Ready to use	20.0 ng/mL	1.3 mL
Standard 4	Transparent	Ready to use	40.0 ng/mL	1.3 mL
Standard 5	Transparent	Ready to use	80.0 ng/mL	1.3 mL
Wash buffer	Brown	Concentrate	10x	100 mL
Conjugate	Red	Concentrate	11x	0.7 mL
Substrate/Chromogen Red chromogen Pro	Brown	Ready to use		10 mL
Stop solution	Yellow	Ready to use		14 mL

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 \times g) + centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner, Art. No. 227261)
- Shaker
- Water bath (50 °C / 122 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12 μ m)
- Graduated pipettes
- Graduated cylinder
- Variable 20 - 200 μ L and 200 - 1000 μ L micropipettes

- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific, Art. No. 95029390 or low binding Greiner bio-one, Art. No. 655901)
- If necessary: 8-channel pipette for 100 µL
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)

5.2 Reagents

- Distilled water (dist. water) or deionized water
- Gluten-free skim milk powder (SMP) (food quality)
- Cocktail (patented) (Art. No. R7006 / R7016) or Cocktail ECO (Art. No. R7080)
- Ethanol solution (**80 %**): i.e. add 120 mL ethanol p.a. to 30 mL dist. water and shake well

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the plate (coated plate and pre-plate, if necessary (see chapter 10.2)). Use separate pipet tips for each standard and each sample extract to avoid cross contamination.

The Cocktail (patented) is harmful to health. It contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 47 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (36 - 47 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation.
- Value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 1.2$) for standard 5.

9. Sample preparation

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol.
- Carry out the sample preparation in a room isolated from the ELISA procedure.
- Check for gliadin contamination of reagents and equipment with the test strips RIDA®QUICK Gliadin (Art. No. R7003 / R7004 / R7005).
- When using the Cocktail (patented), it is recommended to work under a chemical hood, because it contains β -mercaptoethanol.
- β -Mercaptoethanol can disturb the ELISA, therefore dilute the samples at least 1:500.
- The samples should be stored in a cool place, protected against light.

9.1 Extraction with Cocktail (patented) (R7006 / R7016, official AOAC method)

Homogenize (grind thoroughly to powder and mix well or mix well a solution respectively) well a sufficient amount (e.g. 50 g or 50 mL) to ensure taking a representative test portion of sample.

- **Liquid food samples:** use 0.25 mL of the homogenized sample and add 2.5 mL of the Cocktail (patented), close the vial and mix well.
- **Other food samples (e.g. soy and quinoa containing samples):** weigh 0.25 g of the homogenized and add 2.5 mL of the Cocktail (patented), close the vial and mix well.
- **Tannin and polyphenol containing food samples (e.g. chocolate, coffee, cocoa, chestnut flour, buckwheat, millet, and spices):** weigh 0.25 g of the homogenized sample, add 0.25 g of SMP, and add 2.5 mL of the Cocktail (patented), close the vial and mix well.
- **Meat and sausages:** in these matrices gliadin may be not distributed evenly, therefore, weigh 50 g sample and homogenize: weigh 0.25 g of the homogenized sample and add 2.5 mL of the Cocktail (patented), close the vial and mix well.
- **Oat samples:** gliadin may not be distributed evenly, furthermore the samples are difficult to homogenize. Therefore, homogenize 200 g, then carry out the extraction with at least the fourfold amount of reagents: weigh 1 g of the homogenized sample and add 10 mL of the Cocktail (patented), close the vial and mix well.

Please further extract all samples as described in the following:

- Extract for 40 min at 50 °C (122 °F) in a water bath.
- Let the sample cool down shortly (3 - 5 min).
- Add 7.5 mL 80 % ethanol (see chapter 5.2) (with 1 g sample weight (e.g. oat samples): 30 mL 80 % ethanol), close the vial and mix well.
- Shake for 1 h upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter sample or centrifuge for 10 min at > 2,500 x g at room temperature (20 - 25 °C / 68 - 77 °F).
(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The sample extracts (supernatant from centrifugation step or filtrate) can be stored undiluted in a well-sealed container at room temperature (20 - 25 °C / 68 - 77 °F) in the dark (shelf life approx. 4 weeks) until used in the test.

- For use in the test, the extracts (supernatant of centrifugation step or filtrate) must be diluted with buffer (for dilution see chapter 10.2). The diluted sample extracts have a limited shelf life and must be used in the test within 30 minutes.

9.2 Extraction with Cocktail ECO (Art. No. R7080) for processed food

The faster sample preparation using the environmental-friendlier Cocktail ECO (Art. No. R7080) is convenient for the screening of samples. The Cocktail ECO has an extraction efficiency of approx. 70 - 110 % compared to Cocktail (patented). For this fast and environmental-friendlier sample preparation, it is not necessary to work under a chemical hood.

Prepare the necessary amount of Cocktail ECO according to the product information.

Homogenize (grind thoroughly to powder and mix well or mix well a solution respectively) well a sufficient amount (e.g. 50 g or 50 mL) to ensure taking a representative test portion of sample.

- **Liquid food samples:** use 0.25 mL of the homogenized sample and add 2.5 ml of Cocktail ECO, close the vial and mix well.
- **Other food samples (e.g. soy and quinoa containing samples):** weigh 0.25 g of the homogenized and add 2.5 mL of Cocktail ECO, close the vial and mix well.
- **Tannin and polyphenol containing food samples (e.g. chocolate, coffee, cocoa, chestnut flour, buckwheat, millet, and spices):** weigh 0.25 g of the homogenized sample, add 0.25 g of SMP, and add 2.5 mL of Cocktail ECO, close the vial and mix well.
- **Meat and sausages:** in these matrices, gliadin may be not distributed evenly, therefore, weigh 50 g sample and homogenize: weigh 0.25 g of the homogenized sample and add 2.5 mL of Cocktail ECO, close the vial and mix well.
- **Oat samples:** gliadin may not be distributed evenly, furthermore the samples are difficult to homogenize. Therefore, homogenize 200 g, then carry out the extraction with at least the fourfold amount of reagents: weigh 1 g of the homogenized sample and add 10 mL of Cocktail ECO, close the vial and mix well.

Please further extract all samples as described in the following:

- Extract for 10 min at 50 °C (122 °F) in a water bath.
- Let the sample cool down shortly (3 - 5 min).
- Add 7.5 mL 80 % ethanol (see chapter 5.2) (with 1 g sample weight (e.g. oat samples): 30 mL 80 % ethanol), close the vial and mix well.
- Shake for 10 min upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter sample or centrifuge for 5 min at > 2,500 x g at room temperature (20 - 25 °C / 68 - 77 °F).
(Alternatively: transfer 2 ml of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 5 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The sample extracts (supernatant from centrifugation step or filtrate) can be stored undiluted in a well-sealed container at room temperature (20 - 25 °C / 68 - 77 °F) in the dark (shelf life approx. 2 weeks) until used in the test.
- For use in the test, the extracts (supernatant of centrifugation step or filtrate) must be diluted with buffer (for dilution see chapter 10.2). The diluted sample extracts have a limited shelf life and must be used in the test within 30 minutes.

10. Test implementation

10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **buffer** is provided as a 5fold concentrate. Only the amount which actually is needed should be diluted 1:5 (1+4) with dist. water (e.g. 12 mL dist. water + 3 mL concentrate, sufficient for the dilution of 10 samples). The diluted buffer is stable for 1 day. Make sure that the buffer is not contaminated with gliadin.

The **conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) with dist. water (e.g. 2 mL distilled water + 200 µL concentrate, sufficient for 2 microtiter strips). Take care that the water is not contaminated with gliadin.

The **wash buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with dist. water (i.e. (e.g. 900 mL dist. water + 100 mL buffer concentrate). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

Components should be stored at 2 - 8 °C (36 - 47 °F) when no longer required.

10.2 Test procedure

The extracts prepared according to chapter 9.1 and 9.1 must be diluted 1:12.5 (1+11.5) with diluted buffer (see chapter 10.1) before usage (e.g. 920 µL buffer + 80 µL sample). The final dilution factor is 500.

Use the diluted sample extracts immediately (within 30 minutes) in the assay. A longer time period may influence the recovery.

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µL per well) and then exactly 100 µL are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µL of each standard or sample (prepared according to chapter 9 and diluted according to chapter 10.2) in duplicate to the wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid again. Repeat two more times (a total of three wash cycles).
4. Add 100 µL of the diluted conjugate (see chapter 10.1) to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).

5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 μ L diluted wash buffer (see chapter 10.1) and pour out the liquid again. Repeat two more times (a total of three wash cycles).
6. Add 100 μ L of substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 μ L of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

11. Results

A special software, the **RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. It is recommended to use RIDASOFT® Win.NET and cubic spline curve fitting procedure.

For the evaluation it should be clarified, that the quality criteria are fulfilled for the current test run. The course of the standard curve can be taken from the attached analysis certificate.

The assay can be also evaluated when running in single determinations. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not consistent with the AOAC Official Method of Analysis and standards like EN 15633-1 and EN 15842. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

The gliadin concentration in ng/mL (ppb) is read from the RIDASOFT® Win.NET calibration curve and must be further multiplied by the dilution factor of at least 500. This result is then multiplied by 2 in order to obtain the gluten concentration (gliadin represents 50 % of the proteins present in gluten, Codex Definition). The RIDASOFT® Win.NET (version 1.93 or newer) indicates the results in gliadin and gluten.

Calculation example

The absorbance value of a sample corresponds to 10 ng/mL gliadin in the calibration curve. Multiplying by the recommended dilution factor 500 leads to 5,000 ng/mL, corresponding to 5 mg/kg (ppm) gliadin, respectively 0.0005 % gliadin. To calculate the gluten content, it is necessary to multiply by factor 2 which results in 10 mg/kg gluten, respectively 0.001 % gluten. This sample is considered to be gluten-free, because the gluten concentration is below 20 mg/kg.

12. Interpretation of results

Results between LoD and LoQ indicate a low analyte concentration in the sample. Calculated results show a high uncertainty in this area due to the method's high variation below LoQ. Therefore, such results should not be reported with a quantitative value, but qualitative as "< LoQ".

A result below the LoD does not exclude an allergen contamination below the detection limit of the assay or that other allergen components, such as lipids, may be present in a sample. The result should be reported accordingly.

If samples show a higher absorbance than standard 5, the samples can be diluted more strongly and used again in the ELISA. For this purpose, dilute the extract again freshly 1:12.5 with buffer as described in chapter 9.1 or 9.2. Then further dilute the sample extract diluted in this way with the following mixture:

- 2 % Cocktail (patented) or Cocktail ECO (depending on extraction method)
- 6 % of an 80% ethanol solution
- 92 % buffer

This corresponds e.g. to 100 µL Cocktail (patented) or Cocktail ECO, 300 µL 80% ethanol and 4,600 µL buffer. When using this mixture, the composition of the sample medium is maintained according to the extraction in chapter 9.1 or 9.2. In case of a further dilution, the additional dilution factor must be taken into account when calculating the allergen concentration.

Compared to the certificate, higher absorbance values ($A_{450\text{ nm}}$) for the standard curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

Technical limits of the test method are approached outside the designated measurement range. Higher result variation appear typically with repeated testing. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LoD, LoQ, upper limit of measurement range).

An incorrect weigh in of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max $\pm 1\%$.

For detailed results and further information for other food matrices, please refer to the validation report. In addition, data on individual foods may be available from comparative laboratory tests and interlaboratory comparisons.

For the present ELISA, only individual, exemplary foods from different product groups could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

The addition of foreign protein (depending on the test e.g. BSA, gelatine, skim milk powder) during extraction or test procedure may suppress matrix effects.

In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery.

Cross reactivities are side reactions of the used antibody with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

The protein content and protein composition of wheat, rye and barley cultivars may differ. Varying results are thus to be expected for different cultivars.

The main epitope of the R5 antibody is the amino acid sequence QQPFP which is part of many celiac-toxic sequences. The sequence occurs repeatedly in the prolamins from wheat, rye and barley. However, rye and barley contain a higher number of replicates of this sequence, which leads to an overestimation of rye and barley compared to the wheat standard.

Recent research has shown that the real conversion factor from gliadin to gluten is approx. 1.5 ^{[1], [2], [3]}. However, the official conversion factor (Codex Alimentarius) is still 2.

14. Recommendation

Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance, we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control. Gluten-free and gluten-containing (spiked) samples should be used (e.g. Art. No. R7012).
- To do spike experiments to ensure an accurate and correct test procedure. An example of a spike experiment is given in the validation report.
- In case of extremely acidic or basic samples, adjustment of the sample's pH value to neutral (pH 6.5 to 7.5) prior to extraction may be necessary.
- To perform a PCR (e.g. SureFood[®]) for confirmation of the result.

- During the production of foods such as beer or sourdough, proteins are fragmented. In sandwich ELISAs **protein fragments** lead to a reduced recovery, such samples should be analyzed with a competitive ELISA test systems like the **RIDASCREEN® Gliadin competitive** (Art. No. R7021).
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt® / Bolt™) are used.

15. Further application notes

- Sample preparation for processed food with the RIDA® Extraction Solution (colorless) (Art. No. R7098) - **only after validation**
- Sample preparation for raw materials with ethanol
- Sample preparation for polyphenol containing raw materials (e.g. chocolate, coffee, cacao, buckwheat) with fish gelatine and ethanol
- RIDASCREEN®FAST Allergen - Swabbing method for the qualitative analysis of allergens in the production line or for laboratory equipment.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Literature

- [1] Lacorn et al. (2021) Determination of Gliadin as a Measure of Gluten in Food by R5 sandwich ELISA RIDASCREEN® Gliadin Matrix Extension: Collaborative Study 2012.01, J. AOAC Int (publication in progress).
- [2] Schall, E., Scherf, K.A., Bugyi, Z., Hajas, L., Török, K., Koehler, P., Poms, R.E., D'Amico, S., Schoenlechner, R., & Tömösközi, S. (2020) Food Chem. 313, 126049. <https://doi.org/10.1016/j.foodchem.2019.126049>.
- [3] Wieser et al.: Celiac Disease and Gluten. (2014) Elsevier Inc. Amsterdam, ISBN 978-0-12- 420220-7, page 107. [4] Wieser, H. & Koehler, P. (2009) Eur. Food Res. Technol. 229, 9-13.

Version overview

Version number	Chapter and title
2018-02-05	Previous version
2024-07-01	Current version Changes made: <ul style="list-style-type: none">– General linguistic revision– Updated list of related products for gluten / gliadin determination– Completion of additionally required reagents and accessories in chapter 5– Modification of the disposal clause in chapter 6.– Necessary dilution of the extracts for use in the test moved from chapter 9.1 or 9.2 to chapter 10.2– Revision of chapter 11– Addition of chapter 12, 13 and 14– Addition of a further application in chapter 15– Addition of version “Literature”, “Version overview” and “Explanation of symbols”– Addition of marking concerning patents of MORINAGA & Co., Ltd.– Updated disclaimer

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Patent Marking:

If the R-Biopharm product “Cocktail ECO” is used for the processing of food samples, a sulphite-containing extraction agent is used. Food inspection methods using a sulfite-containing extractant as in this product and/or corresponding detection kits are subject to the following patents of MORINAGA & Co., Ltd.: European Patent EP 2 224 239 B1, Australian Patent AU 2008 330 507 B2, United States Patent US 8 859 212 B2, Japanese Patent JP 5 451 854 B2. The patent holder has granted R-Biopharm AG a license to use, and sell products that employ, said protected technology in the above-mentioned territories.

Disclaimer

1. In conformance with the German Civil Code (“BGB”) R-Biopharm AG provides a limited warranty (“Gewährleistung”) against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm’s own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG’s liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG’s sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser’s expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements (“Gewährleistung”) implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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are registered trademarks of R-Biopharm AG.