RIDA® QUICK Gliadin

REF R7003

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Immuno-chromatographic test for the qualitative determination of gluten on surfaces, in cleansing/process water and in food

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AOAC *Official Methods*SM (2015.16)
for processed/non-processed corn products

AOAC *Performance Tested Methods*SM (101702)
for surfaces and cleansing waters

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Please contact for questions and further information:

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RIDA[®], RIDASCREEN[®] and RIDASOFT[®]
are registered trademarks of R-Biopharm AG.
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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Literatur

- [1] Lacorn M. et al.: The Validation of the RIDA®QUICK Gliadin for AOAC Research Institute. J AOAC Int. 2018; 101(5):1548-1557.
- [2] Lacorn M. et al.: Determination of Gluten in Processed and Nonprocessed Corn Products by Qualitative R5 Immunochromatographic Dipstick: Collaborative Study, First Action 2015.16., J AOAC Int. 2016; 99(3):730-737.
- [3] FAO Codex Alimentarius International Food Standards: Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten. CXS 118-1979. Adopted in 1979. Amended in 1983 and 2015. Revised in 2008.
- [4] FAO Codex Alimentarius International Food Standards: Recommended Methods of Analysis and Sampling. CXS 234-1999. Adopted in 1999. Amended in 2021. Pg. 21.

- [5] Valdés I. et al.: Innovative approach to low-level gluten determination in foods using a novel sandwich enzyme-linked immunosorbent assay protocol. Eur J Gastroenterol Hepatol. 2003; 15(5):465-74.
- [6] Immer U., Haas-Lauterbach S.: Gliadin as a measure of gluten in foods containing wheat, rye, and barley-enzyme immunoassay method based on a specific monoclonal antibody to the potentially celiac toxic amino acid prolamin sequences: collaborative study. J AOAC Int. 2012; 95(4):1118-1124.
- [7] Lacorn M. et al.: Determination of Gliadin as a Measure of Gluten in Food by R5 Sandwich ELISA RIDASCREEN® Gliadin Matrix Extension: Collaborative Study 2012.01. J AOAC Int. 2022; 105(2):442-455.
- [8] Garcia E. et al.: Development of a general procedure for complete extraction of gliadins for heat processed and unheated foods. Eur J Gastroenterol Hepatol. 2005 May; 17(5):529-39.

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RIDA®QUICK Gliadin

Brief information

RIDA®QUICK Gliadin (Art. No. R7003) is an immunochromatographic test for the qualitative determination of gluten-contamination on surfaces, in cleansing waters (CIP waters) and in foods.

The test has been approved as AOAC Proficiency Tested Method (PTM 101702) for swabbing and cleansing waters^[1] and as AOAC Official Method of Analysis (OMA 2015.16) for corn based food matrices using ethanol or Cocktail (patented) extraction^[2].

The test kit contains 25 test strips, each of which can be used for one analysis. All reagents required for the swab test are contained in the test kit. Results are evaluated visually.

Time requirement: sampling for swab test..... approx. 1 min
Sample preparation for:
10 CIP water samples* approx. 5 min
10 raw materials approx. 15 min
10 processed foods (R7006) approx. 120 min
10 processed foods (R7080) approx. 35 min
Test implementation (incubation time) 5 min
*CIP = Clean-In-Place (rinse water from cleaning process)

Detection limit: Surfaces approx. 1.6 - 3 µg gluten/100 cm²
Raw material approx. 4.4 mg/kg gluten*
Processed food approx. 6.3 mg/kg gluten**
CIP water (without cleansing agent) approx. 10 ng/mL gluten
CIP water (with cleansing agent) approx. 50-100 ng/mL gluten
*with ethanol extraction **with Cocktail (patented) extraction

Deviation may occur from the above mentioned values depending on investigated matrices and chosen extraction method.

Specificity: The **monoclonal antibody R5** reacts with the gliadin-fractions from wheat and corresponding prolamins from rye and barley.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments (see chapter 11. Limits of the method).

Additional information can be found in the product's validation report.

Related products and accessories for gluten / gliadin determination

RIDA®QUICK Gliadin (single packaged) (Art. No. R7004)

RIDA®QUICK Gliadin (ready to swab) (Art. No. R7005)

RIDASCREEN® Gliadin (Art. No. R7001)

RIDASCREEN®FAST Gliadin (Art. No. R7002)

RIDASCREEN® Gliadin competitive (Art. No. R7021)

RIDASCREEN® Total Gluten (Art. No. R7041)

RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051)

Cocktail (patented) (Art. No. R7006 / R7016)

Cocktail ECO (Art. No. R7080)

RIDA® Extraction solution (colorless) (Art. No. R7098)

Set of 3 processed Gliadin Assay Controls (Art. No. R7012)

SureFood® Allergen Gluten (Art. No. S3606)

SureFood® Allergen 4plex Cereals (Art. No. S7006)

1. Intended use

RIDA®QUICK Gliadin (Art. No. R7003) is an immunochromatographic assay in test strip format for the qualitative determination of gluten

- on surfaces (swab test for hygiene control in production and laboratories),
- in cleansing waters (CIP water),
- in foods after extraction with ethanol, with Cocktail (patented) or with Cocktail ECO.

The test has been developed for the detection of low amounts of gluten and detects gluten contaminations even below 10 mg/kg. Sensitivity of detection varies with investigated matrices and chosen extraction methods.

For the AOAC-RI validation of the assay for testing of surfaces and CIP waters^[1], stainless steel, sealed ceramic, plastic and silicone rubber were validated or waters without and with commercial cleansing agents respectively.

Corn-based raw materials with ethanol extraction and heat-processed corn-based foods with the Cocktail (patented) extraction are part of the AOAC Official Method of Analysis approval^[2].

It can be assumed that the assay is also suitable for the analysis of other foods. However, this must be verified by the user before applying the test kit to the respective food. For detailed results and further information on validation data with other food matrices please refer to our validation report. Further applications are regularly validated in our laboratories, which we publish in our application notes (see chapter 13. Further applications).

2. General information

The use of wheat flour and gluten in foodstuff is extremely common because of their useful effects on e.g. texture, moisture retention and flavor. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley. Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

The Codex Alimentarius Commission has stipulated the limit value for gluten-free food at 20 mg/kg gluten^[3]. This threshold has also been adopted by many national legislations. The prolamin content of gluten (e.g. gliadin in wheat) is set to be 50 %.

The official type I method for gluten determination of the Codex Alimentarius^[4] is an ELISA, which uses the R5 antibody in combination with a special extraction buffer, the Cocktail (patented) (Méndez method)^[5]. The RIDASCREEN® Gliadin ELISA (Art. No. R7001) fulfils this requirement and showed in an additional international collaborative study in 2020 its broad applicability for investigation of foods from different product groups and of different processing levels^[6, 7].

The RIDA®QUICK Gliadin test strips also use the R5 antibody and show a good correlation with the official method the RIDASCREEN® Gliadin ELISA. R-Biopharm AG is the only company that is allowed to use the R5 antibody in immunochromatographic test strips.

3. Test principle

The principle of the test in strip format is an antigen-antibody reaction based on the monoclonal R5-antibody, which is specific for the detection of gliadin from wheat and prolamins from rye and barley. In the lower area, the strip

contains red latex particle-R5-antibody-conjugates (latex conjugates). The test band contains immobilized R5 antibody. The strip is placed in a running buffer containing a prepared sample. Following the chromatographic principle, the liquid (running buffer and sample) moves upwards inside the strip. Latex conjugates are solved by the liquid and pass together with it the test band. If the analyte is present, a sandwich is formed consisting of immobilized R5 antibody, gliadin and red latex-conjugate coloring the test band red. Results are read visually.

4. Reagents provided

Each kit contains sufficient materials for max. 25 determinations. Each test kit contains:

Component	Cap color	Format	Volume
Test strip		Ready to use	25 pieces
Test tube		Ready to use	30 pieces
Disposable pipet			25 pieces
Buffer	Transparent	Ready to use	60 mL
Evaluation card			1 piece

5. Materials required but not provided

5.1 Equipment

- Gloves
- Scale
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Optional: filter paper (pore size 8 - 12 µm)
- Shaker
- Water bath (50 °C / 122 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Graduated pipettes
- Graduated cylinder

5.2 Reagents

Different reagents are recommended for analysis of raw materials and processed food (see chapter 9.5).

- Distilled water (dist. water) or deionized water
- Skimmed milk powder (SMP, food quality) for soy, tannin and polyphenol containing food
- 60 % ethanol solution (mix 150 mL ethanol p.a. with 100 mL dist. water)
- 80 % ethanol solution (mix 120 mL ethanol p.a. with 30 mL dist. water)
- Cocktail (patented) (R7006 / R7016)
- Alternatively: Cocktail ECO (R7080)

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

The Cocktail (patented) contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze the test strips or any test kit components.

Once the test strip container has been opened, store the container at room temperature (20 - 25 °C / 68 - 77 °F).

The test strips are sensitive to moisture and must be absolutely protected against it. Moist test strips can have a negative impact on the result. For this reason, remove the test strips from packaging only after having reached room temperature and immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

8. Indication of instability or deterioration of reagents

- A missing control band or an incomplete test or control band.
- A turbid or flaked buffer.

9. Test implementation

9.1 General information

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol.
- For the analysis, the extraction and the test procedure should be carried out in separate rooms.
- When using the Cocktail (patented), it is recommended to work under a chemical hood, because it contains β -Mercaptoethanol.

9.2 Swab test on surfaces

- Take as many test tubes as surfaces to be analyzed.
- Place 500 μ L of buffer in the test tube (e.g. using the disposable pipette provided).
- Swab the lower end (application area; see also chapter 10. fig. 2) of a dry test strip thoroughly in a cross-hatch technique over a sampling area of 10 x 10 cm (see fig. 1).


Fig. 1: Direct swabbing of surfaces with RIDA®QUICK Gliadin.

- Place the strip into the test tube with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).
- Take out the strip after exactly 5 min (\pm 10 s) and read the result using the evaluation card.

9.3 Analysis of CIP water without cleansing agent

- Take as many test tubes as samples to be analyzed.
- Place 250 μ L of buffer in the test tubes (e.g. using a disposable pipette provided).
- Place 250 μ L of cleansing water in the test tube (e.g. using a disposable pipette provided) and mix gently.
- Place a test strip into the test tube with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).
- Take out the strip after exactly 5 min (\pm 10 s) and read the result using the evaluation card.

9.4 Analysis of CIP water with cleansing agent

- Take as many test tubes as samples to be analyzed.
- Place 500 μ L of buffer in the test tube (e.g. using the disposable pipette provided).
- Place 50 μ L of cleansing water in the test tube (equivalent to 3 drops from the supplied disposable pipette held vertically) and mix gently.
- Place a test strip into the test tube with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).
- Take out the strip after exactly 5 min (\pm 10 s) and read the result using the evaluation card.

9.5 Extraction of foods

For AOAC OMA 2015.16 validation^[2] non-processed corn samples were analyzed using the ethanol extraction and processed corn samples using the Cocktail (patented) extraction.

Depending on the extraction method, the addition of skimmed milk powder (SMP) is necessary for some ingredients to avoid unwanted disturbing reactions.

Food ingredient	Ethanol extraction	Cocktail (patented) / Cocktail ECO extraction
Soy	1 g skimmed milk powder	–
Tannin- and polyphenol containing food (e.g. chocolate, coffee, cacao, chestnut flour, buckwheat, millet, spices)	1 g skimmed milk powder	0.25 g skimmed milk powder

9.5.1 Extraction with ethanol – *recommended for raw materials*

Homogenize well a sufficient amount of sample (grind at least 50 g thoroughly to powder and mix well or mix well 50 mL of a liquid sample respectively).

- **Liquid food samples:** pipet 1 mL of the homogenized sample into a vial (add 1 g of SMP in case of soy, tannin and polyphenol containing samples) and add 9 mL of 60 % ethanol; close the vial and mix well.
- **Soft and hard food samples:** weigh 1 g of the homogenized sample (add 1 g of SMP in case of soy, tannin and polyphenol containing samples) and add 10 mL of 60 % ethanol; close the vial and mix well.
- Shake well for at least 30 sec (e.g. vortex).
- Centrifuge: 10 min / at least 2,500 x g / room temperature (20 - 25 °C / 68 - 77 °F).
- Alternatively: let the sample settle down and / or filtrate.

Principally, ethanol extraction can also be used for fast analysis of other samples than raw materials. **However, this procedure is not part of the AOAC Official Method of Analysis.** Please notice, that the efficiency of ethanol extraction can be reduced especially with heat-processed samples^[5]. Hence, users have to verify suitability of ethanol extraction for the particular sample material. The group of Enrique Mendéz has developed the Cocktail (patented) for improved extraction of gluten from heat-processed samples^[8]. It can generally be used for all kind of foods. As a faster and more environmental-friendly alternative, R-Biopharm has developed the Cocktail ECO, which shows a comparable extraction efficiency.

9.5.2 Extraction with Cocktail (patented) (R7006 / R7016) – *recommended for heat-processed food and in case of unknown treatment*

Homogenize well a sufficient amount of sample (grind at least 50 g thoroughly to powder and mix well or mix well 50 mL of a liquid sample respectively).

- **Liquid food samples:** pipet 0.25 mL of the homogenized sample into a vial (add 0.25 g of SMP in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail (patented); close the vial and mix well.
- **Soft and hard food samples:** weigh 0.25 g of the homogenized sample (add 0.25 g of SMP in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail (patented); close the vial and mix well.
- **Meat and sausages:** uneven gluten distribution is typical for these matrices; therefore, ensure to homogenize a sufficient sample amount: weigh 0.25 g of the homogenized sample (add 0.25 g of SMP in case of tannin and polyphenol containing samples) into a vial and add 2.5 mL of Cocktail (patented); close the vial and mix well.
- **Oat samples:** uneven gluten distribution is typical for oats; furthermore, the samples are difficult to homogenize; therefore, homogenize 200 g and carry out the extraction with at least the fourfold amount of sample and reagents: weigh 1 g of the homogenized sample into a vial and add 10 mL of Cocktail (patented); close the vial and mix well.

Please further extract all Cocktail (patented) samples as described in the following:

- Incubate for 40 min at 50 °C (122 °F) in a water bath.
- Let the sample cool down shortly (1 - 3 min).
- Add 7.5 mL 80 % ethanol (see chapter 5.2) close the vial and mix well (for oat samples, add 30 mL 80 % ethanol accordingly).
- Shake for 1 h upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter sample or centrifuge for 10 min at > 2,500 x g at room temperature (20 - 25 °C) / 68 - 77 °F).

Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge.

- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.

- The undiluted sample extracts (supernatant from centrifugation step or filtrate) can be stored in a well-sealed tube at room temperature (20 - 25 °C / 68 - 77 °F) in the dark (shelf life approx. 4 weeks).

9.5.3 Extraction with Cocktail ECO (R7080) – *recommended for processed food and in case of unknown treatment*

For a faster sample preparation, the more environmental-friendly Cocktail ECO can be used as alternative. This procedure is not part of the AOAC Official Method of Analysis. The Cocktail ECO has an extraction efficiency of approx. 70 - 110 % compared to Cocktail (patented).

Homogenize well a sufficient amount of sample (grind at least 50 g thoroughly to powder and mix well or mix well 50 mL of a liquid sample respectively). Prepare the necessary amount of Cocktail ECO according to the corresponding product information.

- **Liquid food samples:** pipet 0.25 mL of the homogenized sample into a vial (add 0.25 g of SMP in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail ECO; close the vial and mix well.
- **Soft and hard food samples:** weigh 0.25 g of the homogenized sample (add 0.25 g of SMP in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail ECO; close the vial and mix well.
- **Meat and sausages:** uneven gluten distribution is typical for these matrices; therefore, ensure to homogenize a sufficient sample amount: weigh 0.25 g of the homogenized sample (add 0.25 g of SMP in case of tannin and polyphenol containing samples) into a vial and add 2.5 mL of the Cocktail ECO; close the vial and mix well.
- **Oat samples:** uneven gluten distribution is typical for oats; furthermore, the samples are difficult to homogenize; therefore, homogenize 200 g and carry out the extraction with at least the fourfold amount of sample and reagents: weigh 1 g of the homogenized sample into a vial and add 10 mL of Cocktail ECO; close the vial and mix well.

Please further extract all Cocktail ECO samples as described in the following:

- Incubate for 10 min at 50 °C (122 °F) in a water bath.
- Let the sample cool down shortly (1 - 3 min).
- Add 7.5 mL 80 % ethanol (see chapter 5.2) close the vial and mix well (for oat samples, add 30 mL 80 % ethanol accordingly).

- Shake for 10 min upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter sample or centrifuge for 5 min at > 2,500 x g at room temperature (20 - 25 °C) / 68 - 77 °F).
Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 5 min in a microcentrifuge.
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The undiluted sample extracts (supernatant from centrifugation step or filtrate) can be stored in a well-sealed tube at room temperature (20 - 25 °C / 68 - 77 °F) in the dark (shelf life approx. 2 weeks).

9.6 Test procedure for food analysis

- Take as many test tubes as samples to be analyzed.
- Place 500 µL of the buffer into each test tube (e.g. using the disposable pipettes provided).
- Add 50 µL of the extracts prepared according to chapter 9.5.1 - 9.5.3 to the test tubes (e.g. using the provided disposable pipettes, 50 µL correspond to 3 drops from a vertically hold pipette) and mix carefully.
- Place one test strip into the test tube with the arrow pointing down and take care that the strip immerses only with the application area (arrow mark on strip).
- Take out the strip after exactly 5 min (\pm 10 s) and read the result using the evaluation card.

10. Evaluation

The right line in the result area is the red test band (see fig. 2). The sample contains gluten, if the test band becomes visible.

The left line in the result area is the blue control band (see fig. 2). It must be present after each test procedure.


Fig. 2: RIDA®QUICK Gliadin test strip.

Positive result: two complete colored bands

The sample is positive if two complete colored bands (the blue control band and the red test band) are visible within the result area. In the case of swabbing, complete test bands with non-uniform intensity may occur due to an inhomogeneous gluten distribution on the surface or different swabbing procedures. Principally, the test band's red color is stronger with higher analyte concentrations.

Swab test:	> approx. 1.6 - 3 µg gluten/100 cm ²
Raw material:	> approx. 4.4 mg/kg gluten*
Processed food:	> approx. 6.3 mg/kg gluten**
CIP water (without cleansing agent):	> approx. 10 ng/mL gluten
CIP water (with cleansing agent):	> approx. 50 - 100 ng/mL gluten

*ethanol extraction

**Cocktail (patented) extraction

Negative result: only one complete blue band

The sample is negative if only the blue test band is visible within the result area.

Swab test:	< approx. 1.6 - 3 µg gluten/100 cm ²
Raw material:	< approx. 4.4 mg/kg gluten*
Processed food:	< approx. 6.3 mg/kg gluten**
CIP water (without cleansing agent):	< approx. 10 ng/mL gluten
CIP water (with cleansing agent):	< approx. 50 - 100 ng/mL gluten

*ethanol extraction

**Cocktail (patented) extraction

Invalid result: no control band or incomplete bands

If the blue control band is not visible after a run or an incomplete control or test band is visible, the test result is not valid and the test must be repeated with a new test strip. If the result is invalid again, check the expiry date of the test kit and control the content for visible sign of reagent deterioration (see chapter 8.).

11. Limits of the method

The test has been developed for the detection of low amounts of gluten (contamination). No high-dose-hook-effect is observed at high concentrations. However, the red test band may smear at very high gluten concentrations (> 10,000 mg/kg gluten).

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

The limit of detection depends on sample type and extraction efficiency or the properties of the swabbed surface and the kind of contamination respectively.

Inhomogeneous gluten distribution inside a sample and an inadequate sample size may cause a negative result.

Cleansing water containing hypochlorite cannot be analyzed. This cleansing agent destroys gluten very quickly by oxidation. The test strip is not able to detect potentially remaining gluten fragments.

A negative result does not exclude a gluten contamination below the detection limit of the assay, or that other cereal components, such as starch, may be present in a sample. The result should be reported accordingly.

It is recommended to use ethanol extraction with foods, which have not been heated. Extraction efficiency may be reduced with heat-processed samples. Kit users need to verify suitability of an extraction method for particular foods.

Furthermore, proteins may have been altered and / or fragmented in processed food (e.g. heat treatment, dehydration, etc.). This may have an impact on recovery or cross-reactivity.

For evaluation of the cross-reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross-reactivities and exemplary analyzed matrices are described in the validation report.

Cross reactivities are side reactions of the used antibody with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, food processing may induce cross-reactivity in single cases or cause its loss.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive results, but may also reduce or suppress a correct reaction. Such matrix effects (caused by e.g. polyphenols) are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

Addition of foreign protein (e.g. skimmed milk powder) during extraction or test procedure may suppress matrix effects.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and interlaboratory comparisons.

12. Recommendation

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842.
- To carry along test controls for quality control. Gluten-free and gluten-containing (spiked) samples should be used (e.g. R7012).
- To adjust the sample's pH value to neutral (pH 6.5 to 7.5) prior to extraction in case of extremely acidic or basic samples.
- To compare the extraction efficiency of ethanol and Cocktail ECO (R7080) with the Cocktail (patented) (R7006 / R7016).
- To do spike experiments to ensure an accurate and correct test procedure. An example of a spike experiment is given in the validation report.
- To clear unexpected results, use the RIDASCREEN® Gliadin ELISA (R7001) or run a PCR (e.g. SureFood®).
- To use the RIDASCREEN® Gliadin ELISA (R7001) for quantification of a positive result.

For storing the test strip for documentation, cut off the upper part of the strip marked with "Gluten" together with the test bands.

13. Further applications

- Sample preparation for processed food with the RIDA® Extraction solution (colorless) (R7098) - **only after validation**.
- Sample preparation for polyphenol containing raw material using fish gelatine.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Literature

- [1] Lacorn M. et al.: The Validation of the RIDA®QUICK Gliadin for AOAC Research Institute. J AOAC Int. 2018; 101(5):1548-1557.
- [2] Lacorn M. et al.: Determination of Gluten in Processed and Nonprocessed Corn Products by Qualitative R5 Immunochromatographic Dipstick: Collaborative Study, First Action 2015.16., J AOAC Int. 2016; 99(3):730-737.
- [3] FAO Codex Alimentarius International Food Standards: Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten. CXS 118-1979. Adopted in 1979. Amended in 1983 and 2015. Revised in 2008.
- [4] FAO Codex Alimentarius International Food Standards: Recommended Methods of Analysis and Sampling. CXS 234-1999. Adopted in 1999. Amended in 2021. Pg. 21.
- [5] Valdés I. et al.: Innovative approach to low-level gluten determination in foods using a novel sandwich enzyme-linked immunosorbent assay protocol. Eur J Gastroenterol Hepatol. 2003; 15(5):465-74.
- [6] Immer U., Haas-Lauterbach S.: Gliadin as a measure of gluten in foods containing wheat, rye, and barley-enzyme immunoassay method based on a specific monoclonal antibody to the potentially celiac toxic amino acid prolamin sequences: collaborative study. J AOAC Int. 2012; 95(4):1118-1124.
- [7] Lacorn M. et al.: Determination of Gliadin as a Measure of Gluten in Food by R5 Sandwich ELISA RIDASCREEN® Gliadin Matrix Extension: Collaborative Study 2012.01. J AOAC Int. 2022; 105(2):442-455.
- [8] Garcia E. et al.: Development of a general procedure for complete extraction of gliadins for heat processed and unheated foods. Eur J Gastroenterol Hepatol. 2005 May; 17(5):529-39.

Version overview

Version number	Chapter and title
2018-01-16	Previous version
2023-02-23	Current version Changes made: <ul style="list-style-type: none">– General linguistic revision– Modification of the disposal clause in chapter 6– Chapter 8. Indication of instability or deterioration of reagents (new)– Chapter 9.5 specification of usability of different extraction methods– Revision of chapter 10., 11. and 12.– Chapter literature, version overview, explanation of symbols (new)– Disclaimer renewed

Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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