## RIDA®QUICK Gliadin (single packaged)

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Immunochemical test for the qualitative determination of gluten on surfaces, in cleansing/process water and in food

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AOAC *Official Methods*<sup>SM</sup> (2015.16)  
for processed/non-processed corn products

AOAC *Performance Tested Methods*<sup>SM</sup> (101702)  
for surfaces and cleansing waters

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Please contact for questions and further information:

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RIDA<sup>®</sup>, RIDASCREEN<sup>®</sup> and RIDASOFT<sup>®</sup>  
are registered trademarks of R-Biopharm AG.  
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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### Literatur

- [1] Lacorn M. et al.: The Validation of the RIDA®QUICK Gliadin for AOAC Research Institute. J AOAC Int. 2018; 101(5):1548-1557
- [2] Lacorn M. et al.: Determination of Gluten in Processed and Nonprocessed Corn Products by Qualitative R5 Immunochemical Dipstick: Collaborative Study, First Action 2015.16., J AOAC Int. 2016; 99(3):730-737
- [3] FAO Codex Alimentarius International Food Standards: Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten. CXS 118-1979. Adopted in 1979. Amended in 1983 and 2015. Revised in 2008.
- [4] FAO Codex Alimentarius International Food Standards: Recommended Methods of Analysis and Sampling. CXS 234-1999. Adopted in 1999. Amended in 2021. Pg. 21
- [5] Valdés I. et al.: Innovative approach to low-level gluten determination in foods using a novel sandwich enzyme-linked immunosorbent assay protocol. Eur J Gastroenterol Hepatol. 2003; 15(5):465-74
- [6] Immer U., Haas-Lauterbach S.: Gliadin as a measure of gluten in foods containing wheat, rye, and barley-enzyme immunoassay method based on a specific monoclonal antibody to the potentially celiac toxic amino acid prolamin sequences: collaborative study. J AOAC Int. 2012; 95(4):1118-1124
- [7] Lacorn M. et al.: Determination of Gliadin as a Measure of Gluten in Food by R5 Sandwich ELISA RIDASCREENVR Gliadin Matrix Extension: Collaborative Study 2012.01. J AOAC Int. 2022; 105(2):442-455
- [8] Garcia E. et al.: Development of a general procedure for complete extraction of gliadins for heat processed and unheated foods. Eur J Gastroenterol Hepatol. 2005 May; 17(5):529-39

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# RIDA®QUICK Gliadin (single packaged)

## Brief information

RIDA®QUICK Gliadin (single packaged) (Art. No. R7004) is an R5-based immunochromatographic test for the qualitative determination of gluten-contamination on surfaces, in cleansing waters (CIP waters) and in foods.

### The identical test strip RIDA®QUICK Gliadin (Art. No. R7003)

- has AOAC-PTM (101702) for swabbing and cleansing waters<sup>[1]</sup>.
- has been approved as AOAC-OMA 2015.16 for corn based food using Cocktail (patented) or ethanol extraction<sup>[2]</sup>.

The test kit contains 25 test strips (single packaged), each of which can be used for one analysis. All reagents required for the swab test are contained in the test kit. Results are evaluated visually.

Time requirement:    sampling for swab test ..... ca. 1 min  
                                  Sample preparation for:  
                                  10 CIP water samples ..... ca. 5 min  
                                  10 raw materials..... ca. 15 min  
                                  10 processed foods (R7006) ..... ca. 120 min  
                                  10 processed foods (R7080) ..... ca. 35 min  
                                  Test implementation (incubation time)..... 5 min  
CIP = Clean-In-Place (rinse water from cleaning process)

Detection Limit:        Surfaces..... ca. 1.6 - 3 µg gluten/100 cm<sup>2</sup>  
(Matrix-dependent)    Raw material ..... ca. 4.4 mg/kg gluten\*  
                                  Processed food ..... ca. 6.3 mg/kg gluten\*\*  
                                  CIP water (without cleansing agent).....ca. 10 ng/mL gluten  
                                  CIP water (with cleansing agent)... ca. 50-100 ng/mL gluten  
\* with ethanol extraction    \*\* with Cocktail (patented) extraction

Deviation may occur from the above mentioned values depending on investigated matrices and chosen extraction method.

Specificity:                The **monoclonal antibody R5** reacts with the gliadin-fractions from wheat and corresponding prolamins from rye and barley. Further information is contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed/processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments.

## **Related products and accessories for gluten/gliadin determination**

RIDA®QUICK Gliadin (Art. No. R7003)  
RIDA®QUICK Gliadin (ready to swab) (Art. No. R7005)  
RIDASCREEN® Gliadin (Art. No. R7001)  
RIDASCREEN®FAST Gliadin (Art. No. R7002)  
RIDASCREEN® Gliadin competitive (Art. No. R7021)  
RIDASCREEN® Total Gluten (Art. No. R7041)  
RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051)  
Cocktail (patented) (Art. No. R7006 / R7016)  
Cocktail ECO (Art. No. R7080)  
RIDA® Extraction Solution (colorless) (Art. No. R7098)  
Set of 3 processed Gliadin Assay Controls (Art. No. R7012)  
SureFood® ALLERGEN Gluten (Art. No. S3606)  
SureFood® QUANTARD Allergen 40 (Art. No. S3301)  
SureFood® ALLERGEN 4plex Cereals (Art. No. S7006)

### **1. Intended use**

RIDA®QUICK Gliadin (single packaged) (Art. No. R7004) is an R5-based immunochromatographic assay in test strip format for the qualitative determination of gluten

- on surfaces (swab test for hygiene control in production and laboratories)
- in cleansing waters (CIP water).
- in foods after extraction with ethanol, with Cocktail (patented) or with Cocktail ECO.

The test has been developed for the detection of low amounts of gluten and detects gluten contaminations even below 10 mg/kg. Sensitivity of detection varies with investigated matrices and chosen extraction methods.

For the AOAC-PTM validation of the assay for testing of surfaces and CIP waters<sup>[1]</sup>, stainless steel, sealed ceramic, plastic and silicone rubber were validated or waters without and with commercial cleansing agents respectively.

Corn-based raw materials with ethanol extraction and heat-processed corn-based foods with the Cocktail (patented) extraction are part of the AOAC Official Method of Analysis approval<sup>[2]</sup>.

It can be assumed that the assay is also suitable for the analysis of other foods. However, this must be verified by the user before applying the test kit to the respective food. For detailed results and further information on validation data with other food matrices please refer to our validation report. Further applications are regularly validated in our laboratories, which we publish in our Application Notes (see chapter 13).

## 2. General information

The use of wheat flour and gluten in foodstuff is extremely common because of their useful effects on e.g. texture, moisture retention and flavor. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley. Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

The Codex Alimentarius Commission has stipulated in the “Codex Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten” (CODEX STAN 118-1979) the limit value for gluten-free food at 20 mg/kg gluten<sup>[3]</sup>. This threshold has also been adopted by many national legislations. The prolamin content of gluten (e.g. gliadin in wheat) is set to be 50 % (CODEX STAN 118-1979).

The official type I method for gluten determination of the Codex Alimentarius<sup>[4]</sup> is an ELISA, which uses the R5 antibody in combination with a special extraction buffer, the Cocktail (patented) (Mendéz method)<sup>[5]</sup>. The RIDASCREEN® Gliadin ELISA (Art. No. R7001) fulfils this requirement and showed in an additional international collaborative study in 2020 its broad applicability for investigation of foods from different product groups and of different processing levels<sup>[6, 7]</sup>.

The RIDA®QUICK Gliadin test strips also use the R5 antibody and show a good correlation with the official method the RIDASCREEN® Gliadin ELISA.

**Only R-Biopharm AG is authorized to use the R5 antibody in immunochromatographic tests.**

## 3. Test principle

The principle of the test in strip format is an antigen-antibody reaction and bases on the monoclonal R5-antibody, which is specific for the detection of gliadin from wheat and prolamins from rye and barley. In the lower are, the strip contains red latex

particle-R5-antibody-conjugates (latex conjugates). The test band contains immobilized R5 antibody. The strip is placed in a running buffer containing a prepared sample. Following the chromatographic principle, the liquid (running buffer and sample) moves upwards inside the strip. Latex conjugates are solved by the liquid and pass together with it the test band. If the analyte is present, a sandwich is formed consisting of immobilized R5 antibody, gliadin and red latex-conjugate coloring the test band red. Results are read visually.

#### 4. Reagents provided

Each kit contains sufficient materials for max. 25 determinations. Each test kit contains:

| Component       | Cap color   | Format                        | Volume    |
|-----------------|-------------|-------------------------------|-----------|
| Test strip      |             | Ready to use, single packaged | 25 pieces |
| Test tube       |             |                               | 30 pieces |
| Buffer          | Transparent | Ready to use                  | 60 mL     |
| Evaluation card |             |                               | 1 piece   |

#### 5. Materials required but not provided

##### 5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of  $\pm 0.01$  g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x *g*) + centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Optional: filter paper (pore size 8 - 12  $\mu\text{m}$ )
- Shaker
- Water bath (50 °C / 122 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Graduated pipettes
- Graduated cylinder

## 5.2 Reagents

Different reagents are recommended for analysis of raw materials and processed food (see 9.5)

- Distilled water (dist. water) or deionized water
- Skimmed milk powder (SMP) (food quality) for soy, tannin and polyphenol containing food
- 60 % ethanol solution (mix 150 mL ethanol p.a. with 100 mL dist. water)
- 80 % ethanol solution (mix 120 mL ethanol p.a. with 30 mL dist. water)
- Cocktail (patented) (Art. No. R7006 / R7016)
- Alternatively: Cocktail ECO (Art. No. R7080)

## 6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

Only trained laboratory employees should carry out this test. Follow strictly the instruction for use.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

The Cocktail (patented) contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze the test strips or any test kit components.

Once the test strip container has been opened, store the container at room temperature (20 - 25 °C / 68 - 77 °F).

The test strips are sensitive to moisture and must be absolutely protected against it. Moist test strips can have a negative impact on the result. For this reason, remove the test strips from packaging only after having reached room temperature and immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## **8. Indication of instability or deterioration of reagents**

- A missing control band or an incomplete test or control band
- A turbid or flaked buffer

## **9. Test implementation**

9.1 General information

9.2 Swab test on surfaces

9.3 Analysis of CIP water without cleansing agent

9.4 Analysis of CIP water with cleansing agent

9.5 Extraction of foods

9.5.1 Extraction with ethanol

9.5.2 Extraction with Cocktail (patented)

9.5.3 Extraction with Cocktail ECO

9.6 Test procedure for food analysis

9.1 General information

The buffer is ready for use and must be brought to room temperature (20 - 25 °C / 68 - 77 °F) before use. Make sure that the buffer does not become contaminated with gluten during use.

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol.
- For the analysis, the extraction and the test procedure should be carried out in separate rooms.
- When using the Cocktail (patented), it is recommended to work under a chemical hood, because it contains  $\beta$ -Mercaptoethanol.

## 9.2 Swab test on surfaces

- Take as many test tubes as surfaces to be analyzed.
- Place 500  $\mu\text{L}$  of buffer in the test tube.
- Swab the lower end (application area; see also chapter 10 fig. 2) of a dry test strip thoroughly in a cross-hatch technique over a sampling area of 10 x 10 cm (see fig. 1).


Fig. 1: Direct swabbing of surfaces with RIDA®QUICK Gliadin

- Place the strip into the test tube with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).
- Take out the strip after exactly 5 min ( $\pm 10$  s) and read the result using the evaluation card.

## 9.3 Analysis of CIP water without cleansing agent

- Take as many test tubes as samples to be analyzed.
- Place 250  $\mu\text{L}$  of buffer in the test tubes.
- Place 250  $\mu\text{L}$  of cleansing water in the test tube and mix gently.
- Place a test strip into the test tube with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).
- Take out the strip after exactly 5 min ( $\pm 10$  s) and read the result using the evaluation card.

## 9.4 Analysis of CIP water with cleansing agent

- Take as many test tubes as samples to be analyzed.
- Place 500  $\mu\text{L}$  of buffer in the test tube.
- Place 50  $\mu\text{L}$  of cleansing water in the test tube and mix gently.
- Place a test strip into the test tube with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).

- Take out the strip after exactly 5 min ( $\pm$  10 s) and read the result using the evaluation card.

## 9.5 Extraction of foods

For AOAC OMA 2015.16 validation<sup>[2]</sup>, non-processed corn samples were analyzed using the 60 % ethanol extraction and processed corn samples using the Cocktail (patented) extraction.

Depending on the extraction method, the addition of skimmed milk powder (SMP) is necessary for some ingredients to avoid unwanted disturbing reactions.

| Food ingredient   | Ethanol extraction | Cocktail (patented) / Cocktail ECO extraction |
|---|--------------------|---|
| Soy   | 1 g SMP            | -   |
| Tannin- and polyphenol containing food (e.g. chocolate, coffee, cacao, chestnut flour, buckwheat, millet, spices) | 1 g SMP            | 0.25 g SMP                                    |

### 9.5.1 Extraction with ethanol

*(recommended for raw materials)*

Homogenize well a sufficient amount of sample (grind at least 50 g thoroughly to powder and mix well or mix well 50 mL of a liquid sample, respectively).

- **Liquid food samples:** pipet 1 mL of the homogenized sample into a vial (add 1 g of SMP in case of soy, tannin and polyphenol containing samples) and add 9 mL of 60 % ethanol; close the vial and mix well.
- **Soft and hard food samples:** weigh 1 g of the homogenized sample (add 1 g of SMP in case of soy, tannin and polyphenol containing samples) and add 10 mL of 60 % ethanol; close the vial and mix well.
- Shake well for at least 30 sec (e.g. vortex).
- Centrifuge: 10 min / at least 2,500 x *g* / room temperature (20 - 25 °C / 68 - 77 °F).  
(Alternatively: let the sample settle down and/or filtrate).

Principally, ethanol extraction can also be used for fast analysis of other samples than raw materials. **However, this procedure is not part of the AOAC Official Method of Analysis.** Please notice, that efficiency of ethanol extraction can be reduced especially with heat-processed samples [5]. Hence, users have to verify suitability of ethanol extraction for the particular sample material. The group of Enrique Mendéz has developed the Cocktail (patented) for improved extraction of gluten from heat-processed samples<sup>[8]</sup>. It can be used generally for all kind of foods. As a faster and more environmental-friendly alternative, R-Biopharm has developed the Cocktail ECO, which shows a comparable extraction efficiency.

#### 9.5.2 Extraction with Cocktail (patented) (Art. No. R7006 / R7016)

*(recommended for heat-processed food and in case of unknown treatment)*

Homogenize well a sufficient amount of sample (grind at least 50 g thoroughly to powder and mix well or mix well 50 ml of a liquid sample, respectively).

- **Liquid food samples:** pipet 0.25 mL of the homogenized sample into a vial (add 0.25 g of SMP in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail (patented); close the vial and mix well.
- **Soft and hard food samples:** weigh 0.25 g of the homogenized sample (add 0.25 g of SMP in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail (patented); close the vial and mix well.
- **Meat and sausages:** uneven gluten distribution is typical for these matrices; therefore, ensure to homogenize a sufficient sample amount: weigh 0.25 g of the homogenized sample (add 0.25 g of SMP in case of tannin and polyphenol containing samples) into a vial and add 2.5 mL of Cocktail (patented); close the vial and mix well.
- **Oat samples:** uneven gluten distribution is typical for oats; furthermore, the samples are difficult to homogenize; therefore, homogenize 200 g and carry out the extraction with at least the fourfold amount of sample and reagents: weigh 1 g of the homogenized sample into a vial and add 10 mL of Cocktail (patented); close the vial and mix well.

Please further extract all Cocktail (patented) samples as described in the following:

- Incubate for 40 min at 50 °C (122 °F) in a water bath.
- Let the sample cool down shortly (1 - 3 min).
- Add 7.5 mL 80 % ethanol (see chapter 5.2) close the vial and mix well (for oat samples, add 30 mL 80 % ethanol accordingly).

- Shake for 1 h upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter sample or centrifuge for 10 min at > 2,500 x *g* at room temperature (20 - 25 °C) / 68 - 77 °F).  
(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x *g*) for 10 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The undiluted sample extracts (supernatant from centrifugation step or filtrate) can be stored in a well-sealed tube at room temperature (20 - 25 °C / 68 - 77 °F) in the dark (shelf life approx. 4 weeks).

### 9.5.3 Extraction with Cocktail ECO (Art. No. R7080)

*(recommended for processed food and in case of unknown treatment)*

For a faster sample preparation, the more environmental-friendly Cocktail ECO can be used as alternative. This procedure is not part of the AOAC Official Method of Analysis. The Cocktail ECO has an extraction efficiency of approx. 70 - 110 % compared to Cocktail (patented).

Prepare the necessary amount of Cocktail ECO according to the corresponding product information.

Homogenize well a sufficient amount of sample (grind at least 50 g thoroughly to powder and mix well or mix well 50 mL of a liquid sample, respectively).

- **Liquid food samples:** pipet 0.25 mL of the homogenized sample into a vial (add 0.25 g of SMP in case of tannin and polyphenol containing samples) and add 2.5 mL of RIDA® Cocktail ECO; close the vial and mix well.
- **Soft and hard food samples:** weigh 0.25 g of the homogenized sample (add 0.25 g of SMP in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail ECO; close the vial and mix well.
- **Meat and sausages:** uneven gluten distribution is typical for these matrices; therefore, ensure to homogenize a sufficient sample amount: weigh 0.25 g of the homogenized sample (add 0.25 g of SMP in case of tannin and polyphenol containing samples) into a vial and add 2.5 mL of the Cocktail ECO; close the vial and mix well.
- **Oat samples:** uneven gluten distribution is typical for oats; furthermore, the samples are difficult to homogenize; therefore, homogenize 200 g and carry out

the extraction with at least the fourfold amount of sample and reagents: weigh 1 g of the homogenized sample into a vial and add 10 mL of Cocktail ECO; close the vial and mix well.

Please further extract all Cocktail ECO samples as described in the following:

- Incubate for 10 min at 50 °C (122 °F) in a water bath.
- Let the sample cool down shortly (1 - 3 min).
- Add 7.5 mL 80 % ethanol (see chapter 5.2) close the vial and mix well (for oat samples, add 30 ml 80 % ethanol accordingly).
- Shake for 10 min upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter sample or centrifuge for 5 min at > 2,500 x *g* at room temperature (20 - 25 °C) / 68 - 77 °F).  
(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x *g*) for 5 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The undiluted sample extracts (supernatant from centrifugation step or filtrate) can be stored in a well-sealed tube at room temperature (20 - 25 °C / 68 - 77 °F) in the dark (shelf life approx. 2 weeks).

## 9.6 Test procedure for food analysis

- Take as many test tubes as samples to be analyzed.
- Place 500 µL of the buffer into each test tube.
- Add 50 µL of the extracts prepared according to chapter 9.5.1 - 9.5.3 to the test tubes and mix carefully.
- Place one test strip into the test tube with the arrow pointing down and take care that the strip immerses only with the application area (arrow mark on strip).
- Take out the strip after exactly 5 min ( $\pm$  10 s) and read the result using the evaluation card.

## 10. Evaluation

The right line in the result area is the red test band (see fig. 2). The sample contains gluten, if the test band becomes visible.


Fig. 2: Test strip RIDA®QUICK Gliadin

The left line in the result area is the blue control band (see fig. 2). It must be present after each test procedure.

### **Positive result:** two complete colored bands

The sample is positive if two colored bands (the blue control band and the red test band) are visible within the result area. In the case of swabbing, complete test bands with non-uniform intensity may occur due to an inhomogeneous gluten distribution on the surface or different swabbing procedures. Principally, the test band's red color is stronger with higher analyte concentrations.

|                                      |   |
|--------------------------------------|---|
| Swab test:                           | > approx. 1.6 - 3 µg gluten/100 cm <sup>2</sup> |
| Raw material:                        | > approx. 4.4 mg/kg gluten*                     |
| Processed food:                      | > approx. 6.3 mg/kg gluten**                    |
| CIP water (without cleansing agent): | > approx. 10 ng/mL gluten                       |
| CIP water (with cleansing agent):    | > approx. 50 - 100 ng/ml gluten                 |

\* ethanol extraction

\*\* Cocktail (patented) extraction

### **Negative result:** only one complete blue band

The sample is negative if only the blue test band is visible within the result window.

|                                      |   |
|--------------------------------------|---|
| Swab test:                           | < approx. 1.6 - 3 µg gluten/100 cm <sup>2</sup> |
| Raw material:                        | < approx. 4.4 mg/kg gluten*                     |
| Processed food:                      | < approx. 6.3 mg/kg gluten**                    |
| CIP water (without cleansing agent): | < approx. 10 ng/mL gluten                       |
| CIP water (with cleansing agent):    | < approx. 50 - 100 ng/ml gluten                 |

\* ethanol extraction

\*\* Cocktail (patented) extraction

### **Invalid result:** no control band or incomplete bands

If the blue control band is not visible after a run or an incomplete control or test band is visible, the test result is not valid and the test must be repeated with a new test strip. If the result is invalid again, check the expiry date of the test kit and control the content for visible sign of reagent deterioration (see chapter 8).

## 11. Limits of the method

The test has been developed for the detection of low amounts of gluten (contamination). No high-dose-hook-effect is observed at high concentrations. However, the red test band may smear at very high gluten concentrations (> 10,000 mg/kg gluten).

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

The limit of detection depends on sample type and extraction efficiency or the properties of the swabbed surface and the kind of contamination respectively.

Inhomogeneous gluten distribution inside a sample and an inadequate sample size may cause a negative result.

Cleansing water containing hypochlorite cannot be analyzed. This cleansing agent destroys gluten very quickly by oxidation. The test strip is not able to detect potentially remaining gluten fragments.

A negative result does not exclude a gluten contamination below the detection limit of the assay, or that other cereal components, such as starch, may be present in a sample. The result should be reported accordingly.

It is recommended to use ethanol extraction with foods, which have not been heated. Extraction efficiency may be reduced with heat-processed samples. Kit users need to verify suitability of an extraction method for particular foods.

Furthermore, proteins may have been altered and/or fragmented in processed food (e.g. heat treatment, dehydration, etc.). This may have an impact on recovery or cross-reactivity.

For evaluation of the cross-reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross-reactivities and exemplary analyzed matrices are described in the validation report.

Cross reactivities are side reactions of the used antibody with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar

influences (e.g. by heating or drying) as the actual analyte. Therefore, food processing may induce cross-reactivity in single cases or cause its loss.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive results, but may also reduce or suppress a correct reaction. Such matrix effects (caused by e.g. polyphenols) are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

Addition of foreign protein (e.g. skimmed milk powder) during extraction or test procedure may suppress matrix effects.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and interlaboratory comparisons.

## **12. Recommendation**

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842.
- To carry along test controls for quality control. Gluten-free and gluten-containing (spiked) samples should be used (e.g. Art. No. R7012).
- To do spike experiments to ensure an accurate and correct test procedure. An example of a spike experiment is given in the validation report.
- To adjust the sample's pH value to neutral (pH 6.5 to 7.5) prior to extraction in case of extremely acidic or basic samples.
- To compare the extraction efficiency of ethanol and Cocktail ECO (Art. No. R7080) with the Cocktail (patented) (Art. No. R7006 / R7016).
- To clear unexpected results, use the ELISA RIDASCREEN® Gliadin (Art. No. R7001) or run a PCR (e.g. SureFood®).
- To use the ELISA RIDASCREEN® Gliadin (Art. No. R7001) for quantification of a positive result.
- For quick and easy documentation of the test strips, use the "RIDA®SMART APP Allergen" documentation app available free of charge in the Google Play Store (only for Android devices).

For storing the test strip for documentation, cut off the upper part of the strip marked with "Gluten" together with the test bands.

### 13. Further applications

- Sample preparation for processed food with the RIDA® Extraction Solution (colorless) (Art. No. R7098) - **only after validation.**
- Sample preparation for polyphenol containing raw material using fish gelatine.

**For further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

### Literature

- [1] Lacorn M. et al.: Determination of Gluten in Processed and Nonprocessed Corn Products by Qualitative R5 Immuno-chromatographic Dipstick: Collaborative Study, First Action 2015.16., J AOAC Int. 2016; 99(3):730-737
- [2] Lacorn M. et al.: The Validation of the RIDA®QUICK Gliadin for AOAC Research Institute. J AOAC Int. 2018; 101(5):1548-1557
- [3] FAO Codex Alimentarius International Food Standards: Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten. CXS 118-1979. Adopted in 1979. Amended in 1983 and 2015. Revised in 2008.
- [4] FAO Codex Alimentarius International Food Standards: Recommended Methods of Analysis and Sampling. CXS 234-1999. Adopted in 1999. Amended in 2021. Pg. 21
- [5] Valdés I. et al.: Innovative approach to low-level gluten determination in foods using a novel sandwich enzyme-linked immunosorbent assay protocol. Eur J Gastroenterol Hepatol. 2003; 15(5):465-74
- [6] Immer U., Haas-Lauterbach S.: Gliadin as a measure of gluten in foods containing wheat, rye, and barley-enzyme immunoassay method based on a specific monoclonal antibody to the potentially celiac toxic amino acid prolamin sequences: collaborative study. J AOAC Int. 2012; 95(4):1118-1124
- [7] Lacorn M. et al.: Determination of Gliadin as a Measure of Gluten in Food by R5 Sandwich ELISA RIDASCREENVR Gliadin Matrix Extension: Collaborative Study 2012.01. J AOAC Int. 2022; 105(2):442-455
- [8] Garcia E. et al.: Development of a general procedure for complete extraction of gliadins for heat processed and unheated foods. Eur J Gastroenterol Hepatol. 2005 May; 17(5):529-39

## Version overview

| Version number | Chapter and title  |
|----------------|--|
| 2018-01-16     | Previous version   |
| 2023-02-23     | Previous version   |
| 2024-06-27     | Current version<br>Changes made: <ul style="list-style-type: none"><li>– General linguistic and formal revision</li><li>– Addition of related products for gluten/gliadin detection</li><li>– Modification of the disposal clause in chapter 6</li><li>– Chapter 8 “Indication of instability or deterioration of reagents” (new)</li><li>– Addition of buffer handling in chapter 9.1</li><li>– Chapter 9: Specification of usability of different extraction methods (analogous to R7003)</li><li>– Revision of chapter 10, 11 and 12</li><li>– Addition of the availability of the free of charge documentation app in chapter 12</li><li>– Addition of marking concerning patents of MORINAGA &amp; Co., Ltd.</li><li>– New disclaimer</li></ul> |

## Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

## Patent Marking:

If the R-Biopharm product “Cocktail ECO” is used for the processing of food samples, a sulphite-containing extraction agent is used. Food inspection methods using a sulfite-containing extractant as in this product and/or corresponding detection kits are subject to the following patents of MORINAGA & Co., Ltd.: European Patent EP 2 224 239 B1, Australian Patent AU 2008 330 507 B2, United States Patent US 8 859 212 B2, Japanese Patent JP 5 451 854 B2. The patent holder has granted R-Biopharm AG a license to use, and sell products that employ, said protected technology in the above-mentioned territories.

## Disclaimer

1. In conformance with the German Civil Code (“BGB”) R-Biopharm AG provides a limited warranty (“Gewährleistung”) against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.  
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm’s own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
  - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
  - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
  - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
  - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG’s liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG’s sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser’s expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements (“Gewährleistung”) implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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