RIDA® QUICK Gliadin (ready to swab)

REF R7005

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Immuno chromatographic test for the qualitative determination of gliadin and corresponding prolamins on surfaces and in cleansing/process water

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AOAC *Performance Tested Methods*SM (101702)
for surfaces and cleansing waters

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Please contact for questions and further information:

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RIDA[®], RIDASCREEN[®] and RIDASOFT[®]
are registered trademarks of R-Biopharm AG.
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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RIDASCREEN® Gliadin (Art. Nr. R7001)

RIDASCREEN®FAST Gliadin (Art. Nr. R7002)

RIDASCREEN® Gliadin competitive (Art. Nr. R7021)

RIDASCREEN® Total Gluten (Art. Nr. R7041)
RIDASCREEN®FAST Gliadin sensitive (Art. Nr. R7051)
Cocktail (patented) (Art. Nr. R7006 / R7016)
RIDA® Extraction Solution (colorless) (Art. Nr. R7098)
Set of 3 processed Gliadin Assay Controls (Art. Nr. R7012)
SureFood® ALLERGEN Gluten (Art. Nr. S3606)
SureFood® QUANTARD Allergen 40 (Art. Nr. S3301)
SureFood® ALLERGEN 4plex Cereals (Art. Nr. S7006)

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Literatur

- [1] Lacorn M. et al.: The Validation of the RIDA®QUICK Gliadin for AOAC Research Institute. J AOAC Int. 2018; 101(5):1548-1557
- [2] FAO Codex Alimentarius International Food Standards: Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten. CXS 118-1979. Adopted in 1979. Amended in 1983 and 2015. Revised in 2008.
- [3] FAO Codex Alimentarius International Food Standards: Recommended Methods of Analysis and Sampling. CXS 234-1999. Adopted in 1999. Amended in 2021. Pg. 21
- [4] Valdés I. et al.: Innovative approach to low-level gluten determination in foods using a novel sandwich enzyme-linked immunosorbent assay protocol. Eur J Gastroenterol Hepatol. 2003; 15(5):465-74
- [5] Immer U., Haas-Lauterbach S.: Gliadin as a measure of gluten in foods containing wheat, rye, and barley-enzyme immunoassay method based on a specific monoclonal antibody to the potentially celiac toxic amino acid prolamins sequences: collaborative study. J AOAC Int. 2012; 95(4):1118-1124
- [6] Lacorn M. et al.: Determination of Gliadin as a Measure of Gluten in Food by R5 Sandwich ELISA RIDASCREEN® Gliadin Matrix Extension: Collaborative Study 2012.01. J AOAC Int. 2022; 105(2):442-455

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RIDA®QUICK Gliadin

Brief information

RIDA®QUICK Gliadin (ready to swab) (Art. No. R7005) is an R5-based immunochromatographic test for the qualitative determination of gluten-contamination on surfaces and in cleansing waters (CIP waters).

The identical test strip RIDA®QUICK Gliadin (Art. No. R7003) has been approved as AOAC Proficiency Tested Method (PTM 101702) for swabbing and cleansing waters^[1].

The test kit contains 25 test strips (single packaged) and 25 vials filled with 0.5 mL buffer, each of which can be used for one analysis. Results are evaluated visually.

Time requirement: sampling for swab test ca. 1 min
Sample preparation for:
10 CIP water samples ca. 5 min
Test implementation (incubation time)..... 5 min
CIP = Clean-In-Place (rinse water from cleaning process)

Detection Limit: Surfaces..... ca. 1.6 - 3 µg gluten/100 cm²
(Matrix-dependent) CIP water (without cleansing agent)....ca. 10 ng/mL gluten
CIP water (with cleansing agent)...ca. 50 - 100 ng/mL gluten

Specificity: The **monoclonal antibody R5** reacts with the gliadin-fractions from wheat and corresponding prolamins from rye and barley. Further information is contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed/processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments.

Related products and accessories for gluten/gliadin determination

RIDA®QUICK Gliadin (Art. No. R7003)
RIDA®QUICK Gliadin (single packaged) (Art. No. R7004)
RIDASCREEN® Gliadin (Art. No. R7001)
RIDASCREEN®FAST Gliadin (Art. No. R7002)
RIDASCREEN® Gliadin competitive (Art. No. R7021)
RIDASCREEN® Total Gluten (Art. No. R7041)
RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051)
Cocktail (patented) (Art. No. R7006 / R7016)
RIDA® Extraction Solution (colorless) (Art. No. R7098)
Set of 3 processed Gliadin Assay Controls (Art. No. R7012)
SureFood® ALLERGEN Gluten (Art. No. S3606)
SureFood® QUANTARD Allergen 40 (Art. No. S3301)
SureFood® ALLERGEN 4plex Cereals (Art. No. S7006)

1. Intended use

RIDA®QUICK Gliadin (ready to swab) (Art. No. R7005) is an R5-based immunochromatographic assay in test strip format for the qualitative determination of gluten

- on surfaces (swab test for hygiene control in production and laboratories)
- in cleansing waters (CIP water).

The test has been developed for the detection of low amounts of gluten (contamination).

For the AOAC-RI validation of the assay for testing of surfaces and CIP waters^[1], stainless steel, sealed ceramic, plastic and silicone rubber were validated or waters without and with commercial cleansing agents respectively.

2. General information

The use of wheat flour and gluten in foodstuff is extremely common because of their useful effects on e.g. texture, moisture retention and flavor. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley. Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

Especially airborne cereal dust lead to gliadin contamination in the production. The cleaning of surfaces should be carried out with alcoholic solutions. Thereafter, using

swab tests, it is verified if the cleaning of the surfaces was sufficient or the cleaning water is free of residues.

The Codex Alimentarius Commission has stipulated the limit value for gluten-free food at 20 mg/kg gluten^[2]. This threshold has also been adopted by many national legislations. The prolamin content of gluten (e.g. gliadin in wheat) is set to be 50 %.

The official type I method for gluten determination of the Codex Alimentarius^[3] is an ELISA, which uses the R5 antibody in combination with a special extraction buffer, the Cocktail (patented) (Mendéz method)^[4]. The RIDASCREEN® Gliadin ELISA (Art. No. R7001) fulfils this requirement and showed in an additional international collaborative study in 2020 its broad applicability for investigation of foods from different product groups and of different processing levels^[5, 6].

The RIDA®QUICK Gliadin test strips also use the R5 antibody and show a good correlation with the official method the RIDASCREEN® Gliadin ELISA.

Only R-Biopharm AG is authorized to use the R5 antibody in immunochromatographic tests.

3. Test principle

The principle of the test in strip format is an antigen-antibody reaction and bases on the monoclonal R5-antibody, which is specific for the detection of gliadin from wheat and prolamins from rye and barley. In the lower are, the strip contains red latex particle-R5-antibody-conjugates (latex conjugates). The test band contains immobilized R5 antibody. The strip is placed in a running buffer containing a prepared sample. Following the chromatographic principle, the liquid (running buffer and sample) moves upwards inside the strip. Latex conjugates are solved by the liquid and pass together with it the test band. If the analyte is present, a sandwich is formed consisting of immobilized R5 antibody, gliadin and red latex-conjugate coloring the test band red. Results are read visually.

4. Reagents provided

Each kit contains sufficient materials for max. 25 determinations. Each test kit contains:

Component	Cap color	Format	Volume
Test strip		Ready to use, single packaged	25 pieces
Buffer vials	Transparent	Ready to use	25 x 0.5 mL
Evaluation card			1 piece

5. Materials required but not provided

- Gloves

The following is further required for the analysis of cleansing and process water:

- Graduated pipettes

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

Only trained laboratory employees should carry out this test. Follow strictly the instruction for use.

This kit may contain hazardous substances. Please refer to the component safety information in the material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze the test strips or any test kit components.

Once the test strip container has been opened, store the container at room temperature (20 - 25 °C / 68 - 77 °F).

The test strips are sensitive to moisture and must be absolutely protected against it. Moist test strips can have a negative impact on the result. For this reason, remove the test strips from packaging only after having reached room temperature and immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

8. Indication of instability or deterioration of reagents

- A missing control band or an incomplete test or control band
- A turbid or flaked buffer

9. Test implementation

9.1 General information

9.2 Swab test on surfaces

9.3 Analysis of CIP water without cleansing agent

9.4 Analysis of CIP water with cleansing agent

9.1 General information

The buffer is ready for use and must be brought to room temperature (20 - 25 °C / 68 - 77 °F) before use. Make sure that the buffer does not become contaminated with gluten during use.

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol.

9.2 Swab test on surfaces

- Take as many buffer vials as surfaces to be analyzed. Tap the buffer vials onto the table so that the liquid detaches from the lid.
- Swab the lower end (application area; see also chapter 10 fig. 2) of a dry test strip thoroughly in a cross-hatch technique over a sampling area of 10 x 10 cm (see fig. 1).


Fig. 1: Direct swabbing of surfaces with RIDA®QUICK Gliadin

- Place the strip into the test tube with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).

- Take out the strip after exactly 5 min (\pm 10 s) and read the result using the evaluation card.

9.3 Analysis of CIP water without cleansing agent

- Take as many buffer vials as samples to be analyzed.
- Remove 250 μ L of buffer from the vial and discard it.
- Add 250 μ L of cleansing water in the buffer vial and mix gently.
- Place a test strip into the vial with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).
- Take out the strip after exactly 5 min (\pm 10 s) and read the result using the evaluation card.

9.4 Analysis of CIP water with cleansing agent

- Take as many buffer vials as samples to be analyzed.
- Add 50 μ L of cleansing water in the buffer vial and mix gently.
- Place a test strip into the test tube with the arrow pointing down and do not immerse the strip beyond the maximum line (arrow mark on strip).
- Take out the strip after exactly 5 min (\pm 10 s) and read the result using the evaluation card.

10. Evaluation

The right line in the result area is the red test band (see fig. 2). The sample contains gluten, if the test band becomes visible.


Fig. 2: Test strip RIDA®QUICK Gliadin

The left line in the result area is the blue control band (see fig. 2). It must be present after each test procedure.

Positive result: two complete colored bands

The sample is positive if two colored bands (the blue control band and the red test band) are visible within the result area. In the case of swabbing, complete test bands with non-uniform intensity may occur due to an inhomogeneous gluten distribution on the surface or different swabbing procedures. Principally, the test band's red color is stronger with higher analyte concentrations.

Swab test:	> approx. 1.6 - 3 µg gluten/100 cm ²
CIP water (without cleansing agent):	> approx. 10 ng/mL gluten
CIP water (with cleansing agent):	> approx. 50 - 100 ng/ml gluten

Negative result: only one complete blue band

The sample is negative if only the blue test band is visible within the result window.

Swab test:	< approx. 1.6 - 3 µg gluten/100 cm ²
CIP water (without cleansing agent):	< approx. 10 ng/mL gluten
CIP water (with cleansing agent):	< approx. 50 - 100 ng/ml gluten

Invalid result: no control band or incomplete bands

If the blue control band is not visible after a run or an incomplete control or test band is visible, the test result is not valid and the test must be repeated with a new test strip. If the result is invalid again, check the expiry date of the test kit and control the content for visible sign of reagent deterioration (see chapter 8).

11. Limits of the method

The test has been developed for the detection of low amounts of gluten (contamination). No high-dose-hook-effect is observed at high concentrations. However, the red test band may smear at very high gluten concentrations (> 10,000 mg/kg gluten).

Test results may vary depending on the actual test procedure and the laboratory environment.

The limit of detection depends on sample type and extraction efficiency or the properties of the swabbed surface and the kind of contamination respectively.

Cleansing water containing hypochlorite cannot be analyzed. This cleansing agent destroys gluten very quickly by oxidation. The test strip is not able to detect potentially remaining gluten fragments.

A negative result does not exclude a gluten contamination below the detection limit of the assay, or that other cereal components, such as starch, may be present in a sample. The result should be reported accordingly.

For detailed results and further information please refer to the current validation report. In addition, data may be available from comparative laboratory tests and interlaboratory comparisons.

12. Recommendation

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842.
- To carry along test controls for quality control. Gluten-free and gluten-containing (spiked) samples should be used (e.g. Art. No. R7012).
- To do spike experiments to ensure an accurate and correct test procedure. An example of a spike experiment is given in the validation report.
- To adjust the sample's pH value to neutral (pH 6.5 to 7.5) prior to extraction in case of extremely acidic or basic samples.
- To clear unexpected results, use the ELISA RIDASCREEN® Gliadin (Art. No. R7001) or run a PCR (e.g. SureFood®).
- To use the ELISA RIDASCREEN® Gliadin (Art. No. R7001) for quantification of a positive result.
- For quick and easy documentation of the test strips, use the "RIDA®SMART APP Allergen" documentation app available free of charge in the Google Play Store (only for Android devices).

For storing the test strip for documentation, cut off the upper part of the strip marked with "Gluten" together with the test bands.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Literature

- [1] Lacorn M. et al.: Determination of Gluten in Processed and Nonprocessed Corn Products by Qualitative R5 Immuno-chromatographic Dipstick: Collaborative Study, First Action 2015.16., J AOAC Int. 2016; 99(3):730-737
- [2] FAO Codex Alimentarius International Food Standards: Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten. CXS 118-1979. Adopted in 1979. Amended in 1983 and 2015. Revised in 2008.
- [3] FAO Codex Alimentarius International Food Standards: Recommended Methods of Analysis and Sampling. CXS 234-1999. Adopted in 1999. Amended in 2021. Pg. 21
- [4] Valdés I. et al.: Innovative approach to low-level gluten determination in foods using a novel sandwich enzyme-linked immunosorbent assay protocol. Eur J Gastroenterol Hepatol. 2003; 15(5):465-74
- [5] Immer U., Haas-Lauterbach S.: Gliadin as a measure of gluten in foods containing wheat, rye, and barley-enzyme immunoassay method based on a specific monoclonal antibody to the potentially celiac toxic amino acid prolamin sequences: collaborative study. J AOAC Int. 2012; 95(4):1118-1124
- [6] Lacorn M. et al.: Determination of Gliadin as a Measure of Gluten in Food by R5 Sandwich ELISA RIDASCREENVR Gliadin Matrix Extension: Collaborative Study 2012.01. J AOAC Int. 2022; 105(2):442-455

Version overview

Version number	Chapter and title
2018-01-16	Previous version
2023-02-23	Previous version
2024-07-01	Current version Changes made: <ul style="list-style-type: none">– General linguistic revision– Updated list of related products for gluten/gliadin detection– Chapter 5 “Materials required but not provided“ added– Modification of the disposal clause in chapter 6– Chapter 8 “Indication of instability or deterioration of reagents” newly added– Addition of buffer handling in chapter 9.1– Chapter 9: Specification of usability of different extraction methods (analogous to R7003)– Revision of chapter 10, 11 and 12– Addition of the availability of the free of charge documentation app in chapter 12– Chapter „Literature“ newly added– Updated disclaimer

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

1. In conformance with the German Civil Code (“BGB”) R-Biopharm AG provides a limited warranty (“Gewährleistung”) against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
2. ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm’s own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
3. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
4. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
5. R-Biopharm AG’s liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG’s sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser’s expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
6. The foregoing limited warranty is solely intended to fulfill the warranty requirements (“Gewährleistung”) implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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