RIDASCREEN® Gliadin competitive

REF R7021

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Enzyme immunoassay for the quantitative determination
of peptide fragments of gliadins and corresponding proteins

**Approved as AOAC Official Method of Analysis (OMA)
Final Action Status 2015.05**

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Please contact for questions and further information:

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RIDA[®], RIDASCREEN[®] and RIDASOFT[®]
are registered trademarks of R-Biopharm AG.
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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Der R7021 RIDASCREEN® Gliadin competitive ist:

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RIDASCREEN® Gliadin competitive

Brief information

RIDASCREEN® Gliadin competitive (Art. No. R7021) is a competitive enzyme immunoassay for the quantitative determination of prolamine peptide fragments from wheat (gliadin), rye (secalin) and barley (hordein) in foods validated for the method (see chapter 1. Intended use).

Der R5 ELISA RIDASCREEN® Gliadin competitive is:

- official AACCI method (38-55.01) for beer, starch syrup and sourdough
- validated AOAC official method of analysis (AOAC-OMA) Final Action Status (2015.05)
- validated ASBC method (Beer-49)

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction and centrifugation

Time requirement: sample preparation (for 10 samples) approx. 30 min
test implementation (incubation time).....40 min

Standard material: The RIDASCREEN® standard material is a hydrolysate (mixture of wheat, rye and barley; Gessendorf et al. (2009) Anal Bioanal Chem, 395: 1721-1728)

Limit of detection: 2.3 mg/kg (ppm) Gliadin* or to 4.6 mg/kg (ppm) Gluten*
(matrix-dependent) 1.9 - 2.6 mg/kg (ppm) Gliadin
*mean value

Limit of quantification: 5 mg/kg (ppm) Gliadin or to 10 mg/kg (ppm) Gluten

Specificity: The R5 monoclonal antibody recognizes potentially toxic peptide sequences of gliadins from wheat and related prolamins from rye and barley.

Further information are contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments (see chapter 13. Limits of the method).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

Related products for gluten / gliadin determination

RIDASCREEN® Gliadin (Art. No. R7001)

RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051)

RIDASCREEN®FAST Gliadin (Art. No. R7002)

RIDASCREEN® Total Gluten (Art. No. R7041)

RIDASCREEN® Total Gluten Additive TG (Art. No. RA0041)

RIDA®QUICK Gliadin (Art. No. R7003 / R7004 / R7005)

RIDA® Extraction Solution (colorless) (Art. No. R7098)

Cocktail (patented) (Art. No. R7006 / R7016)

Set of 3 processed Gliadin Assay Controls (Art. No. R7012)

SureFood® ALLERGEN Gluten (Art. No. S3606)

SureFood® ALLERGEN 4plex Cereals (Art. No. S7006)

SureFood® QUANTARD Allergen 40 (Art. No. S3301)

1. Intended use

RIDASCREEN® Gliadin competitive (Art. No. R7021) is a competitive enzyme immunoassay test developed for the quantitative determination of prolamine peptide fragments from wheat (gliadin), rye (secalin) and barley (hordein) in fermented or hydrolysed food, which are declared as “gluten-free”. Due to the large number of different food products, the following samples were examined as representative for different food product categories within the scope of the test development: gluten-free beer, rice beer, maize flour, sourdough, and sugar beet syrup. It can be assumed that the assay is also suitable for the

analysis of other foods, however, this must be verified by the user before applying the test kit to that food category.

For detailed results and further information on validation data with other food matrices please refer to the validation report. Further applications are regularly validated in our laboratories, which we make available in our application notes (see chapter 15).

2. General

The use of wheat flour and gluten in foodstuffs is extremely common to improve e.g. texture, moisture retention and flavour. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley. Coeliac disease is a permanent intolerance to gluten that results in damage of the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (CODEX STAN 118-1979) two categories for labeling of food according to the gluten content exist:

- 1) Food products, which contain less than 20 mg/kg (ppm), can be labeled as "**gluten-free**".
- 2) Food products labeled as "**very low gluten**" can have a gluten content above 20 and up to 100 mg/kg (ppm).

The threshold of 20 mg/kg has been adopted by many national legislations in many countries. The prolamin content (e.g. gliadin) of gluten is per definition generally assumed to be 50 % (CODEX STAN 118-1979).

During food processing like e.g. fermentation or hydrolysis intact prolamin molecules are partly or completely degraded to small peptide fragments. However, these can still be dangerous for celiac patients.

Single, small peptide sequences (motifs) cannot be detected by a sandwich ELISA format, because at least two epitopes are necessary for a sandwich ELISA. In contrast, the competitive format can also be used to detect peptide fragments with one epitope.

For the standard material used in the RIDASCREEN® Gliadin competitive (Art. No. R7021) wheat, rye and barley were digested by pepsin and trypsin, the peptide fragments were mixed and lyophilized after protein determination. The protein content was determined according to Dumas. The assay can be related to the prolamin concentration and therefore to the limit values stipulated by the Codex Alimentarius. The result is expressed in mg/kg gliadin.

The hydrolyzed standard was produced by the working group of Prof. Dr. Köhler (German Research Centre for Food Chemistry).

3. Test principle

The principle of the test is the antigen-antibody reaction. The microtiter wells are coated with gliadin. Standards (hydrolysate mixture of wheat, rye and barley), respectively sample solutions and peroxidase labeled anti-gliadin antibodies (conjugate with monoclonal R5-antibodies) are added. Free and immobilized gliadin compete for the gliadin antibody-binding sites (competitive enzyme immunoassay). Any unbound enzyme conjugate is then removed by a washing step. A substrate/chromogen solution is added to the wells and incubated. Bound enzyme conjugate converts the colorless chromogen into a blue end product. The addition of the stop solution leads to a color change from blue to yellow. The absorbance of the solution which is inversely proportional to the gliadin respectively prolamin concentration in the sample is measured photometrically at 450 nm and expressed as ng/mL (ppb) gliadin.

4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Buffer	White	Concentrate	5x	60 mL
Standard 1*	Transparent	Ready to use	0 ng/mL	1.3 mL
Standard 2*	Transparent	Ready to use	10.0 ng/mL	1.3 mL
Standard 3*	Transparent	Ready to use	30.0 ng/mL	1.3 mL
Standard 4*	Transparent	Ready to use	90.0 ng/mL	1.3 mL
Standard 5*	Transparent	Ready to use	270.0 ng/mL	1.3 mL
Wash buffer	Brown	Concentrate	10x	100 mL
Conjugate	Red	Concentrate	11x	0.7 mL
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		10 mL
Stop solution	Yellow	Ready to use		14 mL

*The standards contain a hydrolysate of wheat, rye and barley prolamins. Contrary to the AOAC-OMA Approval 2015.05, the concentrations of the standards are given in gliadin. To convert gliadin into gluten the factor of 2 used (definition Codex Alimentarius).

5. Reagents required but not provided

5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker / rotator
- Fluted filter (pore size 8 - 12 μm)
- Graduated pipettes
- Graduated cylinder
- Variable 20 - 200 μL and 200 - 1000 μL micropipettes
- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100 μL
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)

5.2 Reagents

- Distilled water (dist. water) or deionized water
- For the extraction of samples like syrup, a 60 % ethanol solution is required: transfer e.g. 150 mL ethanol p.a. with 100 mL distilled water and mix well
- For extraction of polyphenol containing samples (e.g. beer, malt, hops): liquid fish gelatine (Serva, Art. No. 22156 or Sigma, Art. No. G-7765, both 45 % solids content)
- 1 M sodium hydroxide solution

6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Do not reuse wells of the microtiter strips (coated microtiter plate and pre-plate, if necessary (see chapter 10.2)). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

Regarding the handling of ethanol, please refer to the safety instructions of the respective manufacturer.

All reagents and materials must be recovered or disposed after use at customer's own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (36 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 1.2$) for the standard 1

9. Preparation of Samples

Wear gloves before starting and during the assay. Airborne grain dust and unclean laboratory equipment can lead to gliadin contamination in the test. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation with 60 % ethanol.
- Carry out the sample preparation in a room isolated from the ELISA procedure.
- Check for gliadin contamination of reagents and equipment with the test strips RIDA[®]QUICK Gliadin (Art. No. R7003, R7004, R7005).

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

For the extraction of samples like syrup, a **60 % ethanol solution** is required. For this purpose, add e.g. 150 mL ethanol p.a. to 100 mL dist. water and mix well.

For the extraction of polyphenol containing samples like beer, malt and hops, an ethanol solution (60 %) containing 10 % liquid fish gelatin is required (**fish gelatin ethanol solution**) must be prepared. For this purpose, fill e.g. 30 mL dist. water in a 100 mL graduated cylinder. Add 10 g of liquid fish gelatin and mix well. Add 60 mL ethanol p.a. and mix well. Adjust the pH value to 8.5 with a pH electrode and 1M NaOH solution and fill up to 100 mL with dist. water. The fish gelatin cannot be dissolve completely in 60 % ethanol but remains as a suspension. Therefore, take the required volume for the sample preparation while stirring. The fish gelatin ethanol solution can be stored for 2 weeks at room temperature (20 - 25 °C / 68 - 77 °F).

9.1 Extraction

For the RIDASCREEN[®] Gliadin competitive (Art. No. R7021) only samples extracted with ethanol should be tested (the Cocktail (patented) should not be used).

Homogenize (grind thoroughly to powder and mix well or mix well a solution respectively) well a sufficient amount (e.g. 50 g or 50 mL) to ensure taking a representative test portion of sample.

- **Solid food (e.g. starch syrup):** transfer 1 g of homogenized sample to a new vial, add 10 mL 60 % ethanol solution (see chapter 9) and mix well.

- **Polyphenol-containing solid food (e.g. malt and hops):** transfer 1 g of homogenized sample, add 10 mL 60 % ethanol solution, containing 10 % liquid fish gelatin (see chapter 9) and mix well.
- **Liquid food (e.g. soy sauce):** mix well 1 mL of the homogenized sample with 9 mL 60 % ethanol solution (see chapter 9).
- **Beer:** mix well 1 mL of the homogenized sample with 9 mL 60 % ethanol solution, containing 10 % liquid fish gelatin (see chapter 9).

Please further extract all samples as described in the following:

- Mix vigorously for at least 30 s (e.g. vortexer).
- Shake well up side down or rotate on a rotator for further 10 min.
- Filter sample or centrifuge for 10 min at $> 2,500 \times g$; if possible at room temperature (20 - 25 °C / 68 - 77 °F).
(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed ($> 10,000 \times g$) for 10 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The extract (supernatant from centrifugation step or filtrate) can be stored undiluted in a well-sealed container at room temperature (20 - 25 °C / 68 - 77 °F) in the dark (shelf life approx. 4 weeks) until used in the test.
- For use in the test, extracts (supernatant from centrifugation step or filtrate) must be diluted with diluted buffer (for dilution see chapter 10.2).

10. Test implementation

10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The yellow stained **buffer** is provided as a 5fold concentrate. On the day of test performance, the required aliquot must be diluted 1:5 (1+4) with dist. water (e.g. 12 mL dist. water + 3 mL concentrate, sufficient for the dilution of 10 samples). The diluted buffer is stable for 1 day. Make sure, that the buffer is not contaminated with gliadin from the sample.

The **conjugate (antibody enzyme conjugate)** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted antibody solution has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the antibody solution should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in dist. water (e.g. 1 mL dist. water + 100 μ L conjugate concentrate, sufficient for 2 microtiter strips). Take care, that the water is not contaminated with gliadin.

The **wash buffer** is provided as a 10fold concentrate. Before use, the buffer concentrate has to be diluted 1:10 (1+9) with dist. water (e.g. 900 mL dist. water + 100 mL buffer concentrate). Prior dilution, dissolve eventually formed crystals by incubating the buffer in a water bath 37 °C (99 °F). The diluted wash buffer is stable at 20 - 25 °C (68 - 77 °F) for approx. 4 weeks.

Components should be stored at 2 - 8 °C (36 - 46 °F) when no longer required.

10.2 Test procedure

The prepared extracts prepared must be diluted 1:50 (1+49) with diluted buffer (see chapter 10.1) before usage (e.g. 980 μ L diluted buffer + 20 μ L sample).

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

It is recommended to pipette the conjugate, substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 μ L of each standard or sample (prepared according to chapter 9 and diluted according to chapter 10.2) in duplicate to the wells.
3. Add 50 μ L of the diluted conjugate (see chapter 10.1) to each well. Mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F).
4. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 μ L diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).

5. Add 100 µL of substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
6. Add 100 µL of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 min after addition of stop solution.

11. Evaluation

A special software, **RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter or cubic spline function. Once an evaluation method has been selected, it should be retained and not be switched between the two functions. Due to its mathematical background, the cubic spline function cannot calculate concentrations outside the measuring range (< standard 2 and > standard 5). The 4-parameter function enables the calculation of values between standard 1 and standard 2 (refer to chapter 13).

For the evaluation it should be clarified, that the quality criteria are fulfilled for the current test run. The course of the standard curve can be taken from the attached analysis certificate.

12. Result Interpretation

The gliadin concentration in ng/mL is read from the standard curve and must be further multiplied by the dilution factor of 500. This factor of 500 is set by default in the RIDASOFT® Win.NET software.

The RIDASOFT® Win.NET software (Version 1.93 or newer) indicates the results in gliadin and gluten. The R5 antibody detects gliadin, secalin and hordein.

Calculation example:

From the standard curve, a concentration of 50 ng/mL (ppb) gliadin is obtained. Multiplying by the recommended dilution factor 500 leads to 25000 ng/kg gliadin, corresponding to 25 mg/kg (ppm) gliadin. To calculate the gluten content, it is necessary to multiply by factor 2, which results in 50 mg/kg (ppm) gluten.

To convert gliadin into gluten the factor of 2 is used (definition Codex Alimentarius). This factor depends strongly on the analyzed sample, especially washed starches may contain a higher proportion of glutelins. In this case, the factor is bigger than 2 and the gluten concentration accordingly higher.

Results between LoD and LoQ may indicate a low content of the analyte tested in the sample. Calculated result show a high uncertainty in this area due to the method's high variation below LoQ. Therefore, such results should not be reported with a quantitative value, but qualitative as "< LoQ". Depending on the matrix, samples that do not contain the analyte can also show a result in this range.

A result below the LoD does not exclude a gluten contamination below the detection limit of the assay, or that other cereal components, such as starch, may be present in a sample. The result should be reported accordingly.

13. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LoD, LoQ, upper limit of measurement range).

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max ± 1 %.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spiking experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

The addition of foreign protein (depending on the test e.g. BSA, gelatine, skim milk powder) during extraction or test procedure may suppress matrix effects.

Cross reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

Heat-treated samples extracted with ethanol may show reduced recovery. For such samples, the extraction should be carried out with Cocktail (patented) followed by analysis in the sandwich ELISA RIDASCREEN® Gliadin (Art. No. R7001). However, Cocktail (patented) extracts cannot be analyzed in the competitive ELISA. Therefore, due to the ethanolic extraction in the RIDASCREEN® Gliadin competitive ELISA, there may be an under-determination in heated samples. In the Sandwich ELISA RIDASCREEN® Gliadin, on the other hand, under-determination may be present due to hydrolysis of the sample. Samples that are both heat-treated and hydrolyzed thus pose a particular analytical challenge. In these cases, analysis in both ELISA is recommended. Due to this problem, the results of both tests should be interpreted cautiously and from a risk perspective, the higher value should be used.

The table below shows the range of variation for the analysis of hydrolyzed samples (further method performance results are published in J. AOAC. Int. 98:1346-1354, 2015):

	Symbol	Beer			Starch syrup		Sourdough	
		gluten-free	15 mg/kg prolamin	50 mg/kg prolamin	gluten-free	naturally contaminated	35 mg/kg prolamin	75 mg/kg prolamin
Total number of laboratories		13	12	11	13	13	13	13
Total number of replicates		26	24	22	26	26	26	26
Overall mean		1.2	13.1	59.7	0.7	5.3	24.2	72.8
Recovery (%)			87	119			69	97
Repeatability SD* (mg/kg)	s(r)	1.2	4.0	18.6	1.0	0.9	5.6	14.2
Reproducibility SD* (mg/kg)	s(R)	1.5	4.8	18.6	1.5	1.8	6.3	20.0
Repeatability rel. SD* (%)	RSD(r)	97.9	30.2	31.2	157.2	16.3	23.1	19.5
Reproducibility rel. SD* (%)	RSD(R)	126.1	36.9	31.2	236.1	34.4	25.9	27.5
Horrat value		8.1	3.4	3.6	13.8	2.8	2.6	3.3

The pepsin / trypsin digest used as a calibrator in this assay does not represent all hydrolysis processes. It is recommended, that users confirm method performance for their specific processes.

The test may not measure or detect all fermented and / or hydrolyzed forms of gluten. Hydrolysis may also produce larger fragments without an epitope for the R5 antibody used, which may pose a risk to celiac patients.

14. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control. Gluten-free and gluten-containing (naturally contaminated or spiked) samples should be used.

- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 to 7.5) to neutral prior to extraction may be necessary.
- To do spike experiments to ensure an accurate and correct test procedure. An example of a spiking experiment is given in the validation report.
- To perform a PCR (e.g. SureFood®) for confirmation of the results.
- To contact sales@r-biopharm.de if automates (e.g. ThunderBolt® / Bolt™) are used.

15. Further application notes

- Analysis of gluten or β -lactoglobulin in supplemental enzymes and unprocessed food containing them.
- More sensitive analysis with RIDASCREEN® Gliadin competitive.

Further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2016-09-21	Release version
2021-07-01	<ul style="list-style-type: none"> – General linguistic revision – Chapter 1, 6, 8, 11, 12, 13, 14
2022-05-06	Current version Changes made: <ul style="list-style-type: none"> – Chapter 10.2. test procedure: Mistake from version 2021-07-01 with "total of four wash cycles" corrected back to "total of three wash cycles"
2024-05-21	Current version Changes made: <ul style="list-style-type: none"> – Update of the list under "Related products for gluten / gliadin detection" – General formal and linguistic adjustments (chapter 3, 5, 6, 8, 11, 12, 13, 14)

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

1. In conformance with the German Civil Code (“BGB”) R-Biopharm AG provides a limited warranty (“Gewährleistung”) against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm’s own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG’s liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG’s sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser’s expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements (“Gewährleistung”) implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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