
r-biopharm®

## **RIDASCREEN® Total Gluten**

**REF R7041**

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## Literatur

[1] Lacorn et al. (2019) Quantification of Wheat, Rye and Barley Gluten in Oat and Oat Products by ELISA RIDASCREEN® Total Gluten: Collaborative Study, First Action 2018.15. *J AOAC Int.*, **102**:5, p. 1535 - 1543.

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# RIDASCREEN® Total Gluten

## Brief information

RIDASCREEN® Total Gluten (Art. No. R7041) is a sandwich enzyme immunoassay for the quantitative determination of intact (non-hydrolyzed) gluten in food validated for the method (see chapter 1).

The ELISA RIDASCREEN® Total Gluten is approved as AOAC Official Method of Analysis – First Action (2018.15)<sup>[1]</sup>.

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization, extraction and centrifugation

Time requirement: sample preparation  
Cocktail (patented) (for 10 samples).....approx. 2 h  
test implementation (incubation time).....50 min

Standard material: The RIDASCREEN® standard material is a total gluten extract from four commercial wheat flours. The results with this standard are traceable to the oat samples from General Mills described in the AOAC SMPR® 2017.021.

Limit of detection: 4 mg/kg (ppm) gluten  
(matrix-dependent)

Limit of quantification: 5 mg/kg (ppm) gluten

Specificity: The used monoclonal antibodies detect gliadins from wheat and corresponding prolamins from rye and barley, high-molecular-weight (HMW) glutenin subunits (GS) from wheat and HMW-secalins from rye as well as low-molecular-weight (LMW)-GS from wheat (see also chapter 2).

There are no known cross reactivities.

Further information are contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments (see chapter 13).

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice manual. It lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA-analysis. The brochure can be retrieved, printed and downloaded from the website <https://food.r-biopharm.com/media/technical-guides/>.

### **Related products for gliadin determination**

RIDASCREEN® Gliadin (Art. No. R7001)  
RIDASCREEN®FAST Gliadin (Art. No. R7002)  
RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051)  
RIDASCREEN® Gliadin competitive (Art. No. R7021)  
RIDASCREEN® Total Gluten Additive TG (Art. No. RA0041)  
RIDA®QUICK Gliadin (Art. No. R7003 / R7004 / R7005)  
RIDA® Extraction Solution (colorless) (Art. No. R7098)  
Cocktail (patented) (Art. No. R7006 / R7016)  
Set of 3 processed Gliadin Assay Controls (Art. No. R7012)  
SureFood® ALLERGEN Gluten (Art. No. S3606)  
SureFood® ALLERGEN 4plex Cereals (Art. No. S7006)  
SureFood® QUANTARD Allergen 40 (Art. No. S3301)

### **1. Intended use**

RIDASCREEN® Total Gluten (Art. No. R7041) is a sandwich enzyme immunoassay for the quantitative determination of intact (non-hydrolyzed) gluten from gluten-containing cereals (wheat, rye and barley) in oat and oat products (oat content > 50 %). The test is AOAC Official Method of Analysis – First Action (2018.15).

Due to the large number of different food products, the following samples were examined representatively for different product groups within the scope of the test development: wheat in oats, rye in oats, barley in oats, oat-based cookies, oat-based porridge, oat flour, oat flakes, oat groats and oat cereals. It can be

assumed, that the test is also suitable for the analysis of other oat foods; this is to be checked by the user before applying the test kit to these products.

For the analysis of oat-free food samples and samples with low oat content (< 50 %) in the RIDASCREEN® Total Gluten ELISA an application note has been prepared (not part of the AOAC approval). For this purpose, the use of an additional reagent is necessary (RIDASCREEN® Total Gluten Additive TG, Art. No. RA0041). Further information on the use of the reagent can be found in the application note "Analysis of food samples with no or low oat content in RIDASCREEN® Total Gluten". It is recommended to analyze samples with unknown oat content also according to this application note.

For detailed results on validation data, please refer to the validation report. Further applications are regularly validated in our laboratories, which we make available in our application notes (see chapter 15).

## 2. General

Traditionally, cereal proteins are subdivided according their solubility (Osborne fractionation):

Soluble in	Fraction
Water	Albumins (non gluten proteins)
0.5 M NaCl	Globulins (non gluten proteins)
60 % Ethanol	Prolamins ( <b>gluten proteins</b> )
Insoluble in 60 % Ethanol	Glutelins ( <b>gluten proteins</b> )

Unfortunately, the subdivision according to solubility is not an accurate criterion, since co-solubilization and co-precipitation occur frequently. Furthermore, additional factors such as temperature can have an influence on the result. Especially, the subdivision into prolamins and glutelins, which form together the gluten proteins, according their solubility is very difficult. Therefore, the modern subdivision of the gluten proteins into protein families according to their electrophoretic mobility and their amino acid sequences is much more accurate.

Due to high sequence homologies between gluten proteins within one cereal species and between different species and due to the combination of four different monoclonal antibodies, the RIDASCREEN® Total Gluten enables the detection of almost all gluten proteins:

Wheat			
Protein family	Fraction	Mean content of gluten protein*	Antibody detection
<b>α/β-gliadins</b>	Mainly prolamins	33 %	R5 antibody
<b>γ-gliadins</b>	Mainly prolamins	27 %	R5 antibody
<b>ω1,2-gliadins</b>	Mainly prolamins	4 %	R5 antibody
<b>ω5-gliadins</b>	Mainly prolamins	3 %	R5 antibody (weak)
<b>LMW-glutenin-subunits</b>	Mainly glutelins	22 %	LMW 1 & 2 antibody
<b>HMW-glutenin-subunits</b>	Mainly glutelins	11 %	HMW antibody

Rye			
Protein family	Fraction	Mean content of gluten protein*	Antibody detection
<b>ω-secalins</b>	Mainly prolamins	18 %	R5 antibody
<b>γ-40k-secalins</b>	Prolamins and glutelins	25 %	R5 antibody
<b>γ-75k-secalins</b>	Prolamins and glutelins	48 %	R5 antibody
<b>HMW secalins</b>	Mainly glutelins	9 %	HMW antibody

Barley			
Protein family	Fraction	Mean content of gluten protein*	Antibody detection
<b>B-hordeins</b>	Mainly prolamins	27%	R5 antibody
<b>C-hordeins</b>	Mainly prolamins	36%	R5 antibody
<b>γ-hordeins</b>	Mainly prolamins	32%	R5 antibody
<b>D-hordeins</b>	Mainly glutelins	5%	No detection

\* Significant differences between different cultivars of the respective cereal and harvest years possible. Data from Wieser et al. Celiac Disease and Gluten (2014) Elsevier Inc. Amsterdam, ISBN 978-0-12-420220-7, page 107.

The use of wheat flour and gluten in foodstuffs is extremely common because of their heat stability and useful effects on e.g. texture, moisture retention and flavor. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley. Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (CODEX STAN 118/1979) two categories for labeling of food according to the gluten content now exist:

- 1) Food products which contain less than 20 mg/kg can be labeled as **"gluten-free"**.
- 2) Food products labeled as **"very low gluten"** can have a gluten content above 20 and up to 100 mg/kg.

The threshold of 20 mg/kg has been adopted by many national legislations in many countries. The prolamin content (e.g. gliadin) of gluten is per definition generally assumed to be 50 % (CODEX STAN 118-1979).

### 3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with antibodies against gluten proteins. By adding the standard or sample solution to the wells, gluten present in the solutions will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Following the washing step, a solution containing antibodies conjugated to peroxidase is added. These conjugates bind to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugates are then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue end product. The addition of the stop solution leads to a color change from blue to yellow. The absorbance of the solution, which is proportional to the gluten concentration in the sample, is measured photometrically at 450 nm and expressed as mg/kg gluten.

### 4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume
<b>Microtiter plate</b>	-	Ready to use	96 wells
<b>Buffer</b>	Transparent	Ready to use	110 mL
<b>Standard 1*</b>	Transparent	Ready to use	0 mg/kg 1.3 mL
<b>Standard 2*</b>	Transparent	Ready to use	5.0 mg/kg 1.3 mL
<b>Standard 3*</b>	Transparent	Ready to use	10.0 mg/kg 1.3 mL
<b>Standard 4*</b>	Transparent	Ready to use	20.0 mg/kg 1.3 mL
<b>Standard 5*</b>	Transparent	Ready to use	40.0 mg/kg 1.3 mL
<b>Standard 6*</b>	Transparent	Ready to use	80.0 mg/kg 1.3 mL
<b>Wash buffer</b>	Brown	<b>Concentrate</b>	<b>10x</b> 100 mL
<b>Conjugate</b>	Red	Ready to use	11 mL
<b>Substrate/Chromogen</b> Red Chromogen Pro	Brown	Ready to use	13 mL
<b>Stop solution</b>	Yellow	Ready to use	14 mL

\* The dilution factor 1,000 from the sample preparation has already been considered when labeling the standards. Therefore, the gluten concentration of samples can directly be read from the standard curve.

## 5. Reagents required but not provided

### 5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of  $\pm 0.01$  g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 ml centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (50 °C / 122 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12  $\mu\text{m}$ )
- Graduated pipettes
- Graduated cylinder
- Variable 20 - 200  $\mu\text{L}$  and 200 - 1000  $\mu\text{L}$  micropipettes
- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100  $\mu\text{L}$
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)

### 5.2 Reagents

- Distilled water (dist. water) or deionized water
- Gluten-free skim milk powder (SMP) (food quality)
- Cocktail (patented) (Art. No. R7006 / R7016)
- Ethanol solution (**80 %**): e.g. add 120 mL ethanol p.a. to 30 mL dist. water and shake well

## 6. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

Do not reuse wells of the plate (coated plate and pre-plate, if necessary, see chapter 10.2). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

The Cocktail (patented) is harmful to health. It contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

Regarding the handling of ethanol, please refer to the safety instructions of the respective manufacturer.

## **7. Storage instructions**

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## **8. Indication of instability or deterioration of reagents**

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 1.2 absorbance units ( $A_{450\text{ nm}} < 1.2$ ) for standard 6

## 9. Sample preparation

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment with 40 - 60 % ethanol or 40 - 60 % 2-propanol.
- Carry out the sample preparation in a room isolated from the ELISA procedure.
- Check for gliadin contamination of reagents and equipment with the test strips RIDA®QUICK Gliadin (R7003 / R7004 / R7005).
- When using the Cocktail (patented), it is recommended to work under a chemical hood, because it contains  $\beta$ -mercaptoethanol.
- $\beta$ -Mercaptoethanol can disturb the ELISA, therefore dilute the samples at least 1:1000.
- Try to homogenize with a clean spatula if samples gelatinizing during preparation.
- The samples should be stored protected against light.

### 9.1 Extraction with Cocktail (patented) (R7006 / R7016) for oat and oat products (official AOAC method)

The gluten distribution in oat samples can be very inhomogeneous. Furthermore, the samples are difficult to homogenize. Therefore, grind and homogenize well at least 200 g.

- Weigh in 1 g of homogenized sample to a new vial, add 10 mL of Cocktail (patented), close the vial and mix well (e.g. vortexer).
- Pay attention to obtain a homogeneous suspension!
- Add 30 mL 80 % ethanol (see chapter 5.2), close the vial and mix well (e.g. vortexer).
- Pay attention to obtain a homogeneous suspension!

*The weighted sample can be further increased for very inhomogeneous samples (> 1 g). In this case, the volume of Cocktail (patented) and 80 % ethanol must be increased accordingly. Choose a sufficient big vial (e.g. Nalgene bottle 250 mL).*

- Extract for 40 min at 50 °C (122 °F) in a water bath.
- Shake for 1 h upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter or centrifuge sample for 10 min at > 2,500 x g at room temperature (20 - 25 °C / 68 - 77 °F).  
(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed > 10,000 x g for 10 min in a microcentrifuge.)
- Transfer the supernatant into a new vial.
- If the supernatant is not free of particles after centrifugation, filter the extract.
- The sample extracts (supernatant from centrifugation step or filtrate) can be stored undiluted in a well-sealed container at room temperature (20 - 25 °C / 68 - 77 °F) in the dark (shelf life approx. 2 weeks) until used in the test.
- For use in the test, the extracts (supernatant of centrifugation step or filtrate) must be diluted with buffer (for dilution see chapter 10.2). The diluted sample extracts have a limited shelf life and must be used in the test within 30 min.
- If further dilutions are required (samples with absorbance values ( $A_{450\text{ nm}}$ ) > standard 6), the following buffer should be used to keep the composition of the buffer constant:
  - Cocktail (patented)            1 %
  - 80 % ethanol solution        3 %
  - Buffer                                96 %
 e.g. 50 µL Cocktail (patented), 150 µL 80 % ethanol solution, 4,800 µL buffer

9.2. Extraction with Cocktail (patented) (R7006 / R7016) and skim milk powder (SMP) for tannin- and polyphenol-containing oat products (e.g. oat products with high content of chocolate, coffee, cocoa, chestnut flour, teff flour, buckwheat, millet or spices)

The gluten distribution in oat samples can be very inhomogeneous. Furthermore, the samples are difficult to homogenize. Therefore, grind and homogenize well at least 200 g.

- Weigh in 1 g of homogenized sample to a new vial, add 1 g SMP and 10 mL of Cocktail (patented), close the vial and mix well (e.g. vortexer).
- Pay attention to obtain a homogeneous suspension!

- Add 30 mL 80 % ethanol (see chapter 5.2), close the vial and mix well (e.g. vortexer).
- Pay attention to obtain a homogeneous suspension!

*The weighted sample can be further increased for very inhomogeneous samples (> 1 g). In this case, the volume of Cocktail (patented) and 80 % ethanol as well as the amount of skim milk powder must be increased accordingly. Choose a sufficient big vial (e.g. Nalgene bottle 250 mL).*

- Extract for 40 min at 50 °C (122 °F) in a water bath.
- Shake for 1 h upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter or centrifuge sample for 10 min at > 2,500 x g at room temperature (20 - 25 °C) / 68 - 77 °F).

(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed > 10,000 x g for 10 min in a microcentrifuge.)

- Transfer the supernatant into a new vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The sample extracts (supernatant from centrifugation step or filtrate) can be stored undiluted in a well-sealed container at room temperature (20 - 25 °C / 68 - 77 °C) (shelf life approx. 2 weeks) until used in the test.
- For use in the test, the extracts (supernatant of centrifugation step or filtrate) must be diluted with buffer (for dilution see chapter 10.2). The diluted sample extracts have a limited shelf life and must be used in the test within 30 min.
- If further dilutions are required (samples with absorbance values ( $A_{450\text{nm}}$ ) > standard 6), the extract should be freshly diluted again 1:25 with buffer for this purpose. Then further dilute the sample extract with the following buffer to keep the composition of the buffer constant:

- Cocktail (patented)            1 %
- 80 % ethanol solution        3 %
- Buffer                                96 %

e.g. 50 µL Cocktail (patented), 150 µL 80 % ethanol solution, 4,800 µL buffer

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (e.g. 900 mL dist. water + 100 mL buffer concentrate). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for four weeks.

Components should be stored at 2 - 8 °C (35 - 46 °F) when no longer required.

### 10.2. Test procedure

**The prepared extracts must be diluted 1:25 (1+24) with buffer before usage (e.g. 960 µL buffer + 40 µL sample). The final dilution factor is 1,000.**

**Use the diluted sample extracts immediately (within 30 min) in the assay. A longer time period may influence the recovery.**

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µL per well) and then exactly 100 µL are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, substrate, chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µL of each standard or sample (prepared according to chapter 9 and diluted according to chapter 10.2) in duplicate to the wells and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).

3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
4. Add 100 µL of the conjugate to each well and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
6. Add 100 µL of substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µL of the stop solution to each well and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

## 11. Results

A special software, **RIDASOFT® Win.NET Food & Feed (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done by use of the 4-parameter function.

For the evaluation it should be clarified, that the quality criteria are fulfilled for the current test run. The course of the standard curve can be taken from the attached certificate of analysis (CoA).

When working in accordance with this extraction, the sample dilution factor is 1,000. Hence, a dilution factor of 1,000 is already taken into account with the standard concentrations (see chapter 4), the gluten concentration can directly be read from the standard curve.

## 12. Interpretation of results

The test result is calculated in mg gluten per kg of food and thus indicates the protein concentration.

Results between LoD and LoQ indicate a low egg concentration in the sample. Calculated result show a high uncertainty in this area due to the method's high variation below LoQ. Therefore, such results should not be reported with a

quantitative value, but qualitative as “< LoQ”. Depending on the matrix, samples that do not contain the analyte can also show a result in this range.

A result below the LoD does not exclude a gluten contamination below the detection limit of the assay, or that other gluten components, such as lipids, may be present in a sample. The result should be reported accordingly.

A further dilution and new detection of samples is recommended for absorbance values ( $A_{450\text{ nm}}$ ) > standard 6. Further dilutions should be made with the respective buffer (see chapter 9.1 and 9.2). In case of a further dilution, the additional dilution factor must be taken into account when calculating the allergen concentration.

Compared to the certificate, higher absorbance values ( $A_{450\text{ nm}}$ ) for the standard curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

### **13. Limits of the method**

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LoD, LoQ, upper limit of measurement range).

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max.  $\pm 1$  %.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results

obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

The addition of foreign protein (depending on the test e.g. BSA, gelatine, skim milk powder) during extraction or test procedure may suppress matrix effects.

Cross reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

## **14. Recommendation**

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842 (e.g. performing duplicate determinations).
- To pre-flush pipette tips with standard or sample extract prior to pipetting.
- To carry along test controls for quality control. Gluten-free and gluten-containing (naturally contaminated or spiked) samples should be used.
- In case of extremely acidic or basic samples, to adjust the sample's pH value to neutral (pH 6.5 to 7.5) prior to extraction may be necessary.

- To do spike experiments to ensure an accurate and correct test procedure. An example of a spiking experiment is given in the validation report.
- When carrying out the test, the specified temperature of 20 - 25 °C (68 - 77 °F) must be observed. Deviations hereof may influence the recovery of the analyte.
- To perform a PCR (e.g. SureFood®) for confirmation of the result.
- During the production of foods such as beer or sourdough, proteins are fragmented. In sandwich ELISAs, **protein fragments** lead to a reduced recovery. Such samples should be analyzed with a competitive ELISA like the **RIDASCREEN® Gliadin competitive** (Art. No. R7021).
- Try to homogenize samples gelatinizing during preparation with a clean spatula.
- To contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de) if automates (e.g. ThunderBolt® / Bolt™) are used.

## 15. Further application notes

- “Analysis of food samples without or with low oat content in RIDASCREEN® Total Gluten“ with RIDASCREEN® Total Gluten Additive TG (Art. No. RA0041)

**Further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Literature

[1] Lacorn et al. (2019) Quantification of Wheat, Rye, and Barley Gluten in Oat and Oat Products by ELISA RIDASCREEN® Total Gluten: Collaborative Study, First Action 2018.15. *J AOAC Int.*, **102**:5, p. 1535 - 1543.

## Version overview

Version number	Chapter and title
2018-11-12	Release version
2023-07-18	Changes made: <ul style="list-style-type: none"><li>– General linguistic revision</li><li>– Inclusion of AOAC-OMA certification</li><li>– Revision of chapter 5.1</li><li>– Modification of the disposal clause in chapter 6</li><li>– Revision of chapter 11 and creation of chapters 12, 13, 14</li><li>– Addition in chapter 11 Evaluation (moved from chapter "Recommendations to ensure high analytical safety")</li><li>– Inclusion of application note for the analysis of oat-free food samples by using Additive TG (RA0041) (chapter 1 and 15)</li></ul>
2024-05-14	Current version Changes made: <ul style="list-style-type: none"><li>– General linguistic revision</li><li>– Update of the content in "Related products for gluten / gliadin determination"</li><li>– Additions to content in chapters 1, 6, 12, 14</li></ul>

## Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM-DD)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM-DD)
	Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.