r-biopharm®

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# RIDASCREEN®FAST Gliadin sensitive

## Brief information

RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051) is a sandwich enzyme immunoassay for the quantitative analysis of prolamins from wheat (gliadin), rye (secalin) and barley (hordein) in food validated for the method (see chapter 1. Intended Use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation  
Cocktail (patented) (for 10 samples) .....2 h  
Cocktail ECO (for 10 samples)..... approx. 35 min  
test implementation (incubation time)..... 30 min

Standard material: The RIDASCREEN® standard material is the standard of the Prolamin Working Group (PWG-Gliadin, <http://www.wgpat.com/handling.html>)

Limit of detection: 0.2 mg/kg (ppm) gliadin or to 0.4 mg/kg (ppm) gluten\*  
(depending on matrix) 0.19 - 0.21 mg/kg (ppm) gliadin  
\*mean value

Limit of quantification: 1.25 mg/kg (ppm) gliadin or to 2.5 mg/kg (ppm) gluten

Specificity: The monoclonal antibody R5 reacts with the gliadin-fractions from wheat and corresponding prolamins from rye and barley.

Further information is contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments (see chapter 13. Limits of the method).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

## **Related product and accessories for gluten / gliadin determination**

RIDASCREEN® Gliadin (Art. No. R7001)

RIDASCREEN® Total Gluten (Art. No. R7041)

RIDA®QUICK Gliadin (Art. No. R7003 / R7004 / R7005)

RIDASCREEN® Gliadin competitive (Art. No. R7021)

RIDASCREEN®FAST Gliadin (Art. No. R7002)

Cocktail (patented) (Art. No. R7006 / R7016)

Cocktail ECO (Art. No. R7080)

RIDA® Extraction Solution (colorless) (Art. No. R7098)

Set of 3 processed Gliadin Assay Controls (Art. No. R7012)

SureFood® Allergen Gluten (Art. No. S3606)

SureFood® Allergen 4plex Cereals (Art. Nr. S7006)

### **1. Intended use**

RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051) is an R5 sandwich enzyme immunoassay based on specific monoclonal antibody to celiac toxic amino acid prolamin sequences to determine gliadin as a measure of gluten in food. Due to the large number of different food products, the following samples were examined representatively for different product groups within the scope of the test development: cookies, rice flour, oat flour, sausages, quinoa, chocolate, millet and various spices. It can be assumed, that the assay is also suitable for the analysis of other foods; this must be checked by the user himself.

For detailed results and further information on validation data with other food matrices please refer to the validation report.

The sample preparation using the Cocktail (patented) (R7006 / R7016) corresponds to the official R5-Mendez method according to the Codex Alimentarius. The faster sample preparation using the environmental-friendlier Cocktail ECO (R7080) is convenient for the screening of samples. The Cocktail ECO has an extraction efficiency of approx. 70 - 110 % compared to Cocktail (patented).

## 2. General information

The use of wheat flour and gluten in foodstuffs is extremely common because of their heat stability and useful effects on e.g. texture, moisture retention and flavour. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley. Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (CODEX STAN 118-1979) two categories for labeling of food according to the gluten content exist:

- 1) Food products which contain less than 20 mg/kg can be labeled as "**gluten-free**".
- 2) Food products labeled as "**very low gluten**" can have a gluten content above 20 and up to 100 mg/kg.

The threshold of 20 mg/kg has been adopted by many national legislations in many countries. The prolamin content of gluten (e.g. gliadin in wheat) is per definition generally assumed to be 50 % (CODEX STAN 118-1979).

## 3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific R5 antibodies against gliadins. By adding the standard or sample solution to the wells, gliadin present in the sample will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Following the washing step, a solution containing R5 antibody conjugated to peroxidase is added. This conjugate is bound to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue end product. A stop solution is added which results in a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance of the solution which is proportional to the gliadin concentration in the sample is measured photometrically at 450 nm and expressed as ng/mL gliadin.

## 4. Reagents provided

Each kit contains sufficient materials for a maximum of 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume
<b>Microtiter plate</b>	-	Ready to use	96 wells
<b>Buffer</b>	Transparent	Ready to use	2 x 100 mL
<b>Standard 1</b>	Transparent	Ready to use	0 ng/mL 1.3 mL
<b>Standard 2</b>	Transparent	Ready to use	2.5 ng/mL 1.3 mL
<b>Standard 3</b>	Transparent	Ready to use	5 ng/mL 1.3 mL
<b>Standard 4</b>	Transparent	Ready to use	10 ng/mL 1.3 mL
<b>Standard 5</b>	Transparent	Ready to use	20 ng/mL 1.3 mL
<b>Wash buffer</b>	Brown	<b>Concentrate</b>	<b>10x</b> 100 mL
<b>Conjugate</b>	Red	Ready to use	10 mL
<b>Substrate/Chromogen</b> Red Chromogen Pro	Brown	Ready to use	10 mL
<b>Stop solution</b>	Yellow	Ready to use	14 mL

## 5. Reagents required but not provided

### 5.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of  $\pm 0.01$  g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials with cap (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (50 °C / 122 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12  $\mu$ m)
- Graduated pipettes
- Graduated cylinder
- Variable 20 - 200  $\mu$ L and 200 - 1000  $\mu$ L micropipettes
- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100  $\mu$ L
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET (Art. No. Z9996FF)

## 5.2 Reagents

- Distilled water (dist. water) or deionized water
- Gluten-free skimmed milk powder (SMP) (food quality)
- Cocktail (patented) (R7006 / R7016)
- Alternatively: Cocktail ECO (R7080)
- Ethanol solution (**80 %**): e.g. add 120 mL ethanol p.a. to 30 mL dist. water and shake well

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

The Cocktail (patented) contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

Do not reuse wells of the microtiter strips (coated microtiter plate and pre-plate if necessary, see chapter 10.2). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- Bluish coloration of the reddish substrate/chromogen prior to test implementation.
- Value of less than 1.2 absorbance units ( $A_{450\text{ nm}} < 1.2$ ) for standard 5.

## 9. Sample preparation

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol.
- Carry out the sample preparation in a room isolated from the ELISA procedure.
- Check for gliadin contamination of reagents and equipment with the test strips RIDA®QUICK Gliadin (R7003 / R7004 / R7005).
- When using the Cocktail (patented), it is recommended to work under a chemical hood, because it contains  $\beta$ -mercaptoethanol.

The samples should be stored in a cool place, protected against light.

### 9.1 Extraction with Cocktail (patented) (R7006 / R7016, corresponds to the Codex Alimentarius method)

Homogenize well a sufficient amount (at least 5 g or 5 mL) of sample (grind it thoroughly to powder and mix well or mix well the solution respectively).

- **Liquid food samples:** pipet 0.25 mL of the homogenized sample into a vial (add 0.25 g of skimmed milk powder in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail (patented); close the vial and mix well.
- **Soft and hard food samples:** weigh 0.25 g of the homogenized sample (add 0.25 g of skimmed milk powder in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail (patented); close the vial and mix well.
- **Meat and sausages:** uneven gluten distribution is typical for these matrices; therefore, ensure to homogenize a sufficient sample amount (at least 50 g): weigh 0.25 g of the homogenized sample (add 0.25 g of skimmed milk

powder in case of tannin and polyphenol containing samples) into a vial and add 2.5 mL of Cocktail (patented); close the vial and mix well.

- **Oat samples:** uneven gluten distribution is typical for oats; furthermore, the samples are difficult to homogenize; therefore, homogenize 200 g and carry out the extraction with at least the fourfold amount of sample and reagents: weigh 1 g of the homogenized sample into a vial and add 10 mL of Cocktail (patented); close the vial and mix well.

**Please further extract all samples as described in the following:**

- Incubate for 40 min at 50 °C (122 °F) in a water bath.
- Let the sample cool down shortly (1 - 3 min).
- Add 7.5 mL 80 % ethanol (see chapter 5.2) (for oat samples: 30 mL 80 % ethanol).
- Close the vial and shake for 1 h upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter sample or centrifuge for 10 min at > 2,500 x g at room temperature (20 - 25 °C / 68 - 77 °F).  
(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The sample extracts (supernatant from centrifugation step or filtrate) can be stored undiluted in a well-sealed container at room temperature (20 - 25 °C / 68 - 77 °C) (shelf life approx. 4 weeks) until used in the test.
- For use in the test, extracts (supernatant of centrifugation step or filtrate) must be diluted with buffer (for dilution see chapter 10.2). The diluted sample extracts have a limited shelf life and must be used in the test within 30 minutes.

## 9.2 Extraction with Cocktail ECO (faster method)

The faster sample preparation using the environmental-friendlier Cocktail ECO (R7080) is convenient for the screening of samples. The Cocktail ECO has an extraction efficiency of approx. 70 - 110 % compared to Cocktail (patented). For this fast and environmental-friendlier sample preparation, it is **not** necessary to work under a chemical hood.

Prepare the necessary amount of Cocktail ECO (R7080) according to the product information.

Homogenize well a sufficient amount (at least 5 g or 5 mL) of sample (grind it thoroughly to powder and mix well or mix well the solution respectively).

- **Liquid food samples:** pipet 0.25 mL of the homogenized sample into a vial (add 0.25 g of skimmed milk powder in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail ECO; close the vial and mix well.
- **Soft and hard food samples:** weigh 0.25 g of the homogenized sample (add 0.25 g of skimmed milk powder in case of tannin and polyphenol containing samples) and add 2.5 mL of Cocktail ECO; close the vial and mix well.
- **Meat and sausages:** uneven gluten distribution is typical for these matrices; therefore, ensure to homogenize a sufficient sample amount (at least 50 g): weigh 0.25 g of the homogenized sample into a vial (add 0.25 g of skimmed milk powder in case of tannin and polyphenol containing samples) and add 2.5 mL of the Cocktail ECO; close the vial and mix well.
- **Oat samples:** uneven gluten distribution is typical for oats; furthermore, the samples are difficult to homogenize; therefore, homogenize 200 g and carry out the extraction with at least the fourfold amount of sample and reagents: weigh 1 g of a homogenized sample into a vial and add 10 mL of Cocktail ECO; close the vial and mix well.

**Please further extract all samples as described in the following:**

- Extract for 10 min at 50 °C (122 °F) in a water bath.
- Let the sample cool down shortly (1 - 3 min).
- Add 7.5 mL 80 % ethanol (see chapter 5.2) (for 1 g samples (e.g. oat samples): 30 mL 80 % ethanol).
- Close the vial and shake for 10 min upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Filter sample or centrifuge for 5 min at > 2,500 x g at room temperature (20 - 25 °C / 68 - 77 °F).  
(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The sample extracts (supernatant from centrifugation step or filtrate) can be stored undiluted in a well-sealed container at room temperature (20 - 25 °C / 68 - 77 °C) (shelf life approx. 2 weeks) until used in the test.

- For use in the test, extracts (supernatant of centrifugation step or filtrate) must be diluted with buffer (for dilution see chapter 10.2). The diluted sample extracts have a limited shelf life and must be used in the test within 30 minutes.

## 10. Test procedure

### 10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **wash buffer** is provided as a 10-fold concentrate. Before use, the buffer has to be diluted 1:10 (1 + 9) with distilled water (e.g. 900 mL dist. water + 100 mL buffer concentrate). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

Components should be stored immediately at 2 - 8 °C (35 - 46 °F) when no longer required.

### 10.2 Test procedure

**The extracts prepared according to chapter 9.1 and 9.2 must be diluted 1:12.5 (1 + 11.5) with buffer before usage (e.g. 920 µL buffer + 80 µL extract). The final sample dilution factor is 500.**

**Use the diluted sample extracts immediately (within 30 minutes) in the assay. A longer time period may influence the recovery.**

Carefully follow the recommended washing procedure to obtain unambiguous results. Do not allow microwells to dry between work steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µL per well). Then, exactly 100 µL are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.

2. Pipette 100 µL of each standard or sample (prepared according to chapter 9. and diluted according chapter 10.2) in duplicate to the wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
4. Add 100 µL of the conjugate to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times in a row) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µL diluted wash buffer (see chapter 10.1) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
6. Add 100 µL of substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Pipette 100 µL of stop solution into each well. Mix gently by shaking the plate manually and measure the extinction at 450 nm. Read within 10 min after addition of stop solution.

## 11. Evaluation

A special software, **RIDASOFT® Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. It is recommended to use RIDASOFT® Win.NET and cubic spline curve fitting procedure.

For the evaluation it should be clarified, that the quality criteria are fulfilled for the current test run. The course of the standard curve can be taken from the attached analysis certificate.

The assay can be also evaluated when running in single determinations. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not consistent with standards like EN 15633-1 and EN 15842. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

The gliadin concentration in ng/mL (ppb) is read from the RIDASOFT® Win.NET calibration curve and must be further multiplied by the dilution factor from the sample preparation of at least 500. This result is then multiplied by 2 in order to obtain the gluten concentration (gliadin represents 50 % of the proteins present in gluten, Codex definition). The RIDASOFT® Win.NET (version 1.93 or newer) indicates the results in gliadin and gluten.

### **Calculation example**

The extinction value of a sample corresponds to 10 ng/mL gliadin in the calibration curve. Multiplying by the recommended dilution factor 500 from the sample preparation leads to 5000 ng/mL, corresponding to 5 mg/kg (ppm) gliadin, respectively 0.0005 % gliadin. To calculate the gluten content, it is necessary to multiply by factor 2 which results in 10 mg/kg gluten, respectively 0.001 % gluten. This sample is considered to be gluten-free, because the gluten concentration is below 20 mg/kg.

## **12. Result interpretation**

Results between LOD and LOQ indicate a low analyte concentration in the sample. Calculated results show a high uncertainty in this area due to the method's high variation below LOQ. Therefore, such results should not be reported with a quantitative value, but qualitative as "< LOQ".

A result below the LOD does not exclude a gluten contamination below the detection limit of the assay, or that other cereal components, such as lipids, may be present in a sample. The result should be reported accordingly.

If samples show a higher absorbance than standard 5, the samples can be further diluted and used again in the ELISA. For this purpose, dilute the extract again freshly 1:12.5 with buffer as described in chapter 9.1 and 9.2. Then, further dilute the sample extract diluted in this way with the following mixture:

- 2 % Cocktail (patented) or Cocktail ECO (depending on extraction method)
- 6 % of an 80 % ethanol solution
- 92 % buffer

This corresponds e.g. to 100 µL Cocktail (patented) / Cocktail ECO, 300 µL 80 % ethanol and 4600 µL buffer. When using this mixture, the composition of the sample medium is maintained according to the extraction in chapter 9.1 or 9.2.

Compared to the certificate, higher absorbance values ( $A_{450\text{ nm}}$ ) for the standard curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

### **13. Limits of the method**

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

The protein content and protein composition of wheat, rye and barley cultivars may differ. Varying results are thus to be expected for different cultivars.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LOD, LOQ, upper limit of measurement range).

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max  $\pm 1$  %.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

The addition of foreign protein (depending on the test e.g. BSA, gelatine, skim milk powder) during extraction or test procedure may suppress matrix effects.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

Cross reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

The main epitope of the R5 antibody is the amino acid sequence QQFPF which is part of many celiac-toxic sequences. The sequence occurs repeatedly in the prolamins from wheat, rye and barley. However, rye and barley contain a higher number of replicates of this sequence, which leads to an overestimation of rye and barley compared to the wheat standard.

Recent research has shown that the real conversion factor from gliadin to gluten is approx. 1.5 [1], [2], [3]. However, the official conversion factor (Codex Alimentarius) is still 2.

## 14. Recommendation

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control. Gluten-free and gluten-containing (spiked) samples should be used.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- To do spike experiments to ensure an accurate and correct test procedure. An example of a spiking experiment is given in the validation report.
- To perform PCR (e.g. SureFood®) for confirmation of the result.
- To contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de) if automates (e.g. ThunderBolt® / Bolt™) are used.
- During the production of foods such as beer or sourdough, proteins are fragmented. In sandwich ELISAs protein fragments lead to a reduced

recovery, such samples should be analyzed with a competitive ELISA test systems like the RIDASCREEN® Gliadin competitive (Art. No. R7021).

## 15. Further application notes

- Sample preparation for processed food with the RIDA® Extraction Solution (colorless) (Art. No. R7098) - **only after validation**
- Sample preparation for raw materials with ethanol
- Sample preparation for polyphenol containing raw materials (e.g. chocolate, coffee, cacao, buckwheat) with fish gelatine and ethanol

**For further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Literature

- [1] Schall, E., Scherf, K. A., Bugyi, Z., Hajas, L., Török, K., Koehler, P., Poms, R. E., D'Amico, S., Schoenlechner, R., & Tömösközi, S. (2020) Food Chem. 313, 126049. <https://doi.org/10.1016/j.foodchem.2019.126049>.
- [2] Wieser et al.: Celiac Disease and Gluten. (2014) Elsevier Inc. Amsterdam, ISBN 978-0-12-420220-7, page 107.
- [3] Wieser, H. & Koehler, P. (2009) Eur. Food Res. Technol. 229, 9 - 13.

## Version overview

Version number	Chapter and title
2017-08-08	Release version
2023-02-09	Current version Changes made: <ul style="list-style-type: none"><li>– General linguistic revision</li><li>– Modification of disposal clause in chapter 6.</li><li>– Additional dilution for use in ELISA moved from chapter 9. to chapter 10.</li><li>– Revision and / or addition of chapter 9, 11., 12., 13. and 14.</li></ul>

## Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.