
r-biopharm®

Cocktail ECO

REF R7080

In vitro Test

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Storage at 2 - 8 °C (36 - 47 °F)

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Cocktail ECO

1. Intended use

The Cocktail ECO (Art. No. R7080) can be used for the sample preparation of raw materials, heat-treated and processed food.

The Cocktail ECO is used in combination with the following test systems:

- RIDASCREEN® Gliadin (Art. No. R7001)
- RIDASCREEN®FAST Gliadin (Art. No. R7002)
- RIDASCREEN®FAST Gliadin sensitive (Art. No. R7051)
- RIDA®QUICK Gliadin (Art. No. R7003)
- RIDA®QUICK Gliadin (single packaged) (Art. No. R7004)

The sample preparation using the **Cocktail (patented)** (Art. No. R7006/R7016) is the official R5-Mendez method according to the Codex Alimentarius and the AOAC.

The faster sample preparation using the more environmental-friendly **Cocktail ECO** is convenient for the screening of samples. The Cocktail ECO has an extraction efficiency of approx. 70 - 110 % compared to Cocktail (patented).

Sample preparation: homogenization, extraction and centrifugation

Time requirement: Cocktail ECO (for 10 samples)..... approx. 35 min

2. Reagents provided

Component	Cap color	Format	Volume
Buffer ECO	Green	Ready to use	2 x 115 mL
Additive ECO	Green	Ready to use	6 g

3. Reagents required but not provided

3.1 Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (50 °C / 122 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12 μm)

3.2 Reagents

- Distilled water (dist. water) or deionized water
- Gluten-free skim milk powder (SMP) (food quality)
- Ethanol solution (80 %): e.g. add 120 mL ethanol p.a. to 30 mL distilled water and shake well

4. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

The Cocktail ECO contains no hazardous substances. It is not necessary to work under a chemical hood. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com or at eifu.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

5. Storage instructions

Store the kit at 2 - 25 °C (35 - 77 °F). Do not freeze any test kit components.

Do not use the product after the expiration date (see test kit label).

6. Preparation and extraction of samples

Before starting and during the assay wear gloves. Airborne cereal dust and dirty laboratory equipment may lead to contamination of the assay. Hence, it is recommended to take the following precautionary measures:

- Clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol before and after each sample preparation,
- Carry out the sample preparation in a room separated from the ELISA procedure,
- Check for gluten contamination of reagents and equipment with the test strips RIDA[®]QUICK Gliadin (Art. No. R7003 / R7004 / R7005).

6.1 Preparation of the Cocktail ECO

To produce the Cocktail ECO it is necessary to add the Additive ECO to the Buffer ECO. Only mix the quantity of Cocktail ECO that is actually needed because the stability of the Cocktail ECO is only one day. The volume needed depends on the weighed portion. The minimal weighed portion is 0.25 g. **Due to the fact that the gluten contamination is often very heterogeneous, a weighed portion of 1 g is recommended.** For a weighed portion of 1 g use 10 mL Cocktail ECO and 30 mL 80 % Ethanol for the following sample preparation.

Number of samples and weigh in	Preparation Cocktail ECO		Needed volume per sample	
	Buffer ECO	Additive ECO	Cocktail ECO	80% Ethanol
1 sample à 0.25 g	3 mL	60 mg	2.5 mL	7.5
5 samples à 0.25 g	15 mL	300 mg	2.5 mL	7.5
10 samples à 0.25 g	27 mL	540 mg	2.5 mL	7.5
20 samples à 0.25 g	54 mL	1,08 g	2.5 mL	7.5
1 sample à 1 g	12 mL	240 mg	10 mL	30
10 samples à 1 g	120 mL	2,4 g	10 mL	30

6.2 Extraction with Cocktail ECO

Homogenize (grind thoroughly to powder and mix well or mix well a solution respectively) well a sufficient amount (e.g. 50 g or 50 mL) to ensure taking a representative test portion of sample.

- **Liquid food samples:** use 0.25 mL of the homogenized sample and add 2.5 mL of the Cocktail ECO, close the vial and mix well.
- **Other food samples (e.g. soy and quinoa containing samples):** weigh 0.25 g of the homogenized sample and add 2.5 mL of the Cocktail ECO, close the vial and mix well.
- **Tannin and polyphenol containing food samples (e.g. chocolate, coffee, cocoa, chestnut flour, buckwheat, millet and spices):** weigh 0.25 g of the homogenized sample, add 0.25 g of skimmed milk powder and add 2.5 mL of the Cocktail ECO, close the vial and mix well.
- **Meat and sausages:** in these matrices, gliadin may not be distributed evenly. Therefore, homogenize 50 g sample: weigh 0.25 g of the homogenized sample and add 2.5 mL of the Cocktail ECO, close the vial and mix well.
- **Oat samples:** gliadin may not be distributed evenly, furthermore the samples are difficult to homogenize. Therefore, homogenize 200 g, then carry out the extraction with at least the fourfold amount of reagents: weigh 1 g of the homogenized sample and add 10 mL of the Cocktail ECO, close the vial and mix well.

Please further extract all samples as described in the following:

- Incubate for 10 min at 50 °C (122 °F).
- Let the sample cool down and then mix it with 7.5 mL 80 % ethanol (for 1 g sample weigh in: 30 mL 80 % ethanol).
- Close the vial and shake for 10 min upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F).
- Centrifuge: 5 min, at least 2,500 x g, at room temperature (20 - 25 °C / 68 - 77 °F) and/or filter the extract.
(Alternatively: transfer 2 mL of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge).
- Transfer the supernatant in a screw top vial.

Remark

The supernatant obtained after the centrifugation or the filtrate can be stored in a tightly closed vial in the dark at room temperature (20 - 25 °C / 68 - 77 °F) up to two weeks.

When using the **ELISA RIDASCREEN® Gliadin** (Art. No. R7001), **RIDASCREEN®FAST Gliadin** (Art. No. R7002) und **RIDASCREEN®FAST Gliadin Sensitive** (Art. No. R7051) (see also the respective test kit insert):

- Dilute the sample 1:12.5 (1+11.5 / 80 µL + 920 µL) with buffer from the ELISA kit: the final dilution factor is 500 (Attention: for R7001 and R7002 the sample dilution buffer needs to be diluted 1:5 before use, see respective test kit insert)
- Use immediately 100 µL per well in the assay

When using the **lateral flow test strip RIDA®QUICK Gliadin** (see also respective test kit insert):

- 500 µL of the sample dilution buffer from the test kit are transferred to a test tube
- 50 µL of the extracted sample are added

7. Limits of the method

An incorrect weight of the sample to be analyzed will have an effect on the measurement result (e.g. a 10 % higher concentration is measured with a weight of + 10 %). A sufficient accuracy is given with a weight deviation of max +/- 1 %.

The addition of foreign protein (depending on the test e.g. BSA, gelatine, skim milk powder) during extraction or test procedure may suppress matrix effects.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

8. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842 (e.g. performing duplicate determinations).
- Pipette tips should be rinsed with standard or sample extract prior to pipetting.
- Carry along test controls for quality control. Gluten-free and gluten-containing (spiked) samples should be used.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 to 7.5) to neutral prior to extraction.
- To do spike experiments to ensure an accurate and correct test procedure. An example of a spike experiment is given in the validation report.
- To perform a PCR (e.g. SureFood®) for confirmation of the results.
- Proteins are fragmented during the production of food such as beer and sourdough. In sandwich ELISA, the recovery is reduced for fragmented proteins; these samples should be analyzed with a competitive ELISA test system such as the RIDASCREEN® Gliadin competitive (Art. No. R7021) therefore.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2017-09-11	Release version (Version 1)
2020-10-09	Version 2 Changes made: The table listed under point 6.1 has been changed. For "20 samples of 0.25 g each" the amount of "Buffer ECO" was changed from 56 ml to 54 ml in order to keep the ratio of Buffer ECO and Additive ECO equal. General formatting
2024-05-28	Current version (Version 3) Changes made: <ul style="list-style-type: none">– General revision– Modification of the disposal clause in chapter 4– Content additions in chapters 3, 4, 7 and 8– Addition of marking concerning patents of MORINAGA & Co., Ltd.

Explanation of symbols

General symbols:

-  Follow the instructions for use
-  Batch number
-  Expiry date (YYYY-MM-DD)
-  Storage temperature
-  Article number
-  Number of test determinations
-  Manufacturing date (YYYY-MM-DD)
-  Manufacturer + address

Patent Marking:

The extraction means in this product contains sulfite. Food inspection methods using a sulfite-containing extractant as in this product and/or corresponding detection kits are subject to the following patents of MORINAGA & Co., Ltd.: European Patent EP 2 224 239 B1, Australian Patent AU 2008 330 507 B2, United States Patent US 8 859 212 B2, Japanese Patent JP 5 451 854 B2. The patent holder has granted R-Biopharm AG a license to use, and sell products that employ, said protected technology in the above-mentioned territories.

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.