r-biopharm®

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RIDASCREEN®EASY Extraction tablets

1. Intended use

The extraction tablets are part of the EASY product line, but may also be enclosed in other test kits. Tablets of this accessory can be used for sample extraction or extract dilution with all test kits, which contain regularly these tablets. For using the tablets, the respective test kit instructions are effective. The tablets can be used independently of the actual lot number.

2. Reagents provided

50 extraction tablets

3. Warnings and precautions for the users

The product / test is only suitable within the scope of its intended use.

This accessory may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Furthermore, apply to warnings and precautions of the corresponding test kit instructions.

4. Storage instructions

Store the tablets at 2 - 25 °C (35 - 77 °F).

To avoid water condensation at the tablets, open the tube for withdrawal of tablets only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Do not use the tablets after the expiration date.

5. Note about extraction tablets

The tablets are identical with all EASY allergen ELISA. Extracts prepared with the tablets can therefore be used in the various EASY Allergen ELISA. The tablets contain all needed chemicals for extraction of foods or for additional dilution of high positive extracts (attention should be paid to the addition of further test-specific additives for certain samples).

For extraction, the tablet is added to each sample according to the respective test kit instructions for use. Then, a defined volume of dist. water is added and the sample is suspended by mixing until the tablet is decomposed. The water should be pre-heated to 60 °C (140 °F). This eases decomposition and guarantees the ideal extraction temperature from the beginning.

If desired, the extraction buffer can be prepared in larger volumes for the extraction of several samples. For this purpose, dissolve the desired number of extraction tablets (1 tablet per extraction) in the corresponding volume of dist. water. Then, add the appropriate amount to the prepared sample (homogenized and weighed).

Further dilution is recommended for highly positive samples that are outside the measuring range of the test. Additional extraction buffer must be prepared for further dilution. For this purpose, dissolve one extraction tablet in 20 mL preheated dist. water (60 °C / 140 °F). The corresponding dilution factor must be taken into account in the evaluation. The additional extraction buffer can also be used to dilute the spiking solution of the test kit.

Pay attention that the readily prepared buffer has only a limited expiry and its activity decreases. Hence, the buffer has to be prepared fresh daily.

The tablets contain insoluble, inert additives, which are needed for tablet pressing. These sediment usually with sample particles during centrifugation. When preparing additional buffer, these additives must be separated by filtration or centrifugation (5 min at 2,500 x g). It is recommended to decant the supernatant into a fresh vial in case of centrifugation.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2024-05-21	Release version

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Patent Marking:

The extraction means in this product contains sulfite. Food inspection methods using a sulfite-containing extractant as in this product and/or corresponding detection kits are subject to the following patents of MORINAGA & Co., Ltd.: European Patent EP 2 224 239 B1, Australian Patent AU 2008 330 507 B2, United States Patent US 8 859 212 B2, Japanese Patent JP 5 451 854 B2. The patent holder has granted R-Biopharm AG a license to use, and sell products that employ, said protected technology in the above-mentioned territories.

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.