Wheat Reference Material

REF RAA7011

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Literatur

- ^[1] Schall, E., Scherf, K.A., Bugyi, Z., Hajas, L., Török, K., Koehler, P., Poms, R.E., D'Amico, S., Schoenlechner, R., & Tömösközi, S. (2020) Food Chem. 313, 126049
<https://doi.org/10.1016/j.foodchem.2019.126049>
- ^[2] Lacorn et al. (2021) Determination of Gliadin as a Measure of Gluten in Food by R5 sandwich ELISA RIDASCREEN® Gliadin Matrix Extension: Collaborative Study 2012.01, J. AOAC Int (Publication open access: <https://doi.org/10.1093/jaoacint/qsab148>).
- ^[3] AOAC Appendix M Guidance on Food Allergen Immunoassay Validation, Official Methods of Analysis of AOAC International, AM-1. <https://doi.org/10.1093/9780197610145.001.0001>

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Wheat Reference Material

1. Intended use

The Wheat Reference Material (Art. No. RAA7011) can be used in gluten analysis. Some of the intended uses include, but are not limited to, application as a quality control sample (e.g. as a starting material for spiking matrices for the preparation of incurred samples or as a spiking solution), calibration material, material for method development, and method validation.

The material can be used in combination with the following R-Biopharm products:

RIDASCREEN® Gliadin (Art. No. R7001)

RIDASCREEN®EASY Gluten (Art. No. RAE7071)

RIDASCREEN® Total Gluten (Art. No. R7041)

RIDA®QUICK Gliadin (Art. No. R7003, R7004, R7005)

RIDA®QUICK Gluten quant. (Art. No. RAL7073)

2. General information

The Wheat Reference Material was developed, characterized, and produced by a consortium of universities and research institutes ^[1]. It is a flour mixture of five different wheat varieties from various origins (Germany, Canada, China, Australia, Hungary). A collaborative study has already been successfully conducted with this material ^[2], demonstrating its suitability for the described purposes.

Quality assurance plays a central role in routine testing in gluten analysis, as the results of these analyses have a direct impact on food safety and consumer protection. An essential part of quality assurance is the use of suitable controls. Positive and negative controls serve to monitor the functionality of the analysis and ensure the validity of the results. Stable and well-characterized controls are essential here. Reference materials are needed for this purpose. Since there is currently no certified reference material, there is an increased need for reliable alternatives.

In addition, new AOAC guidelines emphasize the importance of considering processed samples for the validation of test kits ^[3]. Food processing methods have an impact on the extraction of gluten from food matrices and its analysis using immunological methods. Compared to spiked food samples (gluten is added to the sample immediately before extraction), gluten in so-called

incurred samples undergoes all relevant food processing steps (gluten is added to the food prior to processing).

3. Reagents provided

2 g wheat flour (characterized reference material ^[1])

4. Warnings and precautions for the users

The product is only suitable within the scope of its intended use.

This product is not intended for animal or human consumption.

5. Storage instructions

Store the kit below -8 °C (18 °F).

Before use, bring the material to room temperature (20 - 25 °C / 68 – 77 °F) to minimize the risk of moisture absorption.

After removing the required amount, reseal the aluminum bag containing the reference material and store the material in the zipper bag again at below -8 °C (18 °F). Once opened, the material can be used until the indicated expiration date.

Do not use the test kit after the expiration date (see test kit label).

6. Application

6.1 Quality control sample

The material can be used in various ways as a quality control sample:

1. Extracted as a run control (without matrix)
2. Extracted for spiking a matrix
3. As flour for the production of processed samples (incurred samples)

The recovery may vary depending on whether the sample is spiked immediately prior to extraction or whether the reference material is added to the ingredients and processed together with them before extraction (incurred sample).

When used as a run control or for spiking a matrix, the material is first extracted as recommended for the method to be used.

It may be necessary to perform additional dilutions with the appropriate extraction medium of the extract to reach the measurement range of the method.

Note on incurred samples:

When preparing incurred samples, the wheat flour is added to the food matrices before processing.

One of the biggest challenges is to achieve a homogeneous distribution of the added gluten-containing wheat flour, as the required gluten concentrations are usually very low. It may therefore be advisable to follow one or more of the following recommendations, as also mentioned by the AOAC ^[3]:

- Premix the wheat flour with other ingredients of the sample to be produced that have similar properties to the wheat flour (dry powder) before adding the other ingredients, especially liquid ingredients.
- Contamination is carried out using a dilution series instead of a single high dilution. First, a master mix with a higher wheat flour concentration is prepared (e.g., 1 % wheat flour in rice flour), which is diluted with additional rice flour. Only afterwards, the other ingredients are added for preparation of the incurred matrix.
- If necessary, the reference material can also be suspended in water first and then added as a liquid component to the other ingredients of the incurred sample.

6.2 Method development

Extract the material according to the specified method.

It may be necessary to perform additional dilutions to achieve the appropriate calibration concentrations.

It may also be necessary to dilute in special buffers that ensure the long-term stability of the calibration agents.

6.3 Method and sample validation

For method and sample validation, the reference material can be used both for spiking the test materials and for preparing incurred samples (see 6.1). In this way, performance characteristics such as the recovery of gluten in the sample, the limits of detection and quantification, as well as accuracy and precision, can be determined and evaluated.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Literature

- [1] Schall, E., Scherf, K.A., Bugyi, Z., Hajas, L., Török, K., Koehler, P., Poms, R.E., D'Amico, S., Schoenlechner, R., & Tömösközi, S. (2020) Food Chem. 313, 126049
<https://doi.org/10.1016/j.foodchem.2019.126049>
- [2] Lacorn et al. (2021) Determination of Gliadin as a Measure of Gluten in Food by R5 sandwich ELISA RIDASCREEN® Gliadin Matrix Extension: Collaborative Study 2012.01, J. AOAC Int (Publication open access: <https://doi.org/10.1093/jaoacint/qsab148>).
- [3] AOAC Appendix M Guidance on Food Allergen Immunoassay Validation, Official Methods of Analysis of AOAC International, AM-1. <https://doi.org/10.1093/9780197610145.001.0001>

Version overview

Version number	Chapter and title
2026-07-01	Release version

Explanation of symbols

General symbols:

-  Follow the instructions for use
-  Batch number
-  Expiry date (YYYY-MM-DD)
-  Storage temperature
-  Article number
-  Number of test determinations
-  Manufacturing date (YYYY-MM-DD)
-  Manufacturer + address

Disclaimer

1. R-Biopharm AG provides a limited statutory warranty ("Gewährleistung") as mandated and pursuant to German law against defects in title and material. For products with a defined shelf life, the warranty applies for the stated shelf life. For use cycle limited products, it applies until the stated cycle limit is reached, but not beyond the general warranty period. For all other products, the general warranty period is twelve (12) months from the date the risk transfers.
Such warranty is subject to timely notice of the alleged defect, giving sufficient detail to assess the claim. A product's defectiveness is assessed based on the specifically agreed quality, its suitability for the contractually stated purpose, and whether all agreed-upon accessories, documentation, and instructions were provided.
2. R-Biopharm AG is not liable for issues arising from:
 - a) The failure to read, understand, or follow product instructions and information for use or safety
 - b) The use of unqualified personnel
 - c) The failure to follow industry standards such as Good Laboratory Practices
 - d) The use of unsuitable controls, sample matrices, or procedures
 - e) Other misuse
 - f) unilateral product modifications
 - g) Improper storage by others
 - h) External chemical, electromagnetic, mechanical, or electrolytic influences that exceed the defined parameters.
 - i) Harm caused by external events beyond R-Biopharm's control, including but not limited to burglary, theft, lightning, fire, water, or acts of nature.
3. R-Biopharm AG's liability is limited to foreseeable, contract-typical losses in case of negligent breaches of obligations that are essential to achieving the agreement's purpose and that the other party could reasonably rely on (material contract obligations). R-Biopharm AG is not liable for losses resulting from non-material contract obligations.
4. ALL OTHER WARRANTIES OR GUARANTEES OF ANY KIND, EXPRESS OR IMPLIED, ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, COURSE OF DEALINGS BETWEEN THE PARTIES OR OTHERWISE. R-Biopharm AG assumes no liability for consequential damages, in particular loss of profit, interruptions, or other indirect damages.
5. The foregoing is not intended to limit liability to the extent such liability cannot be excluded or limited under applicable German law or convention, including cases of intent, willful misconduct, or gross negligence, personal injury or death, express warranties or assumptions of procurement risks, strict liability, or other non-waivable legal mandates.
6. The foregoing liability provisions apply equally to R-Biopharm AG, its directors, officers, employees, and agents.

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© 2026 R-Biopharm AG
are registered trademarks of R-Biopharm AG.