RIDA®QUICK Gluten quant.

REF RAL7073

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RIDA®QUICK Gluten quant.

Brief information

RIDA®QUICK Gluten quant. (Art. No. RAL7073) is an R5-based immunochromatographic test for the quantitative determination of gluten from wheat, rye and barley on surfaces, CIP water and for the method validated food samples (see chapter 1).

The test kit contains 15 test strips, each of which can be used for one analysis.

The evaluation is carried out via the RIDA®SMART APP Allergen (see chapter 11).

Time requirement: Sampling for swab testapprox. 1 min
 Sample preparation for:
 10 CIP water samples.....direct use
 10 food samplesapprox. 20 min
 Test implementation (incubation times) 10 min

Calibration material: MoniQA wheat flour (10.6 % gluten).

Limit of detection: Surfaces:0.038 µg/mL gluten
 CIP water:0.072 µg/mL gluten
 Food samples:..... 0.7 mg/kg gluten

Limit of quantification: Surfaces:0.20 µg/mL gluten, RSD_i 22.8 %
 CIP water:0.25 µg/mL gluten, RSD_i 22.7 %
 Food samples:..... 3 mg/kg gluten, RSD_i 14.5 %

Measuring range: Surfaces:0.20 µg/mL - 4 µg/mL gluten
 CIP water:0.25 µg/mL - 4 µg/mL gluten
 Food samples:..... 3 mg/kg - 40 mg/kg gluten

Specificity: The applied **monoclonal antibody R5** reacts with the gliadin-fractions from wheat and corresponding prolamins from rye and barley. The test kit detects the gluten containing cereals einkorn, khorosan-wheat, durum wheat and emmer with recovery rates of approx. 14 % to 39 %. Other gluten containing cereals such as wheat (*triticum aestivum*), rye, barley, spelt and triticale show recovery rates within the AOAC specifications of 50 % to 200 %. No cross-reaction was detected in 107 tested foods. Further information can be found in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e. g. corn flour). In composed/processed food (e. g. maize bread) cross reactivities might be different. Interfering substances (e. g. polyphenols) can be detected by spiking experiments.

Related products for gluten/gliadin determination

RIDASCREEN® Gliadin (Art. No. R7001)
RIDASCREEN®EASY Gluten (Art. No. RAE7071)
RIDASCREEN® Gliadin competitive (Art. No. R7021)
RIDASCREEN® Total Gluten (Art. No. R7041)
RIDA®QUICK Gliadin (Art. No. R7003 / R7004 / R7005)
Cocktail (patented) (Art. No. R7006 / R7016)
RIDA® Extraction Solution (colorless) (Art. No. R7098)
RIDASCREEN®EASY Extraction Tablets (Art. No. RAA0008)
Set of 3 processed Gliadin Assay Controls (Art. No. R7012)
SureFood® ALLERGEN Gluten (Art. No. S3606)
SureFood® QUANTARD Allergen 40 (Art. No. S3301)
SureFood® ALLERGEN 4plex Cereals (Art. No. S7006)

1. Intended use

RIDA®QUICK Gluten quant. (Art. No. RAL7073) is an R5-based immunochromatographic test for the quantitative analysis of contaminations by prolamins from wheat (gliadins), rye (secalins), and barley (hordeins) on surfaces, in cleansing waters (CIP water) and in foods. Due to the large number of different foods, the following samples were examined as representative for food product categories within the scope of test development: baked goods, sauce (pesto), chocolate containing dessert, spices and rice flour.

It can be assumed that the assay is also suitable for the analysis of other foods. However, this must be verified by the user before applying the test kit to the respective food. For detailed results and further information on validation data with other food matrices please refer to our validation report.

2. General information

The use of wheat flour and gluten in foodstuffs is extremely common because of their heat stability and useful effects on e.g. texture, moisture retention and flavor. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley. The prolamin content (e.g. gliadin) of gluten is per definition generally assumed to be 50 % (CODEX STAN 118-1979) ^[1]. Hence, a factor of 2 is used for calculation of gluten from a measured gliadin concentration. However, recent research has shown that the real conversion factor from gliadin to gluten is approx. 1.5 ^{[2], [3], [4]}. The RIDA®QUICK Gluten quant. is not calibrated against gliadin, but against an exemplary wheat flour ^[2] and indicates the result in gluten. Therefore, the calculation factor is omitted. Hereby, lower values are determined being in better agreement with new research findings for the gluten content in comparison to an assay, which calculates the gluten content from a measured gliadin concentration by the factor 2.

Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (CODEX STAN 118-1979) foods for special dietary use for persons intolerant to gluten may contain up to 20 mg/kg gluten. This gluten concentration is usually tolerated by celiac patients.

The threshold of 20 mg/kg has been adopted by national legislations in many countries. Prolamins from wheat are named gliadins. Foods containing < 20

mg/kg gluten can be labelled "gluten-free". The prolamin content of wheat gluten is defined as 50 % (CODEX STAN 118-1979).

The RIDA®QUICK Gluten quant. is calibrated against an exemplary wheat flour and indicates the result in gluten. Only R-Biopharm AG is authorized to use the R5 antibody in immunochromatographic tests.

3. Test principle

The principle of the test in strip format is an antigen-antibody reaction and is based on the monoclonal R5-antibody, which is specific for the detection of gliadin from wheat and prolamins from rye and barley. In the lower area, the strip contains colloidal gold nanoparticle-R5-antibody-conjugates (gold conjugates). The test band contains immobilized R5 antibody. The strip is placed vertically in a running buffer containing the prepared sample. Following the chromatographic principle, the liquid (running buffer and sample) moves upwards inside the strip. Gold conjugates are solved by the liquid and pass together with it the test band. If the analyte is present, a sandwich is formed consisting of immobilized R5 antibody, gliadin and gold conjugate coloring the test band red. The control band should always be visible at the end of the reaction (defined incubation times) to prove the validity of the test. Upon rising gluten concentrations, the hook line becomes weaker. If the sample contains a very high gluten concentration, the hook line is only faintly visible. The intensity of the test band depends on the gluten concentration in the samples. It increases as the gluten concentration increases.

The evaluation is carried out with the RIDA®SMART APP Allergen.

4. Reagents provided

Each kit contains sufficient materials for 15 determinations. Each test kit contains:

Component	Cap color	Format	Volume
Test strip	White	Ready to use	15 pieces
Test tube	-		15 pieces
Extraction tablet	White	Ready to use	15 pieces
Disposable pipet	-		15 pieces
Conjugate 1	Yellow	Ready to use	1.8 mL
Conjugate 2	Blue	Ready to use	1.8 mL
Swabbing buffer	White	Ready to use	5 mL
Swabs	-		16 pieces

5. Materials required, but not provided

- RIDA®QUICK quant. evaluation kit (Art. No. RAL0001); contains test strip holder and RIDA®SMART APP Cover
- Latest version of the RIDA®SMART APP Allergen (available in the Google Play Store; Art. No. ZRSA)
- RIDA®SMART BOX (Art. No. ZRSA-SB)
- Android device (e. g. smartphone or tablet)
- Gloves
- Shaker / Vortexer

For food samples:

- Scale (measurement range at least up to 50 g and precision of ± 0.01 g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- 60 % ethanol (e.g. add 150 mL ethanol (99.5 % denatured) to 100 mL dist. water)
- Centrifugal vials (e.g. 50 mL centrifuge tubes from Greiner Art. No. 227261)
- Fluted filter (pore size 8 - 12 μm)
- Centrifuge (at least 2,500 x g) for corn, tannin and polyphenol containing food samples
- Optional: Variable 20 - 200 μL and 200 - 1000 μL micropipettes

- If required: Gluten-free skim milk powder (SMP) (food quality) for extraction according 9.4
- If required: additional extraction tablets (Art. No. RAA0008)

6. Warnings and precautions

The product / test is only suitable within the scope of its intended use.

This test should only be carried out by trained personnel. Always strictly adhere to the instructions for use of this test.

Do not remove the test strips from the test strip tube until immediately before use in the test. Discard test strips that have not been removed immediately before use, as they may lead to incorrect results.

To reduce the risk of contamination, use the test tubes once.

This kit may contain harmful substances. Please refer to the component safety information in the material safety data sheets (MSDS) for this product, available online at www.r-biopharm.de.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

The extraction tablets are not batch specific. They can also be used in the RIDASCREEN®EASY Allergen ELISA assays. Further, do not interchange individual reagents between kits of different lot numbers.

7. Storage instructions

Store the test at 2 - 8 °C (36 - 46 °F). Do not freeze components of the test kit.

The test strips are sensitive to humidity. Humid test strips may influence the test results negatively. For this reason, keep the strips away from humidity. This applies particularly when the test strip tubes are opened. Therefore, do not remove test strips until immediately before use and after reaching room temperature. The test strip tubes must be closed immediately after the removal of a test strip.

To remove extraction tablets, do not open the container until it has reached room temperature (20 - 25 °C / 68 - 77 °F) to avoid the formation of condensation on the tablets. The extraction tablets can also be stored in the sealed container at room temperature (20 - 25 °C / 68 - 77 °F) until the expiry date printed on the container.

No quality guarantee is accepted after expiry of the kit (see kit label).

8. Indication of instability or deterioration of reagents

- A missing control band

9. Sample preparation

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol.
- Carry out the sample preparation in a room isolated from the LFD procedure.
- Check for gluten contamination of reagents and equipment with the test strips RIDA®QUICK Gliadin (R7003 / R7004 / R7005).

The RIDA®QUICK quant. Gluten Lateral Flow contains a tablet with all needed chemicals for food extraction instead of a liquid extraction buffer concentrate.

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

Perform sample preparation at room temperature (20 - 25 °C / 68 - 77 °F).

9.1 Swabbing

- Take as many test tubes as there are surfaces to be analyzed.
- Add 4 drops each of conjugate 1 (yellow cap) and conjugate 2 (blue cap) to the test tube and mix. Alternatively, the dropper can be removed and 100 µL can be pipetted.

- Moisten the swab with 8 drops of the swabbing buffer and then thoroughly swab the sampling area (e.g. 10 x 10 cm) in a cross-hatch technique (wear gloves!).
- After taking the sample, place the swab into the solution in the test tube from step 1 (mixture of conjugate 1 + 2) and rotate the swab a few times in the liquid. In addition, squeeze the swab by pressing the test tube together. This process should take approx. 15 seconds.
- Pull out the swab while squeezing the test tube to pour out the liquid and collect it in the tube.
- Start the timer immediately and incubate for 5 minutes.

9.2 Cleaning water (CIP water)

- Take as many test tubes as there are samples to be analyzed.
- Add 4 drops each of conjugate 1 (yellow cap) and conjugate 2 (blue cap) to the test tube and mix. Alternatively, the dropper can be removed and 100 µL can be pipetted.
- Add 1 drop – with provided disposable pipette – of CIP water (or 20 µL).
- Mix thoroughly.
- Start the timer immediately and incubate for 5 minutes.

If foam forms, allow the foam to settle first. If a bubble forms on the neck of the tube, wait a few seconds until it bursts on its own or carefully burst it with a clean device (e.g. using a clean pipette tip).

9.3 Food samples

Homogenize (carefully crush, finely grind and mix well or mix the solution well) a sufficient quantity of the food (e.g. 50 g or 50 mL) to ensure that a representative sample is taken.

- Weigh 1 g of the homogenized sample (or 1 mL from a liquid sample) into a 50 mL vial.
- Add 1 extraction tablet.
- Add 10 mL 60 % ethanol (or 9 mL in case of liquid samples).

- Extract at room temperature (20 - 25 °C / 68 - 77 °F) by the following scheme: incubate for 2 minutes, mix vigorously for 20 seconds (e.g. vortexer), incubate for 2 minutes, mix vigorously for 20 seconds.

Caution: The extraction tablet does not dissolve completely; insoluble additives are separated by filtration/centrifugation in the next step.

- Filter the extract into a new, empty tube (filtration should be rapid due to possible evaporation of the ethanol).

Alternatively, the extract can also be centrifuged (5 min at at least 2,500 x *g*). A higher *g*-force (e.g. 5,000 x *g*) results in better sedimentation. This increases the precision of the measurement. If available, the use of a centrifuge with a higher *g*-force is therefore recommended. The duration of the centrifugation remains the same.

- The extract (filtrate or supernatant) can be stored for one day at 20 - 25 °C (68 - 77 °F), two days at 2 - 8 °C or for one month at -20 °C (-4 °F).
- Set up as many test tubes as there are samples to be analyzed.
- Add 4 drops each of conjugate 1 (yellow cap) and conjugate 2 (blue cap) to the test tube and mix. Alternatively, the dropper can be removed and 100 µL can be pipetted.
- Add 1 drop – with provided disposable pipette – of extract (or 20 µL).
- Mix thoroughly.
- Start the timer immediately and incubate for 5 minutes.

9.4 Corn, tannin and polyphenol containing food samples and soy-based milk alternatives

Homogenize (carefully crush, finely grind and mix well or mix the solution well) a sufficient quantity of the food (e.g. 50 g or 50 mL) to ensure that a representative sample is taken.

- Weigh 1 g of homogenized sample (or 1 mL from a liquid sample) into a 50 mL vial.
- Add 1 extraction tablet and 1 g of skim milk powder (SMP).
- Add 10 mL 60 % ethanol.

- Extract at room temperature (20 - 25 °C / 68 - 77 °F) by the following scheme: incubate for 2 minutes, mix vigorously for 20 seconds (e.g. vortexer), incubate for 2 minutes, mix vigorously for 20 seconds.

Caution: The extraction tablet does not dissolve completely; insoluble additives are separated by centrifugation in the next step.

- Due to the additional SMP, there is a risk that the filter will clog and filtration will not be possible. Therefore, prior centrifugation is necessary for this extraction (5 min, 2,500 x g, room temperature (20 - 25 °C / 68 - 77 °F)). Alternatively, allow the sample to settle for 1 - 2 minutes (sedimentation) and only use the resulting supernatant for filtration.
- Filter the supernatant into a new, empty tube (filtration should be rapid due to possible evaporation of the ethanol).
- The extract can be stored for one day at 20 - 25 °C (68 - 77 °F), two days at 2 - 8 °C or for one month at -20 °C (-4 °F).
- Set up as many test tubes as there are samples to be analyzed.
- Add 4 drops each of conjugate 1 (yellow cap) and conjugate 2 (blue cap) to the test tube and mix. Alternatively, the dropper can be removed and 100 µL can be pipetted.
- Add 1 drop – with provided disposable pipette – of filtered extract (or 20 µL).
- Mix thoroughly.
- Start the timer immediately and incubate for 5 minutes.

10. Test implementation

- Place the test strip into the test tube. Make sure the test strip is resting at the bottom of the test tube.

Caution: The test strip must be placed inside the test tube that the front (reaction side) does not touch the test tube's wall!

- Incubate for 5 minutes.
- Remove the test strip from the test tube and remove an excess of liquid (e.g. by dabbing on tissue paper).

Caution: The evaluation must be carried out within 30 seconds after removing the test strip from the sample.

11. Evaluation

The evaluation is carried out with the RIDA®SMART APP Allergen.

Please read the instruction for use of the RIDA®SMART APP Allergen attentively. The test strip evaluation has to be performed by using the RIDA®SMART APP Allergen according to these guidelines.

Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP Allergen on our website www.r-biopharm.de.

The gluten content of each sample is determined automatically in the RIDA®SMART APP Allergen using the 4-parameter function by reference to the stored calibration curve.

The batch specific information regarding the calibration curve is stored in the QR-code on the enclosed certificate of analysis (CoA). When using a new batch, the batch specific QR-code must be scanned and activated by the RIDA®SMART APP Allergen.

Caution: A visual evaluation is not possible!

The accessory kit “RIDA®QUICK quant. evaluation kit” (Art. No. RAL0001) is required for the evaluation. This contains the test strip holder and the RIDA®SMART APP cover. The latter is only required for solely smartphone evaluation (chapter 11.2).

For evaluation according to 11.1 and 11.2, the test strip must be inserted into the test strip holder (direction see fig. 1).


Fig. 1: Test strip alignment for evaluation according to 11.1 and 11.2.

For evaluation according to 11.2, the test strip must be inserted into the test strip holder (direction see fig. 1) and further covered with the RIDA®SMART APP Cover (fig. 2).


Fig. 2: Test strip preparation for evaluation according to 11.2.

The left (highest) line is the control line and must appear after each test run. If the line is missing, the test was not performed properly or the reagents were not functional. In this case, the test must be repeated with a new test strip. If the control line is repeatedly missing, please inform R-Biopharm.

The middle line is the hook line on which the target analyte is immobilized. The intensity of the hook line depends on the prolamin concentration of the sample. The higher the prolamin concentration, the weaker the intensity.

The right (lowest) line is the test line. The intensity of the test line depends on the prolamin concentration in the sample. The higher the prolamin concentration, the higher the intensity.

11.1 Evaluation RIDA®SMART APP Allergen + RIDA®SMART BOX + Android device

- Place the test strip in the test strip holder.
- Place it in the RIDA®SMART BOX and close the drawer.
- Choose the application (Surface Test, CIP Water Test or Food Test) and follow the instructions of the RIDA®SMART APP Allergen.
- The test strip is recorded via the RIDA®SMART BOX.
- The evaluation is carried out with the RIDA®SMART APP Allergen.

Caution: The evaluation must be carried out within 30 seconds after removing the test strip from the sample.

11.2 RIDA®SMART APP Allergen + verified Android Smartphone

(not part of the validation)

As an alternative to the use of the RIDA®SMART BOX, a verified smartphone can be used to record the test strip. Please contact R-Biopharm or visit our website for more information about the verified smartphones.

- Place the test strip in the test strip holder and cover with the RIDA®SMART APP cover.
- Choose the application (Surface Test, CIP Water Test or Food Test) and follow the instructions of the RIDA®SMART APP Allergen.
- The test strip is then recorded with a verified smartphone using the RIDA®SMART APP Allergen.
- The test strip should be within the specified frame when recording. The distance between smartphone and test strip should be approx. 10 cm.
- The evaluation is carried out via the RIDA®SMART APP Allergen.

Caution: The evaluation must be carried out within 30 seconds after removing the test strip from the sample.

Procedure for samples > measuring range:

Additional dilution is recommended for highly positive samples that are outside the measuring range of the test. The corresponding dilution factor must be taken into account during evaluation. If additional extraction buffer is needed for dilution of high positive samples being above the assay's measuring range, one tablet can be solved in 10 mL of 60 % ethanol. Mix thoroughly by the following scheme: incubate for 2 minutes, mix for 20 seconds (e.g. vortexer), incubate for 2 minutes, mix for 20 seconds. Tablets contain insoluble, inert additives, which are needed for tablet pressing. These sediment usually with sample particles during centrifugation. When preparing additional buffer, these additives must be separated by filtration or centrifugation (5 min at 2,500 x g). In case of centrifugation, it is recommended to decant the supernatant into a fresh vial.

It should be noted that the ready-to-use buffer only has limited shelf life. Therefore, the buffer must be prepared daily.

12. Result interpretation

The test is calibrated against gluten. The result therefore indicates the amount of gluten in mg per kg food (**gluten in mg/kg**) for foods and the amount of gluten in µg per mL buffer/sample (**gluten in µg/mL**) for hygiene samples (surfaces and CIP water).

Results between LoD (limit of detection) and LoQ (limit of quantification) indicate a low allergen concentration in the sample and such results should not be reported with a quantitative value, but qualitative as “< LoQ”.

A result below the LoD does not exclude an allergen contamination below the detection limit of the assay, or that other allergen components, such as lipids, may be present in a sample. The result should be reported accordingly.

Results below the limit of quantification (LoQ) are by default reported as < 3 mg/kg (food samples), < 0.20 µg/mL (surfaces) or < 0.25 µg/mL (CIP water). If necessary, for example for the verification of unknown or new matrices, the function “Display results < LoQ” can be activated in the app under settings. This allows numerical values below the LoQ to be output. The use of these values is at your own risk. They should be interpreted with caution and are not suitable for decision-making. Such results can be distorted by background noise or measurement deviations, among other things. The measurement uncertainty increases as the concentration decreases.

Note to surface tests: user cannot always swab a defined area (e.g. the recommended 100 cm²) but rather swab in corners or hard-to-reach spots. The result output of the RIDA®QUICK Gluten quant. gives information what is finally detected in the test tube. Result output in µg/mL simplifies result interpretation by giving information about the concentration in the tube. Customer can then recalculate according to their specific swabbed area to µg/cm². Data acquisition at R-Biopharm AG is always done with 100 cm². For example, 5 µg/mL was measured (in the tube). The swabbed area was 100 cm², resulting in a concentration of 0.05 µg/cm². The volume in the tube is always fixed, due to the volume of conjugate 1 + 2 and the washed swab. This volume does not correspond to 1 mL. However, the RIDA®SMART APP Allergen measures the concentration inside this volume and the result output is given in µg/mL.

13. Limits of the method

For gluten concentrations > 40 mg/kg (food samples) or > 4 µg/mL (hygiene samples), the RIDA®SMART APP Allergen reports the value > 40 mg/kg or > 4 µg/mL, respectively. If the sample is highly contaminated with gluten, the hook line/test line is only faintly visible or not visible at all. In this case, the RIDA®SMART APP Allergen displays the message “Hook effect, calculation not possible. Please dilute”. Please repeat the test run with a diluted sample and a new test strip (see chapter 11).

If the control line is missing, the message “INVALID” appears and the test must be repeated.

Test results may vary depending on the sample matrix, sample preparation, test performance and laboratory conditions. The general quality assurance requirements for laboratories, as listed in the standards EN 15633-1 and 15842, help to recognize such fluctuations (e.g. by performing duplicate determinations). It is therefore recommended to determine the precision of the method by repeated testing, especially for unknown samples.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max. ± 1 %.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

For the present quantitative immunochromatographic test, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

The addition of foreign protein (depending on the test e.g. BSA, gelatine, skim milk powder) during extraction or test procedure may suppress matrix effects.

Cross reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

Food commodities with strong swelling properties such as bovine gelatine, fish gelatine, porcine gelatine, guar gum and xanthan gum have to be diluted with gluten-free rice flour in a ratio of 1:10 (1+9) before extraction. Such foods pass through the membrane more slowly and can clog the pores. This can lead to a lower recovery or to an invalid test result. It should be noted that additional dilution of the sample is associated with a loss of sensitivity of the test. With an additional dilution factor of 10, the limit of quantification is 30 mg/kg gluten.

Strongly yellow-colored foods, such as curcuma or curry, can color the membrane of the test strip and thus impair line detection. The results are displayed as invalid. In these cases, it is recommended to dilute the commodity at least 1:10 (1+9), e.g. with gluten-free rice flour. It should be noted that additional dilution of the sample is associated with a loss of sensitivity of the test. With an additional dilution factor of 10, the limit of quantification is 30 mg/kg gluten.

Strongly acidic or alkaline cleaning agents can have a strong impact on the workflow of the test. For this reason, they must be diluted with water or adjusted to a neutral pH of 6.5 - 7.5 before use. Nevertheless, the protein structure of the gluten can be permanently changed by acid or alkali and may no longer be detected.

The protein content and protein composition of wheat, rye and barley cultivars may differ. Varying results are thus to be expected for different cultivars.

The main epitope of the R5 antibody is the amino acid sequence QQPFP, which is part of many celiac-toxic sequences. The sequence occurs repeatedly in the prolamins from wheat, rye and barley. However, rye and barley contain a higher number of replicates of this sequence, which leads to an overestimation of rye and barley compared to the wheat standard.

The limit of detection (LoD) is an informative value that has been determined during the validation process. This value cannot be determined due to the specified measuring range.

14. Recommendations

In order to ensure a high analytical performance, we recommend:

- Pre-flush pipette tips with sample extract prior to pipetting.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction.
- To do spiking experiments to ensure an accurate and correct test procedure. An example for the spiking procedure is given in the validation report.
- To perform ELISA (e.g. RIDASCREEN®) or PCR (e.g. SureFood®) for confirmation of the result.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Literature

- [1] FAO Codex Alimentarius International Food Standards: Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten. CXS 118-1979. Adopted in 1979. Amended in 1983 and 2015. Revised in 2008.

Version overview

Version number	Chapter and title
2024-09-20	Release version
2025-05-09	Previous version
2026-02-04	Current version Changes made: <ul style="list-style-type: none">– Chapter 2: Additional information on the conversion factor– Test procedure: Additional information on the exact amount of conjugate 1 and 2 when using a pipette– Chapter 9: Note to carry out sample preparation in a room isolated from LFD procedure.– Chapter 12: adjustment of result output in the app; no numerical results < LoQ by default. Additional functionality in the app to display results < LoQ.– Chapter 12: Note on the interpretation of surface tests

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Patent Marking:

The extraction tablet of the test kit contains sulfite. Food inspection methods using a sulfite-containing extractant as used in this test kit and/or corresponding detection kits are subject to the following patents of MORINAGA & Co., Ltd.: European Patent EP 2 224 239 B1, Australian Patent AU 2008 330 507 B2, United States Patent US 8 859 212 B2, Japanese Patent JP 5 451 854 B2. The patent holder has granted R-Biopharm AG a license to use the protected technology in said territories.

Disclaimer

1. In conformance with the German Civil Code ("BGB") R-Biopharm AG provides a limited warranty ("Gewährleistung") against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.
2. ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm's own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
3. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
 - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
 - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
 - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
 - d. any combination thereof.
4. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
5. R-Biopharm AG's liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG's sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser's expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
6. The foregoing limited warranty is solely intended to fulfill the warranty requirements ("Gewährleistung") implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.

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are registered trademarks of R-Biopharm AG.