CONGEN

SureFood® ANIMAL QUANT Beef

Art. No. S1010
2 x 50 rxn

User Manual

September 2025

 Inhalt

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Content

1	General Information	10
1.1	Description	10
1.2	Limit of Detection and Limit of Quantification	10
1.3	DNA-preparation	11
1.4	Kit components and storage	11
1.5	Additionally, required equipment and materials	11
1.6	Precautions for users	11
1.7	Setup.....	12
1.8	Detection channel Set-up	12
2	Qualitative Analysis	13
2.1	Protocol	13
2.1.1	Preparation of the master-mix	13
2.1.2	Preparation of the standard DNA dilutions	14
2.1.3	Preparation of the real-time PCR-mix	14
2.2	Interpretation of results	15
3	Limitations of the method	16
4	Further Information	16
4.1	Product Information	16
4.2	Technical Support	16
4.3	Distribution and Ordering	16

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1 General Information

1.1 Description

The SureFood® ANIMAL QUANT Beef real-time PCR detects the relative beef (*Bos taurus*) DNA amount to the total animal (amniotes) DNA content in meat samples. Therefore, the kit contains two PCR systems, one for detection of a beef-specific gene (**Beef** detection gene) and one for the detection of an animal gene (**Reference** gene).

The real-time PCR assay can be used with established real-time PCR instruments. The internal technical verification of instruments was performed on Roche LightCycler® 480 II, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, Agilent Mx3005P.

1.2 Limit of Detection and Limit of Quantification

The beef specific real-time PCR has a limit of detection of ≤ 5 DNA-copies.

The assay limit of detection depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The limit of quantitation depends on the concentration of the sample DNA used in the analysis. For example, if 50,000 target-sequence copies of the reference gene are present, the relative quantitation limit for beef DNA is 0.04%.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

SureFood® ANIMAL QUANT Beef

(2 x 50 rxn)

Art. No. S1010

September 2025

1.3 DNA-preparation

For DNA-preparation of meat samples the use of SureFood® PREP Basic (Art. No. S1052) is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reference Reaction Mix	2 x 450 µl	Orange
2	Beef Reaction Mix	2 x 450 µl	Yellow
3	Taq Polymerase	1 x 11 µl	Red
4	Dilution Buffer	1 x 1300 µl	White
5	Standard DNA	1 x 25 µl	Dark Blue
6	Positive Control (100% Beef)	1 x 25 µl	Light Blue

Store all reagents at -28 to -16°C and protected from light. The tube of the Taq Polymerase should not be thawed.

1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Basic Art. No. S1052)
- real- time PCR instrument
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

SureFood® ANIMAL QUANT Beef

(2 x 50 rxn)

Art. No. S1010

September 2025

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	10 sec, 95°C	5 sec, 95°C
Annealing	15 sec, 62°C	10 sec, 62°C
Extension (CYCLE)	30 sec, 65°C	15 sec, 65°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	Beef	FAM	+	
	Reference	FAM	+	
Applied Biosystems 7500	Beef	FAM	TAMRA	Check the passive reference option ROX is none.
	Reference	FAM	TAMRA	
Bio-Rad CFX96 / Dx / Opus	Beef	FAM	+	Baseline Settings - Baseline subtracted curve fit - Apply fluorescence drift correction
	Reference	FAM	+	
Qiagen Rotor-Gene Q	Beef	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	Reference	green	+	
Roche LightCycler® 480 II	Beef	465-510	+	For the quantitative analysis use the type "Fit Points"
	Reference	465-510	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples, control reactions and standards) for the specific PCR assay.

Required reactions for the evaluation of the reference detection (see 2.2):

5 reactions for the standard curve

3 reactions for controls* (1x negative control, 2x Positive Control)

For each sample: at least 1 reaction for each sample DNA

Required reactions for the evaluation of beef detection (see 2.2):

5 reactions for the standard curve

3 reactions for controls* (1x negative control, 2x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10% additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

The tube of the Taq Polymerase should not be thawed and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	18.0 µl	198.0 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	18.1 µl	199.1 µl

Mix each master-mix well and centrifuge shortly before use.

***Description of the controls**

- Negative control: only master-mix
- Positive Control: master-mix and within the kit's provided Positive Control

2.1.2 Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different DNA concentrations for the standard curves of the reference gene (**Reference**) and the detection gene (**Beef**). Prepare 5 dilutions of the supplied Standard DNA (**Code 5**) with the supplied Dilution Buffer (**Code 4**). Prepare 5 reaction tubes (labeled S1 to S5) and add 45 µl Dilution Buffer (**Code 4**) each.

The following procedure is recommended:

Standard	Dilutions	Copy number per µl	Final copy number per reaction**
S1	45 µl Dilution Buffer + 5 µl Standard DNA	100,000 copies	200,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	20,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1,000 copies	2,000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	200 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S4	10 copies	20 copies

****Note:** 2 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

2.1.3 Preparation of the real-time PCR-mix

- Pipette 18 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 2 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 2 µl of Positive Control and the standard dilutions into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The calculation for both reaction systems (**Reference, Beef**) has to be made separately. Mark the standards, the controls and the samples for the specific system (**Beef**) and make the evaluation according to the usual analysis program recommended by the real-time PCR instrument manufacturer. Repeat the same procedure for the reference gene system (**Reference**). The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

By using the calculated copy numbers for **Reference** and **Beef** the relative beef content of the sample DNA and the Positive Control can be determined in the following way (example):

Sample Beef	1,350 copies	Positive Control Beef	900 copies
Sample Reference	45,000 copies	Positive Control Reference	1,100 copies

Divide the copy number of the specific system by the copy number of the reference gene system and multiply by 100 to obtain the percentage.

Beef DNA content = **Beef** copy number * 100 / **Reference** copy number

sample **Beef** DNA content = $1,350 * 100 / 45,000$

sample **Beef** DNA content = 3%

For the given example the numbers lead to a **Beef** DNA content of 3.0%, for the sample and 81.8% for the Positive Control with the same calculation.

For a final calculation the use of a correction factor K for the correction of run-to-run fluctuations is necessary. The correction factor is the relation of the true percentage value of the Positive Control and the measured percentage of the Positive Control. The true value of the Positive Control is 100% **Beef** content.

$K = \text{true value} / \text{measured value}$

$K (\text{example}) = 100\% / 81.8\% = 1.2$

The measured DNA content for the sample is multiplied with K to obtain a corrected DNA content.

sample DNA content = measured sample DNA content * K sample (example) = $3.0\% * 1.2 = 3.6\%$

For that example the **Beef DNA content is 3.6%**.

3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- Innards, particularly kidney or liver, show a higher DNA content per weight unit when compared to muscle tissue. This could lead to over-estimation of the species content of the animal from which the innards come from. For meat products with innards content >5 % care must be taken when relative species contents are reported.
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.
- Cross reactivity was observed with DNA extracts from bison (*Bison* sp.) with 100 % and (*Bubalis arnee*) with 0.6%.

4 Further Information

4.1 Product Information

- Microsoft Excel template of calculation and detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

4.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

4.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com