CONGEN

SureFood® PREP Basic

Art. No. S1052
100 extractions

User Manual

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1 General Information

1.1 Description

This kit is intended to be used for the extraction of animal and plant DNA (deoxyribonucleic acid) from raw materials as well as slightly processed food and feed. It is also recommended for the extraction of animal DNA from highly processed food and feed.

1.2 Kit components and storage

Kit Code	Reagent / Material	Amount (per box)
K	Proteinase K	1 Tube
H	PCR grade water	1 x 1.1 ml
L	Lysis Buffer	1 x 20 ml
B	Binding Buffer*	1 x 15 ml
P	Pre-Wash Buffer**	1 x 60 ml
W	Wash Buffer**	1 x 60 ml
E	Elution Buffer	1 x 10 ml
F	Spin Filter, for pre-filtration	1 x 50 Filters
S	Spin Filter, for binding	1 x 50 Filters
R	Receiver Tubes 2,0 ml	2 x 50 Tubes
T	Receiver Tubes 1,5 ml	1 x 50 Tubes

*After adding isopropanol (not supplied with the kit, see 1.3 Additionally required equipment and materials and 2.2 Preparations)

**After adding ethanol (not supplied with the kit, see 1.3 Additionally required equipment and materials and 2.2 Preparations)

The components are equally divided into two boxes.

The Proteinase K should be stored at -28 to -16°C. All other reagents of the kit should be stored dry and at room temperature (14-25°C).

Before every extraction make sure that all components have room temperature. If there are any precipitates within the provided solutions, solve them by warming carefully (up to 30°C).

1.3 Additionally required equipment and materials

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 65°C)
- micro centrifuge (up to 12,000 rpm)
- ethanol for preparation of Wash Buffer (puriss., purity ≥ 96%)
- isopropanol for preparation of Binding Buffer (e.g. Carl Roth 2-Propanol Rotipuran® >99.7; Applichem 2-Propanol for microbiologie; Sigma-Aldrich 2-Propanol)
- disposal bags or waste bin

A flow chart for visualization of the extraction protocol is available at the CONGEN homepage: <http://www.congen.de/en/downloads>.

1.4 Precautions for users

- This test must only be performed by qualified and authorized laboratory personnel.
- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.
- For more details see Material Safety Data Sheet, at www.congen.de/en/eifu/.

2 Protocol

2.1 Principle

1. Preparation of the basic material
2. Lysis of the basic material
3. Pre-filtration and setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

2.2 Preparations

General

Add 1 ml PCR grade water (**Code H**) to the Proteinase K (**Code K**) and mix thoroughly.

Add 11 ml isopropanol to the Binding Buffer (**Code B**) and mix thoroughly.

Add 30 ml ethanol to the Pre-Wash Buffer (**Code P**) and mix thoroughly.

Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

Before each preparation

Preheating the Elution Buffer (**Code E**) - Transfer the needed amount of Elution Buffer (**Code E**) under calculation of a reserve volume into a reaction tube (not provided with the kit).

Equilibrate the Elution Buffer to 65°C. The Elution Buffer is necessary for step 7.

2.3 Extraction protocol

1. Preparation of the basic material

Transfer 50 mg* homogenized sample material into a 2.0 ml reaction tube (not provided with the kit).

For samples with high water content it is necessary to take up to 100 mg* of the basic material.

* The sample amount depends on the specific sample matrix. For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de. Recommendation for processing can be found at: <http://www.congen.de/en/downloads>

2. Lysis of the basic material

Add 400 µl of Lysis Buffer (**Code L**) and 20 µl of Proteinase K (**Code K**) to the reaction tube and mix it briefly.

Incubate on a heating block under continuously shaking for 30 minutes at 65°C.

Note: For strongly swelling samples it can be necessary to add additional volume of Lysis Buffer during the lysis.

3. Pre-filtration and setting of optimal binding conditions

Centrifuge the sample lysate for 1 min at 12,000 rpm.

Place a Spin Filter (**Code F, for pre-filtration**) into a 2.0 ml Receiver Tube (**Code R**). Transfer the liquid supernatant directly onto the Spin Filter.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm. After centrifugation discard the Spin Filter.

Note: For non-processed animal material, the usage of the Spin Filter (**Code F**) is not necessary. Transfer the liquid supernatant directly into a 2.0 ml Receiver Tube (**Code R**).

4. Binding of the nucleic acids on a Spin Filter

Add 200 µl Binding Buffer (**Code B**) to the filtrate and mix.

Place a Spin Filter (**Code S, for binding**) into a new, 2.0 ml Receiver Tube (**Code R**).

Transfer the complete solution onto the Spin Filter. Incubate at room temperature for 1 min.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm. After centrifugation discard the filtrate and place the Spin Filter back into the Receiver Tube.

5. Purification of the bound nucleic acids

Add 550 µl Pre-Wash Buffer (**Code P**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Note: For non-processed animal material the usage of the Pre-Wash Buffer is superfluous.

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Once more add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the Receiver Tube.

6. Drying of the Spin Filter

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

7. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a 1.5 ml Receiver Tube (**Code T**).

Add 100 µl of the preheated (65°C) Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate on a heating block for 3 min at 65°C (no shaking).

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

Note: The DNA can also be eluted with a lower or a higher volume of Elution Buffer (depending on the expected yield of DNA). The minimum volume for the elution is 50 µl.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at +2 to +8°C.

For a storage time of more than 24 hours it should be kept at -20°C.

3 Further Information

3.1 Product Information

- Flow chart (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com

r-biopharm®

4 Safety Information

Lysis Buffer

Warning

H319-H412-
P264-P273-P280-P305+P351+P338-
P337+P313-P501

Pre-Wash Buffer

Warning

H302-H315-H319-
P264-P270-P280-P301+P312-
P302+P352-P305+P351+P338-
P321-P330-P332+P313- P337+P313-
P362+P364-P501

Proteinase K

Danger

H315-H319-H334-H335-
P261-P264-P271-P280-P284-
P302+P352-P304+P340-
P305+P351+P338-P312-P501

H302:	Harmful if swallowed.
H315:	Causes skin irritation
H319:	Causes serious eye irritation.
H334:	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335:	May cause respiratory irritation.
H412:	Harmful to aquatic life with long lasting effects.
P261:	Avoid breathing dust/fume/gas/mist/vapours/spray.
P264:	Wash hands, forearms and face thoroughly after handling.
P270:	Do not eat, drink or smoke when using this product.
P271:	Use only outdoors or in a well-ventilated area.
P273:	Avoid release to the environment.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P284:	Wear respiratory protection.
P301+P312 :	IF SWALLOWED: Call a POISON CENTRE or doctor if you feel unwell.

P302+P352:	IF ON SKIN: Wash with plenty of water.
P304+P340:	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P312:	Call a POISON CENTRE or doctor if you feel unwell.
P321:	Specific treatment (see supplemental first aid instruction on this label).
P330 :	Rinse mouth.
P332+P313 :	If skin irritation occurs: Get medical advice/attention.
P337+P313:	If eye irritation persists: Get medical advice/attention.
P362+P364:	Take off contaminated clothing and wash it before reuse.
P501:	Dispose of contents/container to hazardous or special waste collection point, in

For further information we offer a Material Safety Data Sheet, see at www.congen.de/en/eifu/.
Alternatively please contact your distributor or send an e-mail to sales@r-biopharm.de.