CONGEN

**SureFood® GMO SCREEN 4plex
BAR/PAT/CryIAb/IAC/
CTP2:CP4 EPSPS**

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	DNA-preparation	10
1.4	Kit components and storage	10
1.5	Additionally required equipment and materials	10
1.6	Setup	11
1.7	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	12
2.2	Interpretation of results	13
3	Further Information	14
3.1	Product Information	14
3.2	Technical Support	14
3.3	Distribution and Ordering	14

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1 General Information

1.1 Description

The SureFood® GMO SCREEN 4plex BAR/PAT/CryIAb/IAc/CTP2:CP4 EPSPS is a real-time PCR for the direct, qualitative detection and differentiation of following specific genetically modified DNA - sequences.

- Phosphinothricin-Acetyltransferase gene (BAR) from *Streptomyces hygroscopicus*
- genetically engineered CryIAb-DNA-sequences and CryIAb/Ac-fusion gene sequences
- Phosphinothricin-Acetyltransferase gene (PAT) from *Streptomyces viridochromogenes*
- the transition from CTP2 (Chloroplast-Transpeptide-signal sequence from *Arabidopsis thaliana*) to herbicide tolerance-gene CP4 EPSPS (5-Enolpyruvylshikimat-3-Phosphat Syntheses gene from *Agrobacterium tumefaciens* strain CP4)

The detections are according to the validated methods of the European Commission or the official collection of detection methods of §64 German food law.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The technical validation of instruments was performed on Agilent Mx3005P, Bio-Rad CFX96, Roche LightCycler® 480 II, Applied Biosystems 7500, Qiagen Rotor-Gene Q and R-Biopharm RIDA®CYCLER.

1.2 Limit of Detection

The SureFood® GMO SCREEN 4plex BAR/PAT/CryIAb/IAc/CTP2:CP4 EPSPS real-time PCR has a limit of detection of ≤ 5 DNA copies. This is equivalent to a limit of detection of ≤ 0.01 % based on the genetically modified organisms MON87751 Soya (CryIAb/IAc), DAS-81419-2 Soya (PAT), NK603 Corn (CTP2:CP4 EPSPS) and Bt176 Corn (BAR).

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

Note: Inconsistent mixing ratios* may cause a loss of sensitivity in the low concentration channel in mixed samples especially with high amplicon concentrations (Cp value < 20).

* e.g. 99.9 MON87751 Soya (CryIAb/IAc) and 0.1 % NK603 Corn (CTP2:CP4 EPSPS)

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended. SureFood® PREP Add On is intended to be used for the extraction of DNA from raw materials as well as processed food and feed with sample weight of 2 g. It is used in conjunction with the SureFood® PREP Basic.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053 / SureFood® PREP Add On Art. Nr. S1055)
- real- time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	CryIAb/IAC	FAM	+	
	PAT	HEX	+	
	CTP2:CP4 EPSPS	ROX	+	
	BAR	Cy5	+	
Applied Biosystems 7500	CryIAb/IAC	FAM	None	Check the passive reference option ROX is none.
	PAT	VIC	None	
	CTP2:CP4 EPSPS	ROX	None	
	BAR	Cy5	None	
Bio-Rad CFX96	CryIAb/IAC	FAM	+	
	PAT	VIC/HEX	+	
	CTP2:CP4 EPSPS	ROX	+	
	BAR	Cy5	+	
R-Biopharm RIDA®CYCLER	CryIAb/IAC	green	+	
	PAT	yellow	+	
	CTP2:CP4 EPSPS	orange	+	
	BAR	red	+	
Qiagen Rotor-Gene Q	CryIAb/IAC	green	+	
	PAT	yellow	+	
	CTP2:CP4 EPSPS	orange	+	
	BAR	red	+	
Roche LightCycler® 480 II	CryIAb/IAC	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	PAT	533-580	+	
	CTP2:CP4 EPSPS	533-610	+	
	BAR	618-660	+	
Roche cobas® z 480 Analyzer	CryIAb/IAC	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	PAT	540-580	+	
	CTP2:CP4 EPSPS	540-610	+	
	BAR	610-670	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control and an inhibition control per sample.

For the preparation of the inhibition control, it is recommended to use the SureFood® GMO Plant PLUS (Art. No. S2049), the SureFood® GMO SCREEN 4plex 35S/NOS/FMV+IAC (Art. No. S2126) or the SureFood® GMO SCREEN 35S/NOS/FMV (Art. No. S2026).

Reactions needed for the qualitative CryIAb/IAc, PAT, CTP2:CP4 EPSPS and BAR detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

CryIAb/IAc DNA is detected in the FAM-channel, PAT DNA is detected in the VIC/HEX-channel, CTP2:CP4 EPSPS DNA is detected in the ROX-channel and BAR DNA is detected in the Cy5-channel (see table).

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the external inhibition control of the sample is **positive** with a shift in Cp-value ≤ 2 compared to the negative control.

If the sample DNA in the external inhibition control shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific CryIAb/IAc, PAT, CTP2:CP4 EPSPS or BAR PCR assay.

Result in the respective channel				Result	Interpretation
FAM channel CryIAb/IAc	VIC/HEX channel PAT	ROX channel CTP2:CP4 EPSPS	Cy5 channel BAR	external inhibition control	
positive	negative	negative	negative	positive	CryIAb/IAc DNA detected
negativ	positive	negative	negative	positive	PAT DNA detected
negativ	negative	positive	negative	positive	CTP2:CP4 EPSPS DNA detected
negativ	negative	negative	positive	positive	BAR DNA detected
negativ	negative	negative	negative	negative	invalid

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices
(Download: www.congen.de/en/company/downloads)
- Validation Report upon request

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com