CONGEN

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Reagenz	Menge
GMO Plant Reference Sample	1 x 2,0 g

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1 General Information

1.1 Description

The SureFood® GMO Plant Reference Sample is validated for semi quantitative determination of the plants soya, corn, canola and rice in food and feed samples.

The laboratory reference material SureFood® GMO Plant Reference Sample is based on **wheat flour** material. It contains the plants **soya, corn, canola** and **rice** in a concentration of **0.1%**.

Content of plant substances:

Plant	Material
wheat	wheat flour (<i>Triticum</i> spp.)
soya	soya flour (<i>Glycine max</i>)
corn	corn flour (<i>Zea mays</i>)
canola	canola flour (<i>Brassica napus</i>)
rice	rice flour (<i>Oryza sativa</i>)

The laboratory reference material SureFood® GMO Plant Reference Sample is validated in combination with the DNA extraction method SureFood® PREP Basic (Art. No. S1052) and the SureFood® real-time PCR detection methods.

In case of using alternative methods the sample weight and extraction method in combination with the detection system have to be validated.

1.2 DNA-preparation

SureFood® PREP Basic (Art. No. S1052) or additionally SureFood® PREP Add-On Kit (Art. No. S1055) is recommended for the DNA-preparation. As quantification is performed using an external reference, both the SureFood® GMO Plant Reference Sample and the sample must undergo identical extraction procedures.

Lysis, binding and DNA purification steps must be conducted strictly in accordance with the established protocol.

The SureFood® GMO Plant Reference Sample should be extracted in duplicate in parallel with the samples to be analysed.

1.3 Kit components and storage

Reagent	Amount
GMO Plant Reference Sample	1 x 2.0 g

Store all reagents dry at -28 to -16°C and protected from light.

Note: To avoid moisture absorption the unopened laboratory reference material should thaw at room temperature for at least one hour before use.

1.4 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052; optional SureFood® PREP Add-On Art. No. S1055)
- Real-time PCR methods (SureFood® GMO series)

2 Further Information

2.1 Product Information

- Product-related documents (Download: www.congen.de/en/eifu/)

2.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

2.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com