CONGEN

SureFood® GMO Plant 4plex Corn/Soya/Canola/Cotton

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	DNA-preparation	10
1.4	Kit components and storage	10
1.5	Additionally required equipment and materials	10
1.6	Setup	10
1.7	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	12
2.2	Interpretation of results	13
3	Limitations of the method	14
4	Further Information	14
4.1	Product Information	14
4.2	Technical Support	14
4.3	Distribution and Ordering	14

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1 General Information

1.1 Description

The SureFood® GMO Plant 4plex Corn/Soya/Canola/Cotton is a real-time PCR for the direct, qualitative detection and differentiation of specific DNA sequences of corn (*Zea mays*), soya (*Glycine max*), canola (*Brassica napus*) und cotton (*Gossypium spp.*).

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The technical validation of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Bio-Rad CFX96, Bio-Rad CFX96 Dx, R-Biopharm RIDA®CYCLER, and Roche LightCycler® 480 II.

1.2 Limit of Detection

The SureFood® GMO Plant 4plex Corn/Soya/Canola/Cotton real-time PCR has a limit of detection of ≤ 500 DNA copies. This is equivalent to approx. 0.01% for unprocessed grain.

The limit of detection and quantification of the complete method (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic (Art. No. S1052) and for highly processed food and feed the use of SureFood® PREP Advanced (Art. No. S1053) is recommended. SureFood® PREP Add On (Art. No. S1055) is intended to be used for the extraction of DNA from raw materials as well as processed food and feed with sample weight of 2 g. It is used in conjunction with the SureFood® PREP Basic.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053 / SureFood® PREP Add On Art. No. S1055)
- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD) Cycles	5 min, 95°C 35	1 min, 95°C 35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	Corn	FAM	+	
	Cotton	HEX	+	
	Canola	ROX	+	
	Soya	Cy5	+	
Agilent AriaDx /Mx	Corn	FAM	+	
	Cotton	HEX	+	
	Canola	ROX	+	
	Soya	Cy5	+	
Bio-Rad CFX96/ Dx	Corn	FAM	+	
	Cotton	VIC/HEX	+	
	Canola	ROX	+	
	Soya	Cy5	+	
R-Biopharm RIDA®CYCLER	Corn	green	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter.
	Cotton	yellow	+	
	Canola	orange	+	
	Soya	red	+	
Roche LightCycler® 480 II	Corn	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	Cotton	533-580	+	
	Canola	533-610	+	
	Soya	618-660	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control, Positive Control and inhibition control.

For the preparation of the inhibition control the use of SureFood® GMO Plant PLUS (Art. No. S2049) or SureFood® GMO SCREEN 4plex 35S/NOS/FMV+IAC (Art. No. S2126) is recommended.

Reactions needed for the qualitative corn, soya, canola and cotton detection:

3 reactions for controls (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10% additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Corn DNA is detected in the FAM-channel, cotton DNA is detected in the VIC/HEX-channel, canola DNA is detected in the ROX-channel and soya DNA is detected in the Cy5-channel (see table).

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the external inhibition control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-value ≤ 2 compared to the negative control.

If the sample DNA in the external inhibition control shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific corn, cotton, canola or soya PCR assay.

Result in the respective channel				Result external inhibition-control	Interpretation
FAM channel corn	VIC/HEX channel cotton	ROX channel canola	Cy5 channel soya		
positive	negative	negative	negative	positive	corn DNA detected
negativ	positive	negative	negative	positive	cotton DNA detected
negativ	negative	positive	negative	positive	canola DNA detected
negativ	negative	negative	positive	positive	soya DNA detected
negativ	negative	negative	negative	positive	Negative, corn, cotton, canola, and soya DNA is not detected
negativ	negative	negative	negative	negative	invalid

Note: The results displayed in the table above represent merely an example. Additional combinations are also possible.

3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- A positive test result does not necessarily indicate the presence of a viable organism. It is indicative for the presence of the target DNA.

4 Further Information

4.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

4.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

4.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com

r-biopharm®