**CONGEN**

# **SureFood® GMO ID 4plex Soya III**

Art. No. S2164  
100 rxn

## **User Manual**

**May 2024**

 **Inhalt**

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## **Content**

1	General Information .....	9
1.1	Description .....	9
1.2	Limit of Detection .....	9
1.3	DNA-preparation .....	10
1.4	Kit components and storage .....	10
1.5	Additionally required equipment and materials .....	10
1.6	Setup.....	11
1.7	Detection channel Set-up .....	11
2	Qualitative Analysis .....	12
2.1	Protocol .....	12
2.1.1	Preparation of the master-mix .....	12
2.1.2	Preparation of the real-time PCR-mix .....	12
2.2	Interpretation of results .....	13
3	Limitations of the method .....	14
4	Further Information .....	14
4.1	Product Information .....	14
4.2	Technical Support .....	14
4.3	Distribution and Ordering .....	14

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## **1 General Information**

### **1.1 Description**

The SureFood® GMO ID 4plex Soya III is a real-time PCR for the direct, qualitative detection and differentiation of following specific genetically modified soya DNA sequences.

- FG72 soya (OECD unique identifier MST-FG072-2)
- DAS68416 soya (OECD unique identifier DAS-68416-4)
- GMB151 soya (OECD unique identifier BCS-GM151-6)
- DAS44406 soya (OECD unique identifier DAS-44406-6)

This kit can be used for the screening of genetically modified organisms (GMOs) in food, feed and seeds.

The detections are according to the validated method of the European Commission.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The technical verification of instruments was performed on Roche LightCycler® 480 II, Bio-Rad CFX96 Dx, R-Biopharm RIDA®CYCLER and Agilent AriaMx.

### **1.2 Limit of Detection**

The SureFood® GMO ID 4plex Soya III real-time PCR has a limit of detection of  $\leq 5$  DNA copies. This is equivalent to approx. 0.01 % for unprocessed soybean.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

**Note:** Inconsistent mixing ratios\* may cause a loss of sensitivity in the low concentration channel in mixed samples especially with high amplicon concentrations (Cp value < 20).

\* e.g. 99.9 GMB151 soya and 0.1 % FG72 soya

**1.3 DNA-preparation**

For DNA-preparation of raw material the use of SureFood® PREP Basic (Art. No. S1052), and for highly processed food and feed the use of SureFood® PREP Advanced (Art. No. S1053) is recommended. SureFood® PREP Add On (Art. No. S1055) is intended to be used for the extraction of DNA from raw materials as well as processed food and feed with sample weight of 2 g. It is used in conjunction with the SureFood® PREP Basic.

**1.4 Kit components and storage**

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at –28 bis -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

**1.5 Additionally required equipment and materials**

- DNA-Extraction kit  
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053 / SureFood® PREP Add On Art. No. S1055)
- real- time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

**1.6 Setup**

	<b>Blockcycler &amp; R-Biopharm RIDA®CYCLER</b>	<b>Rotorcycler</b>
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.7 Detection channel Set-up**

<b>Real-time PCR device</b>	<b>Detection</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent AriaDx /Mx</b>	FG72 soya	FAM	+	
	DAS68416 soya	HEX	+	
	GMB151 soya	ROX	+	
	DAS44406 soya	Cy5	+	
<b>Bio-Rad CFX96/Dx/Opus</b>	FG72 soya	FAM	+	
	DAS68416 soya	VIC/HEX	+	
	GMB151 soya	ROX	+	
	DAS44406 soya	Cy5	+	
<b>R-Biopharm RIDA®CYCLER</b>	FG72 soya	green	+	
	DAS68416 soya	yellow	+	
	GMB151 soya	orange	+	
	DAS44406 soya	red	+	
<b>Roche LightCycler® 480 II</b>	FG72 soya	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	DAS68416 soya	533-580	+	
	GMB151 soya	533-610	+	
	DAS44406 soya	618-660	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay. Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control and an inhibition control per sample.

For the preparation of the inhibition control, it is recommended to use the SureFood® GMO Plant PLUS (Art. No. S2049) or the SureFood® GMO SCREEN 4plex 35S/NOS/FMV+IAC (Art. No. S2126).

#### **Reactions needed for the qualitative FG72 soya, DAS68416 soya, GMB151 soya or DAS44406 soya detection:**

3 reactions for controls (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

#### **Example for the calculation and preparation of 10 reactions:**

<b>Components of the master-mix</b>	<b>Amount per reaction</b>	<b>10 reactions (with 10% excess)</b>
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

#### 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

**2.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

FG72 soya DNA is detected in the FAM-channel, DAS68416 soya DNA is detected in the VIC/HEX -channel, GMB151 soya DNA is detected in the ROX-channel and DAS44406 soya DNA is detected in the Cy5-channel (see table).

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the external inhibition control of the sample is **positive** with a shift in Cp-value ≤ 2 compared to the negative control.

If the sample DNA in the external inhibition control shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific FG72 soya, DAS68416 soya, GMB151 soya or DAS44406 soya PCR assay.

Result in the respective channel				Result	Interpretation
FAM channel FG72 soya	VIC/HEX channel DAS68416 soya	ROX channel GMB151 soya	Cy5 channel DAS44406 soya	external inhibition- control	
<b>positive</b>	negative	negative	negative	<b>positive</b>	FG72 soya DNA detected
negativ	<b>positive</b>	negative	negative	<b>positive</b>	DAS68416 soya DNA detected
negativ	negative	<b>positive</b>	negative	<b>positive</b>	GMB151 soya DNA detected
negativ	negative	negative	<b>positive</b>	<b>positive</b>	DAS44406 soya DNA detected
negativ	negative	negative	negative	<b>positive</b>	Negative for FG72, DAS68416, GMB151 and DAS44406 DNA
negativ	negative	negative	negative	negative	invalid

### **3 Limitations of the method**

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.

### **4 Further Information**

#### **4.1 Product Information**

- Detailed information about setup of several real-time PCR devices  
(Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report upon request

#### **4.2 Technical Support**

For further questions please contact your distributor or send an e-mail to [sales@r-biopharm.de](mailto:sales@r-biopharm.de).

#### **4.3 Distribution and Ordering**

R-Biopharm AG  
An der neuen Bergstrasse 17,  
64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 - 81 02-0  
Fax: +49 (0) 61 51 - 81 02-20  
E-Mail: [orders@r-biopharm.de](mailto:orders@r-biopharm.de)  
[www.r-biopharm.com](http://www.r-biopharm.com)