CONGEN

SureFood® QUANTARD Allergen 40

Art. No. S3301
100 rxn

User Manual

September 2026

 Inhalt

1 Allgemeines3

1.1 Beschreibung.....3

1.2 DNA-Präparation4

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1 General Information

1.1 Description

SureFood® QUANTARD Allergen 40 is a laboratory reference material for quantitative determination of 1 mg to 400 mg allergenic substance / kg food sample.

The laboratory reference material SureFood® QUANTARD Allergen 40 is based on **corn flour** material. It contains all potentially allergenic food ingredients according to regulation (EU) No. 1169/2011 (except for **lactose and sulphur dioxide**) in a concentration of **40 mg/kg**.

Content of allergenic substances:

Allergen	Material
Buckwheat	Buckwheat (<i>Fagopyrum esculentum</i>)
Celery	Celery seeds (<i>Apium graveolens</i>)
Crustacean	Black Tiger shrimps (<i>Penaeus monodon</i>), peeled, glazed, blistered, raw, deveined
Edible Nut	Hazelnut kernels (<i>Corylus avellana</i>), without skin Almond (<i>Prunus dulcis</i>), whole, untreated, with skin Walnut kernels (<i>Juglans regia</i>), untreated, with skin Cashew kernels (<i>Anacardium occidentale</i>), untreated, without skin Pecan nuts (<i>Carya illinoensis</i>), untreated, with skin Brazil nut (<i>Bertholletia excelsa</i>), peeled, without skin Pistachio kernels (<i>Pistacia vera</i>), peeled, without skin Macadamia nut (<i>Macadamia ternifolia</i>), roasted in the shell, peeled
Egg	Egg white powder, hen (<i>Gallus gallus</i>)
Fish	Humpback salmon filet (<i>Oncorhynchus gorbuscha</i>), untreated
Gluten containing cereals	Wheat flour type 405 (<i>Triticum</i> spp.)
Lupin	Lupin (<i>Lupinus albus</i>), seed of white lupin, without skin
Milk	Milkpowder, cow (<i>Bos taurus</i>)
Mustard	Mustard seeds, yellow (<i>Sinapis alba</i>)
Molluscs	Sea scallop (<i>Placopecten magellanicus</i>), raw, glazed
Peanut	Peanut kernels (<i>Arachis hypogaea</i>), untreated, without skin
Sesame	Sesam seeds (<i>Sesamum indicum</i>), peeled
Soy	Soybean (<i>Glycine max</i>)

The laboratory reference material SureFood® QUANTARD Allergen 40 is validated in combination with the DNA extraction method SureFood® PREP Advanced (Art. No. S1053), protocol 1, and the SureFood® ALLERGEN real-time PCR detection methods.

Further information about the quantification of allergenic substances is given in the manual of the SureFood® ALLERGEN real-time PCR products.

In case of using alternative methods the sample weight and extraction method in combination with the detection system have to be validated.

1.2 DNA-preparation

SureFood® PREP Advanced (Art. No. S1053) protocol 1 is recommended for the DNA preparation.

As quantification is performed using an external reference, both the SureFood® QUANTARD Allergen 40 standard and the sample must undergo identical extraction procedures. The SureFood® QUANTARD Allergen 40 standard should be extracted in duplicate alongside the test samples to ensure accuracy and reproducibility.

Lysis, binding and DNA purification steps must be conducted strictly in accordance with the established protocol.

1.3 Kit components and storage

Reagent	Amount
QUANTARD Allergen 40	1 x 2.0 g

Store all reagents dry at –28°C to – 16°C and protected from light.

Note: To avoid moisture absorption the unopened laboratory reference material should thaw at room temperature for at least one hour before use.

1.4 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Advanced, Art. No. S1053)
- Real-time PCR methods (SureFood® ALLERGEN series)

2 Further Information

2.1 Protein content (Kjeldahl method)

The protein content of the allergenic food components (100 %) was determined using the Kjeldahl method in accordance with official collection of detection methods of §64 German food law (LFGB). Based on this, the protein content of the allergens in the SureFood® QUANTARD ALLERGEN 40 was calculated.

Allergen	Protein content of allergenic food ingredients in g/100 g	SureFood® QUANTARD ALLERGEN 40 protein content in ppm	Conversion factor from mg allergen/kg food to protein content in ppm
Almond	23.5	9.4	0.24
Brazil Nut	16.7	6.7	0.17
Cashew	22.0	8.8	0.22
Celery	10.1	4.0	0.10
Crustacean	15.3	6.1	0.15
Fish	21.7	8.7	0.22
Hazelnut	16.7	6.7	0.17
Lupin	35.6	14.2	0.36
Macadamia	11.5	4.6	0.12
Molluscs	11.0	4.4	0.11
Mustard	27.2	10.9	0.27
Peanut	24.4	9.8	0.24
Pecan	12.0	4.8	0.12
Pistachio	25.5	10.2	0.26
Sesame	23.7	9.5	0.24
Soy	36.5	14.6	0.37
Walnut	17.9	7.2	0.18
Wheat	14.1	5.6	0.14

2.2 Product Information

- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

2.3 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

2.4 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com