CONGEN

SureFood® ALLERGEN Crustaceans

Art. No. S3612
100 rxn

User Manual

May 2026

 **Inhalt**

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Content

1	General Information	13
1.1	Description	13
1.2	Limit of Detection and Limit of Quantification	13
1.3	DNA preparation	14
1.4	Kit components and storage	14
1.5	Additionally required equipment and materials	14
1.6	Precautions for users	14
1.7	Setup	15
1.8	Detection channel Set-up	15
2	Qualitative Analysis	16
2.1	Protocol	16
2.1.1	Preparation of the master-mix	16
2.1.2	Preparation of the real-time PCR-mix	16
2.2	Interpretation of results	17
3	Quantitative Analysis	18
3.1	Protocol	18
3.1.1	Preparation of the master-mix	18
3.1.2	Preparation of the standard DNA dilutions	19
3.1.3	Preparation of the real-time PCR-mix	19
3.2	Interpretation of results	20
4	Limitations of the method	21
5	Further Information	22
5.1	Product Information	22
5.2	Technical Support	22
5.3	Distribution and Ordering	22

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1 General Information

1.1 Description

The SureFood® ALLERGEN Crustaceans is a real-time PCR for the direct, qualitative and / or quantitative detection of specific crustaceans (*Crustacea*) DNA sequences according to directive (EC) 1169/2011 in food and feed.

The crustacean family contains about 44,000 species. Due to the large number of different species a practical verification of all is impossible. More detailed information is available in the Validation Report of SureFood® ALLERGEN Crustaceans.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected, or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition, spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The internal technical verification of instruments was performed on Agilent AriaDx, Bio-Rad CFX96 Dx, Bio-Rad CFX OPUS 96, Qiagen Rotor-Gene Q, R-Biopharm RIDA®CYCLER and Roche LightCycler® 480 II.

For the quantitative determination the use of the laboratory reference material SureFood® QUANTARD Allergen 40 (Art. No. S3301) containing 40 mg crustaceans / kg food sample and a detection system including a standard series is required.

In general this procedure is applicable for all homogenous food samples. The quantitative determination is inapplicable for swabs, sponges and inhomogeneous fluids/liquids.

The method is verified for quantitative determination of 1 mg to 400 mg allergenic substance / kg food sample using the reference material SureFood® QUANTARD Allergen 40 containing 40 mg allergenic substance / kg food sample.

1.2 Limit of Detection and Limit of Quantification

The SureFood® ALLERGEN Crustaceans PCR has a limit of detection of ≤ 0.4 mg / kg using SureFood® PREP Advanced, protocol 1.

The limit of quantification is 1 mg / kg using SureFood® PREP Advanced, protocol 1.

The limit of detection and quantification were determined using the SureFood® QUANTARD Allergen 40 (matrix: corn flour). The limit of detection and quantification of the complete method (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 DNA preparation

For DNA preparation the use of SureFood® PREP Advanced (Art. No. S1053) is recommended.

If a quantitative determination is performed, the SureFood® QUANTARD Allergen 40 and the sample must be extracted in the same manner. The SureFood® QUANTARD Allergen 40 should be extracted in duplicate alongside the samples to be analyzed.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
4	Dilution Buffer	1 x 1600 µl	White
5	Standard DNA	1 x 200 µl	Dark Blue

Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA extraction kit
(e.g. SureFood® PREP Advanced Art. No. S1053)
- reference material for quantification (SureFood® QUANTARD Allergen 40, Art. No. S3301)
- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

For additional details see cycler operating instructions for real-time PCR device.

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent AriaMx / Dx	Crustaceans	FAM	+	
	IAC	HEX	+	
Bio-Rad CFX96/ Dx/ Opus	Crustaceans	FAM	+	Baseline Settings: <ul style="list-style-type: none"> • Baseline subtracted curve fit • Apply fluorescence drift correction
	IAC	VIC/HEX	+	
Qiagen Rotor-Gene Q	Crustaceans	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
R-Biopharm RIDA®CYCLER	Crustaceans	green	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter.
	IAC	yellow	+	
Roche LightCycler® 480 II	Crustaceans	465-510	+	
	IAC	533-580	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control and Positive Control.

The reaction mix contains an internal amplification control (inhibition control) per reaction.

Reactions needed for the qualitative detection of crustaceans:

3 reactions for controls* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

*** Description of the controls**

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive Control: master-mix and within the kit's provided Positive Control

2.2 Interpretation of results

The evaluation must be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions must show the correct results.

DNA of crustaceans is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive**, if the sample DNA shows amplification in the detection system with a Cp value of ≤ 35 . High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the obtained Cp value is > 35 . The internal amplification control (inhibition control) of the sample must be **positive** (VIC/HEX) with a shift in Cp-Value ≤ 2 compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. It is recommended to repeat the extraction process. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific PCR assay.

Alternatively, a self-selected reference sample at the detection limit can be used to assess the sample. This reference sample is used outside CONGEN's quality and evaluation guidelines.

The following table shows the specification ranges of the kit controls

	Specification range	Internal amplification control (IAC)
Negative control	negative	$25 \leq Cp \leq 34$
Positive Control (FAM – crustaceans)	$25 \leq Cp \leq 33$	not relevant

3 Quantitative Analysis

3.1 Protocol

3.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control, standard curve as positive control and for quantification.

The reaction mix contains an internal amplification control (inhibition control) per reaction.

Reactions needed for the quantitative detection of crustaceans*:

5 reactions for the standard curve

2 reactions for controls (1x negative control, 1x extraction control)

2 reactions for reference material (SureFood® QUANTARD Allergen 40 DNA)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

***Description of the reactions**

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive control: master-mix and within the kit's provided standards

3.1.2 Preparation of the standard DNA dilutions

Use the undiluted Standard DNA (**Code 5**) and dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different concentrations. Prepare a dilution series of 4 steps. Prepare 4 reaction tubes (labelled S2 to S5) and add 45 µl Dilution Buffer each.

The following procedure is recommended:

Standard	Dilution	Copy number per µl	Final copy number per reaction*
S1	Standard DNA undiluted	100,000 copies	500,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	50,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1,000 copies	5,000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	500 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S4	10 copies	50 copies

***Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

3.1.3 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of the SureFood® QUANTARD 40 Allergen DNA and the standard dilutions into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

3.2 Interpretation of results

The evaluation must be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions must show the correct results.

DNA of crustaceans is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive** and can be **quantified**, if the sample DNA shows amplification in the detection system with a Cp value of ≤ 35 . High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the internal amplification control is not needed to obtain a positive result of the Positive control / Standard DNA.

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the obtained Cp value is > 35 . The internal amplification control (inhibition control) of the sample has to be **positive** (VIC/HEX) with a shift in Cp-Value ≤ 2 compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. It is recommended to repeat the extraction process. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific PCR assay.

For quantification mark the standards, controls and samples and make the evaluation according to the analysis program recommended by the real-time PCR instrument manufacturer. The value for the slope of the standard curve must be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

During evaluation, it must be ensured that the sample to be analysed lies within the linear range of the standard curve. If the sample shows amplification before the highest standard, it should be diluted accordingly and measured again.

For the evaluation and calculation of the result, the SureFood® ALLERGEN Quant – Data interpretation sheet (see 5.1 Product information) can be used. Please follow the introductions provided in the template.

As an alternative to the SureFood® QUANTARD Allergen 40, a self-selected reference sample with known allergen content may be used for sample evaluation. The use of such reference sample is outside the quality and assessment guidelines of CONGEN.

The following table shows the specification ranges of the kit controls

	Specification range	Internal amplification control (IAC)
Negative control	negative	$25 \leq Cp \leq 34$
Standard DNA, undiluted (FAM – crustaceans)	$16 \leq Cp \leq 24$	not relevant

For a final calculation of mg allergenic substance / kg food sample following formula is introduced:

$$X[mg/kg] = \frac{10^{\left(\frac{Cp_{sample}}{s}\right)}}{10^{\left(\frac{Cp_{Quantard}}{s}\right)}} * C_{Quantard}$$

s = slope | c = concentration | Cp_{Quantard} = Ø of double assay

Example:

$$X[mg/kg] = \frac{10^{\left(\frac{28.5}{-3.251}\right)}}{10^{\left(\frac{29.2}{-3.251}\right)}} * 40 = 65.7 mg/kg$$

Cp_{sample} = 28.5 | Cp_{Quantard} = Ø 29.2 | s = - 3.251 | c = 40 mg/kg

For this example an amount of **65.7 mg crustaceans / kg food sample** is calculated.

4 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like proteins or lipids for example.
- In highly processed samples, the detection limit may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process– such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.
- The test kit is verified for the quantification of allergen concentrations in the range of 1 to 400 ppm using the SureFood® QUANTARD Allergen 40. Samples can be diluted up to 1:10 without compromising the reliability of the results. Quantification of samples containing more than 4000 ppm (0.4 %) allergen is not recommended.

5 Further Information

5.1 Product Information

- Microsoft Excel template of calculation and detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

5.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

5.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com

r-biopharm®