CONGEN

SureFood® ALLERGEN Pistachio

Art. No. S3614
100 rxn

User Manual

July 2025

 Inhalt

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Content

1	General Information	14
1.1	Description	14
1.2	Limit of Detection and Limit of Quantification	14
1.3	DNA preparation	15
1.4	Kit components and storage	15
1.5	Additionally required equipment and materials	15
1.6	Precautions for users	15
1.7	Setup.....	16
1.8	Detection channel Set-up	16
2	Qualitative Analysis	18
2.1	Protocol	18
2.1.1	Preparation of the master-mix	18
2.1.2	Preparation of the real-time PCR-mix	18
2.2	Interpretation of results	19
3	Quantitative Analysis	20
3.1	Protocol	20
3.1.1	Preparation of the master-mix	20
3.1.2	Preparation of the standard DNA dilutions	20
3.1.3	Preparation of the real-time PCR-mix	21
3.2	Interpretation of results	22
4	Limitations of the method	23
5	Further Information	24
5.1	Product Information	24
5.2	Technical Support	24
5.3	Distribution and Ordering	24

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1 General Information

1.1 Description

The SureFood® ALLERGEN Pistachio is a real-time PCR for the direct, qualitative and / or quantitative detection of specific pistachio (*Pistacia vera*) DNA sequences according to directive (EC) 1169/2011 in food and feed.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected, or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition, spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

For the quantitative determination the use of the laboratory reference material SureFood® QUANTARD Allergen 40 (Art. No. S3301) containing 40 mg pistachio / kg food sample and a detection system including a standard series is required.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The technical verification of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Applied Biosystems 7500, Bio-Rad CFX96, Bio Molecular Systems MIC, Qiagen Rotor-Gene Q, Roche LightCycler® 480 II and Roche cobas® z 480 Analyzer.

In general this procedure is applicable for all homogenous food samples. The quantitative determination is inapplicable for swabs, sponges and inhomogeneous fluids/liquids.

The method is verified for quantitative determination of 1 mg to 400 mg allergenic substance / kg food sample using the reference material SureFood® QUANTARD Allergen 40 containing 40 mg allergenic substance / kg food sample.

1.2 Limit of Detection and Limit of Quantification

The SureFood® ALLERGEN Pistachio PCR has a limit of detection of ≤ 0.4 mg / kg using SureFood® PREP Advanced, protocol 1.

The limit of quantification is 1 mg / kg using SureFood® PREP Advanced, protocol 1.

The limit of detection and quantification were determined using the SureFood® QUANTARD Allergen 40 (matrix: corn flour). The limit of detection and quantification of the complete method (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 DNA preparation

For DNA preparation the use of SureFood® PREP Advanced (Art. No. S1053), protocol 1 is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue
4	Dilution Buffer	1 x 1600 µl	White
5	Standard DNA	1 x 200 µl	Dark Blue

Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA extraction kit (e.g. SureFood® PREP Advanced Art. No. S1053)
- reference material for quantification (SureFood® QUANTARD Allergen 40, Art. No. S3301)
- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler & LightCycler® 480 II
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	Pistachio	FAM	+	
	IAC	HEX	+	
Agilent AriaMx /Dx	Pistachio	FAM	+	
	IAC	HEX	+	
Applied Biosystems 7500	Pistachio	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
Bio-Rad CFX96/ Dx/ Opus	Pistachio	FAM	+	Baseline Settings: • Baseline subtracted curve fit • Apply fluorescence drift correction
	IAC	VIC/HEX	+	
Qiagen Rotor-Gene Q	Pistachio	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
R-Biopharm RIDA®CYCLER/ Bio Molecular Systems MIC	Pistachio	green	+	Ignore cycles before, if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter.
	IAC	yellow	+	

SureFood® ALLERGEN Pistachio
(100 rxn)

Art. No. S3614

July 2025

Real-time PCR device	Detection	Detection channel	Quencher	Note
Roche LightCycler® 480 II	Pistachio	465-510	+	
	IAC	533-580	+	
Roche cobas® z 480 Analyzer	Pistachio	465-510	+	
	IAC	540-580	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Required control reactions for the specific PCR assay: negative control, extraction control, Positive Control.

The reaction mix contains an internal amplification control (Inhibition control) per reaction.

Reactions needed for the qualitative detection of pistachio:

3 reactions for controls* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10% additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

*** Description of the controls**

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive Control: master-mix and within the kit's provided Positive Control

2.2 Interpretation of results

The evaluation must be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions must show the correct results.

Pistachio DNA is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive**, if the sample DNA shows amplification in the detection system with a Cp value of ≤ 35 . High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the obtained Cp value is > 35 . The internal amplification control (inhibition control) of the sample must be **positive** (VIC/HEX) with a shift in Cp-Value ≤ 2 compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific pistachio PCR assay.

Alternatively, a self-selected reference sample at the detection limit can be used to assess the sample. This reference sample is used outside CONGEN's quality and evaluation guidelines.

The following table shows the specification ranges of the kit controls

	Specification range	Internal amplification control (IAC)
negative control	negative	$25 \leq Cp \leq 34$
Positive Control (FAM – pistachio)	$25 \leq Cp \leq 33$	not relevant

3 Quantitative Analysis

3.1 Protocol

3.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay. Recommended control reactions for the specific PCR assay: negative control, extraction control, standard curve as positive control and for quantification.

The reaction mix contains an internal amplification control (Inhibition control) per reaction.

Reactions needed for the quantitative detection of pistachio*:

5 reactions for the standard curve

2 reactions for controls (1x negative control, 1x extraction control)

2 reactions for reference material (SureFood® QUANTARD Allergen 40 DNA)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10% additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

3.1.2 Preparation of the standard DNA dilutions

Use the undiluted Standard DNA (**Code 5**) and dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different concentrations. Prepare a dilution series of 4 steps. Prepare 4 reaction tubes (labelled S2 to S5) and add 45 µl Dilution Buffer each.

The following procedure is recommended:

Standard	Dilution	Copy number per µl	Final copy number per reaction*
S1	Standard DNA undiluted	100,000 copies	500,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	50,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1,000 copies	5,000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	500 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S4	10 copies	50 copies

***Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

3.1.3 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of the SureFood® QUANTARD 40 Allergen DNA and the standard dilutions into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

*** Description of the reactions**

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive control: master-mix and within the kit's provided standards

3.2 Interpretation of results

The evaluation must be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions must show the correct results.

DNA of pistachio is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification control (IAC) is detected.

A sample is stated **positive** and can be **quantified**, if the sample DNA shows amplification in the detection system with a Cp value of ≤ 35 . High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control / Standard DNA.

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the obtained Cp value is > 35 . The internal amplification control (inhibition control) of the sample has to be **positive** (VIC/HEX) with a shift in Cp-Value ≤ 2 compared to the negative control.

If the sample DNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific pistachio PCR assay.

For quantification mark the standards, controls and samples and make the evaluation according to the analysis program recommended by the real-time PCR instrument manufacturer. The value for the slope of the standard curve must be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

The following table shows the specification ranges of the kit controls

	Specification range	Internal amplification control (IAC)
negative control	negative	$25 \leq Cp \leq 34$
Standard DNA, undiluted (FAM – pistachio)	$16 \leq Cp \leq 24$	not relevant

For a final calculation of mg allergenic substance / kg food sample following formula is introduced:

$$X[mg/kg] = \frac{10^{\left(\frac{Cp_{sample}}{s}\right)}}{10^{\left(\frac{Cp_{Quantard}}{s}\right)}} * c_{Quantard}$$

s = slope | c = concentration | Cp_{Quantard} = Ø of double assay

Example:

$$X[mg/kg] = \frac{10^{\left(\frac{28.5}{-3.251}\right)}}{10^{\left(\frac{29.2}{-3.251}\right)}} * 40 = 65.7 mg/kg$$

Cp_{sample} = 28.5 | Cp_{Quantard} = Ø 29.2 | s = - 3.251 | c = 40 mg/kg

For this example an amount of **65.7 mg pistachio / kg food sample** is calculated.

4 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like proteins or lipids for example.
- In highly processed samples, the detection limit may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process– such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.

5 Further Information

5.1 Product Information

- Microsoft Excel template of calculation and detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

5.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

5.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com