CONGEN

SureFood® ANIMAL ID Turkey IAAC

Art. No. S6116
100 rxn

User Manual

March 2025

 **Inhalt**

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	DNA-preparation	10
1.4	Kit components and storage	10
1.5	Additionally, required equipment and materials	10
1.6	Precautions for users	10
1.7	Setup.....	11
1.8	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	12
2.2	Interpretation of results	13
3	Limitations of the method	14
4	Further Information	14
4.1	Product Information	14
4.2	Technical Support	14
4.3	Distribution and Ordering	14

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1 General Information

1.1 Description

The SureFood® ANIMAL ID Turkey IAAC is a real-time PCR for the direct, qualitative detection of a specific turkey (*Meleagris gallopavo*) DNA sequence in food and feed.

Each reaction contains an internal amplification control and an internal detection assay for vertebrate's DNA (IAAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition, spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The technical verification of instruments was performed on Agilent Mx3005P, Applied Biosystems 7500, Bio-Rad CFX96, Qiagen Rotor-Gene Q, and Roche LightCycler® 480 II.

1.2 Limit of Detection

The SureFood® ANIMAL ID Turkey IAAC real-time PCR is developed for the detection of turkey DNA in muscle meat mixture at a relative amount of 0.1%.

The limit of detection of the complete method (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic (Art. No. S1052) is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at -28 to -16°C and protected from light. The tube of the Taq Polymerase should not be thawed.

1.5 Additionally, required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Basic Art. No. S1052)
- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD) Cycles	5 min, 95°C 35	1 min, 95°C 35
Denaturation Annealing/Extension (CYCLE)	15 sec, 95°C 30 sec, 60°C	10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	Turkey	FAM	+	
	IAAC	HEX	+	
Applied Biosystems 7500	Turkey	FAM	None	Check the passive reference option ROX is none.
	IAAC	VIC	None	
Bio-Rad CFX96 / Dx / Opus	Turkey	FAM	+	
	IAAC	VIC/HEX	+	
Qiagen Rotor-Gene Q	Turkey	green	+	Note: Please use only 0.1 ml reaction tube. The gain settings must be set to 5 (factory default) for all channels.
	IAAC	yellow	+	
R-Biopharm RIDA®CYCLER	Turkey	green	+	Ignore cycles before, if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cyclor operating instructions, section 12.1.2 Cycling analysis parameter.
	IAC	yellow	+	
Roche LightCycler® 480 II	Turkey	465-510	+	
	IAAC	533-580	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control, Positive Control.

The reaction mix contains an internal amplification control (Inhibition control) per reaction.

Reactions needed for the qualitative detection of turkey:

3 reactions for controls* (1x negative control, 1x extraction control, 1 Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10% additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

The tube of the Taq Polymerase should not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

*** Description of the controls**

- Negative control: ready for PCR without any addition
- Extraction control: the extraction is performed without the sample – use only Lysis Buffer
- Positive Control: master-mix and within the kit's provided Positive Control

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Turkey DNA is detected in the FAM-channel. In the VIC/HEX-channel it is possible to detect animal DNA in the sample as well as the amplification control (IAC) in a sample with no animal DNA inside.

A sample is stated **positive** for turkey, if the sample DNA shows amplification in the detection system.

A sample is stated **negative** for turkey, if the sample DNA shows no amplification in the detection system. The internal amplification control (inhibition control) of the sample has to be **positive** with a shift in Cp-Value ≤ 2 compared to the negative control in the VIC/HEX-channel.

If the internal control signal in the VIC/HEX-channel of the sample DNA is detected significantly before the signal of the negative control the sample contains animal DNA.

Is the Cp-value of the internal control in the VIC/HEX-channel in the range of the negative control the sample contains no PCR-inhibiting substances but only a low amount or no animal DNA.

If the sample DNA in the VIC/HEX-channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal ($\geq 50\%$) can also show the presence of PCR inhibiting substances. Under these circumstances' DNA isolation and purification of the sample need to be improved. Alternatively, the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific turkey PCR assay.

3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- If the sample contains more than one animal species the DNA mixture can have a competitive influence on the absolute fluorescence. The lower the concentration of the determinant DNA is in a mixture of animal DNAs the lower is the fluorescence level of the amplification curve in the respective channel.
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.

4 Further Information

4.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

4.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

4.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com

r-biopharm®