CONGEN

SureFood® ANIMAL ID 3plex Water Buffalo/Beef+IAAC

Art. No. S6130
100 rxn

User Manual

February 2024

 **Inhalt**

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	DNA-preparation	10
1.4	Kit components and storage	10
1.5	Additionally required equipment and materials	10
1.6	Setup.....	10
1.7	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	12
2.2	Interpretation of results	13
3	Further Information	14
3.1	Product Information	14
3.2	Technical Support	14
3.3	Distribution and Ordering	14

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1 General Information

1.1 Description

The SureFood® ANIMAL ID 3plex Water Buffalo/Beef+IAAC is a real-time PCR for the direct, qualitative detection and differentiation of specific water buffalo (*Bubalis arnee*) and beef (*Bos taurus*) DNA sequences.

Each reaction contains an internal amplification control and an internal detection assay for vertebrates DNA (IAAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of three fluorescence emissions at the channels FAM, VIC/HEX and Cy5 at the same time. The technical verification of instruments was performed on Agilent Mx3005P, Agilent AriaDx, Applied Biosystems 7500, Bio-Rad CFX96, Roche cobas® z 480 Analyzer, Roche LightCycler® 480 II and Qiagen Rotor-Gene Q.

1.2 Limit of Detection

The SureFood® ANIMAL ID 3plex Water Buffalo/Beef+IAAC real-time PCR is developed for the detection of water buffalo and beef DNA in muscle meat mixture at a relative amount of 0.1 %.

The limit of detection of the complete method (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

Note: Inconsistent mixing ratios* may cause a loss of sensitivity in the low concentration channel in mixed samples especially with high amplicon concentrations (Cp value < 20).

* e.g. 99.9 water buffalo and 0.1 % beef

SureFood® ANIMAL ID 3plex Water Buffalo/ Beef+IAAC (100 rxn)

Art. No. S6130

February 2024

1.3 DNA preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic (Art. No. S1052)

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at **-28°C to -16°C and protected from light.**

1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Basic Art. No. S1052)
- real-time PCR instrument with three detection channels (510 nm, 580 nm and 660 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	35	35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

SureFood® ANIMAL ID 3plex Water Buffalo/ Beef+IAAC (100 rxn)

Art. No. S6130

February 2024

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	water buffalo	FAM	+	
	IAAC	HEX	+	
	beef	Cy5	+	
Agilent AriaMx / Dx	water buffalo	FAM	+	
	IAAC	HEX	+	
	beef	Cy5	+	
Applied Biosystems 7500	water buffalo	FAM	None	Check the passive reference option ROX is none.
	IAAC	VIC	None	
	beef	Cy5	None	
Bio-Rad CFX96 / Dx/Opus	water buffalo	FAM	+	
	IAAC	VIC/HEX	+	
	beef	Cy5	+	
Qiagen Rotor-Gene Q	water buffalo	green	+	Note: Please use only 0.1 ml reaction tube. The gain settings must be set to 5 (factory default) for all channels.
	IAAC	yellow	+	
	beef	red	+	
Roche LightCycler® 480 II	water buffalo	465-510	+	
	IAAC	533-580	+	
	beef	618-660	+	
Roche cobas® z 480 Analyzer	water buffalo	465-510	+	
	IAAC	540-580	+	
	beef	610-670	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control. The reaction mix contains an internal amplification control (IAC) per reaction.

Reactions needed for the qualitative detection of Water Buffalo and Beef:

3 reactions for controls (1x negative control, 1x extraction control, 1 positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

The tube of the Taq Polymerase should be kept at -28°C to -16°C, should not be thawed and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Water buffalo DNA is detected in the FAM-channel and beef DNA is detected in the Cy5-channel. In the VIC/HEX-channel it is possible to detect animal DNA in the sample as well as the amplification control (IAC) in a sample with no animal DNA inside. (see also table below)

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel. The internal amplification control (inhibition control) of the sample has to be **positive** with a shift in Cp-Value ≤ 2 compared to the negative control in the VIC/HEX-channel.

If the internal control signal in the VIC/HEX-channel of the sample DNA is detected significantly before the signal of the negative control the sample is stated **positive** for animal DNA.

Is the Cp-value of the internal control in the VIC/HEX-channel in the range of the negative control the sample contains no PCR-inhibiting substances but only a low amount or no animal DNA.

If the sample DNA in the VIC/HEX-channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific water buffalo or beef PCR assay.

result in the respective channel			Interpretation
FAM channel water buffalo	Cy5 channel beef	VIC/HEX channel animal + IAC	
positive	negative	positive	water buffalo DNA detected
negative	positive	positive	beef DNA detected
positive	positive	positive	water buffalo and beef DNA detected
negative	negative	positive, Cp < negative control	vertebrates DNA detected
negative	negative	negative	invalid

Note: If the sample contains more than one animal species the DNA mixture can have a competitive influence on the absolute fluorescence. The lower the concentration of the determinant DNA is in a mixture of animal DNAs the lower is the fluorescence level of the amplification curve.

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/company/downloads)
- Validation Report upon request

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com

r-biopharm®