CONGEN

SureFood® ANIMAL ID 4plex LIVESTOCK Panel

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Content

1	General Information	11
1.1	Description	11
1.2	Limit of Detection	11
1.3	DNA-preparation	12
1.4	Kit components and storage	12
1.5	Additionally required equipment and materials	12
1.6	Precautions for users	12
1.7	Setup.....	13
1.8	Detection channel Set-up	13
2	Qualitative Analysis	15
2.1	Protocol	15
2.1.1	Preparation of the master-mix	15
2.2	Preparation of the real-time PCR-mix.....	15
2.3	Interpretation of results	16
3	Limitations of the method	18
4	Further Information	18
4.1	Product Information	18
4.2	Technical Support.....	18
4.3	Distribution and Ordering	18

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1 General Information

1.1 Description

The SureFood® ANIMAL ID 4plex LIVESTOCK Panel is a real-time PCR for the direct, qualitative detection and differentiation of specific chicken (*Gallus gallus*), turkey (*Meleagris gallopavo*), goose (*Anser anser*), muscovy duck (*Cairina moschata*), American pekin (*Anas platyrhynchos*), pork (*Sus scrofa*), camel (*Camelus dromedarius*), horse (*Equus caballus*), donkey (*Equus asinus*), water buffalo (*Bubalis bubalis*), beef (*Bos taurus*), sheep (*Ovis aries*) and goat (*Capra hircus*) DNA sequences in food and feed.

The test contains an internal amplification control and an internal detection assay for vertebrates DNA (IAAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The internal technical verification of instruments was performed on Agilent AriaDx, Bio-Rad CFX96 Dx, R-Biopharm RIDA®CYCLER, Roche LightCycler® 480 II, Roche Lightcycler® PRO, Applied Biosystems QuantStudio 5 and Qiagen Rotor-Gene Q.

1.2 Limit of Detection

The SureFood® ANIMAL ID 4plex LIVESTOCK Panel real-time PCR is developed for the detection of animal DNA in muscle meat mixture at a relative amount of 0.1 %.

The limit of detection of the complete method (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

Note: Inconsistent mixing ratios may cause a loss of sensitivity in the low concentration channel in mixed samples especially with high amplicon concentrations (Cp value < 20).

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic (Art. No. S1052), SureFast® Mag PREP Food (Art. No. F1060) and for highly processed food and feed the use of SureFood® PREP Advanced (Art. No. S1053) is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053 / SureFast® Mag PREP Food Art. No. F1060)
- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	30	30
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

For additional details see cycler operating instructions for real-time PCR device.

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent AriaMx / Dx	chicken / turkey / goose / duck / pork	FAM	+	
	IAAC	HEX	+	
	camel / horse / donkey / water buffalo	ROX	+	
	beef / sheep / goat	Cy5	+	
Applied Biosystems QuantStudio 5	chicken / turkey / goose / duck / pork	FAM	+	Check the passive reference option ROX is none.
	IAAC	VIC	+	
	camel / horse / donkey / water buffalo	ROX	+	
	beef / sheep / goat	Cy5	+	
Bio-Rad CFX96 / Dx / Opus	chicken / turkey / goose / duck / pork	FAM	+	Baseline Settings: <ul style="list-style-type: none"> • Baseline subtracted curve fit • Apply fluorescence drift correction
	IAAC	VIC/HEX	+	
	camel / horse / donkey / water buffalo	ROX	+	
	beef / sheep / goat	Cy5	+	
R-Biopharm RIDA® CYCLER	chicken / turkey / goose / duck / pork	green	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter.
	IAAC	yellow	+	
	camel / horse / donkey / water buffalo	orange	+	
	beef / sheep / goat	red	+	

SureFood® ANIMAL ID 4plex LIVESTOCK Panel (100 rxn)

Art. No. S6134

June 2026

Real-time PCR device	Detection	Detection channel	Quencher	Note
Qiagen Rotor-Gene Q	chicken / turkey / goose / duck / pork	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAAC	yellow	+	
	camel / horse / donkey / water buffalo	orange	+	
	beef / sheep / goat	red	+	
Roche LightCycler® 480 II	chicken / turkey / goose / duck / pork	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAAC	533-580	+	
	camel / horse / donkey / water buffalo	533-610	+	
	beef / sheep / goat	618-660	+	
Roche Lightcycler® PRO	chicken / turkey / goose / duck / pork	465-510		
	IAAC	533-580		
	camel / horse / donkey / water buffalo	533-580		
	beef / sheep / goat	618-660		

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control and Positive Control.

The reaction mix contains an internal amplification control (Inhibition control) per reaction.

Reactions needed for the qualitative detection:

3 reactions for controls* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19,3 µl	212,3 µl
Taq Polymerase	0,7 µl	7,7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

* Description of the controls

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive Control: master-mix and within the kit's provided Positive Control

2.3 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Chicken, turkey, goose, duck and pork DNA is detected in the FAM-channel; camel, horse, donkey and water buffalo DNA is detected in the ROX-channel and beef, sheep and goat DNA is detected in the Cy5-channel. In the VIC/HEX-channel it is possible to detect animal DNA in the sample as well as the amplification control (IAC) in a sample with no animal DNA inside (see also table below).

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel. The internal amplification control (inhibition control) of the sample has to be **positive** with a shift in Cp-Value ≤ 2 compared to the negative control in the VIC/HEX-channel.

If the internal control signal in the VIC/HEX-channel of the sample DNA is detected significantly before the signal of the negative control the sample is stated **positive** for animal DNA.

Is the Cp-value of the internal control in the VIC/HEX-channel in the range of the negative control the sample contains no PCR-inhibiting substances but only a low amount or no animal DNA.

If the sample DNA in the VIC/HEX-channel shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. It is recommended to repeat the extraction process. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific PCR assay.

result in the respective channel				VIC/HEX channel animal + IAC	Interpretation
FAM channel chicken / turkey / goose / duck / pork	ROX channel camel / horse / donkey / water buffalo	Cy5 channel beef / sheep / goat			
positive	negative	negative	positive	chicken-, turkey-, goose-, duck- and/or pork DNA detected	
negative	positive	negative	positive	camel-, horse-, donkey and/or water buffalo DNA detected	
negative	negative	positive	positive	beef-, sheep- and/or goat DNA detected	
negative	negative	negative	positive, Cp < NTC*	no target DNA, but DNA from other vertebrates is detected	
negative	negative	negative	positive, Cp = NTC*	vertebrates DNA is not detected	
negative	negative	negative	negative	invalid	

Note: The results displayed in the table above represent merely an example. Additional combinations are also possible.

***NTC = negative control**

The following table shows the specification ranges of the kit controls

	Specification range
Positive Control (FAM – chicken / turkey / goose / duck / pork)	$16 \leq Cp \leq 24$
Positive Control (ROX – camel / horse / donkey / water buffalo)	$16 \leq Cp \leq 24$
Positive Control (Cy5 – beef / sheep / goat)	$16 \leq Cp \leq 24$

3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- If the sample contains more than one animal species the DNA mixture can have a competitive influence on the absolute fluorescence. The lower the concentration of the determinant DNA is in a mixture of animal DNAs the lower is the fluorescence level of the amplification curve.
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.

4 Further Information

4.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

4.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

4.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com