CONGEN

SureFood®
QUANT SOFT WHEAT

Art. No. S7010
2 x 50 rxn

User Manual

August 2025

 Inhalt

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Content

1	General Information	10
1.1	Description	10
1.2	Limit of Detection and Limit of Quantification	10
1.3	DNA-preparation	11
1.4	Kit components and storage	11
1.5	Additionally required equipment and materials	11
1.6	Precautions for users	11
1.7	Setup.....	12
1.8	Detection channel Set-up	12
2	Quantitative Analysis	13
2.1	Protocol	13
2.1.1	Preparation of the master-mix	13
2.1.2	Preparation of the standard DNA dilutions	14
2.1.3	Preparation of the real-time PCR-mix	14
2.2	Interpretation of results	15
3	Limitations of the method	16
4	Further Information	16
4.1	Product Information	16
4.2	Technical Support	16
4.3	Distribution and Ordering	16

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1 General Information

1.1 Description

SureFood® QUANT SOFT WHEAT is a real-time PCR for the detection of the relative soft wheat DNA content exclusively in durum wheat samples, such as pasta. Therefore, the kit contains two PCR systems, one for detection of a wheat-specific gene (*Triticum aestivum*) and one for the detection of a wheat gene (*Triticum* spp., reference gene). The Positive Control corresponds to 3 % soft wheat in durum wheat.

The real-time PCR assay can be used with established real-time PCR instruments. The internal technical verification of instruments was performed on Roche LightCycler® 480 II, Bio-Rad CFX96 Dx, Agilent AriaDx, Bio-Rad CFX OPUS 96, R-Biopharm RIDA®CYCLER and Qiagen Rotor-Gene Q.

1.2 Limit of Detection and Limit of Quantification

The SureFood® QUANT SOFT WHEAT specific PCR has a limit of detection of ≤ 5 DNA-copies.

The assay limit of detection depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The limit of quantification depends on the concentration of the sample DNA used in the analysis. For example, if 50,000 target-sequence copies of the reference gene are present, the relative quantification limit for soft wheat DNA is 0.1 %.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

SureFood® QUANT SOFT WHEAT

(2 x 50 rxn)

Art. No. S7010

August 2025

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic (Art. No. S1052), SureFood® PREP Advanced (Art. No. S1053) and SureFast® Mag PREP Food (Art. No. F1060) is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Wheat Reaction Mix	1 x 1050 µl	Orange
2	Soft wheat Reaction Mix	1 x 1050 µl	Yellow
3	Taq Polymerase	1 x 80 µl	Dark Red
4	Dilution Buffer	1 x 1400 µl	White
5	Standard DNA	1 x 45 µl	Dark Blue
6	Positive Control (3% Soft wheat)	1 x 95 µl	Light Blue

Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053 / SureFast® Mag PREP Food Art. No. F1060)
- real- time PCR instrument
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

SureFood® QUANT SOFT WHEAT (2 x 50 rxn)

Art. No. S7010

August 2025

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent AriaDx /Mx	Wheat	FAM	+	
	Soft wheat	FAM	+	
Bio-Rad CFX96/Dx/Opus	Wheat	FAM	+	Baseline Settings: • Baseline subtracted curve fit Apply fluorescence drift correction
	Soft wheat	FAM	+	
R-Biopharm RIDA®CYCLER	Wheat	green	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter.
	Soft wheat	green	+	
Qiagen Rotor-Gene Q	Wheat	green	+	Note: Please use only 0.1 ml reaction tube. The gain settings must be set to 5 (factory default) for all channels.
	Soft wheat	green	+	
Roche LightCycler® 480 II	Wheat	465-510	+	For the quantitative analysis use the type "Fit Points".
	Soft wheat	465-510	+	

2 Quantitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples, control reactions and standards) for the specific PCR assay.

Reactions needed for the wheat detection:

5 reactions for the standard curve

3 reactions for controls* (1x negative control, 2x Positive Control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the soft wheat detection:

5 reactions for the standard curve

3 reactions for controls* (1x negative control, 2x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

*** Description of the reactions**

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive control: master-mix and within the kit's provided Positive Control

2.1.2 Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different DNA concentrations for the standard curves of the reference gene (**wheat**) and the detection gene (**soft wheat**). Prepare 5 dilutions of the supplied Standard DNA (**Code 5**) with the supplied Dilution Buffer (**Code 4**). Prepare 5 reaction tubes (labeled S1 to S5) and add 45 µl Dilution Buffer (**Code 4**) each.

The following procedure is recommended:

Standard	Dilutions	Copy number per µl	Final copy number per reaction**
S1	45 µl Dilution Buffer + 5 µl Standard DNA	100,000 copies	500,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	50,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1,000 copies	5,000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	500 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S4	10 copies	50 copies

****Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

2.1.3 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control and the standard dilutions into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The calculation for both reaction systems (**wheat**, **soft wheat**) has to be made separately. Mark the standards, the controls and the samples for the specific system (**soft wheat**) and make the evaluation according to the usual analysis program recommended by the real-time PCR instrument manufacturer. Repeat the same procedure for the reference gene system (**wheat**). The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

By using the calculated copy numbers for **wheat** and **soft wheat** the relative **soft wheat** content of the sample DNA and the Positive Control can be determined in the following way (example):

Sample soft wheat	1,350 copies	Positive Control soft wheat	400 copies
Sample wheat	45,000 copies	Positive Control wheat	28,000 copies

Divide the copy number of the specific system by the copy number of the reference gene system and multiply by 100 to obtain the percentage.

soft wheat DNA content = **soft wheat** copy number * 100 / **wheat** copy number

sample **soft wheat** DNA content = $1,350 * 100 / 45,000$ sample **soft wheat** DNA content = 3%

For the given example the numbers lead to a **soft wheat** DNA content of 3.0%, for the sample and 1.4% for the Positive Control with the same calculation.

For a final calculation the use of a correction factor K for the correction of run-to-run fluctuations is necessary. The correction factor is the relation of the true percentage value of the Positive Control (1% **soft wheat** content) and the measured percentage of the Positive Control.

$K = \text{true value} / \text{measured value}$

$K (\text{example}) = 1\% / 1.4\% = 0.7$

The measured DNA content for the sample is multiplied with K to obtain a corrected DNA content.

sample DNA content = measured sample DNA content * K sample (example) = $3.0\% * 0.7 = 2.1\%$

For this example the **soft wheat content is 2.1%**.

3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.

4 Further Information

4.1 Product Information

- Microsoft Excel template of calculation and detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

4.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

4.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com