**CONGEN**


# SureFast<sup>®</sup> PREP DNA/RNA Virus

Art. No. F1051  
100 extractions

## User Manual

Efficient DNA preparation of viral or bacterial DNA/RNA

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## Content

1	General Information .....	9
1.1	Description .....	9
1.2	Principle .....	9
1.3	Kit components and storage .....	9
1.4	Additionally required equipment and materials .....	9
2	Protocol .....	10
2.1	Preparations .....	10
2.1.1	General .....	10
2.1.2	Before each preparation .....	10
2.2	Preparation of the basic material .....	10
2.2.1	Liquid samples .....	10
2.2.2	Stool samples .....	10
2.2.3	Swabs and filters .....	10
2.3	Sample preparation for surface swabbing for the detection of SARS-CoV-2 .....	10
2.3.1	Surface swabbing .....	10
2.3.2	Sample preparation for surface swabbing for direct analysis without transport of the swab .....	10
2.3.3	Sample preparation for surface swabbing with transport medium (UTM) after the external transport (48 hours maximal transport duration) .....	11
2.4	Extraction protocol .....	11
3	Further Information .....	12
3.1	Product Information .....	12
3.2	Technical Support .....	12
4	Safety information .....	12

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## 1 General Information

### 1.1 Description

The kit is intended to be used for the isolation of DNA and RNA viruses nucleic acids from serum, urine, plasma, cell culture supernatants, foods (for example wash up fluids from fruits, salads etc.), filters from water samples as well as stool samples and swabs (from work surfaces etc.). SureFast® PREP DNA/RNA virus in combination with SureFast® SARS-CoV-2 PLUS (Art. Nr. F7110) has been validated and certified for stainless steel environmental swabs by AOAC (license no. 022102).

### 1.2 Principle

1. Preparation of the basic material
2. Lysis of the basic material
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

### 1.3 Kit components and storage

Kit Code	Reagent / Material	Amount (per box)	Color
<b>B</b>	Binding Buffer	1 x 30 ml	
<b>E</b>	Elution Buffer	1 x 5 ml	
<b>P</b>	Pre-Wash Buffer*	1 x 40 ml	
<b>W</b>	Wash Buffer*	1 x 60 ml	
<b>ET</b>	Extraction Tubes	1 x 50 Tubes	clear
<b>R</b>	Receiver Tubes 2,0 ml	1 x 50 Tubes	clear
<b>T</b>	Receiver Tubes 1,5 ml	1 x 50 Tubes	clear
<b>S</b>	Spin Filter Set	1 x 50 Filter	pink

\* After adding ethanol (purity  $\geq 96\%$ ; not supplied with the kit)

All reagents of the kit should be stored dry and at room temperature (14-25°C). After opening the stability of the kit is reduced to 12 months.

### 1.4 Additionally required equipment and materials

- Microcentrifuge
- Thermomixer/heating block (to 100°C)
- Pipettes with filter tips
- Reaction tubes 1.5 ml
- Ethanol ( $\geq 96\%$ )
- RNase free PCR grade H<sub>2</sub>O

## 2 Protocol

### 2.1 Preparations

#### 2.1.1 General

Add 20 ml ethanol to the Pre-Wash Buffer (**Code P**) and mix thoroughly.

Add 48 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

#### 2.1.2 Before each preparation

Preheating the Elution Buffer (**Code E**) - Transfer the needed amount under calculation of a reserve volume of Elution Buffer (**Code E**) into a reaction tube (not provided with the kit)

Equilibrate the Elution Buffer to 65°C. The Elution Buffer is necessary for step 6.

### 2.2 Preparation of the basic material

#### 2.2.1 Liquid samples

Mix 200 µl samples fluids with 200 µl RNase free PCR grade H<sub>2</sub>O in a reaction tube (not supplied with the kit). For samples which have a smaller volume than 200 µl please fill up to a total volume of 400 µl with RNase free PCR grade H<sub>2</sub>O. By small yields of nucleic acids use directly 400 µl sample fluid.

Transfer the complete sample in an Extraction Tube (**Code ET**).

#### 2.2.2 Stool samples

Dilute stool samples 10 to 20 fold in RNase free PCR grade H<sub>2</sub>O (for a lower viscosity). Centrifuge the dilution for 2 min at maximum speed. Mix 200 µl supernatant with 200 µl RNase free PCR grade H<sub>2</sub>O. Transfer the complete sample in an Extraction Tube (**Code ET**).

#### 2.2.3 Swabs and filters

Place swabs and filters complete in the Extraction Tube (**Code ET**) (Cut off the shaft of the swab, so that the cap of the extraction tube can be closed). Add 400 µl RNase free PCR grade H<sub>2</sub>O.

### 2.3 Sample preparation for surface swabbing for the detection of SARS-CoV-2

#### 2.3.1 Surface swabbing

The swab should be pre-moisten with 100 µl UTM prior of swabbing the surface. Swab a surface of approximately 25 cm<sup>2</sup> (2 x 2 inch according to AOAC performance tested method (PTM) for stainless steel surfaces) with the wet swab.

For transport to an external laboratory according to the AOAC PTM, put the swab back into the Copan tube with UTM and send it to the designated laboratory. After transport follow with 2.3.3.

#### 2.3.2 Sample preparation for surface swabbing for direct analysis without transport of the swab

Place the swab with the shaft in an Extraction Tube (**Code ET**) and cut the shaft, so that the lid of the Extraction Tube can be closed. Add 400 µl RNase free PCR grade H<sub>2</sub>O.

### 2.3.3 Sample preparation for surface swabbing with transport medium (UTM) after the external transport (48 hours maximal transport duration)

Place the swab with the shaft in an Extraction Tube (**Code ET**) and cut the shaft, so that the lid of the Extraction Tube can be closed. Add 400 µl of the transported UTM. Transportation of the swab samples with UTM should take place at 2-8 °C.

**Note:** After lysis time carefully squeeze out the swab on the wall of the tube and discard the swab.

If the Internal Control RNA is used as an inhibition and extraction control, 20 µl of the ICR must be added after sample preparation in the Extraction tube.

## 2.4 Extraction protocol

### 1. Lysis of the basic material

Add 400 µl of Lysis Buffer (**Code L**) to the reaction tube and mix it briefly and close the reaction tube with caps. Incubate on a heating block under continuously shaking for 10 minutes at 95°C.

### 2. Setting of optimal binding conditions

Add directly 400 µl Binding Buffer (**Code B**) to fluids without particles and mix the sample by vortexing.

**Note:** Does the sample lysate show particles it has to be centrifuged for 1 min at maximum speed. The liquid supernatant has to be transferred into a new reaction tube (not supplied with the kit). Add 400 µl Binding Buffer (**Code B**) to the supernatant and mix the sample by vortexing.

### 3. Binding of the nucleic acids on a Spin Filter

Transfer the complete sample in a Spin Filter Set (**Code S**). Close the cap and centrifuge for 1 minute at 12.000 rpm. Discard the Receiver Tube with the filtrate and place the Spin Filter in a new 2.0 ml Receiver Tube (**Code R**)

### 4. Purification of the bound nucleic acids

Add 550 µl Pre-Wash Buffer (**Code P**) to the Spin Filter and centrifuge at 1 min for 10,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Add 700 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 10,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

### 5. Drying of the Spin Filter

Remove the residual ethanol by final centrifugation for 4 min at maximum speed.

### 6. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 60 µl of the preheated to 65°C Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate for 3 min and centrifuge for 1 min at 10.000 rpm. After centrifugation discard the Spin Filter.

The eluted nucleic acid is ready-to-use for the PCR or RT-PCR. Store the nucleic acids at -20°C or -80°C.

**Note:** The nucleic acids can also be eluted with a lower (but not lower than 40 µl) or a higher volume of Elution Buffer (depends on the expected yield or needed concentration of the nucleic acids).

## 3 Further Information

### 3.1 Product Information

- Flow chart (Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report upon request
- Material Safety Data Sheet upon request

### 3.2 Technical Support

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).

## 4 Safety information

### Binding Buffer

Danger

### Pre-Wash Buffer

Warning

### Extraction Tubes

Warning

H225-319-336 P210-233-(305-351-338)

H302-312-332-412 EUH032 P273-(304+340)

H302-315-319-335-411 P260-262-(305+351+338)-310 EUH208

**H225:** Highly flammable liquid and vapour.

**H319:** Causes serious eye irritation.

**H302:** Harmful if swallowed.

**H312:** Harmful in contact with skin.

**H315:** Causes skin irritation.

**H319:** Causes serious eye irritation.

**H332:** Harmful if inhaled.

**H335:** May cause respiratory irritation.

**H336:** May cause drowsiness or dizziness.

**H411:** Toxic to aquatic life with strong lasting effects

**H412:** Harmful to aquatic life with long lasting effects.

**EUH032:** Contact with acids liberates very toxic gas.

**EUH208:** Contains Proteinase K. May produce an allergic reaction.

**P210:** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

**P233:** Keep container tightly closed.

**P260:** Do not breathe dust / fume / gas / mist / vapours / spray.

**P262:** Do not get in eyes, on skin, or on clothing.

**P273:** Avoid release to the environment.

**P304+340:** IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

**P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information we offer a Material Safety Data Sheet. Please contact your distributor or send an e-mail to [info@congen.de](mailto:info@congen.de).