CONGEN

SureCycle®

Art. No. F4001
260 rxn

User Manual **System Verification**

June 2022

 **Inhalt**

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Content

1	General Information	14
1.1	Intended Use.....	14
1.2	Description	14
1.3	Kit components and storage	15
1.4	Additionally required equipment and materials	15
1.5	Setup.....	15
2	Protocol	17
2.1	Recommended pipetting scheme for system verification.....	17
2.2	Preparation of the master-mix	18
2.3	Preparation of the real-time PCR mix.....	18
3	Interpretation of results	18
3.1	SureCycle® profile 1 Denaturation	18
3.1.1	Denaturation system (VIC channel).....	18
3.1.2	Annealing system (FAM-channel)	20
3.2	SureCycle® Profile 2 Annealing.....	21
3.2.1	Annealing system (FAM-channel)	21
3.2.2	Denaturation system (VIC-channel).....	23
4	Further Information	23
4.1	Product Information	23
4.2	Technical Support	23
4.3	Explanation of Symbols	23

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1 General Information

1.1 Intended Use

The SureCycle®-Kit is intended to use for real-time PCR routine control in all laboratories. Especially laboratories which need a constant control concerning their testing devices (DIN EN ISO 15189 and DIN EN ISO 9001, GMP and GLP) may use the SureCycle®-Kit as a reliable and fast verification method.

The kit is designed to control compatibility of the selected plastic material and to verify the correct temperature control of the thermocycler.

1.2 Description

The real-time Polymerase chain reaction (PCR) is an established method for the specific and sensitive detection of DNA. The area of use spans the clinical infection diagnostics to the analysis of pathogens in food & feed as well as forensics. Each PCR-reaction basically is defined by the primer sequence and the PCR protocol with specific temperature and time intervals. The detection is therefore dependent to the correct temperature system of the thermocycler in use. To prevent false positive or false negative results, the control and verification of the correct temperature system is obligatory.

The SureCycle® -System is a thermosensitive duplex real-time PCR. The temperature sensitivity of the system results from the instability of the primer binding. Binding of the primer to the complementary chain is being destabilized and thereby the annealing- and denaturation reaction become temperature sensitive. Such a system is reacting very sensitively to 'overshooting' of thermocyclers. The test is equally sensitive regarding the inadequate denaturation of double-stranded DNA (dsDNA) due to a too-low denaturation temperature ('undershooting'). Both, 'overshooting' and 'undershooting' of denaturation temperature reduces the efficiency of the PCR.

To control the correct temperature during the annealing- and denaturation reactions, two thermal profiles are used: SureCycle® Profile 1 Denaturation and SureCycle® Profile 2 Annealing.

With Profile 1 an optimal denaturation temperature is programmed, as well as a suboptimal annealing temperature. In this case a positive result is to be expected for the VIC-channel and a negative result in the FAM-channel. Using Profile 2 the opposite situation is evaluated, here the optimal annealing temperature and the suboptimal denaturation temperature is performed. Therefore a positive result in case of the FAM-channel and a negative result in case of the VIC-channel is expected.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for the detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The technical validation of instruments was performed on Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, Bio Molecular Systems MIC, LTF MyGo Pro, Agilent AriaDx and Agilent MX3005P.

1.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	4 x 1400 µl	Yellow
2	Taq Polymerase	1 x 200 µl	Dark Red

Store all reagents at -20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple use on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.4 Additionally required equipment and materials

- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.5 Setup

	SureCycle® Profile 1 Denaturation	SureCycle® Profile 2 Annealing
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	15 sec, 92°C
Annealing/Extension	30 sec, 63°C	30 sec, 60°C
Cooling	10 sec, 40 °C	10 sec, 40 °C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

Real-time PCR device	Detection channel	Quencher	Note
Agilent Mx3005P	FAM	+	-
	HEX	+	
Agilent AriaDx	FAM	+	-
	HEX	+	
Applied Biosystems 7500	FAM	None	Check the passive reference option ROX is none.
	VIC	None	
Bio-Rad CFX96	FAM	+	-
	VIC / HEX	+	
R-Biopharm RIDA®CYCLER	green	+	-
	yellow	+	
LTF MyGo Pro	FAM	+	-
	VIC	+	
Qiagen Rotor-Gene Q	green	+	-
	yellow	+	
Roche LightCycler® 2.0	530	None	The SureCC Color Compensation Kit II (Art. No. F4010) is required. Check the Seek Temperature is 58°C.
	560	None	
Roche LightCycler® 480 II	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	533-580	+	
Roche cobas® z 480 Analyzer	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	540-580	+	

2 Protocol

2.1 Recommended pipetting scheme for system verification

Blockcycler:

For verification of the real-time PCR instrument it is recommended to pipette 4 x 4 wells in the corner positions (highlighted in red) and 4 wells in the center (highlighted in red) per each SureCycle® profile (see figure 1).

	1	2	3	4	5	6	7	8	9	10	11	12
A	26,8	26,6									26,8	26,4
B	26,8	26,8									26,9	26,7
C												
D						26,9	27,0					
E						27,0	27,0					
F												
G	26,7	26,8									26,8	26,7
H	26,6	26,8									26,6	26,7

Figure 1: Blockcycler pipetting scheme in 96-well format

Rotorcycler:

It is recommended to pipette 20 reactions per SureCycle®. 5 reactions should be used as reference reactions (see figure 2, position 1-5, highlighted in red).

In the figure below the Rotor-Gene Q (72-well disc, Qiagen) is represented as an example. The positioning for use with a different rotor size (MIC, LightCycler® 2.0, RIDA®CYCLER) can be adjusted according to the reactions.


Figure 2: Rotorcycler pipetting scheme for e.g. 72-well disc, Rotor-Gene Q (Qiagen)

2.2 Preparation of the master-mix

The total number of reactions needed for the PCR should be calculated according to the recommended plate set-up with 20 reactions per SureCycle® profile (see 2.1).

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 20 reactions:

Components of the master-mix	Amount per reaction	20 reactions (with 10 % excess)
Reaction Mix	19.3 µl	424.6 µl
Taq Polymerase	0.7 µl	15.4 µl
Total volume	20 µl	440 µl

Mix each master-mix well and centrifuge shortly before use.

2.3 Preparation of the real-time PCR mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start run according to the setup.

3 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

A reaction is evaluated as **positive**, if an amplification is shown in the respective channel.

A reaction is evaluated as **negative**, if no amplification is shown in the respective channel.

3.1 SureCycle® profile 1 Denaturation

3.1.1 Denaturation system (VIC channel)

The SureCycle® Profile 1 Denaturation takes place with an optimal denaturation temperature at 95°C and suboptimal Annealing/Extension-temperature of 63°C.

In this case positive results are to be expected for the Denaturation-System (VIC channel). The FAM-channel controls the correct Annealing/Extension. Since the temperature setting is suboptimal, no amplification should occur and hence negative results are to be expected.

Blockcycler

The target value is calculated with the mean value of the reference reactions (Mean value D6/D7/E6/E7, highlighted in red). The target range of control reactions is defined as ± 2 Cp of the reference reactions. The real-time PCR instrument passed the verification, if Cp-values of all control reactions are within the target range.

The real-time PCR instrument did not pass the verification, if one or more Cp-values of control reactions are out of the target range.

In case there is a Cp-value deviation of a control reaction from the target range, a complete inspection of all wells should be performed.

	1	2	3	4	5	6	7	8	9	10	11	12
A	26.8	26.6									26.8	26.4
B	26.8	26.8									26.9	26.7
C												
D						26.9	27.0					
E						27.0	27.0					
F												
G	26.7	26.8									26.8	26.7
H	26.6	26.8									26.6	26.7

Figure 3: Example calculation SureCycle® Profile 1 Denaturation, VIC-channel

Mean value D6/D7/E6/E7: **27.0** Target range ± 2 Cp: **25.0 – 29.0**

Rotorcycler

The target value is calculated with the mean value of the reference reactions (Mean value Position 1-5, highlighted in red). The target range of control reactions is defined as ± 2 Cp of the reference reactions. The real-time PCR instrument passed the verification, if Cp-values of all control reactions are within the target range.

The real-time PCR instrument did not pass the verification, if one or more Cp-values of control reactions are out of the target range.

In the figure below the Rotor-Gene Q (72-well disc, Qiagen) is represented as an example. The positioning for use with a different rotor size can be adjusted according to the reactions.


Figure 4: Example calculation SureCycle® Profile 1 Denaturation, VIC channel

Mean value Pos. 1-5: **27.0** target-range ± 2 Cp: **25.0-29.0**

Table 4: Results Interpretation of SureCycle® Profile 1 Denaturation (VIC-channel)

VIC-channel Denaturation-System		Result	Consequence for further analysis
Reference reaction	Control reaction		
positive	Cp within Target range	correct	Real-time PCR instrument temperature system works correctly
positive	Cp out of target range	incorrect	Denaturation temperature too low → potential false PCR-results
negative	Cp; no target range	incorrect	Denaturation temperature too low → potential false PCR-results

3.1.2 Annealing system (FAM-channel)

All results should be negative.

The real-time PCR instrument passed the verification, if all reactions are negative.

In case control or reference reactions show positive Cp-values, the results have to be interpreted according to table 2.

Table 5: Results Interpretation of SureCycle® Profile 1 Denaturation (FAM-channel)

FAM-channel Annealing-System		Result	Consequence for further analysis
Reference reaction	Control reaction		
negative	negative	correct	Real-time PCR instrument temperature system works correctly
positive	negative	incorrect	Annealing temperature too low → potential false PCR-results
positive	positive	incorrect	Annealing temperature too low → potential false PCR-results

3.2 SureCycle® Profile 2 Annealing

3.2.1 Annealing system (FAM-channel)

The SureCycle® Profile 2 Annealing takes place with a suboptimal denaturation temperature at 92°C and the optimal Annealing/Extension-temperature at 60°C.

In this case negative results are to be expected for the Denaturation-System (VIC-channel). The FAM-channel controls the correct Annealing/Extension. Since the temperature setting is optimal, amplification should occur and hence positive results are to be expected.

Blockcycler

The target value is calculated with the mean value of the reference reactions (Mean value D6/D7/E6/E7, highlighted in red). The target range of control reactions is defined as $\pm 2 \text{ Cp}$ of the reference reactions.

The real-time PCR instrument passed the verification, if Cp-values of all control reactions are within the target range.

The real-time PCR instrument did not pass the verification, if one or more Cp-values of control reactions are out of the target range.

In case there is a Cp-value deviation of a control reaction from the target range, a complete inspection of all wells should be performed.

	1	2	3	4	5	6	7	8	9	10	11	12
A	26.8	26.6									26.8	26.4
B	26.8	26.8									26.9	26.7
C												
D						26.9	27.0					
E						27.0	27.0					
F												
G	26.7	26.8									26.8	26.7
H	26.6	26.8									26.6	26.7

Figure 5: Example calculation SureCycle® Profile 2 Annealing, FAM-channel


Mean value D6/D7/E6/E7: **27.0** Target-range $\pm 2 \text{ Cp}$: **25.0-29.0**

Rotorcycler

The target value is calculated with the mean value of the reference reaction (Mean value Position 1-5, highlighted in red). The target range of control reactions is defined as $\pm 2 C_p$ of the reference reactions. The real-time PCR instrument passed the verification, if C_p -values of all control reactions are within the target range.

The real-time PCR instrument did not pass the verification, if one or more C_p -values of control reactions are out of the target range.

In the figure below the Rotor-Gene Q (72-well disc, Qiagen) is represented as an example. The positioning for use with a different rotor size can be adjusted according to the reactions.

**Figure 6: Example calculation SureCycle® Profile 2
Annealing, FAM-channel**

Mean value Pos. 1-5: **27.0** target-range $\pm 2 C_p$: **25.0-29.0**

Table 6: Result Interpretation of SureCycle® Profile 2 Annealing (FAM-channel)

FAM-channel Annealing-System			
Reference reaction	Control reaction	Result	Consequence for further analysis
positive	C_p within target range	correct	Real-time PCR instrument temperature system works correctly
positive	C_p out of target range	incorrect	Anneling temperature too low → potential false PCR-results
negative	C_p ; no target range	incorrect	Anneling temperature too low → potential false PCR-results

3.2.2 Denaturation system (VIC-channel)

All results should be negative.

The real-time instrument passed the verification, if all reactions are negative.

In case of positive reference reactions the target value is calculated with the mean value of the reference reactions on a Blockcyler D6/D7/E6/E7 and positions 1-5 on a Rotorcyler respectively (highlighted in red).

The target range of control reactions is defined as $\pm 2 C_p$ of the reference reactions.

In case control or reference reactions show positive C_p -values, the results have to be interpreted according to table 4.

Table 4: Result Interpretation of SureCycle® Profile 2 Annealing (VIC-channel)

FAM-channel Denaturation-System		Result	Consequence for further analysis
Reference reaction	Control reaction		
negative	negative	correct	Real-time PCR instrument temperature system works correctly
positive	negative or C_p out of target range	incorrect	Denaturation temperature too high → uncritical, false PCR-results are not to be expected
positive	C_p within target range	incorrect	Denaturation temperature too high → uncritical, false PCR-results are not to be expected

4 Further Information

4.1 Product Information

Detailed information about setup of several real-time PCR devices

(Download: www.congen.de/en/company/downloads)

4.2 Technical Support

For further questions please send an e-mail to info@congen.de.

4.3 Explanation of Symbols

Lot number

Article number

Expiry date

Number of Tests

Storage temperature

Manufacturer