r-biopharm®

## RIDA® QUICK T-2 / HT-2 RQS ECO

**REF** R5304

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# RIDA®QUICK T-2 / HT-2 RQS ECO

## Brief information

RIDA®QUICK T-2 / HT-2 RQS ECO (Art. No. R5304) is an immuno-chromatographic test for the quantitative determination of T-2 / HT-2 toxin in oats, wheat and corn (maize) (see chapter 1. Intended use).

The test kit contains 20 test strips, each of which can be used for one analysis. All reagents required for the assay are contained in the test kit. Results are evaluated with the RIDA®SMART APP software (Art. No. ZRSAM1000) or the RIDA®QUICK SCAN reader.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)... approx. 10 min  
test implementation (incubation time)..... 5 min

Limit of detection: 50 µg/kg (ppb)

Measuring range: 50 - 1,000 µg/kg und 1,000 - 10,000 µg/kg  
Quantitative evaluation with the RIDA®SMART APP

Specificity: The RIDA®QUICK T-2 / HT-2 RQS ECO test detects T-2 / HT-2 toxin in oats, wheat and corn (maize) samples.

The specificity of the RIDA®QUICK T-2 / HT-2 RQS ECO test was determined by analyzing the cross reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

## Related products and accessories for T-2 / HT-2 RQS ECO determination

RIDASCREEN® T-2 Toxin (Art. No. R3801)  
RIDASCREEN® T-2 / HT-2 Toxin (Art. No. R3805)  
RIDASCREEN®FAST T-2 Toxin (Art. No. R5302)

## 1. Intended use

RIDA<sup>®</sup>QUICK T-2 / HT-2 RQS ECO (Art. No. R5304) is a quantitative immunochromatographic test in strip format for the quantitative determination of T-2 / HT-2 toxin in oats, wheat and corn (maize).

## 2. General information

T-2 and HT-2 toxin belong to the trichothecene group of mycotoxins and are formed by fungi of the genus *Fusarium*. T-2 and HT-2 toxin are often found in agricultural commodities, although the incidence and the concentrations found show a broad regional variation. Due to their high cytotoxic and immunosuppressive mode of action T-2 / HT-2 toxins are a threat for human and animal health.

## 3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti- T-2 / HT-2 toxin antibody detects T-2 / HT-2 toxin in the sample. During incubation of the test strip, a band pattern (test line / control line) forms that is used to determine the concentration of T-2 / HT-2 toxin. The intensity of the test line depends on the T-2 / HT-2 toxin concentration of the sample. It increases as the T-2 / HT-2 toxin concentration increases. It must be possible to detect the control line after the reaction has taken place in order to verify the functioning of the test. The control line weakens as the T-2 / HT-2 toxin concentration in the sample increases.

The test strip is evaluated using the RIDA<sup>®</sup>SMART APP software (Art. No. ZRSAM1000) or the RIDA<sup>®</sup>QUICK SCAN reader.


Fig. 1: RIDA<sup>®</sup>QUICK T-2 / HT-2 RQS ECO test strip.

## 4. Reagents provided

Each kit contains sufficient materials for max. 20 determinations. Each test kit contains:

Component	Cap color	Format	Volume
20 x Test strip	-	Ready to use	One determination per test strip (separately packed)
1 x Extraction buffer	Transparent	Concentrate 9x	115 ml
1 x Mobile solvent	White	Ready to use	8 ml
1 x Sample buffer	Transparent	Ready to use	20 ml
1 x RIDA®SMART APP cover	-	Ready to use	Lot-specific evaluation cover for use with RIDA®SMART APP

## 5. Reagents required but not provided

### 5.1 Equipment

- Gloves
- Laboratory grinder
- Balance
- Graduated cylinder
- Optional: Shaker
- Optional: Centrifuge ( $\geq 2000$  g)
- Optional: Filter paper (e.g. Whatman No. 1 or equivalent)
- 100  $\mu$ l, 150  $\mu$ l, 350  $\mu$ l and 900  $\mu$ l micropipettes
- Latest version of the RIDA®SMART APP software (Art. No. ZRSAM1000) activated and installed on a compatible smartphone (or RIDA®QUICK SCAN reader)

### 5.2 Reagents

- Distilled water (dist. Water) or deionized water

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze the test strips or any test kit components.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. For this reason, it is important to protect the strips against moisture. Do not remove the test strips from the test strip packaging until immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Sample preparation

For sample preparation a **ready to use extraction buffer** is needed. For this dilute the extraction buffer (concentrate 9x) 1:9 (1+8) with distilled water. The extraction buffer should be handled under clean working conditions. The ready to use extraction buffer should be stored at 2 - 8 °C (36 - 46 °F) and expires after 8 - 10 weeks.

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

The sample should be stored in a cool place, protected against light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

### 8.1 Detection limit 50 µg/kg

- Weigh 10 g of ground and mixed sample into a suitable container and add 45 ml of ready to use extraction buffer
- Close the tube and shake the sample vigorously for approx. 3 min (manually or with shaker)

- Let the solution settle ( $\leq 5$  min), filtrate or centrifuge (2 min, 2000 g)
- Mix 150  $\mu\text{l}$  of the clear supernatant with 350  $\mu\text{l}$  mobile solvent
- Use 100  $\mu\text{l}$  as sample solution (dilution I) in the test

To analyze samples with a T-2 / HT-2 toxin concentration of  $> 1.000 \mu\text{g}/\text{kg}$  (evaluation with RIDA<sup>®</sup>SMART APP) or  $> 800 \mu\text{g}/\text{kg}$  (evaluation with RIDA<sup>®</sup>QUICK SCAN) the samples solution has to be further diluted. For this please use the sample buffer contained in the test kit and proceed with section 8.2.

## **8.2 Detection limit 1000 $\mu\text{g}/\text{kg}$ (RIDA<sup>®</sup>SMART APP evaluation) or 800 $\mu\text{g}/\text{kg}$ (RIDA<sup>®</sup>QUICK SCAN evaluation)**

- Mix 100  $\mu\text{l}$  of dilution I with 900  $\mu\text{l}$  sample buffer to obtain dilution II
- Use 100  $\mu\text{l}$  of dilution II as sample solution in the test

### Important:

For samples that are evaluated with the RIDA<sup>®</sup>SMART APP please choose the appropriate RIDA<sup>®</sup>SMART APP software application for your sample preparation and matrix (see section 10.1). For samples diluted according to section 8.2 and evaluated with the RIDA<sup>®</sup>QUICK SCAN reader, the obtained result has to be multiplied by the **factor 10** to get the actual T-2 / HT-2 toxin concentration of the sample.

## **9. Test procedure**

Bring all components necessary for sample preparation to room temperature (20 - 25 °C / 68 - 77 °F) before use in the test and perform sample preparation at room temperature.

- Remove a test strip from its foil pouch.
- Apply 100  $\mu\text{l}$  of the sample solution (see section 8. Sample preparation) on the application area of the test strip (see Fig. 1).
- Incubate the test strip at room temperature.
- Evaluate the result after 5 min (+ max. 15 s) incubation time with the RIDA<sup>®</sup>SMART APP or RIDA<sup>®</sup>QUICK SCAN reader.

## **10. Evaluation**

The left line in the reaction area is the test band (test line, see fig. 1). Its appearance and intensity depends on the T-2 / HT-2 toxin concentration of the

sample. The sample is contaminated with T-2 / HT-2 toxin, if the control band (control line) is visible and the test band (test line) is also visible.

The right line in the reaction area is a control band (control line) and must be present after each test procedure. If this line is not visible, the test was not performed correctly or the reagents were not in order. In this case, the test should be repeated with a new test strip. If the right band is missing again, please inform your local distributor.

**Please note:** If the sample is highly contaminated with T-2 / HT-2 toxin, the control line will be only faintly visible. In this case the test result might be reported as “invalid result“(RIDA®SMART APP) or “invalid“ (RIDA®QUICK SCAN reader). In this instance, dilute the sample as described in section 8.2 and repeat the test with a new test strip.

To obtain comparable results, the test strips must always be evaluated after 5 min.

### 10.1 Evaluation with RIDA®SMART APP

First, please read the instructions for use of the RIDA®SMART APP software (Art. No. ZRSAM1000) attentively. The test strip evaluation has to be performed by using the RIDA®SMART APP software according to these guidelines. Please find a short description (Quick Starter Guide) of how to use the RIDA®SMART APP on our website [www.r-biopharm.de](http://www.r-biopharm.de).

For the test strip evaluation the RIDA®SMART APP cover is needed which is included in the test kit. The RIDA®SMART APP cover contains lot-specific information and must be used only with the appropriate lot of the test kit. For further information please check the Quality Assurance Certificate enclosed in the kit. Please make sure the correct RIDA®SMART APP evaluation cover is present prior to starting the test procedure and test strip evaluation with the RIDA®SMART APP software.

RIDA®SMART APP software applications for evaluation:

Matrix	Measuring range	RIDA®SMART APP application
Oats, corn	50 - 1000 µg/kg	Oats_Corn 50-1000ppb
Oats, corn	1.000 - 10000 µg/kg	Oats_Corn 1000-10000ppb
Wheat	50 - 1000 µg/kg	Wheat 50-1000ppb
Wheat	1.000 - 10000 µg/kg	Wheat 1000-10000ppb

## 10.2 Evaluation with RIDA®QUICK SCAN

First, carefully read the user guide for the RIDA®QUICK SCAN and proceed as instructed. Section 3. Test principle contains the description for evaluating the test strips.

To obtain optimal results every test kit contains lot-specific standard curve parameters. These parameters can be found on the Quality Assurance Certificate and have to be used only with this specific lot. The Quality Assurance Certificate contains parameters for several standard curves. Depending on the extraction method used, select the corresponding parameters.

To enter the parameters manually into the RIDA®QUICK SCAN reader use the plain text information on the Quality Assurance Certificate. To enter the data with the barcode scanner use the 2D barcodes. In both cases, the checksum shown on the RIDA®QUICK SCAN screen needs to be identical to the respective checksum on the certificate (number below the 2D barcode).

RIDA®QUICK SCAN software method for test evaluation: **T2/HT2 Quant**

## 11. Recommendation

In order to ensure a high analytical performance we recommend:

- To do spike experiments to ensure an accurate and correct test procedure

## 12. Sensitivity

This test determines T-2 / HT-2 toxin contaminations in the measurement range of  $\geq 50 - 1000 \mu\text{g/kg}$  or  $1000 - 10,000 \mu\text{g/kg}$  in combination with the RIDA®SMART APP software and in the measurement range of  $50 - 800 \mu\text{g/kg}$  or  $500 - 8000 \mu\text{g/kg}$  in combination with the RIDA®QUICK SCAN reader.

## 13. Further application notes

Further application notes are available on request. Please contact your local distributor.

**For further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Version overview

Version number	Chapter and title
2014-09-15	Release version
2020-05-20	General revision
2021-10-12	Current version General revision Changes made: <ul style="list-style-type: none"><li>– Related products and accessories</li><li>– 3. Test principle (Image of new housing)</li></ul>

## Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.