# RIDASCREEN®FAST Ei / Egg Protein

**REF R6402**

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# RIDASCREEN®FAST Ei / Egg Protein

## Brief information

RIDASCREEN®FAST Ei / Egg Protein (Art. No. R6402) is a sandwich enzyme immunoassay for the quantitative analysis of whole egg (-powder) in food validated for the method (see chapter 1. Intended Use).

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization, extraction and centrifugation
Time requirement:	sample preparation (for 10 samples) approx. 20 min test implementation (incubation time).....30 min
Standard material:	NIST Reference Material 8445 (whole egg powder)
Limit of detection: (matrix-dependent)	0.096 mg/kg (ppm) whole egg powder* corresponds to 0.03 mg/kg (ppm) egg white protein 0.052 - 0.157 mg/kg (ppm) Volleipulver *mean value
Limit of quantification:	0.5 mg/kg (ppm) whole egg powder corresponds to 0.13 mg/kg (ppm) egg white protein
Specificity:	The antibodies used specifically detect the egg white proteins ovalbumin and ovomucoid. <b>No</b> cross reactivity was observed to raw or cooked meat of chicken, turkey hen, pork and beef. Further information are contained in the validation report.

Cross reactivities of the antibodies used for this test kit have been determined for the pure food (e.g. corn flour). In composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spiking experiments (see chapter 13. Limits of the method).

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It lists minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

## **Related products for egg determination**

RIDASCREEN®Egg (Art. No. R6411)

RIDASCREEN®FAST Lysozym (Art. No. R6452)

Bioavid Lateral Flow Ei / Egg (Art. No. BL608-10, BL608-25)

Bioavid Lateral Flow Ei / Egg incl. Hook line (Art. No. BLH708-20)

### **1. Intended use**

RIDASCREEN®FAST Ei / Egg Protein (Art. No. R6402) is a sandwich enzyme immunoassay test developed for the quantitative analysis of whole egg (-powder) in food. Due to the large number of different food products, the following samples were examined representative for different food product categories within the scope of the test development: ice cream, noodles, salad dressing and wine. It can be assumed that the assay is also suitable for the analysis of other foods, however, this must be verified by the user before applying the test kit to that food category.

For detailed results and further information on validation data with other food matrices, please refer to the validation report. Further applications are regularly validated in our laboratories, which we make available in our application notes (see chapter 15. Further application notes).

### **2. General**

Egg allergy is very common in children, mostly they lose their egg allergy when getting older. The allergen can be present as an ingredient or as a contamination in raw and cooked products. Egg white contains 9 - 11 % of protein. Four main allergenic proteins come to 80 % of the egg white proteins. These main allergens are ovomucoid (11 %), ovalbumin (54 %), ovotransferrin (12 %) and lysozyme (3.5 %). The allergenic potential of the proteins present in egg yolk is only moderate.

According to the **regulation (EU) No. 1169/2011** from egg and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

### 3. Test principle

The principle of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against egg white proteins. By adding the standard or sample solution to the wells, egg white protein present in the sample will bind to the specific capture antibodies resulting in the formation of an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Following the washing step, a solution containing antibodies conjugated to peroxidase is added. This conjugate binds to the Ab-Ag-complex and an antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugate is then removed in another washing step. A substrate/chromogen solution is added to the wells and incubated. Bound conjugate converts the colorless chromogen into a blue product. A stop solution is added which results in a color change from blue to yellow. The absorbance of the solution which is proportional to the egg white protein concentration in the sample is measured photometrically at 450 nm and expressed as mg/kg whole egg powder.

### 4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
<b>Microtiter plate</b>	-	Ready to use		48 wells
<b>Allergen extraction buffer</b>	Green	<b>Concentrate</b>	<b>10x</b>	100 ml
<b>Standard 1*</b>	Transparent	Ready to use	0 mg/kg	1.3 ml
<b>Standard 2*</b>	Transparent	Ready to use	0.5 mg/kg	1.3 ml
<b>Standard 3*</b>	Transparent	Ready to use	1.5 mg/kg	1.3 ml
<b>Standard 4*</b>	Transparent	Ready to use	4.5 mg/kg	1.3 ml
<b>Standard 5*</b>	Transparent	Ready to use	13.5 mg/kg	1.3 ml
<b>Wash buffer</b>	Brown	<b>Concentrate</b>	<b>10x</b>	100 ml
<b>Conjugate</b>	Red	<b>Concentrate</b>	<b>11x</b>	0.7 ml
<b>Substrate/Chromogen</b> Red Chromogen Pro	Brown	Ready to use		10 ml
<b>Stop solution</b>	Yellow	Ready to use		14 ml

\* The dilution factor 20, which results after sample preparation, has already been considered for the standard concentrations. Therefore, the whole egg powder concentration of samples can directly be read from the standard curve.

## 5. Reagents required but not provided

### 5.1. Equipment

- Gloves
- Scale (measurement range at least up to 50 g and precision of  $\pm 0.01$  g)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge (at least 2,500 x g) + centrifugal vials (e.g. 50 ml centrifuge tubes from Greiner Art. No. 227261)
- Shaker
- Water bath (60 °C / 140 °F; for fluctuation range please refer to the instructions of the water bath manufacturer)
- Fluted filter (pore size 8 - 12  $\mu$ m)
- Graduated pipettes
- Graduated cylinder
- Variable 20 - 200  $\mu$ l and 200 - 1000  $\mu$ l micropipettes
- If necessary: a further microtiter plate (e.g. universal binding, breakable MTP from Thermo Fisher Scientific Art. No. 95029390 or low binding Greiner bio-one Art. No. 655901)
- If necessary: 8-channel pipette for 100  $\mu$ l
- Microtiter plate spectrophotometer (450 nm)
- Optional: RIDASOFT® Win.NET (Art. No. Z9996FF)

### 5.2. Reagents

- Distilled water (dist. water) or deionized water
- If necessary (for samples containing fenugreek, cloves, mustard or celery containing): Casein (e.g. Sigma, Art. No. C5890)

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory personnel. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

Do not reuse wells of the microtiter strips (coated microtiter plate and pre-plate, if necessary (see chapter 10.2. Test procedure). Use separate pipette tips for each standard and each sample extract to avoid cross contamination.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## **7. Storage instructions**

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## **8. Indication of instability or deterioration of reagents**

- Bluish coloration of the reddish substrate/chromogen prior to test implementation
- Value of less than 1.2 absorbance units ( $A_{450\text{ nm}} < 1.2$ ) for standard 5

## **9. Preparation of Samples**

Wear gloves before starting and during the assay. Airborne allergens and dirty laboratory equipment may lead to contamination of the assay. Therefore, please notice the following recommendations:

- Clean surfaces, glass vials, mincers and other equipment before and after each sample preparation.
- Carry out the sample preparation in a room isolated from the ELISA procedure.

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The Allergen extraction buffer is provided as a 10-fold concentrate and must be diluted prior use. Before dilution of the buffer concentrate, dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that, dilute the heated buffer concentrate 1:10 (1+9) with dist. water (i.e. 100 ml buffer concentrate + 900 ml dist. water). The **diluted Allergen extraction buffer (AEB)** is stable at 20 - 25 °C (68 - 77 °F) for approx. 4 weeks or at 2 - 8 °C (35 - 46 °F) for 12 weeks.

Make sure to heat the AEB in time in the 60 °C (140 °F) water bath.

In the following section, the following abbreviation is used:

- AEB: final diluted Allergen extraction buffer

## 9.1. Extraction

**Heat the AEB to 60 °C (140 °F) before sample extraction.**

- Homogenize (grind thoroughly to powder and mix well or mix well a solution respectively) well a sufficient amount (e.g. 50 g or 50 ml) to ensure taking a representative test portion of sample.
- Transfer 1 g (or in case of liquid samples 1 ml) of homogenized sample to a new vial, add 20 ml (or 19 ml in case of liquid samples) pre-heated AEB (see chapter 9.), close the vial and mix well.
- Mix vigorously (e.g. vortexer).
- Extract for 10 min at 60 °C (140 °F) in a water bath.
- Let the sample cool down shortly (e.g. in ice water 3 - 5 min).
- Filter sample or centrifuge for 10 min at > 2,500 x g; if possible at 4 °C (39 °F).  
(Alternatively: transfer 2 ml of the extract into a reaction vial and centrifuge at high speed (> 10,000 x g) for 10 min in a microcentrifuge).
- Transfer the supernatant into a fresh vial.
- If the supernatant is not free of particles after centrifugation, filter the extract additionally.
- The sample extracts (supernatant from centrifugation step or filtrate) can be stored up to one day in a well-closed container at 2 - 8 °C (35 - 46 °F). Unused extracts can be stored at -20 °C (-4 °F) for several months.

## Remark

For fenugreek, cloves, mustard or celery containing samples, add 1 g casein to 1 g or 1 ml sample to avoid disturbing effects. When analyzing samples with unknown composition, casein can be added as a precaution.

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted conjugate has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) in dist. water (e.g. 2 ml dist. water + 200 µl conjugate, sufficient for 2 microtiter strips).

The **wash buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with dist. water (i.e. 900 ml dist. water + 100 ml buffer concentrate). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for 4 weeks.

Components should be stored immediately at 2 - 8 °C (35 - 46 °F) when no longer required.

### 10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time. If more than three strips are needed, a second uncoated plate (see chapter 5.1.) should be used as a pre-plate to avoid a time shift over the microtiter plate. All standards and samples are pipetted into the uncoated plate (at least 150 µl per well) and then exactly 100 µl are quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or sample (prepared according to chapter 9. Sample preparation) in duplicate to the wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see chapter 10.1.) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
4. Add 100 µl of the diluted conjugate (see chapter 10.1.) to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour out the liquid of the wells and tap the microwell holder upside down vigorously (three times) on absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see chapter 10.1.) and pour out the liquid as before. Repeat two more times (a total of three wash cycles).
6. Add 100 µl of substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

## 11. Evaluation

A special software, **RIDASOFT® Win.NET (Art. No. Z9996FF)**, is optional available for evaluation of the RIDASCREEN® enzyme immunoassays. The evaluation should be done using the 4-parameter or cubic spline function. Once an evaluation method has been selected, it should be retained and not be switched between the two functions. Due to its mathematical background, the cubic spline function cannot calculate concentrations outside the measuring range (< standard 2 and > standard 5). The 4-parameter function enables the calculation of values between standard 1 and standard 2 (refer to chapter 13. Limits of the method).

For the evaluation it should be clarified, that the quality criteria are fulfilled for the current test run. The course of the standard curve can be taken from the attached analysis certificate.

The assay can be also evaluated when running in single determinations. This has no influence on the function of the test kit. A special assay evaluation must be written in the RIDASOFT® Win.NET software for this purpose. It is not present by default. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. However, it is not consistent with standards like EN 15633-1 and EN 15842. It should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

When working in accordance with this extraction, the sample dilution factor is 20. Hence, a dilution factor of 20 is already taken into account with the standard concentrations (see chapter 4.\*), the whole egg powder concentration can directly be read from the standard curve.

## 12. Result interpretation

The test is calibrated against NIST Reference Material 8445 (whole egg powder). The result is expressed as mg/kg whole egg powder. The reference material 8445 contains 48 %  $\pm$  1 % total protein. If the result is multiplied by 0.48, then the result can be expressed as mg/kg total protein. If the result is multiplied by 0.263 then mg/kg egg white protein are obtained.

Results between LOD and LOQ indicate a low allergen concentration in the sample. Calculated results show a high uncertainty in this area due to the method's high variation below LOQ. Therefore, such results should not be reported with a quantitative value, but qualitative as "< LOQ".

A result below the LOD does not exclude an allergen contamination below the detection limit of the assay, or that other allergen components, such as lipids, may be present in a sample. The result should be reported accordingly.

A further dilution and new detection of samples is recommended for absorbance values ( $A_{450 \text{ nm}}$ ) > standard 5. In case of a further dilution, the additional dilution factor must be taken into account when calculating the allergen concentration.

Compared to the certificate, higher absorbance values ( $A_{450 \text{ nm}}$ ) for the standard curve, especially for the zero standard, may be a result of insufficient washing or allergen contamination.

### **Calculation example:**

From the standard curve 10 mg/kg whole egg powder are read. If multiplied this result by 0.48 then 4.8 mg/kg whole egg protein is obtained. If multiplied the result by 0.263 then 2.63 mg/kg egg white protein is obtained.

### **13. Limits of the method**

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

Technical limits of the test method are approached outside the designated measurement range resulting in higher variation. This may cause a switch of results between the different areas of the calibration curve especially at the test characteristic boundaries (LOD, LOQ, upper limit of measurement range).

An incorrect weight of the sample to be analyzed will have a 1:1 effect on the measurement result (e.g. a 10 % higher concentration is measured with a weigh in of +10 %). A sufficient accuracy is given with a fluctuation of max  $\pm 1$  %.

For detailed results and further information for other food matrices, please refer to the current validation report. In addition, data on individual foods may be available from comparative laboratory tests and inter-laboratory comparisons.

For the present ELISA, only individual, exemplary foods from different product categories could be validated due to the large number of foods. When analyzing a non-validated matrix, it is recommended to verify the results obtained by means of spike experiments. If necessary, a validation of the sample matrix of interest will need to be performed.

Due to the multitude of food types, matrix effects cannot be excluded. These can lead to false-positive / increased results, but also reduce or suppress a correct reaction. Such matrix effects are independent of the specificity of the antibody used in the test and can be made visible by spiking experiments.

The addition of foreign protein (depending on the test e.g. BSA, gelatine, skim milk powder) during extraction or test procedure may suppress matrix effects.

In processed foods (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented and this may have an impact on the recovery and assay results.

Cross reactivities are side reactions of the antibody used for preparing the test kit with antigen showing similar epitopes as the investigated analyte. These appear especially with antigens from closely related species. In contrast to matrix effects, it is a specific reaction of the antigen with the antibody. The antigen structures are subject to similar influences (e.g. by heating or drying) as the actual analyte. Therefore, cross reactivities may also appear after food processing in single case or are lost.

For evaluation of the cross reactivity only one representative sample was analyzed, other samples may show a different result. All analyzed cross reactivities are described in the validation report.

Egg containing samples that have been heat-treated show a reduced recovery because the proteins denature and are no longer recognized by the antibody used in the ELISA RIDASCREEN®FAST Ei / Egg Protein (Art. No. R6402). The reduction in recovery depends strongly on the temperature and the duration of the heat treatment. If samples are heat treated at high temperatures the recovery can be significantly reduced. For heat-treated samples, we recommend to use our RIDASCREEN® Egg (Art. No. R6411).

## 14. Recommendation

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in standards like EN 15633-1 and EN 15842 (e.g. performing duplicate determinations).
- Pre-flush pipette tips with standard or sample extract prior to pipetting.
- Carry along test controls for quality control and to ensure an accurate and correct test procedure. Egg-free and egg-containing (naturally contaminated or spiked) samples should be used. An example of a spike experiment is given in the validation report.
- In case of extremely acidic or basic samples, adjust the sample's pH value (pH 6.5 - 7.5) to neutral prior to extraction may be necessary.
- To contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de) if automates (e.g. ThunderBolt®) are used.

## 15. Further application notes

- More sensitive allergen analysis for wine (dilution 1:10).
- RIDASCREEN®FAST Ei / Egg Protein (Art. No. R6402) used with automates GEMINI and ChemWell® 2910 Automate.
- RIDASCREEN®FAST Allergen - Swabbing method for the qualitative analysis of allergens in the production line or for laboratory equipment.

**Further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Version overview

Version number	Chapter and title
2015-12-14	Release version
2021-06-16	<ul style="list-style-type: none"><li>– General linguistic revision</li><li>– Chapter 1, 6, 9, 10, 11, 12, 13, 14</li></ul>
2022-05-06	Current version Changes made: <ul style="list-style-type: none"><li>– Chapter 10.2. test procedure: Mistake from version 2021-06-16 with "total of four wash cycles" corrected back to "total of three wash cycles"</li></ul>

## Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.