CONGEN

SureFood[®] GMO QUANT Bt11 Corn

Art. No. S2016
2 x 50 rxn

User Manual

January 2021

 **Inhalt**

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Content

1	General Information	8
1.1	Description	8
1.2	Limit of Detection and Limit of Quantification	8
1.3	DNA-preparation	8
1.4	Kit components and storage	8
1.5	Additionally required equipment and materials	8
1.6	Setup	9
1.7	Detection channel Set-up	9
2	Qualitative Analysis	10
2.1	Protocol	10
2.1.1	Preparation of the master-mix	10
2.1.2	Preparation of the standard DNA dilutions	10
2.1.3	Preparation of the real-time PCR-mix	11
2.2	Interpretation of results	11
3	Further Information	12
3.1	Product Information	12
3.2	Technical Support	12
3.3	Distribution and Ordering	12

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1 General Information

1.1 Description

The test detects the relative Bt11 corn DNA amount to the total corn DNA content. Therefore the kit contains two PCR systems, one for detection of a Bt11 corn specific gene (OECD unique identifier SYN-BTØ11-1) and one specific for corn (reference gene).

The specific Bt11 corn detection is according to the official and validated method of the European Commission.

The real-time PCR assay can be used with established real-time PCR instruments. The technical validation of instruments was performed on Applied Biosystems 7500, BioRad CFX96, Roche LightCycler® LC1.5, Roche LightCycler® LC2.0 and Eppendorf realplex.

1.2 Limit of Detection and Limit of Quantification

The Bt11 corn specific PCR has a limit of detection of ≤ 5 DNA-copies.

The assay limit of detection depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The limit of quantitation depends on the concentration of the sample DNA used in the analysis. For example, if 50,000 target-sequence copies of the reference gene are present, the relative quantitation limit for Bt11 Canola DNA is 0.1 %.

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Ref Reaction Mix	1 x 1100 µl	Orange
2	Bt11 Canola Reaction Mix	1 x 1100 µl	Yellow
3	Taq Polymerase	1 x 11 µl	Red
4	Dilution Buffer	1 x 1300 µl	White
5	Standard DNA	1 x 50 µl	Dark Blue
6	Positive Control (1% Bt11 corn)	1 x 100 µl	Light Blue

Store all reagents at -20°C and protected from light.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053 / SureFood® PREP Add-On Art. No. S1055)
- real-time PCR instrument
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

SureFood® GMO QUANT Bt11 Corn (2 x 50 rxn)

Art. No. S2016

January 2021

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Applied Biosystems 7500	Bt11 corn	FAM	None	Check the passive reference option ROX is none.
	corn	FAM	None	
Bio-Rad CFX96	Bt11 corn	FAM	+	
	corn	FAM	+	
Roche LightCycler® 2.0	Bt11 corn	530/610	+	Check the Seek Temperature is 58°C.
	corn	530/610	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples, control reactions and standards) for the specific PCR assay.

Reactions needed for the corn detection:

5 reactions for the standard curve

3 reactions for controls (1x no-template control, 2x Positive Control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the Bt11 corn detection:

5 reactions for the standard curve

3 reactions for controls (1x no-template control, 2x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different DNA concentrations for the standard curves of the reference gene (**corn**) and the detection gene (**Bt11 corn**). Prepare 5 dilutions of the supplied Standard DNA (**Code 5**) with the supplied Dilution Buffer (**Code 4**). Prepare 5 reaction tubes (labeled S1 to S5) and add 45 µl Dilution Buffer (**Code 4**) each.

The following procedure is recommended:

Standard	Dilutions	Copy number per µl	Final copy number per reaction*
S1	45 µl Dilution Buffer + 5 µl Standard DNA	100,000 copies	500,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	50,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1,000 copies	5,000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	500 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S4	10 copies	50 copies

***Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

2.1.3 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control and the standard dilutions into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The calculation for both reaction systems (**corn**, **Bt11 corn**) has to be made separately. Mark the standards, the controls and the samples for the specific system (**Bt11 corn**) and make the evaluation according to the usual analysis program recommended by the real-time PCR instrument manufacturer. Repeat the same procedure for the reference gene system (**corn**). The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

By using the calculated copy numbers for **corn** and **Bt11 corn** the relative **Bt11 corn** content of the sample DNA and the Positive Control can be determined in the following way (example):

Sample Bt11 corn	1,350 copies	Positive Control Bt11 corn	400 copies
Sample corn	45,000 copies	Positive Control corn	28,000 copies

Divide the copy number of the specific system by the copy number of the reference gene system and multiply by 100 to obtain the percentage.

Bt11 corn DNA content = **Bt11 corn** copy number * 100 / **corn** copy number

sample **Bt11 corn** DNA content = 1,350 * 100 / 45,000 sample **Bt11 corn** DNA content = 3 %

For the given example the numbers lead to a **Bt11 corn** DNA content of 3.0 %, for the sample and 1.4 % for the Positive Control with the same calculation.

For a final calculation the use of a correction factor K for the correction of run-to-run fluctuations is necessary. The correction factor is the relation of the true percentage value of the Positive Control (1 % **Bt11 corn** content) and the measured percentage of the Positive Control.

$K = \text{true value} / \text{measured value}$

$K \text{ (example)} = 1 \% / 1.4 \% = 0.7$

The measured DNA content for the sample is multiplied with K to obtain a corrected DNA content.

sample DNA content = measured sample DNA content * K sample (example) = 3.0 % * 0.7 = 2.1 %

For that example the **Bt11 corn content is 2.1 %**.

3 Further Information

3.1 Product Information

- Microsoft Excel template of calculation and detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/company/downloads)
- Validation Report upon request

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com

r-biopharm®