CONGEN

SureFood[®] GMO Plant Reference Sample

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Content

| | | |
|-----|---|---|
| 1 | General Information | 5 |
| 1.1 | Description | 5 |
| 1.2 | DNA-preparation | 6 |
| 1.3 | Kit components and storage | 6 |
| 1.4 | Additionally required equipment and materials | 6 |
| 2 | Further Information | 6 |
| 2.1 | Technical Support | 6 |
| 2.2 | Distribution and Ordering | 6 |

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| Reagenz | Menge |
|----------------------------|-----------|
| GMO Plant Reference Sample | 1 x 2,0 g |

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1 General Information

1.1 Description

The SureFood® GMO Plant Reference Sample is validated for semiquantitative determination of the plants soya, corn, canola und rice in food and feed samples.

The laboratory reference material SureFood® GMO Plant Reference Sample is based on **wheat flour** material. It contains the plants **soya, corn, canola** and **rice** in a concentration of **0.1 %**.

Content of plant substances:

| Plant | Material |
|--------|--|
| wheat | wheat flour (<i>Triticum spp.</i>) |
| soya | soya flour (<i>Glycine max</i>) |
| corn | corn flour (<i>Zea mays</i>) |
| canola | canola flour (<i>Brassica napus</i>) |
| rice | rice flour (<i>Oryza sativa</i>) |

The laboratory reference material SureFood® GMO Plant Reference Sample is validated in combination with the DNA extraction method SureFood® PREP Basic and the SureFood® real-time PCR detection methods.

In case of using alternative methods the sample weight and extraction method in combination with the detection system have to be validated.

1.2 DNA-preparation

SureFood® PREP Basic (Art. No. S1052) or additionally SureFood® PREP Add-On Kit (Art. No. S1055) is recommended for the DNA- preparation.

1.3 Kit components and storage

| Reagent | Amount |
|----------------------------|-----------|
| GMO Plant Reference Sample | 1 x 2.0 g |

Store all reagents dry at –20°C and protected from light.

Note: To avoid moisture absorption the unopened laboratory reference material should thaw at room temperature for at least one hour before use.

1.4 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFood® PREP Basic, Art. No. S1052; optional SureFood® PREP Add-On Kit, Art. No. S1055)
- Real-time PCR methods (SureFood® GMO series)

2 Further Information

2.1 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

2.2 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com