CONGEN

SureFood[®] GMO ID 4plex Soya I

Art. No. S2161
100 rxn

User Manual

June 2021

 **Inhalt**

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	DNA-preparation	10
1.4	Kit components and storage	10
1.5	Additionally required equipment and materials	10
1.6	Setup	10
1.7	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	12
2.2	Interpretation of results	13
3	Further Information	14
3.1	Product Information	14
3.2	Technical Support	14
3.3	Distribution and Ordering	14

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1 General Information

1.1 Description

The SureFood® GMO ID 4plex Soya I is a real-time PCR for the direct, qualitative detection and differentiation of following specific genetically modified soya DNA sequences.

- FAM-channel: MON87708 Soya (OECD unique identifier MON-87708-9)
- FAM-channel: CV127 Soya (OECD unique identifier BPS-CV127-9)
- VIC-channel: DP305423 Soya (OECD unique identifier DP-305423-1)
- ROX-channel: MON87701 Soya (OECD unique identifier MON87701-2)
- Cy5-channel: MON87769 Soya (OECD unique identifier MON-87769-7)

The detections are according to the validated method of the European Commission.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The technical validation of instruments was performed on Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Applied Biosystems 7500, Bio-Rad CFX96 and Agilent Mx3005P.

1.2 Limit of Detection

The SureFood® GMO ID 4plex Soya I real-time PCR has a limit of detection of ≤ 5 DNA copies. This is equivalent to approx. 0.01 % for unprocessed soybeans.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFood® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

Note: In mixed samples, inconsistent mixing ratios (e.g. 99.9 % MON87708 Soya and 0.1 % MON87769 Soya) may cause a loss of sensitivity in the low concentration channel.

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic (Art. No. S1052), SureFast® Mag PREP Food (Art. No. F1060) and for highly processed food and feed the use of SureFood® PREP Advanced (Art. No. S1053) is recommended. SureFood® PREP Add On (Art. No. S1055) is intended to be used for the extraction of DNA from raw materials as well as processed food and feed with sample weight of 2 g. It is used in conjunction with the SureFood® PREP Basic.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053 / SureFood® PREP Add On Art. No. S1055 / SureFast® Mag PREP Food Art. No. F1060)
- real- time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD) Cycles	5 min, 95°C 45	1 min, 95°C 45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	MON87708 Soya, CV127 Soya	FAM	+	
	DP305423 Soya	HEX	+	
	MON87701 Soya	ROX	+	
	MON87769 Soya	Cy5	+	
Agilent AriaMx / Dx	MON87708 Soya, CV127 Soya	FAM	+	
	DP305423 Soya	HEX	+	
	MON87701 Soya	ROX	+	
	MON87769 Soya	Cy5	+	
Applied Biosystems 7500	MON87708 Soya, CV127 Soya	FAM	None	Check the passive reference option ROX is none.
	DP305423 Soya	VIC	None	
	MON87701 Soya	ROX	None	
	MON87769 Soya	Cy5	None	
Bio-Rad CFX96	MON87708 Soya, CV127 Soya	FAM	+	
	DP305423 Soya	VIC/HEX	+	
	MON87701 Soya	ROX	+	
	MON87769 Soya	Cy5	+	
Roche LightCycler® 480 II	MON87708 Soya, CV127 Soya	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	DP305423 Soya	533-580	+	
	MON87701 Soya	533-610	+	
	MON87769 Soya	618-660	+	
Roche cobas® z 480 Analyzer	MON87708 Soya, CV127 Soya	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	DP305423 Soya	540-580	+	
	MON87701 Soya	540-610	+	
	MON87769 Soya	610-670	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control and an inhibition control per sample.

For the preparation of the inhibition control, it is recommended to use the SureFood® GMO Plant PLUS (Art. No. S2049) or the SureFood® GMO SCREEN 4plex 35S/NOS/FMV+IAC (Art. No. S2126).

Reactions needed for the qualitative MON87708 Soya, CV127 Soya, DP305423 Soya, MON87701 Soya and MON87769 Soya detection:

3 reactions for controls (1x negative control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of positive control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

MON87708 Soya, CV127 Soya DNA is detected in the FAM-channel, DP305423 Soya DNA is detected in the VIC/HEX -channel, MON87701 Soya DNA is detected in the ROX-channel and MON87769 Soya DNA is detected in the Cy5-channel (see table).

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC)

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the external inhibition control of the sample is **positive** with a shift in Cp-value ≤ 2 compared to the negative control.

If the sample DNA in the external inhibition control shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific MON87708 Soya, CV127 Soya, MON87701 Soya, DP305423 Soya or MON87769 Soya PCR assay.

Result in the respective channel				Result	Interpretation
FAM channel MON87708 Soya, CV127 Soya	VIC/HEX channel DP305423 Soya	ROX channel MON87701 Soya	Cy5 channel MON87769 Soya	external inhibition- control	
positive	negative	negative	negative	positive	MON87708 Soya, CV127 Soya DNA detected
negativ	positive	negative	negative	positive	DP305423 Soya DNA detected
negativ	negative	positive	negative	positive	MON87701 Soya DNA detected
negativ	negative	negative	positive	positive	MON87769 Soya DNA detected
negativ	negative	negative	negative	positive	Negative for target DNA
negativ	negative	negative	negative	negative	invalid

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices
(Download: www.congen.de/en/company/downloads)
- Validation Report upon request

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com