CONGEN

SureFood[®] GMO ID 4plex Soya II

Art. No. S2162
100 rxn

User Manual

April 2021

 **Inhalt**

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	DNA-preparation	9
1.4	Kit components and storage	10
1.5	Additionally required equipment and materials	10
1.6	Setup.....	10
1.7	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	12
2.2	Interpretation of results	13
3	Further Information	14
3.1	Product Information.....	14
3.2	Technical Support	14
3.3	Distribution and Ordering	14

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1 General Information

1.1 Description

SureFood® GMO ID 4plex Soya II is a multiplex real-time PCR for the direct, qualitative detection and differentiation of specific genetically modified soya DNA sequences listed as below:

- Roundup Ready Soya (RR Soya) (OECD unique identifier MON-Ø4Ø32-6)
- Roundup Ready 2 Yield Soya (RR2Y Soya) (OECD unique identifier MON-89788-1)
- A2704-12 Soya (OECD unique identifier ACS-GMØØ5-3)
- A5547-127 Soya (OECD unique identifier ACS-GMØØ6-4)

The detections and differentiation are performed according to the validated method of the European Commission respectively official collection of detection methods of §64 German food law.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The technical validation of instruments was performed on Agilent Mx3005P, Applied Biosystems 7500, Bio-Rad CFX96, R-Biopharm RIDA®CYCLER, Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer and Qiagen Rotor-Gene Q.

1.2 Limit of Detection

The SureFood® GMO ID 4plex Soya II real-time PCR has a limit of detection of ≤ 5 DNA copies. This is equivalent to approx. 0.01 % of unprocessed soybean.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

Note: In mixed samples, inconsistent mixing ratios (e.g. 99.9 % RR Soya and 0.1 % A5547-127 Soya) may cause a loss of sensitivity in the low concentration channel.

1.3 DNA-preparation

For DNA-preparation of raw material the use of SureFood® PREP Basic and for highly processed food and feed the use of SureFood® PREP Advanced is recommended. SureFood® PREP Add On is intended to be used for the extraction of DNA from raw materials as well as processed food and feed with sample weight of 2 g. It is used in conjunction with the SureFood® PREP Basic.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1100 µl	Yellow
2	Taq Polymerase	1 x 11 µl	Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at –20°C and protected from light.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053 / SureFood® PREP Add On Art. No. S1055)
- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD)	5 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	RR-Soya	FAM	+	
	RR2Y-Soya	HEX	+	
	A2704-12-Soya	ROX	+	
	A5547-127-Soya	Cy5	+	
Applied Biosystems 7500	RR-Soya	FAM	None	Check the passive reference option ROX is none.
	RR2Y-Soya	VIC	None	
	A2704-12-Soya	ROX	None	
	A5547-127-Soya	Cy5	None	
Bio-Rad CFX96	RR-Soya	FAM	+	
	RR2Y-Soya	VIC	+	
	A2704-12-Soya	ROX	+	
	A5547-127-Soya	Cy5	+	
R-Biopharm RIDA®CYCLER	RR-Soya	green	+	
	RR2Y-Soya	yellow	+	
	A2704-12-Soya	orange	+	
	A5547-127-Soya	red	+	
Qiagen Rotor-Gene Q	RR-Soya	green	+	Note: Please use only 0.1 ml reaction tube. The gain settings must be set to 5 (factory default) for all channels.
	RR2Y-Soya	yellow	+	
	A2704-12-Soya	orange	+	
	A5547-127-Soya	red	+	
Roche LightCycler® 480 II	RR-Soya	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	RR2Y-Soya	533-580	+	
	A2704-12-Soya	533-610	+	
	A5547-127-Soya	618-660	+	
Roche cobas® z 480 Analyzer	RR-Soya	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	RR2Y-Soya	540-580	+	
	A2704-12-Soya	540-610	+	
	A5547-127-Soya	610-670	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, extraction control, positive control and an inhibition control per sample.

For the preparation of the inhibition control the use of the SureFood® GMO Plant Plus (Art. Nr. S2049), SureFood® GMO SCREEN 4plex 35S/NOS/FMV+IAC (Art. No. S2126) respectively the is recommended.

Reactions needed for the qualitative RR-Soya, RR2Y-Soya, A2704-12-Soya and A5547-127-Soya detection:

3 reactions for controls (1x no-template control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

RR-Soya DNA is detected in the FAM-channel, RR2Y-Soya DNA is detected in the VIC/HEX -channel, A2704-12-Soya DNA is detected in the ROX-channel and A5547-127-Soya DNA is detected in the Cy5-channel (see table).

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the external inhibition control of the sample is **positive** with a shift in Cp-value ≤ 2 compared to the negative control.

If the sample DNA in the external inhibition control shows **no amplification** or a shift in Cp-value > 2 compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances DNA isolation and purification of the sample need to be improved. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific RR-Soya, RR2Y-Soya, A2704-12-Soya or A5547-127-Soya PCR assay.

Result in the respective channel				Result external inhibition-control	Interpretation
FAM channel RR-Soya	VIC/HEX channel RR2Y-Soya	ROX channel A2704-12-Soya	Cy5 channel A5547-127-Soya		
positive	negative	negative	negative	positive	RR-Soya DNA detected
negativ	positive	negative	negative	positive	RR2Y-Soya DNA detected
negativ	negative	positive	negative	positive	A2704-12-Soya DNA detected
negativ	negative	negative	positive	positive	A5547-127-Soya DNA detected
negativ	negative	negative	negative	positive	Negative, target DNA is not detected
negativ	negative	negative	negative	negative	invalid

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices
(Download: www.congen.de/en/company/downloads)
- Validation Report upon request

3.2 Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

3.3 Distribution and Ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com